



RESOURCE AND PATIENT MANAGEMENT SYSTEM

Clinical Reporting System (CRS) (BGP)

Technical Manual

Version 6.0
November 2005

Office of Information Technology
Albuquerque, New Mexico

PREFACE

This manual contains the technical manual for the CRS Clinical Reporting System version 6.0, which adds FY 2006 clinical performance measures to existing FY 2002 through FY 2005 measures.

The CRS Clinical Reporting System is an RPMS (Resource and Patient Management System) software application designed for national reporting as well as local and Area monitoring of clinical GPRA and developmental measures. CRS was first released for FY 2002 performance measures (as GPRA+) and is based on a design by the Aberdeen Area (GPRA2000).

The Government Performance and Results Act (GPRA) requires Federal agencies to report annually on how the agency measured up against the performance targets set in its annual Plan. IHS GPRA measures include measures for clinical prevention and treatment, quality of care, infrastructure, and administrative efficiency functions. The CRS Clinical Reporting System is the reporting tool used by the IHS Office of Planning and Evaluation to collect and report clinical performance results annually to the Department of Health and Human Services (DHHS) and to Congress.

Each year, an updated version of CRS software is released to reflect changes in the logic descriptions of the different denominators and numerators. Additional performance measures may also be added. Local facilities can run reports as often as they want to and can also use CRS to transmit data to their Area. The Area Office can use CRS to produce an aggregated Area report for either annual GPRA or Area Director Performance reports.

The CRS Clinical Reporting System will produce reports on demand from local RPMS databases for both GPRA and developmental clinical performance measures that are based on RPMS data. CRS is intended to eliminate the need for manual chart audits for evaluating and reporting clinical measures. Administrative and clinical users will be able to review individual or all measures at any time, and can:

- identify potential data issues in their RPMS, i.e., missing or incorrect data;
- monitor their site's performance against past national performance and upcoming agency goals;
- identify specific areas where the facility is not meeting the measure in order to initiate business process or other changes;
- quickly measure impact of process changes on performance measures;
- identify areas meeting or exceeding measures to provide lessons learned.

To produce reports with comparable data across every facility, the GPRA measure definition was “translated” into programming code with the assistance of clinical subject matter experts. CRS uses pre-defined taxonomies to find data items in PCC to determine if a patient meets the performance measure criteria. Taxonomies contain groups of codes

(e.g., diagnoses or procedures) or site-specific terms. Each performance measure has one or more denominators and numerators defined.

CRS is intended for use by Area and site Quality Improvement staff, Compliance Officers, GPRA Coordinators, clinical staff such as physicians, nurses, nurse practitioners, and other providers, Area Directors, as well as any staff involved with quality assurance initiatives.

TABLE OF CONTENTS

1.0	IMPLEMENTATION AND MAINTENANCE	1
1.1	System Requirements.....	1
2.0	SECURITY KEYS	2
3.0	ROUTINES.....	3
4.0	FILES AND TABLES	21
4.1	File List	21
4.2	File Access	22
4.3	Cross References	23
4.4	Table File	38
5.0	ARCHIVING AND PURGING.....	65
6.0	EXTERNAL RELATIONS	66
6.1	External Calls.....	66
6.2	Callable Routines.....	103
6.3	Exported Options	103
7.0	INTERNAL RELATIONS	108
8.0	HOW TO GENERATE ONLINE DOCUMENTATION	109
8.1	%INDEX.....	109
8.2	Inquire Option	110
8.3	Print Option File	110
8.4	List File Attributes	110
9.0	SAC EXEMPTIONS	112
10.0	GLOSSARY	113
11.0	CONTACT INFORMATION	116

1.0 Implementation and Maintenance

The CLINICAL REPORTING SYSTEM occupies the BGP namespace. Options, security locks/keys, templates, routines, and globals are namespaced BGP.

1.1 System Requirements

- Kernel v8.0 or higher
- FileMan v21 or higher
- IHS Patient Dictionaries (AUPN) v99.1 through patch 15
- PCC Management Reports v3.0 through patch 16
- Taxonomy System v5.1 through patch 8
- XB/ZIB Utilities v3.0 through patch 10
- AUM table update patch 9 (aum*5.1*9)

2.0 Security Keys

Name	Descriptive Name	Description
BGPZAREA	BGP AREA REPORT	This key unlocks the area report options on the main BGPMENU. This key should only be assigned at the Area office level.
BGPZMENU	BGP MAIN MENU - CRS	This key unlocks the main CRS menu. It should be assigned to all users who need to run GPRA/CRS reports.
BGPZ PATIENT LISTS	BGP PATIENT LISTS	This key allows the holder to obtain patient lists. The patient lists contain sensitive patient/clinical information. This key should only be given to those users who have a need for such information.
BGPZ SITE PARAMETERS	BGP SITE PARAMETERS	This key should be given to the those with authority to update the site parameters.
BGPZ TAXONOMY EDIT	BGP TAXONOMY EDIT	This key should be given only to users who will need to update taxonomies.

3.0 Routines

ROUTINE	ROUTINE DESCRIPTION
BGP3BAN	IHS/CMI/LAB - BANNER FOR GPRA 03
BGP3D1	IHS/CMI/LAB - calc indicators
BGP3D10	IHS/CMI/LAB - indicator J-4,J-5,J-6
BGP3D2	IHS/CMI/LAB - indicator 2
BGP3D3	IHS/CMI/LAB - indicator 3
BGP3D4	IHS/CMI/LAB - indicator 3
BGP3D5	IHS/CMI/LAB - indicator 31
BGP3D6	IHS/CMI/LAB - indicator 31
BGP3D7	IHS/CMI/LAB - indicator C
BGP3D8	IHS/CMI/LAB - indicator C
BGP3D9	IHS/CMI/LAB - indicator J
BGP3DA	IHS/CMI/LAB - IHS GPRA - report for area export
BGP3DAR	IHS/CMI/LAB - ihs area GPRA
BGP3DARA	IHS/CMI/LAB - ihs area AA
BGP3DARH	IHS/CMI/LAB - ihs area AA
BGP3DCI	IHS/CMI/LAB - IHS area GPRA
BGP3DG	IHS/CMI/LAB - IHS GPRA - report for area export
BGP3DH	IHS/CMI/LAB - cover page for gpra
BGP3DL	IHS/CMI/LAB - IHS GPRA - report for local use
BGP3DP	IHS/CMI/LAB - IHS gpra print
BGP3DP1	IHS/CMI/LAB - print ind 1
BGP3DP2	IHS/CMI/LAB - print ind 1
BGP3DP3	IHS/CMI/LAB - print ind 1
BGP3DP4	IHS/CMI/LAB - print ind 31
BGP3DP41	IHS/CMI/LAB - print ind 31
BGP3DP5	IHS/CMI/LAB - print ind H
BGP3DP50	IHS/CMI/LAB - print ind H
BGP3DP51	IHS/CMI/LAB - print ind H
BGP3DP6	IHS/CMI/LAB - print ind C-1
BGP3DP60	IHS/CMI/LAB - print ind C-1
BGP3DP7	IHS/CMI/LAB - print ind
BGP3DPH	IHS/CMI/LAB - AREA REPORT HEADER
BGP3DS	IHS/CMI/LAB - IHS gpra print
BGP3DSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP3DSL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP3PARP	IHS/CMI/LAB - IHS gpra print
BGP3PDL	IHS/CMI/LAB - IHS gpra print
BGP3PDL1	IHS/CMI/LAB - print ind 1
BGP3PDL2	IHS/CMI/LAB - print ind 1
BGP3PDL3	IHS/CMI/LAB - print ind 1
BGP3PDL4	IHS/CMI/LAB - print ind 1
BGP3PDL5	IHS/CMI/LAB - print ind 1
BGP3PDL6	IHS/CMI/LAB - print ind 1
BGP3PDL7	IHS/CMI/LAB - print ind 1
BGP3PDL8	IHS/CMI/LAB - print ind 1
BGP3PDL9	IHS/CMI/LAB - print ind 1
BGP3PDLS	IHS/CMI/LAB - IHS gpra print

BGP3PDM1	IHS/CMI/LAB - print ind 1
BGP3PDM2	IHS/CMI/LAB - print ind 1
BGP3TXCH	IHS/CMI/LAB - TAXONOMY CHECK FOR FY03 GPRA REPORT
BGP3ULF	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP3UTL	IHS/CMI/LAB -
BGP3UTL1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP4BAN	IHS/CMI/LAB - BANNER FOR GPRA 04
BGP4D1	IHS/CMI/LAB - calc indicators
BGP4D10	IHS/CMI/LAB - calc indicators
BGP4D2	IHS/CMI/LAB - indicator 1,2,3,4
BGP4D21	IHS/CMI/LAB - indicator 6
BGP4D3	IHS/CMI/LAB - indicator 11
BGP4D31	IHS/CMI/LAB - indicator C
BGP4D4	IHS/CMI/LAB - indicator 3
BGP4D41	IHS/CMI/LAB - indicator 3
BGP4D5	IHS/CMI/LAB - indicator calc
BGP4D6	IHS/CMI/LAB - indicator 31
BGP4D61	IHS/CMI/LAB - indicator 31
BGP4D7	IHS/CMI/LAB - indicator 31
BGP4D71	IHS/CMI/LAB - indicator C
BGP4D8	IHS/CMI/LAB - indicator C
BGP4D9	IHS/CMI/LAB - indicator J
BGP4DA	IHS/CMI/LAB - IHS AREA CLIN 04 REPORT DRIVER
BGP4DAP	IHS/CMI/LAB - IHS GPRA 04 SELECTED REPORT DRIVER ALL PATIENTS
BGP4DAR	IHS/CMI/LAB - ihs area GPRA
BGP4DARA	IHS/CMI/LAB - ihs area AA
BGP4DCHI	IHS/CMI/LAB - calculate HEDIS indicators
BGP4DCI	IHS/CMI/LAB - IHS area GPRA
BGP4DG	IHS/CMI/LAB - IHS GPRA 04 REPORT DRIVER
BGP4DH	IHS/CMI/LAB - cover page for gpra
BGP4DHE	IHS/CMI/LAB - IHS GPRA 04 REPORT DRIVER
BGP4DHE1	IHS/CMI/LAB - indicator HEDIS
BGP4DHEA	IHS/CMI/LAB - ihs area AA
BGP4DHEP	IHS/CMI/LAB - IHS HEDIS print
BGP4DHSL	IHS/CMI/LAB - FY 04 DISPLAY IND LISTS
BGP4DL	IHS/CMI/LAB - IHS GPRA 04 SELECTED REPORT DRIVER
BGP4DP	IHS/CMI/LAB - IHS gpra print
BGP4DP1	IHS/CMI/LAB - print ind 1
BGP4DP11	IHS/CMI/LAB - print ind 1
BGP4DP2	IHS/CMI/LAB - print ind 10
BGP4DP3	IHS/CMI/LAB - print ind 1
BGP4DP4	IHS/CMI/LAB - print ind 31
BGP4DP5	IHS/CMI/LAB - print ind H
BGP4DP50	IHS/CMI/LAB - print ind H
BGP4DP51	IHS/CMI/LAB - print ind H
BGP4DP6	IHS/CMI/LAB - print ind H
BGP4DP60	IHS/CMI/LAB - print ind H
BGP4DPH	IHS/CMI/LAB - AREA REPORT HEADER
BGP4DPP	IHS/CMI/LAB - IHS GPRA 04 SELECTED REPORT DRIVER ALL PATIENTS
BGP4DS	IHS/CMI/LAB - IHS gpra print

BGP4DSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP4DSL	IHS/CMI/LAB - FY 04 DISPLAY IND LISTS
BGP4DSP	IHS/CMI/LAB - IHS summary page
BGP4HEH	IHS/CMI/LAB - cover page for HEDIS
BGP4HEHH	IHS/CMI/LAB - cover page for HEDIS
BGP4HEL	IHS/CMI/LAB - IHS gpra print
BGP4HEL1	IHS/CMI/LAB - print ind 1
BGP4HEP1	IHS/CMI/LAB - print ind 1
BGP4HES	IHS/CMI/LAB - IHS gpra print
BGP4HESL	IHS/CMI/LAB - IHS gpra print
BGP4HUTL	IHS/CMI/LAB -
BGP4LSTF	IHS/CMI/LAB - List 'BG04' files in pub
BGP4PARP	IHS/CMI/LAB - IHS gpra print
BGP4PDL	IHS/CMI/LAB - IHS gpra print
BGP4PDL0	IHS/CMI/LAB - print ind 1
BGP4PDL1	IHS/CMI/LAB - print ind 1
BGP4PDL2	IHS/CMI/LAB - print ind 1
BGP4PDL3	IHS/CMI/LAB - print ind 1
BGP4PDL4	IHS/CMI/LAB - print ind 1
BGP4PDL5	IHS/CMI/LAB - print ind 1
BGP4PDL7	IHS/CMI/LAB - print ind 1
BGP4PDL8	IHS/CMI/LAB - print ind 1
BGP4PDL9	IHS/CMI/LAB - print ind
BGP4PDLH	IHS/CMI/LAB - cover page for gpra del
BGP4PDLS	IHS/CMI/LAB - IHS gpra print
BGP4PDLZ	IHS/CMI/LAB -
BGP4PHEP	IHS/CMI/LAB - IHS gpra print
BGP4TXCH	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 GPRA REPORT
BGP4TXHE	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 HEDIS REPORT
BGP4ULF	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP4UTL	IHS/CMI/LAB -
BGP4UTL1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP50P1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED 20 Dec 2004 9:23 AM
BGP5ASL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5BAN	IHS/CMI/LAB - BANNER FOR CRS 05 06 Jan 2005 9:02 AM
BGP5C1	IHS/CMI/LAB - calc CMS indicators 26 Sep 2004 11:28 AM 04 May 2005 2:38 PM
BGP5C11	IHS/CMI/LAB - calc CMS indicators
BGP5C12	IHS/CMI/LAB - calc CMS indicators 26 Sep 2004 11:28 AM
BGP5C13	IHS/CMI/LAB - calc CMS indicators 26 Sep 2004 11:28 AM
BGP5C2	IHS/CMI/LAB - calc CMS indicators 26 Sep 2004 11:28 AM
BGP5C3	IHS/CMI/LAB - calc CMS indicators 26 Sep 2004 11:28 AM
BGP5CL	IHS/CMI/LAB - IHS GPRA 05 CMS REPORT DRIVER
BGP5CP	IHS/CMI/LAB - IHS gpra print
BGP5CSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5CTL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5CTS	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5CTXC	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT 18 Feb 2005 5:31 PM
BGP5D1	IHS/CMI/LAB - calc indicators 29 Apr 2005 7:38 PM
BGP5D10	IHS/CMI/LAB - calc indicators
BGP5D2	IHS/CMI/LAB - indicator 1,2,3,4

BGP5D21	IHS/CMI/LAB - indicator 6
BGP5D22	IHS/CMI/LAB - indicator I2
BGP5D3	IHS/CMI/LAB - indicator 11
BGP5D31	IHS/CMI/LAB - indicator C
BGP5D32	IHS/CMI/LAB - indicator C
BGP5D33	IHS/CMI/LAB - indicator C
BGP5D34	IHS/CMI/LAB - indicator C
BGP5D35	IHS/CMI/LAB - indicator C
BGP5D4	IHS/CMI/LAB - indicator 3
BGP5D41	IHS/CMI/LAB - indicator 3
BGP5D5	IHS/CMI/LAB - indicator calc
BGP5D6	IHS/CMI/LAB - indicator 31
BGP5D61	IHS/CMI/LAB - indicator 31
BGP5D7	IHS/CMI/LAB - indicator 31
BGP5D71	IHS/CMI/LAB - indicator C
BGP5D72	IHS/CMI/LAB - indicator 31
BGP5D73	IHS/CMI/LAB - indicator 31
BGP5D8	IHS/CMI/LAB - indicator C
BGP5D9	IHS/CMI/LAB - indicator J
BGP5DAP	IHS/CMI/LAB - IHS GPRA 04 SELECTED REPORT DRIVER ALL PATIENTS
BGP5DAR	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP5DCEI	IHS/CMI/LAB - calculate HEDIS indicators
BGP5DCHI	IHS/CMI/LAB - calculate HEDIS indicators
BGP5DCI	IHS/CMI/LAB - IHS area GPRA
BGP5DEL	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER
BGP5DELA	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER
BGP5DELP	IHS/CMI/LAB - IHS HEDIS print
BGP5DESI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5DESL	IHS/CMI/LAB - FY 05 DISPLAY IND LISTS
BGP5DGPA	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP5DGPU	IHS/CMI/LAB - IHS AREA CLIN 05 REPORT DRIVER
BGP5DH	IHS/CMI/LAB - cover page for gpra
BGP5DHE	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER
BGP5DHE1	IHS/CMI/LAB - indicator HEDIS
BGP5DHEA	IHS/CMI/LAB - ihs area AA
BGP5DHEP	IHS/CMI/LAB - IHS HEDIS print
BGP5DHSL	IHS/CMI/LAB - FY 05 DISPLAY IND LISTS
BGP5DL	IHS/CMI/LAB - IHS GPRA 05 SELECTED REPORT DRIVER
BGP5DNG	IHS/CMI/LAB - IHS AREA CLIN 05 REPORT DRIVER
BGP5DP	IHS/CMI/LAB - IHS gpra print
BGP5DP1	IHS/CMI/LAB - print ind 1 21 Mar 2005 12:55 PM
BGP5DP11	IHS/CMI/LAB - print ind 1
BGP5DP2	IHS/CMI/LAB - print ind 10
BGP5DP3	IHS/CMI/LAB - print ind 1
BGP5DP4	IHS/CMI/LAB - print ind 31
BGP5DP5	IHS/CMI/LAB - print ind H
BGP5DP50	IHS/CMI/LAB - print ind H
BGP5DP51	IHS/CMI/LAB - print ind H
BGP5DP6	IHS/CMI/LAB - print ind H
BGP5DP60	IHS/CMI/LAB - print ind H
BGP5DP7	IHS/CMI/LAB - print ind 10

BGP5DPH	IHS/CMI/LAB - AREA REPORT HEADER
BGP5DPP	IHS/CMI/LAB - IHS GPRA 05 SELECTED REPORT DRIVER ALL PATIENTS
BGP5DS	IHS/CMI/LAB - IHS gpra print
BGP5DSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5DSL	IHS/CMI/LAB - FY 05 DISPLAY IND LISTS
BGP5DSP	IHS/CMI/LAB - IHS summary page
BGP5EL10	IHS/CMI/LAB - print ind 1
BGP5EL11	IHS/CMI/LAB - print ind 1
BGP5EL12	IHS/CMI/LAB - print ind 1
BGP5EL13	IHS/CMI/LAB - print ind 1
BGP5EL14	IHS/CMI/LAB - print ind 1
BGP5EL15	IHS/CMI/LAB - print ind 1
BGP5EL16	IHS/CMI/LAB - print ind 1
BGP5EL17	IHS/CMI/LAB - print ind 1
BGP5EL18	IHS/CMI/LAB - print ind 1
BGP5EL19	IHS/CMI/LAB - print ind 1
BGP5EL1A	IHS/CMI/LAB - print ind 1
BGP5EL1B	IHS/CMI/LAB - print ind 1
BGP5EL1C	IHS/CMI/LAB - print ind 1
BGP5EL1D	IHS/CMI/LAB - print ind 1
BGP5EL1E	IHS/CMI/LAB - print ind 1
BGP5EL1F	IHS/CMI/LAB - print ind 1
BGP5EL1G	IHS/CMI/LAB - print ind 1
BGP5EL1H	IHS/CMI/LAB - print ind 1
BGP5EL1I	IHS/CMI/LAB - print ind 1
BGP5EL1J	IHS/CMI/LAB - print ind 1
BGP5EL1K	IHS/CMI/LAB - print ind 1
BGP5EL1L	IHS/CMI/LAB - print ind 1
BGP5EL2	IHS/CMI/LAB - indicator 1,2,3,4 17 Jan 2005 6:49 AM
BGP5EL3	IHS/CMI/LAB - indicator 1,2,3,4 05 Apr 2005 1:44 PM
BGP5EL4	IHS/CMI/LAB - indicator 1,2,3,4 05 Apr 2005 1:44 PM
BGP5ELH	IHS/CMI/LAB - cover page for ELDER
BGP5ELHH	IHS/CMI/LAB - cover page for ELDER
BGP5ELL	IHS/CMI/LAB - IHS gpra print
BGP5ELL1	IHS/CMI/LAB - print ind 1
BGP5ELL2	IHS/CMI/LAB - print ind
BGP5ELL3	IHS/CMI/LAB - print ind
BGP5ELL4	IHS/CMI/LAB - print ind
BGP5ELL5	IHS/CMI/LAB - print ind
BGP5ELL6	IHS/CMI/LAB - print ind
BGP5ELL7	IHS/CMI/LAB - print ind
BGP5ELL8	IHS/CMI/LAB - print ind
BGP5ELL9	IHS/CMI/LAB - print ind
BGP5ELLA	IHS/CMI/LAB - print ind
BGP5ELLB	IHS/CMI/LAB - print ind
BGP5ELLC	IHS/CMI/LAB - print ind
BGP5ELLD	IHS/CMI/LAB - print ind
BGP5ELLE	IHS/CMI/LAB - print ind
BGP5ELLF	IHS/CMI/LAB - print ind
BGP5ELLG	IHS/CMI/LAB - print ind
BGP5ELLH	IHS/CMI/LAB - print ind

BGP5ELLI	IHS/CMI/LAB - print ind
BGP5ELLJ	IHS/CMI/LAB - print ind
BGP5ELLK	IHS/CMI/LAB - print ind
BGP5ELLL	IHS/CMI/LAB - print ind
BGP5ELLM	IHS/CMI/LAB - print ind
BGP5ELLN	IHS/CMI/LAB - print ind 1
BGP5ELLQ	IHS/CMI/LAB - print ind
BGP5ELP1	IHS/CMI/LAB - print ind 1
BGP5ELS	IHS/CMI/LAB - IHS gpra print
BGP5ELSL	IHS/CMI/LAB - IHS gpra print
BGP5ETL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5ETS	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5EUTL	IHS/CMI/LAB -
BGP5GAEL	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP5GAGP	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP5GAHE	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP5GALL	IHS/CMI/LAB - GUI COM REPORT
BGP5GAPU	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP5GCMS	IHS/CMI/LAB - GUI CMS REPORT
BGP5GCOM	IHS/CMI/LAB - GUI COM REPORT
BGP5GELD	IHS/CMI/LAB - GUI COM REPORT
BGP5GGPU	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP5GHED	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP5GLST	IHS/CMI/LAB - GUI CMS REPORT
BGP5GNPL	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP5GNTL	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP5GPP	IHS/CMI/LAB - GUI PP REPORT
BGP5GR	IHS/CMI/LAB - BGPG Visual CRS Reports 12/30/2004 12:29:35 PM
BGP5GRA	IHS/CMI/LAB - BGP Gui Area Reports 5/2/2005 8:38:59 PM
BGP5GT	IHS/CMI/LAB - BGPG Gui CRS Tables 2/2/2005 10:24:22 AM
BGP5GU	IHS/CMI/LAB - BGP Gui Utilities 10/29/2004 3:28:39 PM
BGP5GUPL	IHS/CMI/LAB - GUI CMS REPORT
BGP5HEH	IHS/CMI/LAB - cover page for HEDIS
BGP5HEHH	IHS/CMI/LAB - cover page for HEDIS
BGP5HEL	IHS/CMI/LAB - IHS gpra print
BGP5HEL1	IHS/CMI/LAB - print ind 1
BGP5HEP1	IHS/CMI/LAB - print ind 1
BGP5HES	IHS/CMI/LAB - IHS gpra print
BGP5HESL	IHS/CMI/LAB - IHS gpra print
BGP5HUTL	IHS/CMI/LAB -
BGP5LSTF	IHS/CMI/LAB - List 'BG05' files in pub
BGP5NPL	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP5NPLD	IHS/CMI/LAB - IHS gpra print
BGP5NPLP	IHS/CMI/LAB - IHS gpra print
BGP5PARP	IHS/CMI/LAB - IHS gpra print
BGP5PDL	IHS/CMI/LAB - IHS gpra print
BGP5PDL0	IHS/CMI/LAB - print ind 1
BGP5PDL1	IHS/CMI/LAB - print ind 1
BGP5PDL2	IHS/CMI/LAB - print ind 1
BGP5PDL3	IHS/CMI/LAB - print ind 1
BGP5PDL4	IHS/CMI/LAB - print ind 1
BGP5PDL5	IHS/CMI/LAB - print ind 1

BGP5PDL7	IHS/CMI/LAB - print ind 1
BGP5PDL8	IHS/CMI/LAB - print ind 1
BGP5PDL9	IHS/CMI/LAB - print ind
BGP5PDLH	IHS/CMI/LAB - cover page for gpra del
BGP5PDLS	IHS/CMI/LAB - IHS gpra print
BGP5PHEL	IHS/CMI/LAB - IHS gpra print
BGP5PHEP	IHS/CMI/LAB - IHS gpra print
BGP5SDP	IHS/CMI/LAB - IHS summary page
BGP5TXCH	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP5TXCN	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP5TXEL	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP5TXHE	IHS/CMI/LAB - TAXONOMY CHECK FOR FY05 HEDIS REPORT
BGP5TXL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5TXS	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5TXV	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5TXV1	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5ULF	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP5UTL	IHS/CMI/LAB -
BGP5UTL1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP6AHW	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP6ASL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6AX	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXA	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXB	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXC	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXD	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXE	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXF	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXG	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXGB	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXGC	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXH	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXI	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXJ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXK	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXL	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXM	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXN	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXO	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXP	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXQ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXR	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXS	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXT	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXU	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXV	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXW	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXX	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXY	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXZ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BAN	IHS/CMI/LAB - BANNER FOR CRS 05 06 Jan 2006 9:02 AM
BGP6BX	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005

BGP6BXA	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXB	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBB	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBC	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBD	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBE	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBF	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBG	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXC	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXD	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXE	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXF	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXG	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXH	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXI	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXJ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXK	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXL	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXM	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXN	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXO	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXP	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXQ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXR	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXS	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXT	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXU	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXV	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXW	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXX	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXY	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXZ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6C1	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM 04 May 2006 2:38 PM
BGP6C11	IHS/CMI/LAB - calc CMS measures
BGP6C12	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP6C13	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP6C2	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP6C3	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP6CL	IHS/CMI/LAB - IHS GPRA 05 CMS REPORT DRIVER
BGP6CP	IHS/CMI/LAB - IHS gpra print
BGP6CSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6CTL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6CTS	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6CTXC	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT 18 Feb 2006 5:31 PM
BGP6D1	IHS/CMI/LAB - calc measures 29 Apr 2006 7:38 PM
BGP6D10	IHS/CMI/LAB - calc measures
BGP6D2	IHS/CMI/LAB - measure 1,2,3,4
BGP6D21	IHS/CMI/LAB - measure 6
BGP6D22	IHS/CMI/LAB - measure I2
BGP6D23	IHS/CMI/LAB - measure I2

BGP6D3	IHS/CMI/LAB - measure 11
BGP6D31	IHS/CMI/LAB - measure C
BGP6D32	IHS/CMI/LAB - measure C
BGP6D33	IHS/CMI/LAB - measure C
BGP6D34	IHS/CMI/LAB - measure C
BGP6D35	IHS/CMI/LAB - measure C
BGP6D4	IHS/CMI/LAB - measure 3
BGP6D41	IHS/CMI/LAB - measure 3
BGP6D42	IHS/CMI/LAB - measure 11
BGP6D5	IHS/CMI/LAB - measure calc
BGP6D51	IHS/CMI/LAB - measure I2
BGP6D52	IHS/CMI/LAB - measure 31
BGP6D6	IHS/CMI/LAB - measure 31
BGP6D61	IHS/CMI/LAB - measure 31
BGP6D7	IHS/CMI/LAB - measure 31
BGP6D71	IHS/CMI/LAB - measure C
BGP6D72	IHS/CMI/LAB - measure 31
BGP6D73	IHS/CMI/LAB - measure 31
BGP6D8	IHS/CMI/LAB - measure C
BGP6D9	IHS/CMI/LAB - measure J
BGP6DAP	IHS/CMI/LAB - IHS GPRA 04 SELECTED REPORT DRIVER ALL PATIENTS
BGP6DAR	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP6DCEI	IHS/CMI/LAB - calculate HEDIS measures
BGP6DCHA	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP6DCHI	IHS/CMI/LAB - calculate HEDIS measures
BGP6DCHW	IHS/CMI/LAB - calc measures 29 Apr 2006 7:38 PM
BGP6DCI	IHS/CMI/LAB - IHS area GPRA
BGP6DCL	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP6DCLD	IHS/CMI/LAB - IHS gpra print
BGP6DCLP	IHS/CMI/LAB - IHS gpra print
BGP6DEL	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER
BGP6DELA	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER
BGP6DELP	IHS/CMI/LAB - IHS HEDIS print
BGP6DESI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6DESL	IHS/CMI/LAB - FY 06 DISPLAY IND LISTS
BGP6DGPA	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP6DGPU	IHS/CMI/LAB - IHS AREA CLIN 05 REPORT DRIVER
BGP6DH	IHS/CMI/LAB - cover page for gpra
BGP6DHE	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER
BGP6DHE1	IHS/CMI/LAB - measure HEDIS
BGP6DHEA	IHS/CMI/LAB - ihs area AA
BGP6DHEP	IHS/CMI/LAB - IHS HEDIS print
BGP6DHSL	IHS/CMI/LAB - FY 06 DISPLAY IND LISTS
BGP6DL	IHS/CMI/LAB - IHS GPRA 05 SELECTED REPORT DRIVER
BGP6DLT	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP6DMT	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP6DNG	IHS/CMI/LAB - IHS AREA CLIN 05 REPORT DRIVER
BGP6DP	IHS/CMI/LAB - IHS gpra print
BGP6DP1	IHS/CMI/LAB - print ind 1 21 Mar 2006 12:55 PM
BGP6DP11	IHS/CMI/LAB - print ind 1
BGP6DP12	IHS/CMI/LAB - print ind 19 AGE DIST

BGP6DP2	IHS/CMI/LAB - print ind 10
BGP6DP3	IHS/CMI/LAB - print ind 1
BGP6DP4	IHS/CMI/LAB - print ind 31
BGP6DP5	IHS/CMI/LAB - print ind H
BGP6DP50	IHS/CMI/LAB - print ind H
BGP6DP51	IHS/CMI/LAB - print ind H
BGP6DP6	IHS/CMI/LAB - print ind H
BGP6DP60	IHS/CMI/LAB - print ind H
BGP6DP7	IHS/CMI/LAB - print ind 10
BGP6DPH	IHS/CMI/LAB - AREA REPORT HEADER
BGP6DPP	IHS/CMI/LAB - IHS GPRA 05 SELECTED REPORT DRIVER ALL PATIENTS
BGP6DS	IHS/CMI/LAB - IHS gpra print
BGP6DSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6DSL	IHS/CMI/LAB - FY 06 DISPLAY IND LISTS
BGP6DSP	IHS/CMI/LAB - IHS summary page
BGP6DU	IHS/CMI/LAB - gpra utility calls
BGP6EL10	IHS/CMI/LAB - print ind 1
BGP6EL11	IHS/CMI/LAB - print ind 1
BGP6EL12	IHS/CMI/LAB - print ind 1
BGP6EL13	IHS/CMI/LAB - print ind 1
BGP6EL14	IHS/CMI/LAB - print ind 1
BGP6EL15	IHS/CMI/LAB - print ind 1
BGP6EL16	IHS/CMI/LAB - print ind 1
BGP6EL17	IHS/CMI/LAB - print ind 1
BGP6EL18	IHS/CMI/LAB - print ind 1
BGP6EL19	IHS/CMI/LAB - print ind 1
BGP6EL1A	IHS/CMI/LAB - print ind 1
BGP6EL1B	IHS/CMI/LAB - print ind 1
BGP6EL1C	IHS/CMI/LAB - print ind 1
BGP6EL1D	IHS/CMI/LAB - print ind 1
BGP6EL1E	IHS/CMI/LAB - print ind 1
BGP6EL1F	IHS/CMI/LAB - print ind 1
BGP6EL1G	IHS/CMI/LAB - print ind 1
BGP6EL1H	IHS/CMI/LAB - print ind 1
BGP6EL1I	IHS/CMI/LAB - print ind 1
BGP6EL1J	IHS/CMI/LAB - print ind 1
BGP6EL1K	IHS/CMI/LAB - print ind 1
BGP6EL1L	IHS/CMI/LAB - print ind 1
BGP6EL1M	IHS/CMI/LAB - print ind 1
BGP6EL2	IHS/CMI/LAB - measure 1,2,3,4 17 Jan 2006 6:49 AM
BGP6EL3	IHS/CMI/LAB - measure 1,2,3,4 05 Apr 2006 1:44 PM
BGP6EL4	IHS/CMI/LAB - measure 1,2,3,4 05 Apr 2006 1:44 PM
BGP6ELH	IHS/CMI/LAB - cover page for ELDER
BGP6ELHH	IHS/CMI/LAB - cover page for ELDER
BGP6ELL	IHS/CMI/LAB - IHS gpra print
BGP6ELL1	IHS/CMI/LAB - print ind 1
BGP6ELL2	IHS/CMI/LAB - print ind
BGP6ELL3	IHS/CMI/LAB - print ind
BGP6ELL4	IHS/CMI/LAB - print ind
BGP6ELL5	IHS/CMI/LAB - print ind
BGP6ELL6	IHS/CMI/LAB - print ind

BGP6ELL7	IHS/CMI/LAB - print ind
BGP6ELL8	IHS/CMI/LAB - print ind
BGP6ELL9	IHS/CMI/LAB - print ind
BGP6ELLA	IHS/CMI/LAB - print ind
BGP6ELLB	IHS/CMI/LAB - print ind
BGP6ELLC	IHS/CMI/LAB - print ind
BGP6ELLD	IHS/CMI/LAB - print ind
BGP6ELLE	IHS/CMI/LAB - print ind
BGP6ELLF	IHS/CMI/LAB - print ind
BGP6ELLG	IHS/CMI/LAB - print ind
BGP6ELLH	IHS/CMI/LAB - print ind
BGP6ELLI	IHS/CMI/LAB - print ind
BGP6ELLJ	IHS/CMI/LAB - print ind
BGP6ELLK	IHS/CMI/LAB - print ind
BGP6ELLL	IHS/CMI/LAB - print ind
BGP6ELLM	IHS/CMI/LAB - print ind
BGP6ELLN	IHS/CMI/LAB - print ind 1
BGP6ELLO	IHS/CMI/LAB - print ind
BGP6ELLQ	IHS/CMI/LAB - print ind
BGP6ELP1	IHS/CMI/LAB - print ind 1
BGP6ELS	IHS/CMI/LAB - IHS gpra print
BGP6ELSL	IHS/CMI/LAB - IHS gpra print
BGP6ETL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6EUTL	IHS/CMI/LAB -
BGP6EX	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXA	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXB	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXBB	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXC	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXD	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXE	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXF	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXG	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXH	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXHB	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXHC	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXI	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXJ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXK	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXL	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXM	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXN	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXO	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXP	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXQ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXR	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXS	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXT	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXU	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXV	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXW	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXX	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005

BGP6EXY	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXZ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FX	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXA	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXB	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXC	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXD	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXE	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXF	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXG	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6GACW	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP6GAEL	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP6GAGP	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP6GAHE	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP6GALL	IHS/CMI/LAB - GUI COM REPORT
BGP6GAPU	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP6GCMP	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GCMS	IHS/CMI/LAB - GUI CMS REPORT
BGP6GCOM	IHS/CMI/LAB - GUI COM REPORT
BGP6GELD	IHS/CMI/LAB - GUI COM REPORT
BGP6GGPU	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GHED	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GLST	IHS/CMI/LAB - GUI CMS REPORT
BGP6GLTX	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GMTX	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GNPL	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GNTL	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GPP	IHS/CMI/LAB - GUI PP REPORT
BGP6GR	IHS/CMI/LAB - BGPG Visual CRS Reports 12/30/2004 12:29:35 PM
BGP6GRA	IHS/CMI/LAB - BGP Gui Area Reports 5/2/2005 8:38:59 PM
BGP6GT	IHS/CMI/LAB - BGPG Gui CRS Tables 2/2/2005 10:24:22 AM
BGP6GU	IHS/CMI/LAB - BGP Gui Utilities 10/29/2004 3:28:39 PM
BGP6GUPL	IHS/CMI/LAB - GUI CMS REPORT
BGP6HEH	IHS/CMI/LAB - cover page for HEDIS
BGP6HEHH	IHS/CMI/LAB - cover page for HEDIS
BGP6HEL	IHS/CMI/LAB - IHS gpra print
BGP6HEL1	IHS/CMI/LAB - print ind 1
BGP6HEP1	IHS/CMI/LAB - print ind 1
BGP6HES	IHS/CMI/LAB - IHS gpra print
BGP6HESL	IHS/CMI/LAB - IHS gpra print
BGP6HUTL	IHS/CMI/LAB -
BGP6HX	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXA	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXB	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXC	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXD	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXE	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXF	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXG	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXH	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXI	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXJ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005

BGP6HXK	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXL	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXM	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXN	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXO	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXP	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXQ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXR	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXS	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXT	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXU	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXV	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXW	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXX	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXY	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXZ	IHS/CMI/LAB-CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6LSTF	IHS/CMI/LAB - List 'BG06' files in pub
BGP6NPL	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP6NPLD	IHS/CMI/LAB - IHS gpra print
BGP6NPLP	IHS/CMI/LAB - IHS gpra print
BGP6PARP	IHS/CMI/LAB - IHS gpra print
BGP6PDL	IHS/CMI/LAB - IHS gpra print
BGP6PDL0	IHS/CMI/LAB - print ind 1
BGP6PDL1	IHS/CMI/LAB - print ind 1
BGP6PDL2	IHS/CMI/LAB - print ind 1
BGP6PDL3	IHS/CMI/LAB - print ind 1
BGP6PDL4	IHS/CMI/LAB - print ind 1
BGP6PDL5	IHS/CMI/LAB - print ind 1
BGP6PDL7	IHS/CMI/LAB - print ind 1
BGP6PDL8	IHS/CMI/LAB - print ind 1
BGP6PDL9	IHS/CMI/LAB - print ind
BGP6PDLA	IHS/CMI/LAB - print ind
BGP6PDLH	IHS/CMI/LAB - cover page for gpra del
BGP6PDLS	IHS/CMI/LAB - IHS gpra print
BGP6PHEL	IHS/CMI/LAB - IHS gpra print
BGP6PHEP	IHS/CMI/LAB - IHS gpra print
BGP6SDP	IHS/CMI/LAB - IHS summary page
BGP6ULF	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP6UTL	IHS/CMI/LAB -
BGP6UTL1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP6XTCH	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP6XTCN	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP6XTEL	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP6XTHE	IHS/CMI/LAB - TAXONOMY CHECK FOR FY06 HEDIS REPORT
BGP6XTL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6XTS	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6XTV	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6XTV1	IHS/CMI/LAB - DISPLAY IND LISTS
BGPAP	IHS/CMI/LAB - IHS gpra print
BGPAP1	IHS/CMI/LAB - print ind 1
BGPAP12	IHS/CMI/LAB - print ind 12
BGPAP13	IHS/CMI/LAB - print ind 13

BGPAP14	IHS/CMI/LAB - print ind 14
BGPAP1B	IHS/CMI/LAB - print ind 1B
BGPAP2	IHS/CMI/LAB - print ind 2
BGPAP22	IHS/CMI/LAB - print ind 22
BGPAP221	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGPAP22A	IHS/CMI/LAB - print ind 22
BGPAP22B	IHS/CMI/LAB - print ind 22
BGPAP22C	IHS/CMI/LAB - print ind 22
BGPAP23	IHS/CMI/LAB - print ind 23
BGPAP24	IHS/CMI/LAB - print ind 24
BGPAP3	IHS/CMI/LAB - print ind 3
BGPAP30	IHS/CMI/LAB - IHS gpra print
BGPAP30A	IHS/CMI/LAB - IHS gpra print
BGPAP30B	IHS/CMI/LAB - IHS gpra print
BGPAP30C	IHS/CMI/LAB - IHS gpra print
BGPAP4	IHS/CMI/LAB -print ind 4
BGPAP5	IHS/CMI/LAB - print ind 5
BGPAP6	IHS/CMI/LAB - print ind 6
BGPAP7	IHS/CMI/LAB - print ind 7
BGPAP8	IHS/CMI/LAB - print ind 8
BGPAPA	IHS/CMI/LAB - IHS gpra print
BGPAPB	IHS/CMI/LAB - IHS gpra print
BGPAPC	IHS/CMI/LAB - IHS gpra print
BGPAPD	IHS/CMI/LAB - IHS gpra print
BGPAPH	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGPD	IHS/CMI/LAB - IHS GPRA - report for local use
BGPD0	IHS/CMI/LAB - DISPLAY IND LISTS
BGPD1	IHS/CMI/LAB - IHS area GPRA
BGPD12	IHS/CMI/LAB - indicator 12
BGPD13	IHS/CMI/LAB - indicator 13
BGPD14	IHS/CMI/LAB - indicator 14
BGPD1B	IHS/CMI/LAB - ind 1b
BGPD2	IHS/CMI/LAB - indicator 2
BGPD22	IHS/CMI/LAB - indicator 18
BGPD23	IHS/CMI/LAB - indicator 23
BGPD24	IHS/CMI/LAB - indicator 24
BGPD29	IHS/CMI/LAB - indicator 29
BGPD3	IHS/CMI/LAB - indicator 3
BGPD30	IHS/CMI/LAB - indicator 30
BGPD31	IHS/CMI/LAB - indicator 31
BGPD4	IHS/CMI/LAB - indicator 4
BGPD5	IHS/CMI/LAB - indicator 5
BGPD6	IHS/CMI/LAB - indicator 6
BGPD7	IHS/CMI/LAB - indicator 7
BGPD8	IHS/CMI/LAB - indicator 8
BGPDA	IHS/CMI/LAB - indicator A
BGPDAREA	IHS/CMI/LAB - ihs area GPRA
BGPDARP	IHS/CMI/LAB - IHS area GPRA
BGPDB	IHS/CMI/LAB - indicator B
BGPDC	IHS/CMI/LAB - indicator C
BGPDD	IHS/CMI/LAB - indicator D
BGPDFTA	IHS/CMI/LAB - IHS area GPRA

BGPDH	IHS/CMI/LAB - cover page for gpra
BGPDL	IHS/CMI/LAB - DISPLAY IND LISTS
BGPDP	IHS/CMI/LAB - IHS gpra print
BGPDP1	IHS/CMI/LAB - print ind 1
BGPDP12	IHS/CMI/LAB - print ind 12
BGPDP13	IHS/CMI/LAB - print ind 13
BGPDP14	IHS/CMI/LAB - print ind 14
BGPDP1B	IHS/CMI/LAB - print ind 1B
BGPDP2	IHS/CMI/LAB - print ind 2
BGPDP22	IHS/CMI/LAB - print ind 22
BGPDP221	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGPDP22A	IHS/CMI/LAB - print ind 22
BGPDP22B	IHS/CMI/LAB - print ind 22
BGPDP22C	IHS/CMI/LAB - print ind 22
BGPDP23	IHS/CMI/LAB - print ind 23
BGPDP24	IHS/CMI/LAB - print ind 24
BGPDP29	IHS/CMI/LAB - IHS gpra print
BGPDP3	IHS/CMI/LAB - print ind 3
BGPDP30	IHS/CMI/LAB - IHS gpra print
BGPDP30A	IHS/CMI/LAB - IHS gpra print
BGPDP30B	IHS/CMI/LAB - IHS gpra print
BGPDP30C	IHS/CMI/LAB - IHS gpra print
BGPDP31	IHS/CMI/LAB - print ind 31
BGPDP32	IHS/CMI/LAB - print ind 32
BGPDP4	IHS/CMI/LAB -print ind 4
BGPDP5	IHS/CMI/LAB - print ind 5
BGPDP6	IHS/CMI/LAB - print ind 6
BGPDP7	IHS/CMI/LAB - print ind 7
BGPDP8	IHS/CMI/LAB - print ind 8
BGPDPA	IHS/CMI/LAB - IHS gpra print
BGPDPARP	IHS/CMI/LAB - IHS gpra print
BGPDPB	IHS/CMI/LAB - IHS gpra print
BGPDPC	IHS/CMI/LAB - IHS gpra print
BGPDPD	IHS/CMI/LAB - IHS gpra print
BGPDPH	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGPDS	IHS/CMI/LAB - IHS gpra print
BGPDT	IHS/CMI/LAB - 2001 DIABETES AUDIT
BGPDU	IHS/CMI/LAB - gpra utility calls
BGPDUP	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGPGR	IHS/CMI/LAB - BGPG Visual CRS Reports 12/30/2004 12:29:35 PM
BGPGRA	IHS/CMI/LAB - BGP Gui Area Reports 5/2/2005 8:38:59 PM
BGPGBT	IHS/CMI/LAB - BGPG Gui CRS Tables 2/2/2005 10:24:22 AM
BGPGU	IHS/CMI/LAB - BGP Gui Utilities 10/29/2004 3:28:39 PM
BGPPOST	IHS/CMI/LAB - NO DESCRIPTION PROVIDED 28 Jan 2005 1:34 PM
BGPPOST1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED 28 Jan 2005 1:34 PM
BGPTX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXEB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005

BGPTXEC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXED	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXEE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXEF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXEG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXGB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXPB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005

BGPVTXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 26, 2005
BGPWTXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 24, 2005
BGPWTXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 24, 2005
BGPWTXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 26, 2005

4.0 Files and Tables

4.1 File List

File	Global	Filename
90240.01	^BGPD(BGP IHS GPRA DATA 02
90241.01	^BGPCTRL(BGP CONTROL FILE
90241.02	^BGPSITE(BGP SITE PARAMETERS
90242	^BGPINDC(BGP INDIVIDUAL INDICATOR DEFINITIONS
90242.01	^BGPIND(BGP INDICATORS
90243	^BGPGPDC(BGP 03 DATA CURRENT
90243.01	^BGPGPDP(BGP 03 DATA PREVIOUS
90243.02	^BGPGPDB(BGP 03 DATA BASELINE
90244.01	^BGPINDF(BGP 04 INDICATORS
90244.02	^BGPINDFC(BGP 04 INDIVIDUAL INDICATORS
90244.03	^BGPGPDCF(BGP 04 DATA CURRENT
90244.04	^BGPGPDPF(BGP 04 DATA PREVIOUS
90244.05	^BGPGPDBF(BGP 04 DATA BASELINE
90370.01	^BGPHEIF(BGP 04 HEDIS INDICATORS
90370.02	^BGPHEIIF(BGP 04 HEDIS INDIVIDUAL INDICATORS
90370.03	^BGPHEDCF(BGP 04 HEDIS DATA CURRENT
90370.04	^BGPHEDPF(BGP 04 HEDIS DATA PREVIOUS
90370.05	^BGPHEDBF(BGP 04 HEDIS DATA BASELINE
90371.01	^BGPINDV(BGP 05 INDICATORS
90371.02	^BGPINDVC(BGP 05 INDIVIDUAL INDICATORS
90371.03	^BGPGPDCV(BGP 05 DATA CURRENT
90371.04	^BGPGPDPV(BGP 05 DATA PREVIOUS
90371.05	^BGPGPDBV(BGP 05 DATA BASELINE
90371.06	^BGPVNPL(BGP 05 NATIONAL PATIENT LISTS
90371.08	^BGP TAXV(BGP 05 TAXONOMIES
90372.01	^BGPHEIV(BGP 05 HEDIS INDICATORS
90372.02	^BGPHEIIV(BGP 05 HEDIS INDIVIDUAL INDICATORS
90372.03	^BGPHEDCV(BGP 05 HEDIS DATA CURRENT
90372.04	^BGPHEDPV(BGP 05 HEDIS DATA PREVIOUS
90372.05	^BGPHEDBV(BGP 05 HEDIS DATA BASELINE
90372.08	^BPGUI(BGP 05 GUI REPORT OUTPUT
90373.01	^BGP CMSIF(BGP 05 CMS INDICATORS
90373.02	^BGP ELIIV(BGP 05 ELDER INDIVIDUAL INDICATORS
90373.03	^BGP ELD CV(BGP 05 ELDER DATA CURRENT
90373.04	^BGP ELD PV(BGP 05 ELDER DATA PREVIOUS
90373.05	^BGP ELD BV(BGP 05 ELDER DATA BASELINE
90374.01	^BGPINDS(BGP 06 INDICATORS
90374.02	^BGPINDSC(BGP 06 INDIVIDUAL INDICATORS
90374.03	^BGP GPDCS(BGP 06 DATA CURRENT

90374.04	^BGPGPDP(BGP 06 DATA PREVIOUS
90374.05	^BGPGPDBS(BGP 06 DATA BASELINE
90374.06	^BGPSNPL(BGP 06 NATIONAL PATIENT LISTS
90374.08	^BGPTAXS(BGP 06 TAXONOMIES
90375.01	^BGPHEIS(BGP 06 HEDIS INDICATORS
90375.02	^BGPHEIIS(BGP 06 HEDIS INDIVIDUAL INDICATORS
90375.03	^BGPHEDCS(BGP 06 HEDIS DATA CURRENT
90375.04	^BGPHEDPS(BGP 06 HEDIS DATA PREVIOUS
90375.05	^BGPHEDBS(BGP 06 HEDIS DATA BASELINE
90375.08	^BPGUI(BGP 06 GUI REPORT OUTPUT
90376.01	^BGP CMSIS(BGP 06 CMS INDICATORS
90376.02	^BGP ELIIS(BGP 06 ELDER INDIVIDUAL INDICATORS
90376.03	^BGP ELDCS(BGP 06 ELDER DATA CURRENT
90376.04	^BGP ELDPS(BGP 06 ELDER DATA PREVIOUS
90376.05	^BGP ELDDBS(BGP 06 ELDER DATA BASELINE
90376.06	^BGP ELIS(BGP 06 ELDER INDICATORS

4.2 File Access

FILE (#)	GLOBAL	RD	WR	LYG	DD	DEL
90240.01	^BPGD	@	@	@	@	@
90241.01	^BGPCTRL	@	@	@	@	@
90241.02	^BGP SITE	M	M	M	@	@
90242	^BGPINDC	@	@	@	@	@
90242.01	^BGPIND	@	@	@	@	@
90243	^BGP GPDC	@	@	@	@	@
90243.01	^BGP GPDP	@	@	@	@	@
90243.02	^BGP GPDB	@	@	@	@	@
90244.01	^BGP INDF(0)	@	@	@	@	@
90244.02	^BGP INDFC(0)	@	@	@	@	@
90244.03	^BGP GPDCF(0)	@	@	@	@	@
90244.04	^BGP GPDPF(0)	@	@	@	@	@
90244.05	^BGP GPDBF(0)	@	@	@	@	@
90370.01	^BGP HEIF(0)	@	@	@	@	@
90370.02	^BGP HEIIF(0)	@	@	@	@	@
90370.03	^BGP HEDCF(0)	@	@	@	@	@
90370.04	^BGP HEDPF(0)	@	@	@	@	@
90370.05	^BGP HEDBF(0)	@	@	@	@	@
90370.01	^BGP HEIF(0)	@	@	@	@	@
90370.02	^BGP HEIIF(0)	@	@	@	@	@
90370.03	^BGP HEDCF(0)	@	@	@	@	@
90370.04	^BGP HEDPF	@	@	@	@	@
90370.05	^BGP HEDBF	@	@	@	@	@
90371.01	^BGP INDV(0)	@	@	@	@	@
90371.02	^BGP INDVC(0)	@	@	@	@	@

FILE (#)	GLOBAL	RD	WR	LYG	DD	DEL
90371.03	^BGPGPDCV(0)	@	@	@	@	@
90371.04	^BGPGPDPV(0)	@	@	@	@	@
90371.05	^BGPGPDBV(0)	@	@	@	@	@
90371.06	^BGPVNPL(0)	@	@	@	@	@
90371.08	^BGPTAXV(@	@	@	@	@
90372.01	^BGPHEIV(0)	@	@	@	@	@
90372.02	^BGPHEIIV(0)	@	@	@	@	@
90372.03	^BGPHEDCV(0)	@	@	@	@	@
90372.04	^BGPHEDPV	@	@	@	@	@
90372.05	^BGPHEDBV	@	@	@	@	@
90372.08	^BGPGUI(@	@	@	@	@
90373.01	^BGPcmsif(@	@	@	@	@
90373.02	^BGPeliiv(@	@	@	@	@
90373.03	^BGPelDCV(@	@	@	@	@
90373.04	^BGPeldPV(@	@	@	@	@
90373.05	^BGPelDBV(@	@	@	@	@
90374.01	^BGPINDS(0)	@	@	@	@	@
90374.02	^BGPINDSC(0)	@	@	@	@	@
90374.03	^BGPgpDCS(0)	@	@	@	@	@
90374.04	^BGPgpDPs(0)	@	@	@	@	@
90374.05	^BGPgpDBS(0)	@	@	@	@	@
90374.06	^BGPsnPL(0)	@	@	@	@	@
90374.08	^BGPTAXS(@	@	@	@	@
90375.01	^BGPHEIS(0)	@	@	@	@	@
90375.02	^BGPHEIIS(0)	@	@	@	@	@
90375.03	^BGPHEDCS(0)	@	@	@	@	@
90375.04	^BGPHEDPS	@	@	@	@	@
90375.05	^BGPHEDBS	@	@	@	@	@
90375.08	^BPGuis(@	@	@	@	@
90376.01	^BGPcmsis(@	@	@	@	@
90376.02	^BGPeliis(@	@	@	@	@
90376.03	^BGPelDCS(@	@	@	@	@
90376.04	^BGPelDPs(@	@	@	@	@
90376.05	^BGPelDBS(@	@	@	@	@
90376.06	^BGPelIS(@	@	@	@	@

4.3 Cross References

90240.01 (BGP IHS GPRA DATA 02)

.01 BEGINNING DATE

B Regular type cross reference

18000 PHN VISIT PRIMARY DX (90240.018)

.01 PHN VISIT PRIMARY DX

B Regular type cross reference

24000 CURRENT INJURIES (90240.24)
.01 CURRENT INJURIES
 B Regular type cross reference
24000 E CODES (90240.2424)
.01 E CODES
 B Regular type cross reference
25000 INJURIES 98 PREV (90240.25)
.01 INJURIES 98 PREV
 B Regular type cross reference
25000 E CODE TALLY (90240.2525)
.01 E CODE TALLY
 B Regular type cross reference
26000 INJURIES 98 (90240.26)
.01 INJURIES 98
 B Regular type cross reference
26000 E CODES (90240.2626)
.01 E CODES
 B Regular type cross reference
28000 COMMUNITIES (90240.28)
.01 COMMUNITIES
 B Regular type cross reference
48000 PHN VISIT PRIMARY DX PR (90240.048)
.01 PHN VISIT PRIMARY DX PR
 B Regular type cross reference
88000 PHN PRIMARY DIAGNOSIS BL (90240.088)
.01 PHN PRIMARY DIAGNOSIS BL
 B Regular type cross reference

90241.01 (BGP CONTROL FILE)

.01 FY
 B Regular type cross reference
1101 CLINICS FOR CLINICAL USER (90241.0111)
.01 CLINIC CODE
 B Regular type cross reference
1201 CLINIC CODES FOR 2ND VISIT (90241.0112)
.01 CLINIC CODES FOR 2ND VISIT
 B Regular type cross reference

90241.02 (BGP SITE PARAMETERS)

.01 LOCATION
 B Regular type cross reference

90242 (BGP INDIVIDUAL INDICATOR DEFINITIONS)

.01 INDICATOR
 B Regular type cross reference
.04 INDICATOR ID
 C Regular type cross reference

AA MUMPS type cross reference
.05 GPRA
AGPRA Regular type cross reference
.06 AREA APP
AREA Regular type cross reference
.07 EISS
AEISS Regular type cross reference
.08 DENOM FIELD NUMBER
ADF Regular type cross reference
.09 NUMERATOR FIELD NUMBER
ANF Regular type cross reference
.12 DM INDICATOR?
ADM Regular type cross reference
.13 CARDIO INDICATOR?
ACARD Regular type cross reference
.14 EISS NUMERATOR PIECE
APO MUMPS type cross reference

90242.01 (BGP INDICATORS)

.01 FY
B Regular type cross reference
AO MUMPS type cross reference
.06 ORDER TO DISPLAY
AOTOO MUMPS type cross reference

90243 BGP 03 DATA CURRENT

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90243.09999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90243.099999)
.01 PATIENTS EISS
B Regular type cross reference

90243.01 (BGP 03 Data Previous)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90243.10999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90243.109999)
.01 PATIENTS EISS
B Regular type cross reference

90243.02 (BGP 03 Data Baseline)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90243.11999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90243.119999)
.01 PATIENTS EISS
B Regular type cross reference

90244.01 BGP 04 INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.02 ORDER ON SELECTION/REPORT
AO Regular type cross reference
.03 TITLE
C Regular type cross reference
.07 GPRA INDICATOR
GPRA Regular type cross reference
.08 AREA INDICATOR
AREA Regular type cross reference
.09 HEDIS INDICATOR
HEDIS Regular type cross reference
.11 LOCAL REPORT
LOCAL Regular type cross reference

90244.02 BGP 04 INDIVIDUAL INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.04 INDICATOR ID
C Regular type cross reference
AA Mumps cross reference
.05 GPRA INDICATOR
AGPRA Regular type cross reference
.06 AREA APP
B Regular type cross reference
.07 EISS
AEISS Regular type cross reference
.08 DENOMINATOR FIELD
ADF Regular type cross reference
.09 NUMERATOR FIELD
ANF Regular type cross reference
.12 DM INDICATOR?
ADM Regular type cross reference
.13 CARDIO INDICATOR?
ACARD Regular type cross reference
.14 EISS PIECE

APO Mumps type cross reference
1202 WOMEN'S HEALTH
AWH Regular type cross reference
1203 PREVENTION RELATED?
APR Regular type cross reference
1204 ORDER
ORDER Mumps type cross reference
AB Mumps type cross reference
OR Regular type cross reference

90244.03 (BGP 04 Data Current)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90244.119999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISSIONS (90243.129999)
.01 PATIENTS EISSIONS
B Regular type cross reference

90244.04 (BGP 04 Data Previous)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90244.13999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISSIONS (90243.139999)
.01 PATIENTS EISSIONS
B Regular type cross reference

90244.05 (BGP 04 Data Baseline)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90244.14999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISSIONS (90244.1499999)
.01 PATIENTS EISSIONS
B Regular type cross reference

90370.01 BGP 04 HEDIS INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.02 ORDER ON SELECTION/REPORT
AO Regular type cross reference
.03 TITLE
C Regular type cross reference

.07 GPRA INDICATOR
GPRA Regular type cross reference
.08 AREA INDICATOR
AREA Regular type cross reference
.09 HEDIS INDICATOR
HEDIS Regular type cross reference
.11 LOCAL REPORT
LOCAL Regular type cross reference

90370.02 BGP 04 HEDIS INDIVIDUAL INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.04 INDICATOR ID
C Regular type cross reference
AA Mumps cross reference
.05 GPRA INDICATOR
AGPRA Regular type cross reference
.06 AREA APP
B Regular type cross reference
.07 EISS
AEIIS Regular type cross reference
.08 DENOMINATOR FIELD
ADF Regular type cross reference
.09 NUMERATOR FIELD
ANF Regular type cross reference
.12 DM INDICATOR?
ADM Regular type cross reference
.13 CARDIO INDICATOR?
ACARD Regular type cross reference
.14 EISS PIECE
APO Mumps type cross reference
1205 WOMEN'S HEALTH
AWH Regular type cross reference
1206 PREVENTION RELATED?
APR Regular type cross reference
1207 ORDER
ORDER Mumps type cross reference
AB Mumps type cross reference
OR Regular type cross reference

90370.03 (BGP 04 HEDIS Data Current)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90370.119999)
.01 COMMUNITIES
B Regular type cross reference

99999 PATIENTS EISS (90243.129999)
.01 PATIENTS EISS
B Regular type cross reference

90370.04 (BGP 04 HEDIS Data Previous)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90370.13999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90243.139999)
.01 PATIENTS EISS
B Regular type cross reference

90370.05 (BGP 04 HEDIS Data Baseline)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90370.14999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90370.1499999)
.01 PATIENTS EISS
B Regular type cross reference

90371.01 BGP 05 INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.02 ORDER ON SELECTION/REPORT
AO Regular type cross reference
.03 TITLE
C Regular type cross reference
.07 GPRA INDICATOR
GPRA Regular type cross reference
.08 AREA INDICATOR
AREA Regular type cross reference
.09 HEDIS INDICATOR
HEDIS Regular type cross reference
.11 LOCAL REPORT
LOCAL Regular type cross reference

90371.02 BGP 05 INDIVIDUAL INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.04 INDICATOR ID
C Regular type cross reference
AA Mumps cross reference

.05 GPRA INDICATOR
AGPRA Regular type cross reference
.06 AREA APP
B Regular type cross reference
.07 EISS
AEISS Regular type cross reference
.08 DENOMINATOR FIELD
ADF Regular type cross reference
.09 NUMERATOR FIELD
ANF Regular type cross reference
.12 DM INDICATOR?
ADM Regular type cross reference
.13 CARDIO INDICATOR?
ACARD Regular type cross reference
.14 EISS PIECE
APO Mumps type cross reference
1208 WOMEN'S HEALTH
AWH Regular type cross reference
1209 PREVENTION RELATED?
APR Regular type cross reference
1210 ORDER
ORDER Mumps type cross reference
AB Mumps type cross reference
OR Regular type cross reference

90371.03 (BGP 05 Data Current)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90371.119999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90243.129999)
.01 PATIENTS EISS
B Regular type cross reference

90371.04 (BGP 05 Data Previous)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90371.13999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90243.139999)
.01 PATIENTS EISS
B Regular type cross reference

90371.05 (BGP 05 Data Baseline)

.01 BEGINNING DATE

B Regular type cross reference
9999 COMMUNITIES (90371.14999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90371.1499999)
.01 PATIENTS EISS
B Regular type cross reference

90371.08 BGP 05 TAXONOMIES

.01 TAXONOMY NAME
B Regular type cross reference
1200 REPORTS
B Regular type cross reference

90372.01 BGP 05 HEDIS INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.02 ORDER ON SELECTION/REPORT
AO Regular type cross reference
.03 TITLE
C Regular type cross reference
.07 GPRA INDICATOR
GPRA Regular type cross reference
.08 AREA INDICATOR
AREA Regular type cross reference
.09 HEDIS INDICATOR
HEDIS Regular type cross reference
.11 LOCAL REPORT
LOCAL Regular type cross reference

90372.02 BGP 05 HEDIS INDIVIDUAL INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.04 INDICATOR ID
C Regular type cross reference
AA Mumps cross reference
.05 GPRA INDICATOR
AGPRA Regular type cross reference
.06 AREA APP
B Regular type cross reference
.07 EISS
AEISS Regular type cross reference
.08 DENOMINATOR FIELD
ADF Regular type cross reference
.09 NUMERATOR FIELD
ANF Regular type cross reference

.12 DM INDICATOR?
ADM Regular type cross reference
.13 CARDIO INDICATOR?
ACARD Regular type cross reference
.14 EISS PIECE
APO Mumps type cross reference
1211 WOMEN'S HEALTH
AWH Regular type cross reference
1212 PREVENTION RELATED?
APR Regular type cross reference
1213 ORDER
ORDER Mumps type cross reference
AB Mumps type cross reference
OR Regular type cross reference

90372.03 (BGP 05 HEDIS Data Current)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90372.119999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90243.129999)
.01 PATIENTS EISS
B Regular type cross reference

90372.04 (BGP 05 HEDIS Data Previous)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90372.13999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90243.139999)
.01 PATIENTS EISS
B Regular type cross reference

90372.05 (BGP 05 HEDIS Data Baseline)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90372.14999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90372.149999)
.01 PATIENTS EISS

B Regular type cross reference

90372.08 BGP 05 GUI REPORT OUTPUT

.01 NAME
 B Regular type cross reference
.02 USER WHO REQUESTED REPORT
 AUSR Regular type cross reference
.03 START TIME
 AC Regular type cross reference

90373.01 BGP 05 CMS INDICATORS

.01 INDICATOR NUMBER
 B Regular type cross reference
.02 ORDER ON SELECTION/REPORT
 AO Regular type cross reference
.03 TITLE
 C Regular type cross reference

90373.02 BGP 05 ELDER INDIVIDUAL INDICATORS

.01 INDICATOR
 B Regular type cross reference
.04 INDICATOR ID
 C Regular type cross reference
1204 ORDER
 ABC MUMPS
 AB MUMPS
 OR Regular type cross reference
 AP MUMPS

90373.03 BGP 05 ELDER DATA CURRENT

.01 BEGINNING DATE
 B Regular type cross reference

90373.04 BGP 05 ELDER DATA PREVIOUS

.01 BEGINNING DATE
 B Regular type cross reference

90373.05 BGP 05 ELDER DATA BASELINE

.01 BEGINNING DATE
 B Regular type cross reference

90373.06 BGP 05 ELDER INDICATORS

.01 INDICATOR NUMBER
 B Regular type cross reference
.02 ORDER ON SELECTION/REPORT
 AO Regular type cross reference
.03 TITLE
 C Regular type cross reference
.07 NATIONAL REPORT

GPRA Regular type cross reference
.11 LOCAL REPORT
LOCAL Regular type cross reference
.01 DENOMINATOR DEFINITIONS
B Regular type cross reference
.01 NUMERATOR DEFINITIONS
B Regular type cross reference

90374.01 BGP 06 INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.02 ORDER ON SELECTION/REPORT
AO Regular type cross reference
.03 TITLE
C Regular type cross reference
.07 NATIONAL REPORT
GPRA Regular type cross reference
.11 LOCAL REPORT
LOCAL Regular type cross reference
.01 DENOMINATOR DEFINITIONS
B Regular type cross reference
.01 NUMERATOR DEFINITIONS
B Regular type cross reference

90374.02 BGP 06 INDIVIDUAL INDICATORS

.01 INDICATOR
B Regular type cross reference
.04 INDICATOR ID
C Regular type cross reference
.05 NATIONAL
AGPRA Regular type cross reference
.12 DM INDICATOR?
ADM Regular type cross reference
.13 CARDIO INDICATOR?
ACARD Regular type cross reference
1202 WOMEN'S HEALTH
AWH Regular type cross reference
1203 ELDER RELATED
AEL Regular type cross reference
1204 ORDER
ABC MUMPS
AB MUMPS
OR Regular type cross reference

90374.03 BGP 06 DATA CURRENT

.01 BEGINNING DATE
B Regular type cross reference

90374.04 BGP 06 DATA PREVIOUS

- .01 BEGINNING DATE
 - B Regular type cross reference

90374.05 BGP 06 DATA BASELINE

- .01 BEGINNING DATE
 - B Regular type cross reference

90374.06 BGP 06 NATIONAL PATIENT LISTS

- .01 POINTER
 - B Regular type cross reference

90374.08 BGP 06 TAXONOMIES

- .01 TAXONOMY NAME
 - B Regular type cross reference
- .01 REPORTS
 - B Regular type cross reference
- .01 COMMUNITIES
 - B Regular type cross reference
- .01 PATIENTS EISS
 - B Regular type cross reference
- .01 COMMUNITIES
 - B Regular type cross reference
- .01 PATIENTS EISS
 - B Regular type cross reference
- .01 COMMUNITIES
 - B Regular type cross reference
- .01 PATIENTS EISS
 - B Regular type cross reference

90375.01 BGP 06 HEDIS INDICATORS

- .01 INDICATOR NUMBER
 - B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT
 - AO Regular type cross reference
- .03 TITLE
 - C Regular type cross reference
- .09 HEDIS INDICATOR
 - HEDIS Regular type cross reference
- .01 DENOMINATOR DEFINITIONS
 - B Regular type cross reference
- .01 NUMERATOR DEFINITIONS
 - B Regular type cross reference

90375.02 BGP 06 HEDIS INDIVIDUAL INDICATORS

.01 INDICATOR
 B Regular type cross reference
.04 INDICATOR ID
 C Regular type cross reference
 AA MUMPS
.08 DENOM FIELD NUMBER
 ADF Regular type cross reference
.09 NUMERATOR FIELD NUMBER
 ANF Regular type cross reference
1204 ORDER
 AB MUMPS
 OR Regular type cross reference

90375.03 BGP 06 HEDIS DATA CURRENT

.01 BEGINNING DATE
 B Regular type cross reference

90375.04 BGP 06 HEDIS DATA PREVIOUS

.01 BEGINNING DATE
 B Regular type cross reference

90375.05 BGP 06 HEDIS DATA BASELINE

.01 BEGINNING DATE
 B Regular type cross reference

90375.08 BGP 06 GUI REPORT OUTPUT

.01 NAME
 B Regular type cross reference
.03 START TIME
 AC Regular type cross reference
 AUSR MUMPS
.01 COMMUNITIES
 B Regular type cross reference
.01 COMMUNITIES
 B Regular type cross reference
.01 COMMUNITIES
 B Regular type cross reference

90376.01 BGP 06 CMS INDICATORS

.01 INDICATOR NUMBER
 B Regular type cross reference
.02 ORDER ON SELECTION/REPORT
 AO Regular type cross reference
.03 TITLE
 C Regular type cross reference

90376.02 BGP 06 ELDER INDIVIDUAL INDICATORS

.01 INDICATOR
 B Regular type cross reference
.04 INDICATOR ID
 C Regular type cross reference
1204 ORDER
 ABC MUMPS
 AB MUMPS
 OR Regular type cross reference
 AP MUMPS

90376.03 BGP 06 ELDER DATA CURRENT

.01 BEGINNING DATE
 B Regular type cross reference

90376.04 BGP 06 ELDER DATA PREVIOUS

.01 BEGINNING DATE
 B Regular type cross reference

90376.05 BGP 06 ELDER DATA BASELINE

.01 BEGINNING DATE
 B Regular type cross reference

90376.06 BGP 06 ELDER INDICATORS

.01 INDICATOR NUMBER
 B Regular type cross reference
.02 ORDER ON SELECTION/REPORT
 AO Regular type cross reference
.03 TITLE
 C Regular type cross reference
.07 NATIONAL REPORT
 GPRA Regular type cross reference
.11 LOCAL REPORT
 LOCAL Regular type cross reference
.01 DENOMINATOR DEFINITIONS
 B Regular type cross reference
.01 NUMERATOR DEFINITIONS
 B Regular type cross reference
.01 COMMUNITIES
 B Regular type cross reference
.01 PATIENTS EISS
 B Regular type cross reference
.01 COMMUNITIES
 B Regular type cross reference
.01 PATIENTS EISS
 B Regular type cross reference

.01 COMMUNITIES
B Regular type cross reference
.01 PATIENTS EISS
B Regular type cross reference

4.4 Table File

NOTE: The following 16 files have over 1,000 fields per file. These files are used to hold calculated data values for the GPRA reports. The table files are not listed here to conserve space in this document. If you want a listing use the ^XBFLD or FM data dictionary listing utility to list the data dictionary in detail

File: 90240.01 BGP IHS GPRA DATA 02

Global: ^BGPD(

FILE: 90243 BGP 03 DATA CURRENT

GLOBAL: ^BGPGPDC(

FILE: 90243.01 BGP 03 DATA PREVIOUS

GLOBAL: ^BGPGPDP(

FILE: 90243.02 BGP 03 DATA BASELINE

GLOBAL: ^BGPGPDB(

FILE: 90244.03 BGP 04 DATA CURRENT

GLOBAL: ^BGPGPDCF(

FILE: 90244.04 BGP 04 DATA PREVIOUS

GLOBAL: ^BGPGPDPF(

FILE: 90244.05 BGP 04 DATA BASELINE

GLOBAL: ^BGPHEDBF(

FILE: 90370.03 BGP 04 HEDIS DATA CURRENT

GLOBAL: ^BGPHEDCF(

FILE: 90370.04 BGP 04 HEDIS DATA PREVIOUS

GLOBAL: ^BGPHEDPF(

FILE: 90370.05 BGP 04 HEDIS DATA BASELINE

GLOBAL: ^BGPHEDBV(

FILE: 90371.03 BGP 05 DATA CURRENT

GLOBAL: ^BGPGPDCV(

FILE: 90371.04 BGP 05 DATA PREVIOUS

GLOBAL: ^BGPGPDPV(

FILE: 90371.05 BGP 05 DATA BASELINE

GLOBAL: ^BGPHEDBV(

FILE: 90372.03 BGP 05 HEDIS DATA CURRENT

GLOBAL: ^BGPHEDCV(

FILE: 90372.04 BGP 05 HEDIS DATA PREVIOUS

GLOBAL: ^BGPHEDPV(

FILE: 90372.05 BGP 05 HEDIS DATA BASELINE

GLOBAL: ^BGPHEDBV(

FILE: 90373.03 BGP 05 ELDER DATA CURRENT

GLOBAL: ^BGPELDCV(

FILE: 90373.04 BGP 05 ELDER DATA PREVIOUS

GLOBAL: ^BGPELDPV(

FILE: 90373.05 BGP 05 ELDER DATA BASELINE

GLOBAL: ^BGPELDBV(

FILE: 90374.03 BGP 06 DATA CURRENT

GLOBAL: ^BGPGPDCV(

FILE: 90374.04 BGP 06 DATA PREVIOUS

GLOBAL: ^BGPGPDPV(

FILE: 90374.05 BGP 06 DATA BASELINE

GLOBAL: ^BGPHEDBV(

FILE: 90375.03 BGP 06 HEDIS DATA CURRENT

GLOBAL: ^BGPHEDCV(

FILE: 90375.04 BGP 06 HEDIS DATA PREVIOUS

GLOBAL: ^BGPHEDPV(

FILE: 90375.05 BGP 06 HEDIS DATA BASELINE

GLOBAL: ^BGPHEDBV(

FILE: 90376.03 BGP 06 ELDER DATA CURRENT

GLOBAL: ^BGPELDCV(

FILE: 90376.04 BGP 06 ELDER DATA PREVIOUS

GLOBAL: ^BGPELDPV(

FILE: 90376.05 BGP 06 ELDER DATA BASELINE

GLOBAL: ^BGPELDBV(

File: 90241.01 BGP CONTROL FILE
Global: ^BGPCTRL(

Field #	Field Name	Subscript	Piece	Type
.01	FY	D0,0	1	F
1101	Clinics for Clinical User (90241.0111)			
.01	Clinic Code	D0,11,D1,0	1	F
.019	Clinic Name	Computed		
1201	Clinic Codes for 2ND Visit (90241.0112)			
.01	Clinic Codes for 2ND Visit	D0,12,D1,0	1	F

FILE: 90241.02 BGP SITE PARAMETERS
GLOBAL: ^BGPSITE(

Field #	Field Name	Subscript	Piece	Type
.01	LOCATION	D0,0	1	P
.02	HOME LOCATION	"	2	P
.03	AREA EXPORT?	"	3	S
.04	EISS EXPORT?	"	4	S

FILE: 90242 BGP INDIVIDUAL INDICATOR DEFINITIONS
GLOBAL: ^BGPINDC(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.02	INDICATOR SECOND LEVEL	D0,0	2	F
.03	INDICATOR THIRD LEVEL	"	3	F
.04	INDICATOR ID	"	4	F
.05	GPRA	"	5	S
.06	AREA APP	"	6	S
.07	EISS	"	7	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.11	EISS DENOM PIECE	"	11	N
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS NUMERATOR PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.21	LOCAL DENOM TITLE 3	"	21	F
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1100	DENOMINATOR (90242.11)			
.01	DENOMINATOR	D0,11,D1,0	1	W
1101	GPRA DENOM LINE TITLE	D0,11	1	F
1102	GPRA DENOM TITLE 2	"	2	F
1103	GPRA DENOM TITLE 3	"	3	F
1104	AREA DENOM TITLE 1	"	4	F
1105	AREA DENOM TITLE 2	"	5	F
1106	AREA DENOM TITLE 3	"	6	F
1200	NUMERATRO (90242.12)			
.01	NUMERATRO	D0,12,D1,0	1	W
1300	LOGIC (90242.13)			
.01	LOGIC	D0,13,D1,0	1	W

FILE: 90242.01 BGP INDICATORS
GLOBAL: ^BGPIND(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	FY	D0,0	1	F
.02	INDICATOR NUMBER	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	GPRA INDICATOR	"	7	S
.08	AREA INDICATOR	"	8	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	LOGIC USED (90242.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	DESCRIPTION (90242.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	LOCAL REPORT PRINT HEADER (90242.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90242.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
5400	GPRA REPORT PRINT HEADER			

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
	(90242.0154)			
.01	GPRA REPORT PRINT HEADER	D0,54,D1,0	1	W
5500	AREA DIRECTORS REPORT HEADER (90242.0155)			
.01	AREA DIRECTORS REPORT HEADER	D0,55,D1,0	1	W
7100	REPORT LIST 1ST PAGE (90242.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90242.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90244.01 BGP 04 INDICATORS
GLOBAL: ^BGPINDF(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	GPRA INDICATOR	"	7	S
.08	AREA INDICATOR	"	8	S
.09	HEDIS INDICATOR	"	9	S
.11	LOCAL REPORT	"	11	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90244.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90244.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90244.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90244.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90244.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
1	DENOM DESCRIPTION (90244.01611)			

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	DENOM DESCRIPTION	D0,61,D1,1,D2, 0	1	W
6200	NUMERATOR DEFINITIONS (90244.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
1	NUMERATOR DESCRIPTION (90244.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D2, 0	1	W
7100	REPORT LIST 1ST PAGE (90244.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90244.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90244.02 BGP 04 INDIVIDUAL INDICATORS
GLOBAL: ^BGPINDFC(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.001	NUMBER			N
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.02	INDICATOR SECOND LEVEL	D0,0	2	F
.03	INDICATOR THIRD LEVEL	"	3	F
.04	INDICATOR ID	"	4	F
.05	GPRA	"	5	S
.06	AREA APP	"	6	S
.07	EISS	"	7	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.11	EISS DENOM PIECE	"	11	N
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS NUMERATOR PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1101	GPRA DENOM LINE TITLE	D0,11	1	F

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
1102	GPRA DENOM TITLE 2	"	2	F
1103	GPRA DENOM TITLE 3	"	3	F
1104	AREA DENOM TITLE 1	"	4	F
1105	AREA DENOM TITLE 2	"	5	F
1106	AREA DENOM TITLE 3	"	6	F
1201	HEDIS?	D0,12	1	S
1202	WOMEN'S HEALTH	"	2	S
1203	PREVENTION RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	NATIONAL PREVIOUS	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N

FILE: 90370.01 BGP 04 HEDIS INDICATORS**GLOBAL: ^BGPHEIF(**

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.09	HEDIS INDICATOR	"	9	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90244.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90244.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90244.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE			

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
	(90244.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90244.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
1	DENOM DESCRIPTION (90244.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D2, 0	1	W
6200	NUMERATOR DEFINITIONS (90244.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
1	NUMERATOR DESCRIPTION (90244.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D2, 0	1	W
7100	REPORT LIST 1ST PAGE (90244.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90244.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90370.02 BGP 04 HEDIS INDIVIDUAL INDICATORS
GLOBAL: ^BGPHEIIIF(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.001	NUMBER			N
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.02	INDICATOR SECOND LEVEL	D0,0	2	F
.03	INDICATOR THIRD LEVEL	"	3	F
.04	INDICATOR ID	"	4	F
.07	EISS	"	7	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.11	EISS DENOM PIECE	"	11	N
.13	CARDIO INDICATOR?	"	13	S
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1101	GPRA DENOM LINE TITLE	D0,11	1	F
1102	GPRA DENOM TITLE 2	"	2	F
1103	GPRA DENOM TITLE 3	"	3	F
1104	AREA DENOM TITLE 1	"	4	F
1105	AREA DENOM TITLE 2	"	5	F
1106	AREA DENOM TITLE 3	"	6	F
1201	HEDIS?	D0,12	1	S
1202	WOMEN'S HEALTH	"	2	S
1203	PREVENTION RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	NATIONAL PREVIOUS	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N

FILE: 90371.01 BGP 05 INDICATORS**GLOBAL: ^BGPINDV(**

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	NATIONAL REPORT	"	7	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90371.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION			

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
	(90371.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90371.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90371.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90371.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	DENOM DESCRIPTION (90371.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D2,0	1	W
6200	NUMERATOR DEFINITIONS (90371.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	NUMERATOR DESCRIPTION (90371.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D2,0	1	W
7100	REPORT LIST 1ST PAGE (90371.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90371.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90371.02 BGP 05 INDIVIDUAL INDICATORS
GLOBAL: ^BGPINDVC(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.04	INDICATOR ID	D0,0	4	F
.05	NATIONAL	"	5	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
.22	ADD BANG	"	22	S
.23	NATIONAL LIST	"	23	S
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1202	WOMEN'S HEALTH	D0,12	2	S
1203	ELDER RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	GOAL 2010	"	3	F
1404	TITLE	"	4	F

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N
1407	TITLE 2	"	7	F

FILE: 90371.06 BGP 05 NATIONAL PATIENT LISTS
GLOBAL: ^BGPVNPL(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	POINTER	D0,0	1	P
.02	INDIVIDUAL INDICATOR ID	"	2	F
.03	LIST TITLE	"	3	F
1101	LIST TITLE FOR REPORT (90371.61101)			
.01	LIST TITLE FOR REPORT	D0,11,D1,0	1	W
1201	SET LIST VALUE	D0,12		K

FILE: 90371.08 BGP 05 TAXONOMIES
GLOBAL: ^BGPTAXV(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	TAXONOMY NAME	D0,0	1	F
.02	TYPE	"	2	S
.03	LAST CHECKED/EXPORTED	"	3	D
.04	EDITABLE	"	4	S
1100	DESCRIPTION (90371.0811)	"	5	W
1200	REPORTS (90371.0812)	"	1	
.01	REPORTS	"	1	S

FILE: 90372.01 BGP 05 HEDIS INDICATORS
GLOBAL: ^BGPHEIV(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.09	HEDIS INDICATOR	"	9	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90372.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90372.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90372.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90372.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90372.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
1	DENOM DESCRIPTION (90372.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D2,0	1	W
6200	NUMERATOR DEFINITIONS (90372.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
1	NUMERATOR DESCRIPTION (90372.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D2,0	1	W
7100	REPORT LIST 1ST PAGE (90372.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90372.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90372.02 BGP 05 HEDIS INDIVIDUAL INDICATORS
GLOBAL: ^BGPHEIIV(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.001	NUMBER			N
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.02	CORRESPONDING GPRA	D0,0	2	F
.03	CORRESPONDING GPRA ORDER	"	3	F
.04	INDICATOR ID	"	4	F
.07	EISS	"	7	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.11	EISS DENOM PIECE	"	11	N
.13	CARDIO INDICATOR?	"	13	S
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1101	GPRA DENOM LINE TITLE	D0,11	1	F
1102	GPRA DENOM TITLE 2	"	2	F
1103	GPRA DENOM TITLE 3	"	3	F
1104	AREA DENOM TITLE 1	"	4	F
1105	AREA DENOM TITLE 2	"	5	F
1106	AREA DENOM TITLE 3	"	6	F
1201	HEDIS?	D0,12	1	S
1202	WOMEN'S HEALTH	"	2	S
1203	PREVENTION RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	NATIONAL PREVIOUS	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N

FILE: 90373.01 BGP 05 CMS INDICATORS

GLOBAL: ^BGPCMSV(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.06	ORDER TO DISPLAY	"	6	N
1	M CODE TO EXECUTE	D0,1		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
5100	INDICATOR DESCRIPTION (90372.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	LIST 2 HEADER (90372.0152)			
.01	LIST 2 HEADER	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90372.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
5400	LIST 1 HEADER (90372.0154)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
5500	LIST 1 TITLE (90372.0155)			
.01	LIST 1 TITLE	D0,72,D1,0	1	W

FILE: 90372.08 BGP 05 GUI REPORT OUTPUT**GLOBAL: ^BGPGUI(**

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	NAME	D0,0	1	F
.02	USER WHO REQUESTED REPORT	"	2	P
.03	START TIME	"	3	D
.04	END TIME	"	4	D
.05	TYPE OF REPORT	"	5	F
.06	STATUS	"	6	S
.07	TYPE OF OUTPUT	"	7	S
1100	PRINT FORMATTED OUTPUT			
.01	OUTPUT	D0,11,D1,0	1	W
1200	DELIMITED FORMATTED OUTPUT	D0,12,d1,0	1	W

FILE: 90373.06 BGP 05 ELDER INDICATORS**GLOBAL: ^BGPELIV(**

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	NATIONAL REPORT	"	7	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90371.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90371.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90371.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90371.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90371.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	DENOM DESCRIPTION (90371.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D2,0	1	W
6200	NUMERATOR DEFINITIONS (90371.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	NUMERATOR DESCRIPTION (90371.01621)			

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D2,0	1	W
7100	REPORT LIST 1ST PAGE (90371.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90371.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1, 0	1	W

FILE: 90373.02 BGP 05 ELDER INDIVIDUAL INDICATORS
GLOBAL: ^BGPELIIC(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.04	INDICATOR ID	D0,0	4	F
.05	NATIONAL	"	5	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
.22	ADD BANG	"	22	S
.23	NATIONAL LIST	"	23	S
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1202	WOMEN'S HEALTH	D0,12	2	S
1203	ELDER RELATED	"	3	S

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	GOAL 2010	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N
1407	TITLE 2	"	7	F

FILE: 90374.01 BGP 06 INDICATORS
GLOBAL: ^BGPINDS(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	NATIONAL REPORT	"	7	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90371.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90371.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90371.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
5300	REPORT PRINT 2ND PAGE (90371.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90371.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	DENOM DESCRIPTION (90371.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D2,0	1	W
6200	NUMERATOR DEFINITIONS (90371.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	NUMERATOR DESCRIPTION (90371.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D2,0	1	W
7100	REPORT LIST 1ST PAGE (90371.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90371.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90374.02 BGP 06 INDIVIDUAL INDICATORS
GLOBAL: ^BGPINDSC(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR	D0,0	1	P
.04	INDICATOR ID	D0,0	4	F
.05	NATIONAL	"	5	S
.06	GRANT1 PIECE	D0,0	6	N
.07	GRANT2 PIECE	D0,0	7	N
.08	DENOM FIELD	"	8	N

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
	NUMBER			
.09	NUMERATOR FIELD NUMBER	"	9	N
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
.22	ADD BANG	"	22	S
.23	NATIONAL LIST	"	23	S
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
3	M CODE FOR COMP PL	D0,3		K
1201	COLUMN TEXT FOR CMP PL	D0,12	1	F
1202	WOMEN'S HEALTH	D0,12	2	S
1203	ELDER RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	GOAL 2010	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N
1407	TITLE 2	"	7	F

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
1408	GOAL 06	"	8	F
1409	SECOND LINE GOAL	"	9	F
1410	NATIONAL CURRENT 2	"	10	F
1411	SECOND LINE 2010 GOAL	"	11	F

FILE: 90375.06 BGP 06 NATIONAL PATIENT LISTS
GLOBAL: ^BGPSNPL(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	POINTER	D0,0	1	P
.02	INDIVIDUAL INDICATOR ID	"	2	F
.03	LIST TITLE	"	3	F
1101	LIST TITLE FOR REPORT (90371.61101)			
.01	LIST TITLE FOR REPORT	D0,11,D1,0	1	W
1201	SET LIST VALUE	D0,12		K

FILE: 90375.08 BGP 06 TAXONOMIES
GLOBAL: ^BGPTAXS(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	TAXONOMY NAME	D0,0	1	F
.02	TYPE	"	2	S
.03	LAST CHECKED/EXPORTED	"	3	D
.04	EDITABLE	"	4	S
1100	DESCRIPTION (90371.0811)	"	5	W
1200	REPORTS (90371.0812)	"	1	
.01	REPORTS	"	1	S

FILE: 90375.01 BGP 05 HEDIS INDICATORS
GLOBAL: ^BGPHEIS(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.09	HEDIS INDICATOR	"	9	S

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90372.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90372.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90372.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90372.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90372.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
1	DENOM DESCRIPTION (90372.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D2,0	1	W
6200	NUMERATOR DEFINITIONS (90372.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
1	NUMERATOR DESCRIPTION (90372.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D2,0	1	W
7100	REPORT LIST 1ST PAGE (90372.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90372.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90375.02 BGP 06 HEDIS INDIVIDUAL INDICATORS
GLOBAL: ^BGPHEIIS(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.001	NUMBER			N
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.02	CORRESPONDING GPRA	D0,0	2	F
.03	CORRESPONDING GPRA ORDER	"	3	F
.04	INDICATOR ID	"	4	F
.07	EISS	"	7	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.11	EISS DENOM PIECE	"	11	N
.13	CARDIO INDICATOR?	"	13	S
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1101	GPRA DENOM LINE TITLE	D0,11	1	F
1102	GPRA DENOM TITLE 2	"	2	F
1103	GPRA DENOM TITLE 3	"	3	F
1104	AREA DENOM TITLE 1	"	4	F
1105	AREA DENOM TITLE 2	"	5	F
1106	AREA DENOM TITLE 3	"	6	F
1201	HEDIS?	D0,12	1	S
1202	WOMEN'S HEALTH	"	2	S
1203	PREVENTION RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	NATIONAL PREVIOUS	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N

FILE: 90376.01 BGP 06 CMS INDICATORS
GLOBAL: ^BGPCMSS(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.06	ORDER TO DISPLAY	"	6	N
1	M CODE TO EXECUTE	D0,1		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
5100	INDICATOR DESCRIPTION (90372.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	LIST 2 HEADER (90372.0152)			
.01	LIST 2 HEADER	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90372.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
5400	LIST 1 HEADER (90372.0154)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
5500	LIST 1 TITLE (90372.0155)			
.01	LIST 1 TITLE	D0,72,D1,0	1	W

FILE: 90375.08 BGP 06 GUI REPORT OUTPUT
GLOBAL: ^BGP_GUI(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	NAME	D0,0	1	F
.02	USER WHO REQUESTED REPORT	"	2	P
.03	START TIME	"	3	D
.04	END TIME	"	4	D
.05	TYPE OF REPORT	"	5	F
.06	STATUS	"	6	S
.07	TYPE OF OUTPUT	"	7	S
1100	PRINT FORMATTED OUTPUT			
.01	OUTPUT	D0,11,D1,0	1	W
1200	DELIMITED FORMATTED OUTPUT	D0,12,d1,0	1	W

FILE: 90376.06 BGP 06 ELDER INDICATORS
GLOBAL: ^BGPELIS(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	NATIONAL REPORT	"	7	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90371.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90371.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90371.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90371.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90371.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	DENOM DESCRIPTION (90371.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D2,0	1	W
6200	NUMERATOR DEFINITIONS (90371.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	NUMERATOR DESCRIPTION (90371.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D2,0	1	W
7100	REPORT LIST 1ST PAGE (90371.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90371.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90376.02 BGP 06 ELDER INDIVIDUAL INDICATORS
GLOBAL: ^BGPELIIS(

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.04	INDICATOR ID	D0,0	4	F
.05	NATIONAL	"	5	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
.22	ADD BANG	"	22	S
.23	NATIONAL LIST	"	23	S
1	M CODE NUMERATOR	D0,1		K

FIELD #	FIELD NAME	SUBSCRIPT	PIECE	TYPE
2	M CODE DENOMINATOR	D0,2		K
1202	WOMEN'S HEALTH	D0,12	2	S
1203	ELDER RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	GOAL 2010	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N
1407	TITLE 2	"	7	F

5.0 Archiving and Purging

There is no archiving and purging in this package.

6.0 External Relations

6.1 External Calls

This package calls the following documented entry points:

ROUTINE	IS INVOKED BY
^%DT	BGP3UTL1,BGP4UTL1,BGP5UTL1,BGP6UTL1, dd90240.01 dd90240.0197, dd90243, dd90243.01, dd90243.02 dd90244.03, dd90244.04, dd90244.05, dd90370.03 dd90370.04, dd90370.05, dd90371.03, dd90371.04 dd90371.05, dd90371.08, dd90372.03, dd90372.04
^%DT	dd90372.05, dd90372.08, dd90373.03, dd90373.04 dd90373.05, dd90374.03, dd90374.04, dd90374.05 dd90374.08, dd90375.03, dd90375.04, dd90375.05 dd90375.08, dd90376.03, dd90376.04, dd90376.05
COMMA^%DTC	BGP3DP,BGP3DP1,BGP3DP2,BGP3DP3,BGP3DP4,BGP3DP41,BGP 3DP5 BGP3DP50,BGP3DP51,BGP3DP6,BGP3DP60,BGP3DP7,BGP3DPH BGP3PDL,BGP3PDL1,BGP3PDL2,BGP3PDL3,BGP3PDL4,BGP3PDL5 BGP3PDL6,BGP3PDL7,BGP3PDL8,BGP3PDL9,BGP3PDM1,BGP3PD M2 BGP4DHEP,BGP4DP,BGP4DP1 BGP4DP11,BGP4DP2,BGP4DP3,BGP4DP4,BGP4DP5,BGP4DP50 BGP4DP51,BGP4DP6,BGP4DP60,BGP4DPH,BGP4HEL,BGP4HEL1 BGP4HEP1,BGP4PDL,BGP4PDL0,BGP4PDL1,BGP4PDL2,BGP4PDL 3 BGP4PDL4,BGP4PDL5,BGP4PDL7,BGP4PDL8,BGP4PDL9,BGP5CP BGP5DELP,BGP5DHEP,BGP5DP BGP5DP1,BGP5DP11,BGP5DP2,BGP5DP3,BGP5DP4,BGP5DP5 BGP5DP50,BGP5DP51,BGP5DP6,BGP5DP60,BGP5DP7,BGP5DPH BGP5EL10,BGP5EL11,BGP5EL12,BGP5EL13,BGP5EL14,BGP5EL15 BGP5EL16,BGP5EL17,BGP5EL18,BGP5EL19,BGP5EL1A,BGP5EL1B BGP5EL1C,BGP5EL1D,BGP5EL1E BGP5EL1F,BGP5EL1G,BGP5EL1H,BGP5EL1I,BGP5EL1J,BGP5EL1K BGP5EL1L,BGP5ELL,BGP5ELL1,BGP5ELL2,BGP5ELL3,BGP5ELL4 BGP5ELL5,BGP5ELL6,BGP5ELL7,BGP5ELL8,BGP5ELL9,BGP5ELLA BGP5ELLB,BGP5ELLC,BGP5ELLD,BGP5ELLE,BGP5ELLF,BGP5ELL G BGP5ELLH,BGP5ELLI BGP5ELLJ,BGP5ELLK,BGP5ELLL,BGP5ELLM,BGP5ELLN,BGP5ELL Q BGP5ELP1,BGP5HEL,BGP5HEL1,BGP5HEP1,BGP5PDL,BGP5PDL0 BGP5PDL1,BGP5PDL2,BGP5PDL3,BGP5PDL4,BGP5PDL5,BGP5PDL 7 BGP5PDL8,BGP5PDL9,BGP6CP,BGP6DELP,BGP6DHEP,BGP6DP BGP6DP1,BGP6DP11,BGP6DP12 BGP6DP2,BGP6DP3,BGP6DP4,BGP6DP5,BGP6DP50,BGP6DP51 BGP6DP6,BGP6DP60,BGP6DP7,BGP6DPH,BGP6EL10,BGP6EL11 BGP6EL12,BGP6EL13,BGP6EL14,BGP6EL15,BGP6EL16,BGP6EL17

BGP6EL18,BGP6EL19,BGP6EL1A,BGP6EL1B,BGP6EL1C,BGP6EL1D
BGP6EL1E,BGP6EL1F,BGP6EL1G
BGP6EL1H,BGP6EL1I,BGP6EL1J,BGP6EL1K,BGP6EL1L,BGP6EL1M
BGP6ELL,BGP6ELL1,BGP6ELL2,BGP6ELL3,BGP6ELL4,BGP6ELL5
BGP6ELL6,BGP6ELL7,BGP6ELL8,BGP6ELL9,BGP6ELLA,BGP6ELLB
BGP6ELLC,BGP6ELLD,BGP6ELLE,BGP6ELLF,BGP6ELLG,BGP6ELLH
BGP6ELLI,BGP6ELLJ
BGP6ELLK,BGP6ELLL,BGP6ELLM,BGP6ELLN,BGP6ELLO,BGP6ELLQ
BGP6ELP1,BGP6HEL,BGP6HEL1,BGP6HEP1,BGP6PDL,BGP6PDL0
BGP6PDL1,BGP6PDL2,BGP6PDL3,BGP6PDL4,BGP6PDL5,BGP6PDL7
BGP6PDL8,BGP6PDL9,BGP6PDLA,BGPAP,BGPAP1,BGPAP12,BGPAP13
BGPAP14,BGPAP1B
BGPAP2,BGPAP22,BGPAP22A,BGPAP22B,BGPAP22C,BGPAP23
BGPAP24,BGPAP3,BGPAP30,BGPAP30A,BGPAP30B,BGPAP30C,BGPAP4
BGPAP5,BGPAP6,BGPAP7,BGPAP8,BGPAPA,BGPAPB,BGPAPC,BGPAPD
BGPAPH,BGPDP,BGPDP1,BGPDP12,BGPDP13,BGPDP14,BGPDP1B
BGPDP2,BGPDP22,BGPDP22A
BGPDP22B,BGPDP22C,BGPDP23,BGPDP24,BGPDP29,BGPDP3,BGPDP30
BGPDP30A,BGPDP30B,BGPDP30C,BGPDP31,BGPDP32,BGPDP4,BGPDP5
BGPDP6,BGPDP7,BGPDP8,BGPDP8A,BGPDPARP,BGPDPB,BGPDPC
BGPDPD,BGPDPH
BGP6GRA,BGP6GU
BGP3DA,BGP3DAR,BGP3DARA,BGP3DARH,BGP3DG,BGP3DL,BGP4DA
BGP4DAP,BGP4DAR,BGP4DARA,BGP4DG,BGP4DHE,BGP4DHEA,BGP4DL
BGP4DPP,BGP5CL,BGP5DAP,BGP5DAR,BGP5DEL,BGP5DELA,BGP5DGA
BGP5DGPU,BGP5DHE,BGP5DHEA,BGP5DL,BGP5DNG,BGP5DPP,BGP5NPL
BGP6AHW
BGP6CL,BGP6DAP,BGP6DAR,BGP6DCL,BGP6DEL,BGP6DELA,BGP6DGA
BGP6DGPU,BGP6DHE,BGP6DHEA,BGP6DL,BGP6DLT,BGP6DMT,BGP6DNG
BGP6DPP,BGP6NPL,BGP6DAREA,BGP6DARP,BGP6DFTA
BGP3TXCH,BGP3ULF,BGP4LSTF,BGP4TXCH,BGP4TXHE,BGP4ULFN
BGP50P1,BGP5CTXC,BGP5GUPL,BGP5LSTF,BGP5TXCH,BGP5TXC
BGP5TXEL,BGP5TXHE,BGP5ULF,BGP6CTXC,BGP6DLT,BGP6DMT
BGP6GUPL,BGP6LSTF,BGP6ULF,BGP6XTCH,BGP6XTCN,BGP6XTE
BGP6XTHE,BGP6DARP,BGP6DFTA

	BGPDUP
^%ZISC	BGP3DA,BGP3DAR,BGP3DARA,BGP3DARH,BGP3DG,BGP3DL,BGP3DP BGP3PARP,BGP3ULF,BGP4DA,BGP4DAP,BGP4DAR,BGP4DARA,BGP4DG BGP4DHE,BGP4DHEA,BGP4DHEP,BGP4DL,BGP4DP,BGP4DPP,BGP4PARP BGP4PHEP,BGP4ULF,BGP4UTL,BGP5CL,BGP5DAP,BGP5DAR,BGP5DEL BGP5DELA BGP5DELP,BGP5DGPA,BGP5DGPU,BGP5DHE,BGP5DHEA,BGP5DHEP BGP5DL,BGP5DNG,BGP5DP,BGP5DPP,BGP5GUPL,BGP5NPL,BGP5PARP BGP5PHEL,BGP5PHEP,BGP5ULF,BGP5UTL,BGP6AHW,BGP6CL,BGP6DAP BGP6DAR,BGP6DCHA,BGP6DCL,BGP6DEL,BGP6DELA,BGP6DELP BGP6DGPA,BGP6DGPU BGP6DHE,BGP6DHEA,BGP6DHEP,BGP6DL,BGP6DLT,BGP6DMT,BGP6DNG BGP6DP,BGP6DPP,BGP6GUPL,BGP6NPL,BGP6PARP,BGP6PHEL BGP6PHEP,BGP6ULF,BGP6UTL,BGPD,BGPDAREA,BGPDARP,BGPDFTA BGPDUP
\$\$LIST^%ZISH	BGP4LSTF,BGP5GLST,BGP5GU,BGP5LSTF,BGP6GLST,BGP6GU BGP6LSTF,BGPGU
\$\$OPEN^%ZISH	BGP3ULF,BGP4ULF,BGP4UTL,BGP5GUPL,BGP5ULF,BGP5UTL BGP6GUPL,BGP6ULF,BGP6UTL,BGPDUP
\$\$STATUS^%ZISH	BGP3ULF,BGP4ULF,BGP5GUPL,BGP5ULF,BGP6GUPL,BGP6ULF,BGP6DUP
^%ZTLOAD	BGP3DA,BGP3DAR,BGP3DARA,BGP3DARH,BGP3DG,BGP3DL,BGP4DA BGP4DAP,BGP4DAR,BGP4DARA,BGP4DG,BGP4DHE,BGP4DHEA,BGP4DL BGP4DPP,BGP5CL,BGP5DAP,BGP5DAR,BGP5DEL,BGP5DELA,BGP5DGPA BGP5DGPU,BGP5DHE,BGP5DHEA,BGP5DL,BGP5DNG,BGP5DPP BGP5GAEL,BGP5GAGP BGP5GAHE,BGP5GALL,BGP5GAPU,BGP5GCMS,BGP5GCOM,BGP5GELD BGP5GGPU,BGP5GHED,BGP5GLST,BGP5GNPL,BGP5GNTL,BGP5GPP BGP5NPL,BGP6AHW,BGP6CL,BGP6DAP,BGP6DAR,BGP6DCL,BGP6DEL BGP6DELA,BGP6DGPA,BGP6DGPU,BGP6DHE,BGP6DHEA,BGP6DL BGP6DLT,BGP6DMT,BGP6DNG
^%ZTLOAD	BGP6DPP,BGP6GACW,BGP6GAEL,BGP6GAGP,BGP6GAHE,BGP6GALL BGP6GAPU,BGP6GCMP,BGP6GCMS,BGP6GCOM,BGP6GELD,BGP6GGPU BGP6GHED,BGP6GLST,BGP6GLTX,BGP6GMTX,BGP6GNPL,BGP6GNTL BGP6GPP,BGP6NPL,BGPD,BGPDAREA,BGPDARP,BGPDFTA
^APCDCHKJ	BGP5C1,BGP6C1

\$\$START1^APCLDF BGP3D1,BGP3D10,BGP3D2,BGP3D3,BGP3D4,BGP3D5,BGP3D6,BGP3D7
 BGP3D8,BGP3D9,BGP4D1,BGP4D10,BGP4D2,BGP4D21,BGP4D3
 BGP4D31,BGP4D4,BGP4D41,BGP4D5,BGP4D6,BGP4D61,BGP4D7
 BGP4D71,BGP4D8,BGP4D9,BGP5C1,BGP5C11,BGP5C13,BGP5C2
 BGP5C3,BGP5D1,BGP5D10
 BGP5D2,BGP5D21,BGP5D22,BGP5D3,BGP5D31,BGP5D33,BGP5D34
 BGP5D35,BGP5D4,BGP5D41,BGP5D5,BGP5D6,BGP5D61,BGP5D7
 BGP5D71,BGP5D72,BGP5D8,BGP5D9,BGP5DHE1,BGP5EL2,BGP5EL3
 BGP5EL4,BGP6C1,BGP6C11,BGP6C13,BGP6C2,BGP6C3,BGP6D1
 BGP6D10,BGP6D2,BGP6D21
 BGP6D22,BGP6D23,BGP6D3,BGP6D31,BGP6D33,BGP6D34,BGP6D35
 BGP6D4,BGP6D41,BGP6D5,BGP6D51,BGP6D52,BGP6D6,BGP6D61
 BGP6D7,BGP6D71,BGP6D72,BGP6D8,BGP6D9,BGP6DCHW,BGP6DHE1
 BGP6EL2,BGP6EL3,BGP6EL4,BGPD1,BGPD12,BGPD13,BGPD14
 BGPD1B,BGPD2,BGPD22
 BGPD24,BGPD29,BGPD3,BGPD31,BGPD4,BGPD5,BGPD6,BGPD7,BGPD8A
 BGPDB,BGPDC,BGPDD
 ^APCLTAX0 |opt
 \$\$CLINIC^APCLV BGP3D1,BGP3D2,BGP3D3,BGP3D4,BGP3D6,BGP3D9,BGP4D1,BGP4D10
 BGP4D2,BGP4D21,BGP4D41,BGP4D7,BGP4D71,BGP5D1,BGP5D10
 BGP5D2,BGP5D21,BGP5D22,BGP5D41,BGP5D7,BGP5D71,BGP5DH
 BGP5EL2,BGP5EL4,BGP6D1,BGP6D10,BGP6D2,BGP6D21,BGP6D22
 BGP6D23,BGP6D41
 BGP6D51,BGP6D7,BGP6D71,BGP6DHE1,BGP6EL2,BGP6EL4,BGPD14
 BGPD2,BGPD22,BGPD3,BGPD8,BGPDD
 \$\$DSCHDATE^APCLV BGP5D22,BGP5D72,BGP5EL3,BGP6D22,BGP6D51,BGP6D72,BGP6EL3
 \$\$PRIMPOV^APCLV BGP3D2,BGP3D3,BGP3D8,BGP4D21,BGP4D31,BGP4D71,BGP5C1
 BGP5C2,BGP5C3,BGP5D21,BGP5D22,BGP5D31,BGP5D71,BGP5DH
 BGP6C1,BGP6C2,BGP6C3,BGP6D21,BGP6D22,BGP6D31,BGP6D51
 BGP6D71,BGP6DHE1,BGPD14,BGPD8,BGPDD
 \$\$PRIMPROV^APCLV BGP3D1,BGP3D2,BGP4D1,BGP4D10,BGP4D21,BGP4D71,BGP5D1
 BGP5D10,BGP5D21,BGP5D71,BGP5DHE1,BGP6D1,BGP6D10,BGP6D21
 BGP6D51,BGP6D71,BGP6DHE1,BGPD14,BGPD2,BGPDD
 \$\$ICD^ATXCHK BGP3D1,BGP3D6,BGP3D7,BGP3D8,BGP3D9,BGP4D1,BGP4D10
 BGP4D31,BGP4D5,BGP4D7,BGP4D8,BGP4D9,BGP5C1,BGP5C11
 BGP5C13,BGP5C2,BGP5C3,BGP5D1,BGP5D10,BGP5D22,BGP5D31
 BGP5D4,BGP5D5,BGP5D61,BGP5D7,BGP5D72,BGP5D8,BGP5D9
 BGP5DHE1,BGP5EL3,BGP5UTL1
 BGP6C1,BGP6C11,BGP6C13,BGP6C2,BGP6C3,BGP6D1,BGP6D10
 BGP6D22,BGP6D31,BGP6D4,BGP6D5,BGP6D51,BGP6D61,BGP6D7

BGP6D72,BGP6D8,BGP6D9,BGP6DHE1,BGP6DU,BGP6EL3,BGP6U
TL1
BGPD2,BGPDU
BULL^ATXSTX2 BGP6AXA,BGP6AXB,BGP6AXC,BGP6AXD,BGP6AXE,BGP6AXF,BGP
6AXG
BGP6AXH,BGP6AXI,BGP6AXJ,BGP6AXK,BGP6AXL,BGP6AXM,BGP6
AXN
BGP6AXO,BGP6AXP,BGP6AXQ,BGP6AXR,BGP6AXS,BGP6AXT,BG
P6AXU
BGP6AXV,BGP6AXW,BGP6AXX,BGP6AXY,BGP6AXZ,BGP6BXA,BG
P6BXB
BGP6BXC
BGP6BXD,BGP6BXE,BGP6BXF,BGP6BXG,BGP6BXH,BGP6BXI,BGP
6BXJ
BGP6BXK,BGP6BXL,BGP6BXM,BGP6BXN,BGP6BXO,BGP6BXP,BG
P6BXQ
BGP6BXR,BGP6BXS,BGP6BXT,BGP6BXU,BGP6BXV,BGP6BXW,BG
P6BXX
BGP6BXY,BGP6BXZ,BGP6EXA,BGP6EXB,BGP6EXC,BGP6EXD,BGP
6EXE
BGP6EXF
BGP6EXG,BGP6EXH,BGP6EXI,BGP6EXJ,BGP6EXK,BGP6EXL,BGP6
EXM
BGP6EXN,BGP6EXO,BGP6EXP,BGP6EXQ,BGP6EXR,BGP6EXS,BG
P6EXT
BGP6EXU,BGP6EXV,BGP6EXW,BGP6EXX,BGP6EXY,BGP6EXZ,BG
P6FXA
BGP6FXB,BGP6FXC,BGP6FXD,BGP6FXE,BGP6FXF,BGP6FXG,BGP
6HXA
BGP6HXB
BGP6HXC,BGP6HxD,BGP6HxE,BGP6HxF,BGP6HXG,BGP6HXH,BG
P6HXI
BGP6HXJ,BGP6HXK,BGP6HXL,BGP6HXM,BGP6HXN,BGP6HXO,BG
P6HXP
BGP6HXQ,BGP6HXR,BGP6HXS,BGP6HXT,BGP6HXU,BGP6HXV,BG
P6HXW
BGP6HXX,BGP6HXY,BGP6HXZ,BGPTXA,BGPTXB,BGPTXC,BGPTX
D
BGPTXE,BGPTXF
BGPTXG,BGPTXH,BGPTXI,BGPTXJ,BGPTXK,BGPTXL,BGPTXM,BG
PTXN
BGPTXO,BGPTXP,BGPTXQ,BGPTXR,BGPTXS,BGPTXT,BGPTXU,B
GPTXV
BGPTXW,BGPTXX,BGPTXY,BGPTXZ,BGPUTXA,BGPUTXB,BGPUTX
C
BGPUTXD,BGPUTXE,BGPUTXF,BGPUTXG,BGPUTXH,BGPUTXI,BG
PUTXJ
BGPUTXK,BGPUTXL
BGPUTXM,BGPUTXN,BGPUTXO,BGPUTXP,BGPUTXQ,BGPUTXR,B
GPUTXS
BGPUTXT,BGPUTXU,BGPUTXV,BGPUTXW,BGPUTXX,BGPUTXY,B
GPUTXZ
BGPVTXA,BGPVTXB,BGPVTXC,BGPVTXD,BGPVTXE,BGPVTXF,BG
PVTXG

BGPVTXH,BGPVTXI,BGPVTXJ,BGPVTXK,BGPVTXL,BGPVTXM,BGP
VTXN
BGPVTXO
BGPVTXP,BGPVTXQ,BGPVTXR,BGPVTXS,BGPVTXT,BGPVTXU,BG
PVTXV
BGPVTXW,BGPVTXX,BGPVTXY,BGPVTXZ,BGPWTXA,BGPWTXB,BG
PWTXC
BGPWTXD,BGPWTXE,BGPWTXF,BGPWTXG,BGPWTXH,BGPWTXI,
BGPWTXJ
BGPWTXK,BGPWTXL,BGPWTXM,BGPWTXN,BGPWTXO,BGPWTXP,
BGPWTXQ
BGPWTXR
BGPWTXS,BGPWTXT,BGPWTXU
KILL^ATXSTX2 BGP6AXA,BGP6AXB,BGP6AXC,BGP6AXD,BGP6AXE,BGP6AXF,BGP
6AXG
BGP6AXH,BGP6AXI,BGP6AXJ,BGP6AXK,BGP6AXL,BGP6AXM,BGP6
AXN
BGP6AXO,BGP6AXP,BGP6AXQ,BGP6AXR,BGP6AXS,BGP6AXT,BG
P6AXU
BGP6AXV,BGP6AXW,BGP6AXX,BGP6AXY,BGP6AXZ,BGP6BXA,BG
P6BXB
BGP6BXC
BGP6BXD,BGP6BXE,BGP6BXF,BGP6BXG,BGP6BXH,BGP6BXI,BGP
6BXJ
BGP6BXK,BGP6BXL,BGP6BXM,BGP6BXN,BGP6BXO,BGP6BXP,BG
P6BXQ
BGP6BXR,BGP6BXS,BGP6BXT,BGP6BXU,BGP6BXV,BGP6BXW,BG
P6BXX
BGP6BXY,BGP6BXZ,BGP6EXA,BGP6EXB,BGP6EXC,BGP6EXD,BGP
6EXE
BGP6EXF
BGP6EXG,BGP6EXH,BGP6EXI,BGP6EXJ,BGP6EXK,BGP6EXL,BGP6
EXM
BGP6EXN,BGP6EXO,BGP6EXP,BGP6EXQ,BGP6EXR,BGP6EXS,BG
P6EXT
BGP6EXU,BGP6EXV,BGP6EXW,BGP6EXX,BGP6EXY,BGP6EXZ,BG
P6FXA
BGP6FXB,BGP6FCX,BGP6FXD,BGP6FXE,BGP6FXF,BGP6FXG,BGP
6HXA
BGP6HXB
BGP6HXC,BGP6HxD,BGP6HxE,BGP6HXF,BGP6HXG,BGP6HXH,BG
P6HXI
BGP6HXJ,BGP6HXK,BGP6HXL,BGP6HXM,BGP6HXN,BGP6HXO,BG
P6HXP
BGP6HXQ,BGP6HXR,BGP6HXS,BGP6HXT,BGP6HXU,BGP6HXV,BG
P6HXW
BGP6HXX,BGP6HXY,BGP6HXZ,BGPTXA,BGPTXB,BGPTXC,BGPTX
D
BGPTXE,BGPTXF
KILL^ATXSTX2 BGPTXG,BGPTXH,BGPTXI,BGPTXJ,BGPTXK,BGPTXL,BGPTXM,BG
PTXN
BGPTXO,BGPTXP,BGPTXQ,BGPTXR,BGPTXS,BGPTXT,BGPTXU,BG
PTXV
BGPTXW,BGPTXX,BGPTXY,BGPTXZ,BGPUTXA,BGPUTXB,BGPUTX
C

BGPUTXD,BGPUTXE,BGPUTXF,BGPUTXG,BGPUTXH,BGPUTXI,BG
PUTXJ
BGPUTXK,BGPUTXL
BGPUTXM,BGPUTXN,BGPUTXO,BGPUTXP,BGPUTXQ,BGPUTXR,B
PUTXS
BGPUTXT,BGPUTXU,BGPUTXV,BGPUTXW,BGPUTXX,BGPUTXY,B
PUTXZ
BGPVTXA,BGPVTXB,BGPVTXC,BGPVTXD,BGPVTXE,BGPVTXF,BG
PVTXG
BGPVTXH,BGPVTXI,BGPVTXJ,BGPVTXK,BGPVTXL,BGPVTXM,BGP
VTXN
BGPVTXO
BGPVTXP,BGPVTXQ,BGPVTXR,BGPVTXS,BGPVTXT,BGPVTXU,BG
PVTXV
BGPVTXW,BGPVTXX,BGPVTXY,BGPVTXZ,BGPWTXA,BGPWTXB,B
GPWTXC
BGPWTXD,BGPWTXE,BGPWTXF,BGPWTXG,BGPWTXH,BGPWTXI,
BGPWTXJ
BGPWTXK,BGPWTXL,BGPWTXM,BGPWTXN,BGPWTXO,BGPWTXP,
BGPWTXQ
BGPWTXR
BGPWTXS,BGPWTXT,BGPWTXU
TAX^ATXSTX2
BGP6AXA,BGP6AXB,BGP6AXC,BGP6AXD,BGP6AXE,BGP6AXF,BGP
6AXG
BGP6AXH,BGP6AXI,BGP6AXJ,BGP6AXK,BGP6AXL,BGP6AXM,BGP6
AXN
BGP6AXO,BGP6AXP,BGP6AXQ,BGP6AXR,BGP6AXS,BGP6AXT,BG
P6AXU
BGP6AXV,BGP6AXW,BGP6AXX,BGP6AXY,BGP6AXZ,BGP6BXA,BG
P6BXB
BGP6BC
BGP6BXD,BGP6BXE,BGP6BXF,BGP6BXG,BGP6BXH,BGP6BXI,BGP
6BXJ
BGP6BXK,BGP6BXL,BGP6BXM,BGP6BXN,BGP6BXO,BGP6BXP,BG
P6BXQ
BGP6BXR,BGP6BXS,BGP6BXT,BGP6BXU,BGP6BXV,BGP6BXW,BG
P6BXX
BGP6BXY,BGP6BXZ,BGP6EXA,BGP6EXB,BGP6EXC,BGP6EXD,BGP
6EXE
BGP6EXF
BGP6EXG,BGP6EXH,BGP6EXI,BGP6EXJ,BGP6EXK,BGP6EXL,BGP6
EXM
BGP6EXN,BGP6EXO,BGP6EXP,BGP6EXQ,BGP6EXR,BGP6EXS,BG
P6EXT
BGP6EXU,BGP6EXV,BGP6EXW,BGP6EXX,BGP6EXY,BGP6EXZ,BG
P6FXA
BGP6FXB,BGP6FC,BGP6FXD,BGP6FXE,BGP6FXF,BGP6FXG,BGP
6HXA
BGP6HXB
BGP6HXC,BGP6HxD,BGP6HxE,BGP6HxF,BGP6HXG,BGP6HXH,BG
P6HXI
BGP6HXJ,BGP6HXK,BGP6HXL,BGP6HXM,BGP6HXN,BGP6HXO,BG
P6HXP
BGP6HXQ,BGP6HXR,BGP6HXS,BGP6HXT,BGP6HXU,BGP6HXV,BG
P6HXW

	BGP6HXX,BGP6HXY,BGP6HXZ,BGPTXA,BGPTXB,BGPTXC,BGPTX D BGPTXE,BGPTXF BGPTXG,BGPTXH,BGPTXI,BGPTXJ,BGPTXK,BGPTXL,BGPTXM,BG PTXN BGPTXO,BGPTXP,BGPTXQ,BGPTXR,BGPTXS,BGPTXT,BGPTXU,B GPTXV BGPTXW,BGPTXX,BGPTXY,BGPTXZ,BGPUTXA,BGPUTXB,BGPUTX C BGPUTXD,BGPUTXE,BGPUTXF,BGPUTXG,BGPUTXH,BGPUTXI,BG PUTXJ BGPUTXK,BGPUTXL BGPUTXM,BGPUTXN,BGPUTXO,BGPUTXP,BGPUTXQ,BGPUTXR,B GPUTXS BGPUTXT,BGPUTXU,BGPUTXV,BGPUTXW,BGPUTXX,BGPUTXY,B GPUTXZ BGPVTXA,BGPVTXB,BGPVTXC,BGPVTXD,BGPVTXE,BGPVTXF,BG PVTXG BGPVTXH,BGPVTXI,BGPVTXJ,BGPVTXK,BGPVTXL,BGPVTXM,BGP VTXN BGPVTXO BGPVTXP,BGPVTXQ,BGPVTXR,BGPVTXS,BGPVTXT,BGPVTXU,BG PVTXV BGPVTXW,BGPVTXX,BGPVTXY,BGPVTXZ,BGPWTXA,BGPWTXB,B GPWTXC BGPWTXD,BGPWTXE,BGPWTXF,BGPWTXG,BGPWTXH,BGPWTXI, BGPWTXJ BGPWTXK,BGPWTXL,BGPWTXM,BGPWTXN,BGPWTXO,BGPWTXP, BGPWTXQ BGPWTXR TAX^ATXSTX2 \$\$AGE^AUPNPAT \$\$BEN^AUPNPAT \$\$COMMRES^AUPNPA T \$\$DOB^AUPNPAT \$\$DOD^AUPNPAT	
	BGPWTXS,BGPWTXT,BGPWTXU BGP3D1,BGP3D3,BGP4D1,BGP4D10,BGP4D7,BGP5C1,BGP5C2,BG P5C3 BGP5D1,BGP5D10,BGP5D6,BGP5D7,BGP6C1,BGP6C2,BGP6C3,BG P6D1 BGP6D10,BGP6D23,BGP6D6,BGP6D7,BGPD1,BGPD2,BGPD22 BGP3D1,BGP4D1,BGP4D10,BGP5D1,BGP5D10,BGP6D1,BGP6D10 BGPD1 BGP6DCHW BGP3D1,BGP3D2,BGP3D3,BGP3D4,BGP3DCI,BGP4D1,BGP4D10 BGP4D2,BGP4D3,BGP4D4,BGP5C1,BGP5C11,BGP5C13,BGP5C2 BGP5D1,BGP5D10,BGP5D2,BGP5D22,BGP5D3,BGP5D32,BGP5D33 BGP5D34,BGP5D35,BGP5D4,BGP5D61,BGP5D71,BGP5D72,BGP5D 9 BGP5EL2,BGP5EL3,BGP6C1 BGP6C11,BGP6C13,BGP6C2,BGP6D1,BGP6D10,BGP6D2,BGP6D22 BGP6D3,BGP6D32,BGP6D33,BGP6D34,BGP6D35,BGP6D4,BGP6D6 1 BGP6D71,BGP6D72,BGP6D9,BGP6DHE1,BGP6EL2,BGP6EL3,BGPD 1 BGPD24 BGP3D1,BGP4D1,BGP4D10,BGP5D1,BGP5D10,BGP6D1,BGP6D10	

	BGPD1
\$\$HRN^AUPNPAT	BGP3DS,BGP3PDLS,BGP4DS,BGP4HES,BGP4HESL,BGP4PDLS,BG P5C1
	BGP5C2,BGP5C3,BGP5DS,BGP5ELS,BGP5ELSL,BGP5HES,BGP5H ESL
	BGP5NPLD,BGP5NPLP,BGP5PDLS,BGP6C1,BGP6C2,BGP6C3,BGP 6DCLD
	BGP6DCLP,BGP6DS,BGP6ELS,BGP6ELSL,BGP6HES,BGP6HESL
	BGP6NPLD,BGP6NPLP
	BGP6PDLS,BGPDS
\$\$TRIBE^AUPNPAT	BGP6DCHW
KILL^AUPNPAT	BGP3DA,BGP3DAR,BGP3DARA,BGP3DARH,BGP3DG,BGP3DL,BGP 4DA
	BGP4DAP,BGP4DAR,BGP4DARA,BGP4DG,BGP4DHE,BGP4DHEA,B GP4DL
	BGP4DPP,BGP5CL,BGP5DAP,BGP5DAR,BGP5DEL,BGP5DELA,BGP 5DGPA
	BGP5DGPU,BGP5DHE,BGP5DHEA,BGP5DL,BGP5DNG,BGP5DPP
	BGP5GAEL,BGP5GAGP
	BGP5GAHE,BGP5GALL,BGP5GAPU,BGP5GCMS,BGP5GCOM,BGP5 GELD
	BGP5GGPU,BGP5GHED,BGP5GLST,BGP5GNPL,BGP5GNTL,BGP5 GPP
	BGP5NPL,BGP6AHW,BGP6CL,BGP6DAP,BGP6DAR,BGP6DCHA,BG P6DCL
	BGP6DEL,BGP6DELA,BGP6DGPA,BGP6DGPU,BGP6DHE,BGP6DHE A
	BGP6DL,BGP6DNG,BGP6DPP
	BGP6GACW,BGP6GAEL,BGP6GAGP,BGP6GAHE,BGP6GALL,BGP6 GAPU
	BGP6GCMP,BGP6GCMS,BGP6GCOM,BGP6GELD,BGP6GGPU,BGP 6GHED
	BGP6GLST,BGP6GLTX,BGP6GMTX,BGP6GNPL,BGP6GNTL,BGP6G PP
	BGP6NPL,BGPD,BGPDAREA,BGPDARP,BGPDFTA
^BGP3BAN	opt
^BGP3D1	BGP3DA,BGP3DG,BGP3DL
\$\$LASTHF^BGP3D6	BGP3D8
\$\$PED^BGP3D6	BGP3D9
\$\$FIRSTIHD^BGP3D9	BGP3D1
\$\$V2IHD^BGP3D9	BGP3D1
^BGP3DA	opt
^BGP3DAR	opt
^BGP3DARA	opt
CALCIND^BGP3DCI	BGP3D1
^BGP3DG	opt
^BGP3DH	BGP3DP
AREACP^BGP3DH	BGP3PARP
^BGP3DL	opt
^BGP3DP	BGP3DA,BGP3DG,BGP3DL
PRINT1^BGP3DP	BGP3PARP
PI^BGP3DP1	BGP3DP3,BGP3DP4,BGP3DP5,BGP3DP6
I1AGE^BGP3DP41	BGP3DP4

I1AGEP^BGP3DP50	BGP3DP5,BGP3DP51
I1AGE^BGP3DP51	BGP3DP5
I1AGE^BGP3DP60	BGP3DP6
H1^BGP3DPH	BGP3DP1,BGP3DP4,BGP3DP5,BGP3DP6,BGP3DP7
H2^BGP3DPH	BGP3DP1,BGP3DP6,BGP3DP7
H3^BGP3DPH	BGP3DP1
HEADER^BGP3DPH	BGP3DH,BGP3DP,BGP3DP1,BGP3DP3,BGP3DP4,BGP3DP41,BGP3DP5 BGP3DP50,BGP3DP51,BGP3DP6,BGP3DP60,BGP3DP7
^BGP3DS	BGP3DP
EN^BGP3DSI	BGP3DL
EN^BGP3DSL	BGP3DL
PT^BGP3DSL	BGP3DA,BGP3DAR,BGP3DARA,BGP3DARH,BGP3DG,BGP3DL
RT^BGP3DSL	BGP3DL
PRINT^BGP3PARP	BGP3DAR,BGP3DARA,BGP3DARH
^BGP3PDL	BGP3DP,BGP3PARP
H1^BGP3PDL1	BGP3PDL3,BGP3PDL4,BGP3PDL6,BGP3PDL7
PI^BGP3PDL1	BGP3PDL3,BGP3PDL4,BGP3PDL6,BGP3PDL7
I1AGE^BGP3PDL5	BGP3PDL4
I1AGEP^BGP3PDL8	BGP3PDL7,BGP3PDL9
I1AGE^BGP3PDL9	BGP3PDL7
^BGP3PDLS	BGP3PDL
I1AGE^BGP3PDM1	BGP3PDL6
^BGP3TXCH	opt
TAXCHK^BGP3TXCH	BGP3DA,BGP3DG,BGP3DL
^BGP3ULF	opt
\$\$DATE^BGP3UTL	BGP3D10,BGP3D2,BGP3D3,BGP3D4,BGP3D5,BGP3D6,BGP3D7,BG P3D8 BGP3D9,BGP4D31
EISSGS^BGP3UTL	BGP3DA
GETIEN^BGP3UTL	BGP3ULF
GS^BGP3UTL	BGP3DA,BGP3DG
REPORT^BGP3UTL	BGP3DA,BGP3DG,BGP3DL
STMP^BGP3UTL	BGP3DCI BGP3D2,BGP3D3,BGP3D4,BGP3D6,BGP3D8
\$\$REFUSAL^BGP3UTL1	
^BGP4BAN	opt
^BGP4D1	BGP4DA,BGP4DAP,BGP4DG,BGP4DHE,BGP4DL
^BGP4D10	BGP4DPP
\$\$HDL^BGP4D2	BGP4D41,BGP4D9
\$\$HGBA1C^BGP4D2	BGP4DHE1
\$\$LDL^BGP4D2	BGP4D41,BGP4D9,BGP4DHE1
\$\$LIPID^BGP4D2	BGP4D41,BGP4D9
\$\$TRIG^BGP4D2	BGP4D41,BGP4D9
\$\$EYE^BGP4D21	BGP4DHE1
\$\$POSUR^BGP4D21	BGP4DHE1
\$\$PNEU^BGP4D31	BGP4D3
\$\$HYSTER^BGP4D4	BGP4D3
\$\$MEANBP^BGP4D41	BGP4D9
\$\$CHL^BGP4D5	BGP4DHE1
\$\$BMI^BGP4D6	BGP4D71

\$\$OB^BGP4D6	BGP4D71
\$\$OW^BGP4D6	BGP4D71
\$\$FOB^BGP4D61	BGP4DHE1
\$\$IBH^BGP4D61	BGP4DHE1
\$\$RECT^BGP4D61	BGP4DHE1
\$\$DENT^BGP4D7	BGP4D71
\$\$DX^BGP4D7	BGP4D71
\$\$TOBACCO^BGP4D7	BGP4D71
^BGP4DA	opt
^BGP4DAP	opt
^BGP4DAR	opt
^BGP4DARA	opt
CALCIND^BGP4DCHI	BGP4D1
CALCIND^BGP4DCI	BGP4D1,BGP4D10
^BGP4DG	opt
^BGP4DH	BGP4DP
AREACP^BGP4DH	BGP4PARP
^BGP4DHE	opt
^BGP4DHEA	opt
^BGP4DHEP	BGP4DHE
PRINT1^BGP4DHEP	BGP4PHEP
EN^BGP4DHSL	BGP4DHE
PT^BGP4DHSL	BGP4DHE,BGP4DHEA
RT^BGP4DHSL	BGP4DHE
^BGP4DL	opt
^BGP4DP	BGP4DA,BGP4DAP,BGP4DG,BGP4DL,BGP4DPP
PRINT1^BGP4DP	BGP4PARP
SETEXCEL^BGP4DP	BGP4PARP
PI^BGP4DP1	BGP4DP3,BGP4DP4,BGP4DP5,BGP4DP6
SETEXA^BGP4DP1	BGP4DP2,BGP4DP3,BGP4PDL1,BGP4PDL2,BGP4PDL3
SETN^BGP4DP1	BGP4DP2,BGP4DP3,BGP4PDL0,BGP4PDL1,BGP4PDL2,BGP4PDL3
I1AGE^BGP4DP11	BGP4DP1
I1AGEP^BGP4DP50	BGP4DP5,BGP4DP51
I1AGEP^BGP4DP60	BGP4DP6
H1^BGP4DPH	BGP4DP1,BGP4DP4,BGP4DP5,BGP4DP6,BGP4HEP1
H2^BGP4DPH	BGP4DP1,BGP4HEP1
H3^BGP4DPH	BGP4DP11
HEADER^BGP4DPH	BGP4DHEP,BGP4DP,BGP4DP1,BGP4DP11,BGP4DP3,BGP4DP4,BG P4DP5 BGP4DP50,BGP4DP51,BGP4DP6,BGP4DP60,BGP4DS,BGP4DSP BGP4HEP1,BGP4HES
^BGP4DPP	opt
^BGP4DS	BGP4DP
EN^BGP4DSI	BGP4DAP,BGP4DL,BGP4DPP
EN^BGP4DSL	BGP4DAP,BGP4DL,BGP4DPP
PT^BGP4DSL	BGP4DA,BGP4DAP,BGP4DAR,BGP4DARA,BGP4DG,BGP4DL,BGP4 DPP
RT^BGP4DSL	BGP4DAP,BGP4DL,BGP4DPP
^BGP4DSP	BGP4DP
^BGP4HEH	BGP4DHEP
AREACP^BGP4HEH	BGP4PHEP

^BGP4HEHH	BGP4HEL
^BGP4HEL	BGP4DHEP,BGP4PHEP
^BGP4HES	BGP4DHEP
^BGP4HESL	BGP4HEL
GETIEN^BGP4HUTL	BGP4ULF
GS^BGP4HUTL	BGP4DHE
REPORT^BGP4HUTL	BGP4DHE
STMP^BGP4HUTL	BGP4DCII
^BGP4LSTF	opt
PRINT^BGP4PARP	BGP4DAR,BGP4DARA
^BGP4PDL	BGP4DP,BGP4PARP
HEADER^BGP4PDL	BGP4PDLZ
H1^BGP4PDL1	BGP4PDL0,BGP4PDL3,BGP4PDL4,BGP4PDL7
PI^BGP4PDL1	BGP4PDL0,BGP4PDL3,BGP4PDL4,BGP4PDL7
I1AGEP^BGP4PDL5	BGP4PDL0
I1AGEP^BGP4PDL8	BGP4PDL7
I1AGE^BGP4PDL9	BGP4PDL1
^BGP4PDLH	BGP4PDL
^BGP4PDLS	BGP4PDL
PRINT^BGP4PHEP	BGP4DHEA
^BGP4TXCH	opt
TAXCHK^BGP4TXCH	BGP4DA,BGP4DAP,BGP4DG,BGP4DL,BGP4DPP
TAXCHK^BGP4TXHE	BGP4DHE
^BGP4ULF	opt
\$\$D^BGP4UTL	BGP4DA,BGP4DAR,BGP4DARA,BGP4DG
\$\$DATE^BGP4UTL	BGP4D2,BGP4D21,BGP4D3,BGP4D31,BGP4D4,BGP4D41,BGP4D5 BGP4D61,BGP4D7,BGP4D71,BGP4D9,BGP4DHE1,BGP4DP
\$\$LZERO^BGP4UTL	BGP4DAR,BGP4DARA
\$\$RZERO^BGP4UTL	BGP4DA,BGP4DAR,BGP4DARA,BGP4DG
EXCELGS^BGP4UTL	BGP4DA,BGP4DAR,BGP4DARA,BGP4DG
GETIEN^BGP4UTL	BGP4ULF
GS^BGP4UTL	BGP4DA,BGP4DG
JRNL^BGP4UTL	BGP4D1,BGP4D10
REPORT^BGP4UTL	BGP4DA,BGP4DAP,BGP4DG,BGP4DL,BGP4DPP
STMP^BGP4UTL	BGP4DCI BGP4D21,BGP4D3,BGP4D31,BGP4D4,BGP4D5,BGP4D61,BGP4D8
\$\$REFUSAL^BGP4UTL1	
EN^BGP5ASL	BGP5DAR,BGP5DELA,BGP5DGPA,BGP5DHEA
GET^BGP5ASL	BGP5GU,BGPGU
^BGP5BAN	opt
^BGP5C1	BGP5CL,BGP5GCMS
GETMEDS^BGP5C1	BGP5C11,BGP5C13,BGP5C3
ACEIALG1^BGP5C11	BGP5C2
ACEICON1^BGP5C11	BGP5C2
ACEIRX1^BGP5C11	BGP5C2
EN^BGP5C11	BGP5C1
^BGP5C12	BGP5C11
EN^BGP5C12	BGP5C2,BGP5C3
ARBALG1^BGP5C13	BGP5C11
ARBCON1^BGP5C13	BGP5C11
ARBRX1^BGP5C13	BGP5C11

^BGP5CL	opt
^BGP5CP	BGP5CL
EP^BGP5CSI	BGP5CL
EP^BGP5CTL	BGP5CTS,BGP5ETS
^BGP5CTS	opt
^BGP5CTXC	opt
TAXCHK^BGP5CTXC	BGP5CL
^BGP5D1	BGP5DAP,BGP5DEL,BGP5DELA,BGP5DGPU,BGP5DHE,BGP5DL,B GP5DNG
	BGP5GALL,BGP5GCOM,BGP5GELD,BGP5GGPU,BGP5GHED,BGP5 GNPL
	BGP5GNTL,BGP5NPL
^BGP5D10	BGP5DPP,BGP5PPP
\$\$HGBA1C^BGP5D2	BGP5D22,BGP5DHE1
\$\$LDL^BGP5D2	BGP5D22,BGP5D41,BGP5D72,BGP5D73,BGP5DHE1,BGP5EL3
\$\$MEANBP^BGP5D2	BGP5D22
\$\$DENTSRV^BGP5D21	BGP5DHE1,BGP5EL2
\$\$DEP^BGP5D21	BGP5D41,BGP5EL3
\$\$DEPSCR^BGP5D21	BGP5D41,BGP5EL3
\$\$DNKA^BGP5D21	BGP5DHE1
\$\$EYE^BGP5D21	BGP5D22,BGP5EL2
\$\$GFR^BGP5D21	BGP5D22,BGP5EL2
\$\$POSUR^BGP5D21	BGP5D22,BGP5EL2
\$\$FLU^BGP5D3	BGP5EL3
\$\$HOSPA^BGP5D31	BGP5EL3
\$\$LAST^BGP5D31	BGP5EL3
\$\$PNEU^BGP5D31	BGP5D3,BGP5EL3
\$\$V2ASTH^BGP5D31	BGP5EL3
\$\$IMMREF^BGP5D32	BGP5D33,BGP5D34,BGP5D35
GETIMMS^BGP5D32	BGP5D33,BGP5D34,BGP5D35
\$\$MMR^BGP5D33	BGP5D32
\$\$OPV^BGP5D33	BGP5D32
\$\$DTAP^BGP5D34	BGP5D32
\$\$HEP^BGP5D35	BGP5D32
\$\$HIB^BGP5D35	BGP5D32
\$\$VAR^BGP5D35	BGP5D32
\$\$MAM^BGP5D4	BGP5EL3
\$\$MAS^BGP5D4	BGP5EL3
\$\$MEANBP^BGP5D41	BGP5D9,BGP5EL3
\$\$CHL^BGP5D5	BGP5DHE1
\$\$DVDX^BGP5D5	BGP5EL3
\$\$DVEX^BGP5D5	BGP5EL3
\$\$DVPED^BGP5D5	BGP5EL3
\$\$REFDV^BGP5D5	BGP5EL3
\$\$BMI^BGP5D6	BGP5D41,BGP5D71,BGP5EL3
\$\$OB^BGP5D6	BGP5D71,BGP5EL3
\$\$OW^BGP5D6	BGP5D71,BGP5EL3
\$\$CRC^BGP5D61	BGP5DHE1,BGP5EL3
\$\$FOB^BGP5D61	BGP5DHE1,BGP5EL3
\$\$IBH^BGP5D61	BGP5DHE1,BGP5EL3
\$\$RECT^BGP5D61	BGP5DHE1,BGP5EL3

\$\$REF^BGP5D61	BGP5DHE1,BGP5EL3
\$\$DENT^BGP5D7	BGP5D41,BGP5D71,BGP5EL3
\$\$DX^BGP5D7	BGP5D71,BGP5EL3
\$\$PED^BGP5D7	BGP5D41,BGP5EL3
\$\$PREG^BGP5D7	BGP5D8
\$\$TOBACCO^BGP5D7	BGP5D41,BGP5D71,BGP5EL3
\$\$MEDNUTR^BGP5D71	BGP5D41
\$\$OTHREL^BGP5D71	BGP5D41
\$\$SPECEX^BGP5D71	BGP5D41
\$\$SPECNUTR^BGP5D71	BGP5D41
\$\$CHOL^BGP5D73	BGP5EL3
\$\$FIRSTIHD^BGP5D9	BGP5D21,BGP5D41,BGP5D73,BGP5DHE1
\$\$V2IHD^BGP5D9	BGP5D21,BGP5D41,BGP5D73,BGP5DHE1
^BGP5DAP	opt
^BGP5DAR	opt
CALCIND^BGP5DCEI	BGP5D1
CALCIND^BGP5DCHI	BGP5D1
CALCIND^BGP5DCI	BGP5D1,BGP5D10
^BGP5DEL	opt
^BGP5DELA	opt
^BGP5DELP	BGP5DEL,BGP5DELA
PRINT1^BGP5DELP	BGP5PHEL
EN^BGP5DESI	BGP5DEL
EN^BGP5DESL	BGP5DEL,BGP5DELA
PT^BGP5DESL	BGP5DEL,BGP5DELA
RT^BGP5DESL	BGP5DEL,BGP5DELA
^BGP5DGPA	opt
^BGP5DGPU	opt
^BGP5DH	BGP5DP
AREACP^BGP5DH	BGP5PARP
^BGP5DHE	opt
^BGP5DHEA	opt
^BGP5DHEP	BGP5DHE
PRINT1^BGP5DHEP	BGP5PHEP
EN^BGP5DHSL	BGP5DHE
PT^BGP5DHSL	BGP5DHE,BGP5DHEA
RT^BGP5DHSL	BGP5DHE
^BGP5DL	opt
^BGP5DNG	opt
^BGP5DP	BGP5DAP,BGP5DGPU,BGP5DL,BGP5DNG,BGP5DPP,BGP5NPL
PRINT1^BGP5DP	BGP5PARP
SETEXCEL^BGP5DP	BGP5PARP
PI^BGP5DP1	BGP5DP3,BGP5DP4,BGP5DP5,BGP5DP6
PI1^BGP5DP1	BGP5DP4,BGP5DP5
SETEXA^BGP5DP1	BGP5DP2,BGP5DP3,BGP5DP7,BGP5PDL1,BGP5PDL2,BGP5PDL3
SETN^BGP5DP1	BGP5DP2,BGP5DP3,BGP5DP7,BGP5PDL0,BGP5PDL1,BGP5PDL2
	BGP5PDL3
I1AGE^BGP5DP11	BGP5DP1
I1AGEP^BGP5DP50	BGP5DP5,BGP5DP51
I1AGEP^BGP5DP60	BGP5DP6

H1^BGP5DPH	BGP5DP1,BGP5DP4,BGP5DP5,BGP5DP6,BGP5ELP1,BGP5HEP1
H2^BGP5DPH	BGP5DP1,BGP5ELP1,BGP5HEP1
H3^BGP5DPH	BGP5DP11
HEADER^BGP5DPH	BGP5DELP,BGP5DHEP,BGP5DP,BGP5DP1,BGP5DP11,BGP5DP3 BGP5DP4,BGP5DP5,BGP5DP50,BGP5DP51,BGP5DP6,BGP5DP60 BGP5DS,BGP5DSP,BGP5EL10,BGP5EL11,BGP5EL12,BGP5EL13 BGP5EL14,BGP5EL15,BGP5EL16,BGP5EL17,BGP5EL18,BGP5EL19 BGP5EL1A,BGP5EL1B,BGP5EL1C BGP5EL1D,BGP5EL1E,BGP5EL1F,BGP5EL1G,BGP5EL1H,BGP5EL1I BGP5EL1J,BGP5EL1K,BGP5EL1L,BGP5ELP1,BGP5ELS,BGP5HEP1 BGP5HES,BGP5SDP
^BGP5DPP	opt
^BGP5DS	BGP5DP
EN^BGP5DSI	BGP5DAP,BGP5DL,BGP5DPP,BGP5NPL
EN^BGP5DSL	BGP5DAP,BGP5DL,BGP5DPP
PT^BGP5DSL	BGP5DAP,BGP5DAR,BGP5DGPA,BGP5DGPU,BGP5DL,BGP5DNG,BGP5DPP BGP5NPL
RT^BGP5DSL	BGP5DAP,BGP5DL,BGP5DPP,BGP5NPL
^BGP5DSP	BGP5DP
I1AGE^BGP5EL10	BGP5ELP1
I1AGE^BGP5EL11	BGP5ELP1
I1AGE^BGP5EL12	BGP5ELP1
I1AGE^BGP5EL13	BGP5ELP1
I1AGE^BGP5EL14	BGP5ELP1
I1AGE^BGP5EL15	BGP5ELP1
I1AGE^BGP5EL16	BGP5ELP1
I1AGE^BGP5EL17	BGP5ELP1
I1AGE^BGP5EL18	BGP5ELP1
I1AGE^BGP5EL19	BGP5ELP1
I1AGE^BGP5EL1A	BGP5ELP1
I1AGE^BGP5EL1B	BGP5ELP1
I1AGE^BGP5EL1C	BGP5ELP1
I1AGE^BGP5EL1D	BGP5ELP1
I1AGE^BGP5EL1E	BGP5ELP1
I1AGE^BGP5EL1F	BGP5ELP1
I1AGE^BGP5EL1G	BGP5ELP1
I1AGE^BGP5EL1H	BGP5ELP1
I1AGE^BGP5EL1I	BGP5ELP1
I1AGE^BGP5EL1J	BGP5ELP1
I1AGE^BGP5EL1K	BGP5ELP1
I22^BGP5EL1L	BGP5ELP1
\$\$FRACTURE^BGP5EL3	BGP5D22
\$\$FUNCTION^BGP5EL4	BGP5EL3
\$\$PHNV^BGP5EL4	BGP5EL3
\$\$TXBMD^BGP5EL4	BGP5D22,BGP5EL3
^BGP5ELH	BGP5DELP
AREACP^BGP5ELH	BGP5PHEL
^BGP5ELHH	BGP5ELL
^BGP5ELL	BGP5DELP,BGP5PHEL

SETN^BGP5ELL1	BGP5ELLN
I1AGE^BGP5ELL2	BGP5ELL1
I1AGE^BGP5ELL3	BGP5ELL1
I1AGE^BGP5ELL4	BGP5ELL1
I1AGE^BGP5ELL5	BGP5ELL1
I1AGE^BGP5ELL6	BGP5ELL1
I1AGE^BGP5ELL7	BGP5ELL1
I1AGE^BGP5ELL8	BGP5ELL1
I1AGE^BGP5ELL9	BGP5ELL1
I1AGE^BGP5ELLA	BGP5ELL1
I1AGE^BGP5ELLB	BGP5ELL1
I1AGE^BGP5ELLC	BGP5ELL1
I1AGE^BGP5ELLD	BGP5ELL1
I1AGE^BGP5ELLE	BGP5ELL1
I1AGE^BGP5ELLF	BGP5ELL1
I1AGE^BGP5ELLG	BGP5ELL1
I1AGE^BGP5ELLH	BGP5ELL1
I1AGE^BGP5ELLI	BGP5ELL1
I1AGE^BGP5ELLJ	BGP5ELL1
I1AGE^BGP5ELLK	BGP5ELL1
I1AGE^BGP5ELL	BGP5ELL1
I1AGE^BGP5ELLM	BGP5ELL1
I22^BGP5ELLN	BGP5ELL1
I1AGEP^BGP5ELLQ	BGP5ELL3
SETN^BGP5ELP1	BGP5EL1L
^BGP5ELS	BGP5DELP
^BGP5ELSL	BGP5ELL
GS^BGP5EUTL	BGP5DEL,BGP5DELA,BGP5GELD
REPORT^BGP5EUTL	BGP5DEL,BGP5GELD
STMP^BGP5EUTL	BGP5DCEI
EP^BGP5GAEL	BGP5GRA,BGPGRA
EP^BGP5GAGP	BGP5GRA,BGPGRA
EP^BGP5GAHE	BGP5GRA,BGPGRA
EP^BGP5GALL	BGP5GR,BGPGR
EP^BGP5GAPU	BGP5GRA,BGPGRA
EP^BGP5GCMS	BGP5GR,BGPGR
EP^BGP5GCOM	BGP5GR,BGPGR
EP^BGP5GELD	BGP5GR,BGPGR
EP^BGP5GGPU	BGP5GR,BGPGR
EP^BGP5GHED	BGP5GR,BGPGR
EP^BGP5GLST	BGP5GU,BGPGU
EP^BGP5GNPL	BGP5GR,BGPGR
EP^BGP5GNTL	BGP5GR,BGPGR
EP^BGP5GPP	BGP5GR,BGPGR
EP^BGP5GUPL	BGP5GU,BGPGU
^BGP5HEH	BGP5DHEP
AREACP^BGP5HEH	BGP5PHEP
^BGP5HEHH	BGP5HEL
^BGP5HEL	BGP5DHEP,BGP5PHEP
^BGP5HES	BGP5DHEP

^BGP5HESL	BGP5HEL
GETIEN^BGP5HUTL	BGP5GUPL,BGP5ULF
GS^BGP5HUTL	BGP5DHE,BGP5GHED
REPORT^BGP5HUTL	BGP5DHE,BGP5GHED
STMP^BGP5HUTL	BGP5DCHI
^BGP5LSTF	opt
^BGP5NPL	opt
NPL1^BGP5NPLD	BGP5PDL
NPL1^BGP5NPLP	BGP5DP
PRINT^BGP5PARP	BGP5DAR,BGP5DGPA
^BGP5PDL	BGP5DP,BGP5PARP
\$\$SB^BGP5PDL1	BGP5D6,BGP5EL3
H1^BGP5PDL1	BGP5PDL0,BGP5PDL3,BGP5PDL4,BGP5PDL7
PI^BGP5PDL1	BGP5PDL0,BGP5PDL3,BGP5PDL4,BGP5PDL7
PI1^BGP5PDL1	BGP5PDL4,BGP5PDL7
I1AGEP^BGP5PDL5	BGP5PDL0
I1AGEP^BGP5PDL8	BGP5PDL7
I1AGE^BGP5PDL9	BGP5PDL1
^BGP5PDLH	BGP5PDL
^BGP5PDLS	BGP5PDL
PRINT^BGP5PHEL	BGP5DELA
PRINT^BGP5PHEP	BGP5DHEA
^BGP5SDP	BGP5DSP
^BGP5TXCH	opt
TAXCHK^BGP5TXCH	BGP5DAP,BGP5DL,BGP5DPP
^BGP5TXCN	opt
TAXCHK^BGP5TXCN	BGP5DGPU,BGP5DNG,BGP5NPL
^BGP5TXEL	opt
TAXCHK^BGP5TXEL	BGP5DEL
^BGP5TXHE	opt
TAXCHK^BGP5TXHE	BGP5DHE
EP^BGP5TXL	BGP5TXS
^BGP5TXV	opt
EP^BGP5TXV1	BGP5CTS,BGP5TXV
^BGP5ULF	opt
\$\$D^BGP5UTL	BGP5DAR,BGP5DNG,BGP5GAGP
\$\$DATE^BGP5UTL	BGP5ASL,BGP5C1,BGP5C11,BGP5C12,BGP5C13,BGP5C2,BGP5C3 BGP5D2,BGP5D21,BGP5D22,BGP5D3,BGP5D31,BGP5D4,BGP5D41 BGP5D5,BGP5D6,BGP5D61,BGP5D7,BGP5D71,BGP5D72,BGP5D73 BGP5D8,BGP5D9,BGP5DGPA,BGP5DHE1,BGP5DP,BGP5EL2,BGP5EL3 BGP5EL4,BGP6DHE1
\$\$LZERO^BGP5UTL	BGP5DAR,BGP5GAGP
\$\$RZERO^BGP5UTL	BGP5DAR,BGP5DNG,BGP5GAGP
EXCELGS^BGP5UTL	BGP5DAR,BGP5DNG,BGP5GAGP
GETIEN^BGP5UTL	BGP5GUPL,BGP5ULF
GS^BGP5UTL	BGP5DGPU,BGP5DNG,BGP5GGPU,BGP5GNTL
JRNL^BGP5UTL	BGP5D1,BGP5D10
REPORT^BGP5UTL	BGP5DAP,BGP5DGPU,BGP5DL,BGP5DNG,BGP5DPP,BGP5GALL BGP5GCOM,BGP5GGPU,BGP5GNPL,BGP5GNTL,BGP5GPP,BGP5NPL
STMP^BGP5UTL	BGP5DCI

\$\$FIRSTPRC^BGP5UTL1	BGP5EL3
\$\$LASTDX^BGP5UTL1	BGP5D72
\$\$LASTPRC^BGP5UTL1	BGP5D72,BGP5DHE1,BGP5EL4
\$\$LASTPRC\^BGP5UTL1	BGP5D3,BGP5D61
\$\$RADREF^BGP5UTL1	BGP5D61
\$\$REFUSAL^BGP5UTL1	BGP5D21,BGP5D3,BGP5D31,BGP5D4,BGP5D5,BGP5D61,BGP5D8 BGP5DHE1
EN^BGP6ASL	BGP6DAR,BGP6DCHA,BGP6DELA,BGP6DGPA,BGP6DHEA
GET^BGP6ASL	BGP6GU
^BGP6AX	BGPPOST
^BGP6AXA	BGP6AX
^BGP6AXB	BGP6AX
^BGP6AXC	BGP6AX
^BGP6AXD	BGP6AX
^BGP6AXE	BGP6AX
^BGP6AXF	BGP6AX
^BGP6AXG	BGP6AX
^BGP6AXGB	BGP6AXG
^BGP6AXGC	BGP6AXG
^BGP6AXH	BGP6AX
^BGP6AXI	BGP6AX
^BGP6AXJ	BGP6AX
^BGP6AXK	BGP6AX
^BGP6AXL	BGP6AX
^BGP6AXM	BGP6AX
^BGP6AXN	BGP6AX
^BGP6AXO	BGP6AX
^BGP6AXP	BGP6AX
^BGP6AXQ	BGP6AX
 ^BGP6AXR	BGP6AX
^BGP6AXS	BGP6AX
^BGP6AXT	BGP6AX
^BGP6AXU	BGP6AX
^BGP6AXV	BGP6AX
^BGP6AXW	BGP6AX
^BGP6AXX	BGP6AX
^BGP6AXY	BGP6AX
^BGP6AXZ	BGP6AX
^BGP6BAN	opt
^BGP6BX	BGPPOST
^BGP6BXA	BGP6BX
^BGP6BXB	BGP6BX
^BGP6BXBB	BGP6BXB
^BGP6BXBC	BGP6BXB
^BGP6BXBD	BGP6BXB

^BGP6BXBE	BGP6BXB
^BGP6BXBF	BGP6BXB
^BGP6BXBG	BGP6BXB
^BGP6BXC	BGP6BX
^BGP6BXD	BGP6BX
^BGP6BXE	BGP6BX
^BGP6BXF	BGP6BX
^BGP6BXG	BGP6BX
^BGP6BXH	BGP6BX
^BGP6BXI	BGP6BX
^BGP6BXJ	BGP6BX
^BGP6BXK	BGP6BX
^BGP6BXL	BGP6BX
^BGP6BXM	BGP6BX
^BGP6BXN	BGP6BX
^BGP6BXO	BGP6BX
^BGP6BXP	BGP6BX
^BGP6BXQ	BGP6BX
^BGP6BXR	BGP6BX
^BGP6BXS	BGP6BX
^BGP6BXT	BGP6BX
^BGP6BXU	BGP6BX
^BGP6BXV	BGP6BX
^BGP6BXW	BGP6BX
^BGP6BXX	BGP6BX
^BGP6BXY	BGP6BX
^BGP6BXZ	BGP6BX
^BGP6C1	BGP6CL,BGP6GCMS
GETMEDS^BGP6C1	BGP6C11,BGP6C13,BGP6C3
ACEIALG1^BGP6C11	BGP6C2
ACEICON1^BGP6C11	BGP6C2
ACEIRX1^BGP6C11	BGP6C2
EN^BGP6C11	BGP6C1
^BGP6C12	BGP6C11
EN^BGP6C12	BGP6C2,BGP6C3
ARBALG1^BGP6C13	BGP6C11
ARBCON1^BGP6C13	BGP6C11
ARBRX1^BGP6C13	BGP6C11
^BGP6CL	opt
^BGP6CP	BGP6CL
EP^BGP6CSI	BGP6CL
EP^BGP6CTL	BGP6CTS
^BGP6CTS	opt
INIT^BGP6CTS	BGP6GU
^BGP6CTXC	opt
TAXCHK^BGP6CTXC	BGP6CL
^BGP6D1	BGP6AHW,BGP6DAP,BGP6DCL,BGP6DEL,BGP6DELA,BGP6DGPU BGP6DHE,BGP6DL,BGP6DNG,BGP6GALL,BGP6GCMP,BGP6GCOM BGP6GELD,BGP6GGPU,BGP6GHED,BGP6GNPL,BGP6GNTL,BGP6NPL
^BGP6D10	BGP6DPP,BGP6GPP

\$\$HGBA1C^BGP6D2	BGP6D22,BGP6DHE1
\$\$LDL^BGP6D2	BGP6D22,BGP6D23,BGP6D41,BGP6D72,BGP6D73,BGP6DHE1,BGP6EL4
\$\$MEANBP^BGP6D2	BGP6D22
\$\$DENTSRV^BGP6D21	BGP6DHE1,BGP6EL2
\$\$DEP^BGP6D21	BGP6D23,BGP6D41,BGP6EL3
\$\$DEPREF^BGP6D21	BGP6D23,BGP6EL3
\$\$DEPSCR^BGP6D21	BGP6D23,BGP6D41,BGP6EL3
\$\$DNKA^BGP6D21	BGP6DHE1
\$\$EYE^BGP6D21	BGP6D22,BGP6EL2
\$\$GFR^BGP6D21	BGP6D22,BGP6EL2
\$\$POSUR^BGP6D21	BGP6D22,BGP6D23,BGP6EL2
\$\$FLU^BGP6D3	BGP6EL3
\$\$HOSPA^BGP6D31	BGP6EL3
\$\$LAST^BGP6D31	BGP6EL3
\$\$PNEU^BGP6D31	BGP6D3,BGP6EL3
\$\$V2ASTH^BGP6D31	BGP6EL3
\$\$IMMREF^BGP6D32	BGP6D33,BGP6D34,BGP6D35
GETIMMS^BGP6D32	BGP6D33,BGP6D34,BGP6D35
\$\$MMR^BGP6D33	BGP6D32
\$\$OPV^BGP6D33	BGP6D32
\$\$DTAP^BGP6D34	BGP6D32
\$\$HEP^BGP6D35	BGP6D32
\$\$HIB^BGP6D35	BGP6D32
\$\$VAR^BGP6D35	BGP6D32
\$\$MAM^BGP6D4	BGP6EL3
\$\$MAS^BGP6D4	BGP6EL3
\$\$LIFE^BGP6D41	BGP6D23
\$\$MEANBP^BGP6D41	BGP6D9,BGP6EL3
\$\$OSTEOSCR^BGP6D42	BGP6EL3
\$\$CHL^BGP6D5	BGP6DHE1
\$\$DVDX^BGP6D5	BGP6EL3
\$\$DVEX^BGP6D5	BGP6EL3
\$\$DVPED^BGP6D5	BGP6EL3
\$\$REFDV^BGP6D5	BGP6EL3
\$\$EAPT^BGP6D52	BGP6D51
\$\$BMI^BGP6D6	BGP6D23,BGP6D41,BGP6D71,BGP6EL3
\$\$OB^BGP6D6	BGP6D71,BGP6EL3
\$\$OW^BGP6D6	BGP6D71,BGP6EL3
\$\$REF^BGP6D6	BGP6D23,BGP6EL3
\$\$CRC^BGP6D61	BGP6DHE1,BGP6EL3
\$\$FOB^BGP6D61	BGP6DHE1,BGP6EL3
\$\$IBH^BGP6D61	BGP6DHE1,BGP6EL3
\$\$RECT^BGP6D61	BGP6DHE1,BGP6EL3
\$\$REF^BGP6D61	BGP6DHE1,BGP6EL3
\$\$DENT^BGP6D7	BGP6D23,BGP6D41,BGP6D71,BGP6DHE1,BGP6EL3
\$\$DX^BGP6D7	BGP6D23,BGP6D71,BGP6DHE1,BGP6EL3
\$\$PED^BGP6D7	BGP6D23,BGP6D41,BGP6EL3
\$\$PREG^BGP6D7	BGP6D8
\$\$TOBACCO^BGP6D7	BGP6D23,BGP6D41,BGP6D71,BGP6DHE1,BGP6EL3
	BGP6D41

\$\$MEDNUTR^BGP6D71	
\$\$OTHREL^BGP6D71	BGP6D41
\$\$SPECEX^BGP6D71	BGP6D41
\$\$SPECNUTR^BGP6D71	BGP6D41
\$\$CHOL^BGP6D73	BGP6EL4
\$\$FIRSTIHD^BGP6D9	BGP6D21,BGP6D41,BGP6D73,BGP6DHE1
\$\$V2IHD^BGP6D9	BGP6D21,BGP6D41,BGP6D73,BGP6DHE1
^BGP6DAP	opt
^BGP6DAR	opt
CALCIND^BGP6DCEI	BGP6D1
^BGP6DCHA	opt
\$\$D^BGP6DCHA	BGP6GACW
GS^BGP6DCHA	BGP6GACW
CALCIND^BGP6DCHI	BGP6D1
GEN^BGP6DCHW	BGP6D1
CALCIND^BGP6DCI	BGP6D1,BGP6D10
^BGP6DCL	opt
CPPL1^BGP6DCLD	BGP6PDL
CPPL1^BGP6DCLP	BGP6DP
^BGP6DEL	opt
^BGP6DELA	opt
^BGP6DELP	BGP6DEL,BGP6DELA
PRINT1^BGP6DELP	BGP6PHEL
EN^BGP6DESI	BGP6DEL
EN^BGP6DESL	BGP6DEL,BGP6DELA
PT^BGP6DESL	BGP6DEL,BGP6DELA
RT^BGP6DESL	BGP6DEL,BGP6DELA
^BGP6DGPA	opt
^BGP6DGPU	opt
^BGP6DH	BGP6DP
AREACP^BGP6DH	BGP6PARP
EOP^BGP6DH	BGP6GUPL,BGP6LSTF,BGP6ULF
^BGP6DHE	opt
^BGP6DHEA	opt
^BGP6DHEP	BGP6DHE
PRINT1^BGP6DHEP	BGP6PHEP
EN^BGP6DHSL	BGP6DHE
PT^BGP6DHSL	BGP6DHE,BGP6DHEA
RT^BGP6DHSL	BGP6DHE
^BGP6DL	opt
^BGP6DLT	opt
^BGP6DMT	opt
^BGP6DNG	opt
^BGP6DP	BGP6AHW,BGP6DAP,BGP6DCL,BGP6DGPU,BGP6DL,BGP6DNG,BGP6DPP BGP6NPL
PRINT1^BGP6DP	BGP6PARP
SETEXCEL^BGP6DP	BGP6PARP
PI^BGP6DP1	BGP6DP3,BGP6DP4,BGP6DP5,BGP6DP6
PI1^BGP6DP1	BGP6DP4,BGP6DP5
SETEXA^BGP6DP1	BGP6DP2,BGP6DP3,BGP6DP7,BGP6PDL1,BGP6PDL2,BGP6PDL3

SETN^BGP6DP1	BGP6DP2,BGP6DP3,BGP6DP7,BGP6PDL0,BGP6PDL1,BGP6PDL2 BGP6PDL3
I1AGE^BGP6DP11	BGP6DP1
I1AGE^BGP6DP12	BGP6DP1
I1AGEP^BGP6DP50	BGP6DP5,BGP6DP51
I1AGEP^BGP6DP60	BGP6DP6
H1^BGP6DPH	BGP6DP1,BGP6DP4,BGP6DP5,BGP6DP6,BGP6ELP1,BGP6HEP1
H2^BGP6DPH	BGP6DP1,BGP6ELP1,BGP6HEP1
H3^BGP6DPH	BGP6DP11
H4^BGP6DPH	BGP6DP12
HEADER^BGP6DPH	BGP6DELP,BGP6DHEP,BGP6DP,BGP6DP1,BGP6DP11,BGP6DP12 BGP6DP3,BGP6DP4,BGP6DP5,BGP6DP50,BGP6DP51,BGP6DP6 BGP6DP60,BGP6DS,BGP6DSP,BGP6EL10,BGP6EL11,BGP6EL12 BGP6EL13,BGP6EL14,BGP6EL15,BGP6EL16,BGP6EL17,BGP6EL18 BGP6EL19,BGP6EL1A,BGP6EL1B BGP6EL1C,BGP6EL1D,BGP6EL1E,BGP6EL1F,BGP6EL1G,BGP6EL1H BGP6EL1I,BGP6EL1J,BGP6EL1K,BGP6EL1L,BGP6EL1M,BGP6ELP1 BGP6ELS,BGP6HEP1,BGP6HES,BGP6SDP
^BGP6DPP	opt
^BGP6DS	BGP6DP
EN^BGP6DSI	BGP6AHW,BGP6DAP,BGP6DL,BGP6DPP,BGP6NPL
EN^BGP6DSL	BGP6DAP,BGP6DL,BGP6DPP
PT^BGP6DSL	BGP6AHW,BGP6DAP,BGP6DAR,BGP6DCL,BGP6DGPA,BGP6DGPU,BGP6DL BGP6DNG,BGP6DPP,BGP6NPL
RT^BGP6DSL	BGP6AHW,BGP6DAP,BGP6DCL,BGP6DL,BGP6DPP,BGP6NPL
^BGP6DSP	BGP6DP
\$\$CPT^BGP6DU	BGP6D22,BGP6D3,BGP6D4,BGP6D42,BGP6D5,BGP6D61,BGP6D72 BGP6D8,BGP6DHE1,BGP6EL4
\$\$CPTI^BGP6DU	BGP6D2,BGP6D21,BGP6D23,BGP6D31,BGP6D61,BGP6D71,BGP6D72 BGP6D73,BGP6D8,BGP6D9,BGP6DHE1,BGP6EL2
\$\$PLCODE^BGP6DU	BGP6D34
\$\$PLTAX^BGP6DU	BGP6D33,BGP6D34,BGP6D35
\$\$RAD^BGP6DU	BGP6D42,BGP6D61,BGP6EL4
\$\$WH^BGP6DU	BGP6D3,BGP6D4
I1AGE^BGP6EL10	BGP6ELP1
I1AGE^BGP6EL11	BGP6ELP1
I1AGE^BGP6EL12	BGP6ELP1
I1AGE^BGP6EL13	BGP6ELP1
I1AGE^BGP6EL14	BGP6ELP1
I1AGE^BGP6EL15	BGP6ELP1
I1AGE^BGP6EL16	BGP6ELP1
I1AGE^BGP6EL17	BGP6ELP1
I1AGE^BGP6EL18	BGP6ELP1
I1AGE^BGP6EL19	BGP6ELP1
I1AGE^BGP6EL1A	BGP6ELP1
I1AGE^BGP6EL1B	BGP6ELP1
I1AGE^BGP6EL1C	BGP6ELP1
I1AGE^BGP6EL1D	BGP6ELP1
I1AGE^BGP6EL1E	BGP6ELP1
I1AGE^BGP6EL1F	BGP6ELP1

I1AGE^BGP6EL1G	BGP6ELP1
I1AGE^BGP6EL1H	BGP6ELP1
I1AGE^BGP6EL1I	BGP6ELP1
I1AGE^BGP6EL1J	BGP6ELP1
I1AGE^BGP6EL1K	BGP6ELP1
IELDPHA^BGP6EL1L	BGP6ELP1
I1AGE^BGP6EL1M	BGP6ELP1
\$\$FRACTURE^BGP6EL3	BGP6D22
\$\$FUNCTION^BGP6EL4	BGP6EL3
\$\$PHNV^BGP6EL4	BGP6EL3
\$\$TXBMD^BGP6EL4	BGP6D22,BGP6EL3
^BGP6ELH	BGP6DELP
AREACP^BGP6ELH	BGP6PHEL
^BGP6ELHH	BGP6ELL
^BGP6ELL	BGP6DELP,BGP6PHEL
SETN^BGP6ELL1	BGP6ELLN
I1AGE^BGP6ELL2	BGP6ELL1
I1AGE^BGP6ELL3	BGP6ELL1
I1AGE^BGP6ELL4	BGP6ELL1
I1AGE^BGP6ELL5	BGP6ELL1
I1AGE^BGP6ELL6	BGP6ELL1
I1AGE^BGP6ELL7	BGP6ELL1
I1AGE^BGP6ELL8	BGP6ELL1
I1AGE^BGP6ELL9	BGP6ELL1
I1AGE^BGP6ELLA	BGP6ELL1
I1AGE^BGP6ELLB	BGP6ELL1
I1AGE^BGP6ELLC	BGP6ELL1
I1AGE^BGP6ELLD	BGP6ELL1
I1AGE^BGP6ELLE	BGP6ELL1
I1AGE^BGP6ELLF	BGP6ELL1
I1AGE^BGP6ELLG	BGP6ELL1
I1AGE^BGP6ELLH	BGP6ELL1
I1AGE^BGP6ELLI	BGP6ELL1
I1AGE^BGP6ELLJ	BGP6ELL1
I1AGE^BGP6ELLK	BGP6ELL1
I1AGE^BGP6ELL	BGP6ELL1
I1AGE^BGP6ELLM	BGP6ELL1
IELDPHA^BGP6ELLN	BGP6ELL1
I1AGE^BGP6ELLO	BGP6ELL1
I1AGEP^BGP6ELLQ	BGP6ELL3
SETN^BGP6ELP1	BGP6EL1L
^BGP6ELS	BGP6DELP
^BGP6ELSL	BGP6ELL
GS^BGP6EUTL	BGP6DEL,BGP6DELA,BGP6GELD
REPORT^BGP6EUTL	BGP6DEL,BGP6GELD
STMP^BGP6EUTL	BGP6DCEI
^BGP6EX	BGPPOST
^BGP6EXA	BGP6EX
^BGP6EXB	BGP6EX

^BGP6EXBB	BGP6EXB
^BGP6EXC	BGP6EX
^BGP6EXD	BGP6EX
^BGP6EXE	BGP6EX
^BGP6EXF	BGP6EX
^BGP6EXG	BGP6EX
^BGP6EXH	BGP6EX
^BGP6EXHB	BGP6EXH
^BGP6EXHC	BGP6EXH
^BGP6EXI	BGP6EX
^BGP6EXJ	BGP6EX
^BGP6EXK	BGP6EX
^BGP6EXL	BGP6EX
^BGP6EXM	BGP6EX
^BGP6EXN	BGP6EX
^BGP6EXO	BGP6EX
^BGP6EXP	BGP6EX
^BGP6EXQ	BGP6EX
^BGP6EXR	BGP6EX
^BGP6EXS	BGP6EX
^BGP6EXT	BGP6EX
^BGP6EXU	BGP6EX
^BGP6EXV	BGP6EX
^BGP6EXW	BGP6EX
^BGP6EXX	BGP6EX
^BGP6EXY	BGP6EX
^BGP6EXZ	BGP6EX
^BGP6FX	BGPPOST
^BGP6FXA	BGP6FX
^BGP6FXB	BGP6FX
^BGP6FXC	BGP6FX
^BGP6FXD	BGP6FX
^BGP6FXE	BGP6FX
^BGP6FXF	BGP6FX
^BGP6FXG	BGP6FX
EP^BGP6GACW	BGP6GRA
EP^BGP6GAEL	BGP6GRA
EP^BGP6GAGP	BGP6GRA
EP^BGP6GAHE	BGP6GRA
EP^BGP6GALL	BGP6GR
EP^BGP6GAPU	BGP6GRA
EP^BGP6GCMP	BGP6GR
EP^BGP6GCMS	BGP6GR
EP^BGP6GCOM	BGP6GR
EP^BGP6GELD	BGP6GR
EP^BGP6GGPU	BGP6GR
EP^BGP6GHED	BGP6GR
EP^BGP6GLST	BGP6GU
EP^BGP6GLTX	BGP6GU
EP^BGP6GMTX	BGP6GU

EP^BGP6GNPL	BGP6GR
EP^BGP6GNTL	BGP6GR
EP^BGP6GPP	BGP6GR
EP^BGP6GUPL	BGP6GU
^BGP6HEH	BGP6DHEP
AREACP^BGP6HEH	BGP6PHEP
^BGP6HEHH	BGP6HEL
^BGP6HEL	BGP6DHEP,BGP6PHEP
^BGP6HES	BGP6DHEP
^BGP6HESL	BGP6HEL
GETIEN^BGP6HUTL	BGP6GUPL,BGP6ULF
GS^BGP6HUTL	BGP6DHE,BGP6GHED
REPORT^BGP6HUTL	BGP6DHE,BGP6GHED
STMP^BGP6HUTL	BGP6DCHI
^BGP6HX	BGPPOST
^BGP6HXA	BGP6HX
^BGP6HXB	BGP6HX
^BGP6HXC	BGP6HX
^BGP6HXD	BGP6HX
^BGP6HXE	BGP6HX
^BGP6HXF	BGP6HX
^BGP6HXG	BGP6HX
^BGP6HXH	BGP6HX
^BGP6HXI	BGP6HX
^BGP6HXJ	BGP6HX
^BGP6HXK	BGP6HX
^BGP6HXL	BGP6HX
^BGP6HXM	BGP6HX
^BGP6HZN	BGP6HX
^BGP6HXO	BGP6HX
^BGP6HXP	BGP6HX
^BGP6HXQ	BGP6HX
^BGP6HXR	BGP6HX
^BGP6HXS	BGP6HX
^BGP6HXT	BGP6HX
^BGP6HXU	BGP6HX
^BGP6HXV	BGP6HX
^BGP6HXW	BGP6HX
^BGP6HXX	BGP6HX
^BGP6HXY	BGP6HX
^BGP6HXZ	BGP6HX
^BGP6LSTF	lopt
^BGP6NPL	lopt
NPL1^BGP6NPLD	BGP6PDL
NPL1^BGP6NPLP	BGP6DP
PRINT^BGP6PARP	BGP6DAR,BGP6DGPA
^BGP6PDL	BGP6DP,BGP6PARP
\$\$SB^BGP6PDL1	BGP6D23,BGP6D6,BGP6EL3
H1^BGP6PDL1	BGP6PDL0,BGP6PDL3,BGP6PDL4,BGP6PDL7
PI^BGP6PDL1	BGP6PDL0,BGP6PDL3,BGP6PDL4,BGP6PDL7

PI1^BGP6PDL1	BGP6PDL4,BGP6PDL7
I1AGEP^BGP6PDL5	BGP6PDL0
I1AGEP^BGP6PDL8	BGP6PDL7
I1AGE^BGP6PDL9	BGP6PDL1
I1AGE^BGP6PDLA	BGP6PDL1
^BGP6PDLH	BGP6PDL
^BGP6PDLS	BGP6PDL
PRINT^BGP6PHEL	BGP6DELA
PRINT^BGP6PHEP	BGP6DHEA
^BGP6SDP	BGP6DSP
^BGP6ULF	opt
\$\$D^BGP6UTL	BGP6DAR,BGP6DNG,BGP6GAGP
\$\$DATE^BGP6UTL	BGP6ASL,BGP6C1,BGP6C11,BGP6C12,BGP6C13,BGP6C2,BGP6C3 BGP6D2,BGP6D21,BGP6D22,BGP6D23,BGP6D3,BGP6D31,BGP6D4 BGP6D41,BGP6D42,BGP6D5,BGP6D51,BGP6D52,BGP6D6,BGP6D61 BGP6D7,BGP6D71,BGP6D72,BGP6D73,BGP6D8,BGP6D9,BGP6DGPA BGP6DHE1,BGP6DP,BGP6EL2 BGP6EL3,BGP6EL4
\$\$LZERO^BGP6UTL	BGP6DAR,BGP6GAGP
\$\$RZERO^BGP6UTL	BGP6DAR,BGP6DNG,BGP6GACW,BGP6GAGP
EXCELGS^BGP6UTL	BGP6DAR,BGP6DNG,BGP6GAGP
GETIEN^BGP6UTL	BGP6GUPL,BGP6ULF
GS^BGP6UTL	BGP6DGPU,BGP6DNG,BGP6GGPU,BGP6GNTL
JRNL^BGP6UTL	BGP6D1,BGP6D10
REPORT^BGP6UTL	BGP6AHW,BGP6DAP,BGP6DCL,BGP6DGPU,BGP6DL,BGP6DNG,BGP6DPP BGP6GALL,BGP6GCMP,BGP6GCOM,BGP6GGPU,BGP6GNPL,BGP6GNTL BGP6GPP,BGP6NPL
STMP^BGP6UTL	BGP6DCI
XTMP^BGP6UTL	BGP6C1,BGP6D1
\$\$FIRSTPRC^BGP6UTL1	BGP6EL3
\$\$LASTDX^BGP6UTL1	BGP6D23,BGP6D42,BGP6D72,BGP6EL3
\$\$LASTPRC^BGP6UTL1	BGP6D72,BGP6DHE1,BGP6EL4
\$\$LASTPRCI^BGP6UTL1	BGP6D3,BGP6D42,BGP6D61
\$\$RADREF^BGP6UTL1	BGP6D42,BGP6D61
\$\$REFUSAL^BGP6UTL1	BGP6D21,BGP6D3,BGP6D31,BGP6D4,BGP6D5,BGP6D6,BGP6D61 BGP6D8,BGP6DHE1
^BGP6XTCH	opt
TAXCHK^BGP6XTCH	BGP6DAP,BGP6DL,BGP6DPP
^BGP6XTCN	opt
TAXCHK^BGP6XTCN	BGP6AHW,BGP6DCL,BGP6DGPU,BGP6DNG,BGP6NPL
^BGP6XTEL	opt
TAXCHK^BGP6XTEL	BGP6DEL
^BGP6XTHE	opt
TAXCHK^BGP6XTHE	BGP6DHE
EP^BGP6XTL	BGP6XTS
^BGP6XTV	opt
INIT^BGP6XTV	BGP6GU

EP^BGP6XTV1	BGP6CTS,BGP6XTV
^BGPAP	BGPDPARP
I1^BGPAP1	BGPAP
I12^BGPAP12	BGPAP
I13^BGPAP13	BGPAP
I14^BGPAP14	BGPAP
I1B^BGPAP1B	BGPAP
I2A^BGPAP2	BGPAP
I2B^BGPAP2	BGPAP
I2C^BGPAP2	BGPAP
I22^BGPAP22	BGPAP
I23^BGPAP23	BGPAP
I24^BGPAP24	BGPAP
I3A^BGPAP3	BGPAP
I3B^BGPAP3	BGPAP
I3C^BGPAP3	BGPAP
I4A^BGPAP4	BGPAP
I4B^BGPAP4	BGPAP
I4C^BGPAP4	BGPAP
I5A^BGPAP5	BGPAP
I5B^BGPAP5	BGPAP
I5C^BGPAP5	BGPAP
I6^BGPAP6	BGPAP
I6A^BGPAP6	BGPAP
I7^BGPAP7	BGPAP
I8^BGPAP8	BGPAP
IA^BGPAP8A	BGPAP
IB^BGPAP8B	BGPAP
IC^BGPAP8C	BGPAP
ID^BGPAPD	BGPAP
H2^BGPAPH	BGPAP13,BGPAP22B,BGPAP22C,BGPAP30,BGPAP30A,BGPAP30B BGPAP30C
HEADER^BGPAPH	BGPAP1,BGPAP12,BGPAP13,BGPAP14,BGPAP1B,BGPAP2,BGPAP22 BGPAP22A,BGPAP22B,BGPAP22C,BGPAP23,BGPAP24,BGPAP3 BGPAP30,BGPAP30A,BGPAP30B,BGPAP30C,BGPAP4,BGPAP5,BGPAP6 BGPAP7,BGPAP8,BGPAP8A,BGPAP8B,BGPAP8C,BGPAPD,BGPAPH
^BGPD	lopt
EN^BGPD0	BGPD
^BGPD1	BGPD,BGPDARP,BGPDTA
\$\$DM^BGPD1	BGPD24,BGPDA,BGPDD
I12^BGPD12	BGPD1
I13^BGPD13	BGPD1
I14^BGPD14	BGPD1
I1B^BGPD1B	BGPD1
I2A^BGPD2	BGPD1
I2B^BGPD2	BGPD1
I2C^BGPD2	BGPD1
I22^BGPD22	BGPD1
I24^BGPD24	BGPD1
I29^BGPD29	BGPD1

I3A^BGPD3	BGPD1
I3B^BGPD3	BGPD1
I3C^BGPD3	BGPD1
I30^BGPD30	BGPD1
I4A^BGPD4	BGPD1
I4B^BGPD4	BGPD1
I4C^BGPD4	BGPD1
I5A^BGPD5	BGPD1
I5B^BGPD5	BGPD1
I5C^BGPD5	BGPD1
I6^BGPD6	BGPD1
I6A^BGPD6	BGPD1
I7^BGPD7	BGPD1
I8^BGPD8	BGPD1
IA^BGPDA	BGPD1
^BGPDAREA	opt
^BGPDARP	opt
IB^BGPDB	BGPD1
IC^BGPDC	BGPD1
ID^BGPDD	BGPD1
^BGPDFTA	opt
^BGPDH	BGPDP
EOP^BGPDH	BGP3ULF,BGP4LSTF,BGP4ULF,BGP5GUPL,BGP5LSTF,BGP5ULF BGPDUPL
EN^BGPDL	BGPD,BGPDARP,BGPDFTA
^BGPDP	BGPD,BGPDARP,BGPDFTA
PRINT1^BGPDP	BGPDPARP
I1^BGPDP1	BGPDP
I12^BGPDP12	BGPDP
I13^BGPDP13	BGPDP
I14^BGPDP14	BGPDP
I1B^BGPDP1B	BGPDP
I2A^BGPDP2	BGPDP
I2B^BGPDP2	BGPDP
I2C^BGPDP2	BGPDP
I22^BGPDP22	BGPDP
POV^BGPDP221	BGPDP22
AAGE^BGPDP22A	BGPDP22
MAGE^BGPDP22B	BGPDP22
FAGE^BGPDP22C	BGPDP22
I23^BGPDP23	BGPDP
I24^BGPDP24	BGPDP
I29^BGPDP29	BGPDP
I3A^BGPDP3	BGPDP
I3B^BGPDP3	BGPDP
I3C^BGPDP3	BGPDP
I30^BGPDP30	BGPDP
^BGPDP30A	BGPDP30
^BGPDP30B	BGPDP30
^BGPDP30C	BGPDP30

I4A^BGPDP4	BGPDP
I4B^BGPDP4	BGPDP
I4C^BGPDP4	BGPDP
I5A^BGPDP5	BGPDP
I5B^BGPDP5	BGPDP
I5C^BGPDP5	BGPDP
I6^BGPDP6	BGPDP
I6A^BGPDP6	BGPDP
I7^BGPDP7	BGPDP
I8^BGPDP8	BGPDP
IA^BGPDPA	BGPDP
PRINT^BGPDPARP	BGPDAREA
IB^BGPDPB	BGPDP
IC^BGPDPD	BGPDP
ID^BGPDPD	BGPDP
AREACP^BGPDPH	BGPDPARP
H1^BGPDPH	BGPDP1,BGPDP12,BGPDP13,BGPDP14,BGPDP1B,BGPDP2,BGPDP22 BGPDP22B,BGPDP22C,BGPDP23,BGPDP24,BGPDP29,BGPDP3,BGPDP30 BGPDP30A,BGPDP30B,BGPDP30C,BGPDP31,BGPDP32,BGPDP4,BGPDP5 BGPDP6,BGPDP7,BGPDP8,BGPDPB,BGPDPB,BGPDPB,BGPDPD BGPDP1,BGPDP12,BGPDP13,BGPDP14,BGPDP1B,BGPDP2,BGPDP22 BGPDP22B,BGPDP22C,BGPDP23,BGPDP24,BGPDP29,BGPDP3,BGPDP30 BGPDP30A,BGPDP30B,BGPDP30C,BGPDP31,BGPDP32,BGPDP4,BGPDP5 BGPDP6,BGPDP7,BGPDP8,BGPDPB,BGPDPB,BGPDPB,BGPDPD
H2^BGPDPH	BGPDP1,BGPDP12,BGPDP13,BGPDP14,BGPDP1B,BGPDP2,BGPDP22 BGPDP22B,BGPDP22C,BGPDP23,BGPDP24,BGPDP29,BGPDP3,BGPDP30 BGPDP30A,BGPDP30B,BGPDP30C,BGPDP31,BGPDP32,BGPDP4,BGPDP5 BGPDP6,BGPDP7,BGPDP8,BGPDPB,BGPDPB,BGPDPB,BGPDPD
H3^BGPDPH	BGPDP1
H4^BGPDPH	BGPDP1B
H6^BGPDPH	BGPDPB
HEADER^BGPDPH	BGP3DPH,BGPAP1,BGPAP12,BGPAP13,BGPAP14,BGPAP1B,BGPAP2 BGPAP22,BGPAP23,BGPAP24,BGPAP3,BGPAP4,BGPAP5,BGPAP6 BGPAP7,BGPAP8,BGPAPB,BGPAPD,BGPDP1,BGPDP12,BGPDP13 BGPDP14,BGPDP1B,BGPDP2,BGPDP22,BGPDP22A,BGPDP22B BGPDP22C,BGPDP23,BGPDP24 BGPDP29,BGPDP3,BGPDP30,BGPDP30A,BGPDP30B,BGPDP30C BGPDP31,BGPDP32,BGPDP4,BGPDP5,BGPDP6,BGPDP7,BGPDP8 BGPDPB,BGPDPARP,BGPDPB,BGPDPB,BGPDPD,BGPDPH
^BGPDS	BGPDP
TAXCHK^BGPDT	BGPD,BGPDARP,BGPDFTA, opt
^BGPDU	BGP3D6,BGP4D61,BGPDB
\$\$CPT^BGPDU	BGP3D3,BGP3D4,BGP3D8,BGP4D3,BGP4D4,BGP4D5,BGP4D61,BGP4D8 BGP5D22,BGP5D3,BGP5D4,BGP5D5,BGP5D61,BGP5D72,BGP5D8 BGP5DHE1,BGP5EL4,BGPD24,BGPD6,BGPD7
\$\$CPTI^BGPDU	BGP3D2,BGP3D4,BGP3D6,BGP3D7,BGP3D9,BGP4D2,BGP4D21 BGP4D31,BGP4D61,BGP4D7,BGP4D71,BGP4D8,BGP5D2,BGP5D21 BGP5D31,BGP5D61,BGP5D7,BGP5D71,BGP5D72,BGP5D73,BGP5D8 BGP5D9,BGP5EL2,BGPD4,BGPDB
\$\$PLCODE^BGPDU	BGP5D34
\$\$PLTAX^BGPDU	BGP5D33,BGP5D34,BGP5D35
\$\$RAD^BGPDU	BGP5D61,BGP5EL4
\$\$WH^BGPDU	BGP3D3,BGP4D3,BGP4D4,BGP5D3,BGP5D4,BGP6,BGP7
BANNER^BGPDU	opt

^BGPDUP	opt
CATSTR^BGPGR	BGP5GRA,BGP6GRA,BGPGRA
CLTAX^BGPPOST1	BGPPOST
DRUGS^BGPPOST1	BGPPOST
LAB^BGPPOST1	BGPPOST
MHTAX^BGPPOST1	BGPPOST
PRVTAX^BGPPOST1	BGPPOST
^BGPTXA	BGPTX
^BGPTXB	BGPTX
^BGPTXC	BGPTX
^BGPTXD	BGPTX
^BGPTXE	BGPTX
^BGPTXEB	BGPTXE
^BGPTXEC	BGPTXE
^BGPTXED	BGPTXE
^BGPTXEE	BGPTXE
^BGPTXEF	BGPTXE
^BGPTXEG	BGPTXE
^BGPTXF	BGPTX
^BGPTXG	BGPTX
^BGPTXGB	BGPTXG
^BGPTXH	BGP50P1,BGPTX
^BGPTXI	BGPTX
^BGPTXJ	BGPTX
^BGPTXK	BGPTX
^BGPTXL	BGPTX
^BGPTXM	BGPTX
^BGPTXN	BGPTX
^BGPTXO	BGPTX
^BGPTXP	BGPTX
^BGPTXPB	BGPTXP
^BGPTXQ	BGPTX
^BGPTXR	BGPTX
^BGPTXS	BGPTX
^BGPTXT	BGPTX
^BGPTXU	BGPTX
^BGPTXV	BGPTX
^BGPTXW	BGPTX
^BGPTXX	BGPTX
^BGPTXY	BGPTX
^BGPTXZ	BGPTX
^BGPUTXA	BGPUTX
^BGPUTXAB	BGPUTXA
^BGPUTXAC	BGPUTXA
^BGPUTXAD	BGPUTXA
^BGPUTXAE	BGPUTXA
^BGPUTXAF	BGPUTXA
^BGPUTXAG	BGPUTXA
^BGPUTXAH	BGPUTXA
^BGPUTXAI	BGPUTXA

^BGPUTXAJ	BGPUTXA
^BGPUTXAK	BGPUTXA
^BGPUTXAL	BGPUTXA
^BGPUTXAM	BGPUTXA
^BGPUTXAN	BGPUTXA
^BGPUTXAO	BGPUTXA
^BGPUTXAP	BGPUTXA
^BGPUTXAQ	BGPUTXA
^BGPUTXAR	BGPUTXA
^BGPUTXB	BGPUTX
^BGPUTXC	BGPUTX
^BGPUTXD	BGPUTX
^BGPUTXE	BGPUTX
^BGPUTXF	BGPUTX
^BGPUTXG	BGPUTX
^BGPUTXH	BGPUTX
^BGPUTXI	BGPUTX
^BGPUTXJ	BGPUTX
^BGPUTXK	BGPUTX
^BGPUTXL	BGPUTX
^BGPUTXM	BGPUTX
^BGPUTXN	BGPUTX
^BGPUTXO	BGPUTX
^BGPUTXP	BGPUTX
^BGPUTXQ	BGPUTX
^BGPUTXR	BGPUTX
^BGPUTXS	BGPUTX
^BGPUTXT	BGPUTX
^BGPUTXU	BGPUTX
^BGPUTXV	BGPUTX
^BGPUTXW	BGPUTX
^BGPUTXX	BGPUTX
^BGPUTXY	BGPUTX
^BGPUTXZ	BGPUTX
^BGPUTXZB	BGPUTXZ
^BGPUTXZC	BGPUTXZ
^BGPUTXZD	BGPUTXZ
^BGPUTXZE	BGPUTXZ
^BGPVTXA	BGPVTX
^BGPVTXB	BGPVTX
^BGPVTXC	BGPVTX
^BGPVTXD	BGPVTX
^BGPVTXE	BGPVTX
^BGPVTXF	BGPVTX
^BGPVTXG	BGPVTX
^BGPVTXH	BGPVTX
^BGPVTXI	BGPVTX
^BGPVTXIB	BGPVTXI

^BGPVTXIC	BGPVTXI
^BGPVTXJ	BGPVTX
^BGPVTXK	BGPVTX
^BGPVTXL	BGPVTX
^BGPVTXM	BGPVTX
^BGPVTXN	BGPVTX
^BGPVTXO	BGPVTX
^BGPVTXP	BGPVTX
^BGPVTXQ	BGPVTX
^BGPVTXR	BGPVTX
^BGPVTXS	BGPVTX
^BGPVTXT	BGPVTX
^BGPVTXU	BGPVTX
^BGPVTXV	BGPVTX
^BGPVTXW	BGPVTX
^BGPVTXX	BGPVTX
^BGPVTXY	BGPVTX
^BGPVTXZ	BGPVTX
^BGPWTXA	BGPWTX
^BGPWTXB	BGPWTX
^BGPWTXBB	BGPWTXB
^BGPWTXBC	BGPWTXB
^BGPWTXBD	BGPWTXB
^BGPWTXBE	BGPWTXB
^BGPWTXBF	BGPWTXB
^BGPWTXBG	BGPWTXB
^BGPWTXC	BGPWTX
^BGPWTXD	BGPWTX
^BGPWTXE	BGPWTX
^BGPWTXF	BGPWTX
^BGPWTXG	BGPWTX
^BGPWTXH	BGPWTX
^BGPWTXI	BGPWTX
^BGPWTXJ	BGPWTX
^BGPWTXK	BGPWTX
^BGPWTXL	BGPWTX
^BGPWTXM	BGPWTX
^BGPWTXN	BGPWTX
^BGPWTXO	BGPWTX
^BGPWTXP	BGPWTX
^BGPWTXQ	BGPWTX
^BGPWTXR	BGPWTX
^BGPWTXS	BGPWTX
^BGPWTXT	BGPWTX
^BGPWTXU	BGPWTX
IMMFORC^BIRPC	BGPD23
^DIC	BGP3DA,BGP3DG,BGP3DL,BGP3UTL,BGP4DA,BGP4DG,BGP4DHE BGP4DL,BGP4DPP,BGP4HUTL,BGP4UTL,BGP5CL,BGP5CTL,BGP5DEL BGP5DGPU,BGP5DHE,BGP5DL,BGP5DNG,BGP5DPP,BGP5ETL,BGP5EUTL BGP5HUTL,BGP5NPL,BGP5TXL,BGP5UTL,BGP6AHW,BGP6CL,BGP6CTL

	BGP6DCL,BGP6DEL BGP6DGPU,BGP6DHE,BGP6DL,BGP6DNG,BGP6DPP,BGP6ETL,BGP6EUTL BGP6HUTL,BGP6NPL,BGP6UTL,BGP6XTL,BGPD,BGPDARP,BGPDFTA BGPOST,BGPOST1, dd90241.0111, dd90241.02
IX^DIC	BGPOST
FILE^DICN	BGP3ULF,BGP4ULF,BGP5CTL,BGP5ETL,BGP5GAEL,BGP5GAGP BGP5GAHE,BGP5GALL,BGP5GAPU,BGP5GCMS,BGP5GCOM,BGP5GELD BGP5GGPU,BGP5GHED,BGP5GLST,BGP5GNPL,BGP5GNTL,BGP5GPP BGP5GUPL,BGP5TXL,BGP5ULF,BGP6CTL,BGP6ETL,BGP6GAEL BGP6GAGP,BGP6GAHE,BGP6GALL BGP6GAPU,BGP6GCMP,BGP6GCMS,BGP6GCOM,BGP6GELD,BGP6GGPU BGP6GHED,BGP6GLST,BGP6GLTX,BGP6GMTX,BGP6GNPL,BGP6GNTL BGP6GPP,BGP6GUPL,BGP6ULF,BGP6XTL,BGPDUP,BGPOST
DT^DICRW	BGP50P1
^DIE	BGP5CTL,BGP5ETL,BGP5GAEL,BGP5GAGP,BGP5GAHE,BGP5GALL BGP5GAPU,BGP5GCMS,BGP5GCOM,BGP5GELD,BGP5GGPU,BGP5GHED BGP5GLST,BGP5GNPL,BGP5GNTL,BGP5GPP,BGP5TXL,BGP6CTL BGP6ETL,BGP6GACW,BGP6GAEL,BGP6GAGP,BGP6GAHE,BGP6GALL BGP6GAPU,BGP6GCMP,BGP6GCMS BGP6GCOM,BGP6GELD,BGP6GGPU,BGP6GHED,BGP6GLST,BGP6GLTX BGP6GMTX,BGP6GNPL,BGP6GNTL,BGP6GPP,BGP6XTL,BGPOST
FILE^DIE	BGPOST1
UPDATE^DIE	BGP5GU,BGP6GU,BGPGU
^DIK	BGP5GU,BGP6GU,BGPGU BGP3DA,BGP3DG,BGP3DL,BGP3ULF,BGP4DA,BGP4DAP,BGP4DG BGP4DHE,BGP4DL,BGP4DPP,BGP4ULF,BGP50P1,BGP5DAP,BGP5DEL BGP5DGPU,BGP5DHE,BGP5DL,BGP5DNG,BGP5DPP,BGP5GU,BGP5GUPL BGP5NPL,BGP5ULF,BGP6AHW,BGP6DAP,BGP6DCL,BGP6DEL,BGP6DGPU BGP6DHE,BGP6DL BGP6DNG,BGP6DPP,BGP6GU,BGP6GUPL,BGP6NPL,BGP6ULF,BGPD BGPDARP,BGPDFTA,BGPDUP,BGPGU,BGPOST,BGPOST1
IX1^DIK	BGP3ULF,BGP4ULF,BGP5GUPL,BGP5ULF,BGP6GUPL,BGP6ULF,BGPDUP BGPOST,BGPOST1
IXALL^DIK	BGPDI1
^DIM	dd90242, dd90242.01, dd90244.01, dd90244.02, dd90370.01 dd90370.02, dd90371.01, dd90371.02, dd90371.06 dd90372.01, dd90372.02, dd90373.01, dd90373.02 dd90373.06, dd90374.01, dd90374.02, dd90374.06 dd90375.01, dd90375.02, dd90376.01 dd90376.02, dd90376.06
\$\$GET1^DIQ	BGP5GU,BGP6DCHW,BGP6GU,BGPGU
^DIR	BGP3DA,BGP3DAR,BGP3DARA,BGP3DARH,BGP3DG,BGP3DH,BGP3DL BGP3DP,BGP3DPH,BGP3DS,BGP3DSI,BGP3DSL,BGP3PDLS,BGP3TXCH BGP3ULF,BGP4DA,BGP4DAP,BGP4DAR,BGP4DARA,BGP4DG,BGP4DH BGP4DHE,BGP4DHEA,BGP4DHEP,BGP4DHSL,BGP4DL,BGP4DP,BGP4DPH BGP4DPP,BGP4DSI,BGP4DSL,BGP4HEH,BGP4HEHH,BGP4HES,BGP4HESL BGP4LSTF,BGP4PDLS,BGP4TXCH,BGP4TXHE,BGP4ULF,BGP5ASL BGP5CL,BGP5CP,BGP5CSI,BGP5CTL,BGP5CTS,BGP5CTXC,BGP5DAP BGP5DAR,BGP5DEL,BGP5DELA,BGP5DELP,BGP5DESI,BGP5DESL

	BGP5DGPA,BGP5DGPU,BGP5DH BGP5DHE,BGP5DHEA,BGP5DHEP,BGP5DHSL,BGP5DL,BGP5DNG,BGP5DP BGP5DPH,BGP5DPP,BGP5DS,BGP5DSI,BGP5DSL,BGP5ELH,BGP5ELHH BGP5ELS,BGP5ELSL,BGP5ETL,BGP5ETS,BGP5HEH,BGP5HEHH BGP5HES,BGP5HESL,BGP5LSTF,BGP5NPL,BGP5NPLP,BGP5PDLS BGP5TXCH,BGP5TXCN BGP5TXEL,BGP5TXHE,BGP5TXL,BGP5TXS,BGP5TXV,BGP5ULF BGP6AHW,BGP6ASL,BGP6CL,BGP6CP,BGP6CSI,BGP6CTL,BGP6CTS BGP6CTXC,BGP6DAP,BGP6DAR,BGP6DCHA,BGP6DCL,BGP6DCLP BGP6DEL,BGP6DELA,BGP6DELP,BGP6DESI,BGP6DESL,BGP6DGPA BGP6DGPU,BGP6DH,BGP6DHE BGP6DHEA,BGP6DHEP,BGP6DHSL,BGP6DL,BGP6DLT,BGP6DMT BGP6DNG,BGP6DP,BGP6DPH,BGP6DPP,BGP6DS,BGP6DSI,BGP6DSL BGP6ELH,BGP6ELHH,BGP6ELS,BGP6ELSL,BGP6ETL,BGP6HEH BGP6HEHH,BGP6HES,BGP6HESL,BGP6LSTF,BGP6NPL,BGP6NPLP BGP6PDLS,BGP6ULF,BGP6XTCH BGP6XTCN,BGP6XTEL,BGP6XTHE,BGP6XTL,BGP6XTS,BGP6XTV,BGPAP BGPAPH,BGPD,BGPD0,BGPDAREA,BGPDARP,BGPDFTA,BGPDH,BGPDL BGPDP,BGPDPH,BGPDS,BGPDT,BGPDUP
^DIR	BGP50P1 BGP3DSI,BGP3DSL,BGP4DHSL,BGP4DSI,BGP4DSL,BGP5ASL,BGP5CSI BGP5CTL,BGP5CTS,BGP5DESI,BGP5DESL,BGP5DHSL,BGP5DSI BGP5DSL,BGP5ETL,BGP5ETS,BGP5TXL,BGP5TXS,BGP5TXV,BGP5TXV1 BGP6ASL,BGP6CSI,BGP6CTL,BGP6CTS,BGP6DESI,BGP6DESL BGP6DHSL,BGP6DSI BGP6DSL,BGP6ETL,BGP6XTL,BGP6XTS,BGP6XTV,BGP6XTV1,BGPD0 BGPDL
TERM^VALM0	BGP3DSI,BGP3DSL,BGP4DHSL,BGP4DSI,BGP4DSL,BGP5ASL,BGP5CSI BGP5CTL,BGP5CTS,BGP5DESI,BGP5DESL,BGP5DHSL,BGP5DSI BGP5DSL,BGP5ETL,BGP5ETS,BGP5TXL,BGP5TXS,BGP5TXV,BGP5TXV1 BGP6ASL,BGP6CSI,BGP6CTL,BGP6CTS,BGP6DESI,BGP6DESL BGP6DHSL,BGP6DSI BGP6DSL,BGP6ETL,BGP6XTL,BGP6XTS,BGP6XTV,BGP6XTV1,BGPD0 BGPDL
CLEAR^VALM1	BGP3DSI,BGP3DSL,BGP4DHSL,BGP4DSI,BGP4DSL,BGP5ASL,BGP5CSI BGP5CTL,BGP5CTS,BGP5DESI,BGP5DESL,BGP5DHSL,BGP5DSI BGP5DSL,BGP5ETL,BGP5ETS,BGP5TXL,BGP5TXS,BGP5TXV,BGP5TXV1 BGP6ASL,BGP6CSI,BGP6CTL,BGP6CTS,BGP6DESI,BGP6DESL BGP6DHSL,BGP6DSI BGP6DSL,BGP6ETL,BGP6XTL,BGP6XTS,BGP6XTV,BGP6XTV1,BGPD0 BGPDL
FULL^VALM1	BGP3DSI,BGP3DSL,BGP4DHSL,BGP4DSI,BGP4DSL,BGP5ASL,BGP5CSI BGP5CTL,BGP5CTS,BGP5DESI,BGP5DESL,BGP5DHSL,BGP5DSI BGP5DSL,BGP5ETL,BGP5ETS,BGP5TXL,BGP5TXS,BGP5TXV,BGP5TXV1 BGP6ASL,BGP6CSI,BGP6CTL,BGP6CTS,BGP6DESI,BGP6DESL BGP6DHSL,BGP6DSI BGP6DSL,BGP6ETL,BGP6XTL,BGP6XTS,BGP6XTV,BGP6XTV1,BGPD0 BGPDL
^XBDBQUE	BGP3DA,BGP3DG,BGP3DL,BGP3TXCH,BGP4DA,BGP4DAP,BGP4DG BGP4DHE,BGP4DL,BGP4DPP,BGP4TXCH,BGP4TXHE,BGP5CTXC

BGP5DAP,BGP5DEL,BGP5DELA,BGP5DGPU,BGP5DHE,BGP5DL,BGP5DNG
 BGP5DPP,BGP5NPL,BGP5TXCH,BGP5TCN,BGP5TXEL,BGP5TXHE
 BGP6AHW,BGP6CTXC,BGP6DAP
 BGP6DCL,BGP6DEL,BGP6DELA,BGP6DGPU,BGP6DHE,BGP6DL,BGP6DNG
 BGP6DPP,BGP6NPL,BGP6XTCH,BGP6XTCN,BGP6XTEL,BGP6XTHE,BGPD
 BGPDAREA,BGPDT
\$\$VAL^XBDIQ1
 BGP3D6,BGP3D9,BGP3DAR,BGP3DARA,BGP3DARH,BGP4D5,BGP4D7
 BGP4DAR,BGP4DARA,BGP4DH,BGP4DHEA,BGP4HEH,BGP4HEHH
 BGP4PDLH,BGP5C1,BGP5C11,BGP5C13,BGP5C2,BGP5C3,BGP5CTL
 BGP5CTS,BGP5D5,BGP5D7,BGP5D72,BGP5DH,BGP5ELH,BGP5ELHH
 BGP5ETL,BGP5ETS,BGP5HEH
 BGP5HEHH,BGP5PDLH,BGP5TXV,BGP5TXV1,BGP6C1,BGP6C11
 BGP6C13,BGP6C2,BGP6C3,BGP6CTL,BGP6CTS,BGP6D5,BGP6D7
 BGP6D72,BGP6DCHW,BGP6DH,BGP6ELH,BGP6ELHH,BGP6ETL,BGP6HEH
 BGP6HEHH,BGP6PDLH,BGP6XTV,BGP6XTV1,BGPD30,BGPDAREA
\$\$VALI^XBDIQ1
^XBFMK
 BGP3D3,BGP4D7,BGP5D7,BGP5EL4,BGP6D7,BGP6EL4,BGPD22
 BGP3DA,BGP3DAR,BGP3DARA,BGP3DARH,BGP3DG,BGP3DL,BGP3ULF
 BGP3UTL,BGP4DA,BGP4DAP,BGP4DAR,BGP4DARA,BGP4DG,BGP4DHE
 BGP4DHEA,BGP4DL,BGP4DPP,BGP4HUTL,BGP4ULF,BGP4UTL,BGP5CL
 BGP5CTL,BGP5DAP,BGP5DAR,BGP5DEL,BGP5DELA,BGP5DGPA
 BGP5DGPU,BGP5DHE
 BGP5DHEA,BGP5DL,BGP5DNG,BGP5DPP,BGP5ETL,BGP5EUTL
 BGP5GAEL,BGP5GAGP,BGP5GAHE,BGP5GALL,BGP5GAPU,BGP5GCMS
 BGP5GCOM,BGP5GELD,BGP5GGPU,BGP5GHED,BGP5GLST,BGP5GNPL
 BGP5GNTL,BGP5GPP,BGP5GUPL,BGP5HUTL,BGP5NPL,BGP5TXL
 BGP5ULF,BGP5UTL,BGP6AHW
 BGP6CL,BGP6CTL,BGP6DAP,BGP6DAR,BGP6DCHA,BGP6DCL,BGP6DEL
 BGP6DELA,BGP6DGPA,BGP6DGPU,BGP6DHE,BGP6DHEA,BGP6DL
 BGP6DLT,BGP6DMT,BGP6DNG,BGP6DPP,BGP6ETL,BGP6EUTL
 BGP6GACW,BGP6GAEL,BGP6GAGP,BGP6GAHE,BGP6GALL,BGP6GAPU
 BGP6GCMP,BGP6GCMS,BGP6GCOM
 BGP6GELD,BGP6GGPU,BGP6GHED,BGP6GLST,BGP6GLTX,BGP6GMTX
 BGP6GNPL,BGP6GNTL,BGP6GPP,BGP6GUPL,BGP6HUTL,BGP6NPL
 BGP6ULF,BGP6UTL,BGP6XTL,BGPD,BGPDAREA,BGPDARP,BGPDFTA
 BGPDUP,BGPOST,BGPOST1
^XBGSAVE
 BGP3PDL,BGP3UTL,BGP4HEL,BGP4HUTL,BGP4PDL,BGP4UTL,BGP5ELL
 BGP5EUTL,BGP5HEL,BGP5HUTL,BGP5PDL,BGP5UTL,BGP6DCHA
 BGP6ELL,BGP6EUTL,BGP6HEL,BGP6HUTL,BGP6PDL,BGP6UTL
 BGPDARP,BGPDFTA
^XBKVAR
 BGP5GAEL,BGP5GAGP,BGP5GAHE,BGP5GALL,BGP5GAPU,BGP5GCMS
 BGP5GCOM,BGP5GELD,BGP5GGPU,BGP5GHED,BGP5GLST,BGP5GNPL
 BGP5GNTL,BGP5GPP,BGP5GUPL,BGP6GACW,BGP6GAEL,BGP6GAGP
 BGP6GAHE,BGP6GALL,BGP6GAPU,BGP6GCMP,BGP6GCMS,BGP6GCOM
 BGP6GELD,BGP6GGPU
 BGP6GHED,BGP6GLST,BGP6GLTX,BGP6GMTX,BGP6GNPL,BGP6GNTL
 BGP6GPP,BGP6GUPL
GUIR^XBLM
 BGP5GAEL,BGP5GAGP,BGP5GAHE,BGP5GALL,BGP5GAPU,BGP5GCMS
 BGP5GCOM,BGP5GELD,BGP5GGPU,BGP5GHED,BGP5GNPL,BGP5GNTL
 BGP5GPP,BGP5GU,BGP6GAEL,BGP6GAGP,BGP6GAHE,BGP6GALL
 BGP6GAPU,BGP6GCMP,BGP6GCMS,BGP6GCOM,BGP6GELD,BGP6GGPU

EN^XBVK	BGP6GHED,BGP6GLTX BGP6GMTX,BGP6GNPL,BGP6GNTL,BGP6GPP,BGP6GU,BGPGU BGP3DA,BGP3DAR,BGP3DARA,BGP3DARH,BGP3DG,BGP3DL,BGP3ULF BGP4DA,BGP4DAP,BGP4DAR,BGP4DARA,BGP4DG,BGP4DHE,BGP4DHEA BGP4DL,BGP4DPP,BGP4LSTF,BGP4ULF,BGP5CL,BGP5CTL,BGP5CTS BGP5DAP,BGP5DAR,BGP5DEL,BGP5DELA,BGP5DGPA,BGP5DGPU BGP5DHE,BGP5DHEA BGP5DL,BGP5DNG,BGP5DPP,BGP5ETL,BGP5ETS,BGP5GAEL,BGP5GAGP BGP5GAHE,BGP5GALL,BGP5GAPU,BGP5GCMS,BGP5GCOM,BGP5GELD BGP5GGPU,BGP5GHED,BGP5GLST,BGP5GNPL,BGP5GNTL,BGP5GPP BGP5GR,BGP5GRA,BGP5GU,BGP5GUPL,BGP5LSTF,BGP5NPL,BGP5TXL BGP5TXS,BGP5TXV BGP5TXV1,BGP6AHW,BGP6CL,BGP6CTL,BGP6CTS,BGP6DAP,BGP6DAR BGP6DCHA,BGP6DCL,BGP6DEL,BGP6DELA,BGP6DGPA,BGP6DGPU BGP6DHE,BGP6DHEA,BGP6DL,BGP6DLT,BGP6DMT,BGP6DNG,BGP6DPP BGP6ETL,BGP6GACW,BGP6GAEL,BGP6GAGP,BGP6GAHE,BGP6GALL BGP6GAPU,BGP6GCMP
EN^XBVK	BGP6GCMS,BGP6GCOM,BGP6GELD,BGP6GGPU,BGP6GHED,BGP6GLST BGP6GLTX,BGP6GMTX,BGP6GNPL,BGP6GNTL,BGP6GPP,BGP6GR BGP6GRA,BGP6GU,BGP6GUPL,BGP6LSTF,BGP6NPL,BGP6XTL,BGP6XTS BGP6XTV,BGP6XTV1,BGPD,BGPDAREA,BGPDARP,BGPDFTA,BGPDUP BGPGR,BGPGRA,BGPGU BGPOST
\$\$DT^XLFDT	BGP5GAEL,BGP5GAGP,BGP5GAHE,BGP5GALL,BGP5GAPU,BGP5GCMS BGP5GCOM,BGP5GELD,BGP5GGPU,BGP5GHED,BGP5GLST,BGP5GNPL BGP5GNTL,BGP5GPP,BGP5GUPL,BGP6GACW,BGP6GAEL,BGP6GAGP BGP6GAHE,BGP6GALL,BGP6GAPU,BGP6GCMP,BGP6GCMS,BGP6GCOM BGP6GELD,BGP6GGPU BGP6GHED,BGP6GLST,BGP6GLTX,BGP6GMTX,BGP6GNPL,BGP6GNTL BGP6GPP,BGP6GUPL
\$\$FMADD^XLFDT	BGP3BAN,BGP3D1,BGP3D10,BGP3D3,BGP3D4,BGP3D5,BGP3D6 BGP3D7,BGP3D8,BGP3D9,BGP4BAN,BGP4D1,BGP4D10,BGP4D21 BGP4D3,BGP4D31,BGP4D4,BGP4D41,BGP4D5,BGP4D6,BGP4D61 BGP4D7,BGP4D8,BGP4D9,BGP5ASL,BGP5BAN,BGP5C1,BGP5C11 BGP5C12,BGP5C13,BGP5C3,BGP5CL BGP5D1,BGP5D10,BGP5D22,BGP5D3,BGP5D31,BGP5D4,BGP5D41 BGP5D5,BGP5D6,BGP5D61,BGP5D7,BGP5D71,BGP5D72,BGP5D73 BGP5D8,BGP5D9,BGP5DEL,BGP5DELA,BGP5DGPA,BGP5DGPU BGP5DHE1,BGP5EL3,BGP5GAEL,BGP5GAHE,BGP5GAPU,BGP5GCMS BGP5GELD,BGP5GGPU,BGP5GNPL BGP5NPL,BGP5UTL1,BGP6AHW,BGP6ASL,BGP6BAN,BGP6C1,BGP6C11 BGP6C12,BGP6C13,BGP6C3,BGP6CL,BGP6D1,BGP6D10,BGP6D21 BGP6D22,BGP6D23,BGP6D3,BGP6D31,BGP6D4,BGP6D41,BGP6D42 BGP6D5,BGP6D51,BGP6D52,BGP6D6,BGP6D61,BGP6D7,BGP6D71 BGP6D72,BGP6D73 BGP6D8,BGP6D9,BGP6DCL,BGP6DEL,BGP6DELA,BGP6DGPA,BGP6DGPU BGP6DHE1,BGP6DU,BGP6EL3,BGP6EL4,BGP6GAEL,BGP6GAHE BGP6GAPU,BGP6GCMP,BGP6GCMS,BGP6GELD,BGP6GGPU,BGP6GNPL BGP6NPL,BGP6UTL,BGP6UTL1,BGPD,BGPD1,BGPD12,BGPD14,BGPD1B BGPD2,BGPD23

\$\$FMDIFF^XLFDT BGPD24,BGPD29,BGPD3,BGPD4,BGPD5,BGPD6,BGPD7,BGPD8,BGPDA
BGPDB,BGPDC,BGPDD,BGPDU
BGP3D1,BGP3D3,BGP3D9,BGP4D1,BGP4D10,BGP4D7,BGP4D9,BGP5D1
BGP5D10,BGP5D32,BGP5D33,BGP5D34,BGP5D35,BGP5D7,BGP5D72
BGP5D9,BGP6D1,BGP6D10,BGP6D32,BGP6D33,BGP6D34,BGP6D35
BGP6D52,BGP6D7,BGP6D72,BGP6D9,BGPD2,BGPD22
BGP3D1,BGP3D10,BGP3D2,BGP3D3,BGP3D4,BGP3D5,BGP3D6,BGP3D7
BGP3D8,BGP3D9,BGP3DA,BGP3DAR,BGP3DARA,BGP3DARH,BGP3DG
BGP3DH,BGP3DL,BGP3DPH,BGP3DS,BGP3PDL,BGP3PDLS,BGP4D1
BGP4D10,BGP4D2,BGP4D21,BGP4D3,BGP4D31,BGP4D4,BGP4D41
BGP4D5,BGP4D6
BGP4D61,BGP4D7,BGP4D71,BGP4D8,BGP4D9,BGP4DA,BGP4DAP
BGP4DAR,BGP4DARA,BGP4DG,BGP4DH,BGP4DHE,BGP4DHEA,BGP4DL
BGP4DPH,BGP4DPP,BGP4DS,BGP4HEH,BGP4HEHH,BGP4HEL,BGP4HES
BGP4HESL,BGP4PDL,BGP4PDHL,BGP4PDLS,BGP5C1,BGP5C11
BGP5C13,BGP5C2,BGP5C3
BGP5CL,BGP5CP,BGP5D1,BGP5D10,BGP5D2,BGP5D21,BGP5D22
BGP5D3,BGP5D31,BGP5D4,BGP5D41,BGP5D5,BGP5D6,BGP5D61
BGP5D7,BGP5D71,BGP5D72,BGP5D8,BGP5D9,BGP5DAP,BGP5DAR
BGP5DEL,BGP5DELA,BGP5DGPA,BGP5DGPU,BGP5DH,BGP5DHE
BGP5DHE1,BGP5DHEA,BGP5D
BGP5DNG,BGP5DPH,BGP5DPP,BGP5DS,BGP5EL2,BGP5EL3,BGP5EL4
BGP5ELH,BGP5ELHH,BGP5ELL,BGP5ELS,BGP5ELSL,BGP5GU,BGP5HEH
BGP5HEHH,BGP5HEL,BGP5HES,BGP5HESL,BGP5NPL,BGP5NPLD
BGP5NPLP,BGP5PDL,BGP5PDHL,BGP5PDLS,BGP6AHW,BGP6C1
BGP6C11,BGP6C13,BGP6C2
BGP6C3,BGP6CL,BGP6CP,BGP6D1,BGP6D10,BGP6D2,BGP6D21
BGP6D22,BGP6D23,BGP6D3,BGP6D31,BGP6D4,BGP6D41,BGP6D5
BGP6D51,BGP6D52,BGP6D6,BGP6D61,BGP6D7,BGP6D71,BGP6D72
BGP6D8,BGP6D9,BGP6DAP,BGP6DAR,BGP6DCHA,BGP6DCHW,BGP6DCL
BGP6DCLD,BGP6DCLP
BGP6DEL,BGP6DELA,BGP6DGPA,BGP6DGPU,BGP6DH,BGP6DHE
BGP6DHE1,BGP6DHEA,BGP6DL,BGP6DLT,BGP6DMT,BGP6DNG,BGP6DPH
BGP6DPP,BGP6DS,BGP6DU,BGP6EL2,BGP6EL3,BGP6EL4,BGP6ELH
BGP6ELHH,BGP6ELL,BGP6ELS,BGP6ELSL,BGP6GU,BGP6HEH
BGP6HEHH,BGP6HEL,BGP6HES
BGP6HESL,BGP6NPL,BGP6NPLD,BGP6NPLP,BGP6PDL,BGP6PDHL
BGP6PDLS,BGPAPH,BGPD,BGPD1,BGPD12,BGPD13,BGPD14,BGPD1B
BGPD2,BGPD22,BGPD24,BGPD29,BGPD3,BGPD31,BGPD4,BGPD5
BGPD6,BGPD7,BGPDA,BGPDAREA,BGPDARP,BGPDB,BGPDD,BGPDFTA
BGPDH,BGPDPH,BGPDS
BGPDU,BGPGU
\$\$NOW^XLFDT BGP3UTL,BGP4DA,BGP4DAR,BGP4DARA,BGP4DG,BGP5DAR,BGP5DNG
BGP5GAEL,BGP5GAGP,BGP5GAHE,BGP5GALL,BGP5GAPU,BGP5GCMS
BGP5GCOM,BGP5GELD,BGP5GGPU,BGP5GHED,BGP5GLST,BGP5GNPL
BGP5GNTL,BGP5GPP,BGP5GUPL,BGP6DAR,BGP6DCHA,BGP6DNG
BGP6GACW,BGP6GAEL
BGP6GAGP,BGP6GAHE,BGP6GALL,BGP6GAPU,BGP6GCMP,BGP6GCMS
BGP6GCOM,BGP6GELD,BGP6GGPU,BGP6GHED,BGP6GLST,BGP6GLTX
BGP6GMTX,BGP6GNPL,BGP6GNTL,BGP6GPP,BGP6GUPL

\$\$CJ^XLFSTR	BGPPOST
\$\$REPEAT^XLFSTR	BGP6DLT,BGP6DMT
\$\$STRIP^XLFSTR	BGP4DP1,BGP4HEL1,BGP4HEP1,BGP4PDL3,BGP5DP1,BGP5ELL1 BGP5ELLN,BGP5ELP1,BGP5HEL1,BGP5HEP1,BGP5PDL3,BGP6D23 BGP6DCHW,BGP6DP1,BGP6EL2,BGP6ELL1,BGP6ELLN,BGP6ELP1 BGP6HEL1,BGP6HEP1,BGP6PDL3
\$\$UP^XLFSTR	BGP3D2,BGP3ULF,BGP4D2,BGP4ULF,BGP5C1,BGP5C11,BGP5C12 BGP5C13,BGP5C2,BGP5C3,BGP5D2,BGP5D72,BGP5EL2,BGP5ULF BGP6C1,BGP6C11,BGP6C12,BGP6C13,BGP6C2,BGP6C3,BGP6D2 BGP6D72,BGP6EL2,BGP6ULF,BGPD2,BGPD4
^XMD	BGP50P1
MES^XPDUML	BGPPOST
DISP^XQORM1	BGP3DSI,BGP3DSL,BGP4DHSL,BGP4DSI,BGP4DSL,BGP5ASL,BGP5CSI BGP5CTL,BGP5CTS,BGP5DESI,BGP5DESL,BGP5DHSL,BGP5DSI BGP5DSL,BGP5ETL,BGP5ETS,BGP5TXL,BGP5TXS,BGP5TXV,BGP5TXV1 BGP6ASL,BGP6CSI,BGP6CTL,BGP6CTS,BGP6DESI,BGP6DESL BGP6DHSL,BGP6DSI
DISP^XQORM1	BGP6DSL,BGP6ETL,BGP6XTL,BGP6XTS,BGP6XTV,BGP6XTV1,BGPD0 BGPDL
\$\$NOJOURN^ZIBGCHAR	BGP3D1,BGP4HUTL,BGP4UTL,BGP5EUTL,BGP5HUTL,BGP5UTL BGP6EUTL,BGP6HUTL,BGP6UTL,BGPD1

6.2 Callable Routines

There are no published entry points in this package.

6.3 Exported Options

Option Name	Description
BGP 02 MENU	GPRA+ FY02
BGP 03 AREA AREA	Run AREA Annual Clinical Report (Area Only)
BGP 03 AREA GPRA	Run AREA GPRA Report (to be used at Area only)
BGP 03 AREA MENU	Area Options
BGP 03 EDIT SITE	
PARAMETERS	Site Parameters
BGP 03 MAIN MENU	GPRA+ FY03
BGP 03 MENU REPORTS	Reports
BGP 03 MENU SETUP	System Setup
BGP 03 RPT AREA PERF	Run Area Annual Performance Report for Area Export
BGP 03 RPT GPRA AREA	
EXPT	Run GPRA Report for Area Export
BGP 03 RPT LOCAL USE	Run Indicator Reports for Local Use
BGP 03 TAXONOMY	
CHECK	Check for Taxonomies Required by the GPRA Report
BGP 03 TAXONOMY	
SETUP	Taxonomy Setup
BGP 03 UPLOAD TO	
AREA	Upload Report Files from Site
BGP 04 AREA AREA	Run AREA Area Director Report (Area Only)

Option Name	Description
BGP 04 AREA DIR REPORT	Area Director's Performance Report
BGP 04 AREA GPRA	Run AREA GPRA Report
BGP 04 AREA HEDIS	Run Area HEDIS Report
BGP 04 AREA MENU	Area Office Options
BGP 04 HEDIS AREA REPORT	Run AREA HEDIS Report
BGP 04 HEDIS REPORT	HEDIS Performance Report
BGP 04 LIST FILES	List files in a directory
BGP 04 MAIN MENU	GPRA+ FY04
BGP 04 REPORT GPRA	GPRA Performance Report
BGP 04 REPORTS MENU	Reports
BGP 04 SEL ALL PATS	Selected Indicators with No Community Defined
BGP 04 SEL PATIENT PANEL	Selected Indicators w/Patient Panel Population
BGP 04 SELECTED IND REPORT	Selected Indicators w/Community Taxonomy Defined
BGP 04 SITE PARAMETERS	Site Parameters
BGP 04 SYSTEM SETUP	System Setup
BGP 04 TAXONOMY CHECK	Check for Taxonomies Required by FY 04 GPRA Report
BGP 04 TAXONOMY SETUP	Taxonomy Setup
BGP 04 UPLOAD FILES	Upload Report Files from Site
BGP 05 AREA ELDER REPORT	Run AREA Elder Care Report
BGP 05 AREA GPRA	Run AREA National GPRA Report
BGP 05 AREA GPU	Run AREA GPRA Performance Report
BGP 05 AREA HEDIS	Run Area HEDIS Report
BGP 05 AREA MENU	Area Options
BGP 05 CMS PERF RPT	CMS Performance Report
BGP 05 CMS TAXONOMY CHECK	Taxonomy Check-CMS Report
BGP 05 CMS TAXONOMY SETUP	Taxonomy Setup-CMS Report
BGP 05 CRS TAXONOMY SETUP	Taxonomy Setup-All CRS Reports
BGP 05 ELDER REPORT	Elder Care Report
BGP 05 ELDER TAXONOMY CHECK	Taxonomy Check-Elder Care Report
BGP 05 GPU GPRA PERFORMANCE	GPRA Performance Report
BGP 05 HEDIS AREA REPORT	Run AREA HEDIS Report
BGP 05 HEDIS REPORT	HEDIS Performance Report
BGP 05 HEDIS TAXONOMY CHECK	Taxonomy Check-HEDIS Report
BGP 05 HOSPITAL REPORTS	Hospital Reports
BGP 05 LIST FILES	List files in a directory
BGP 05 LOCAL REPORTS	Reports for Local Use: IHS Clinical Indicators

Option Name	Description
MENU	
BGP 05 LOCAL	
TAXONOMY CHECK	Taxonomy Check-Local CRS Reports
BGP 05 MAIN MENU	CRS 2005
BGP 05 MENU NATIONAL	National GPRA Reports
BGP 05 NATIONAL GPRA	
REPORT	National GPRA Report
BGP 05 NATIONAL PAT	
LISTS	National GPRA Report Patient List
BGP 05 NGR TAXONOMY	Taxonomy Check-National GPRA/GPRA Performance
CHECK	Rpts
BGP 05 NGR TAXONOMY	Taxonomy Setup-National GPRA/GPRA Performance
SETUP	Rpts
BGP 05 OTHER CMS	
REPORTS	Other CMS Reports
BGP 05 OTHER NTL RPTS	
MENU	Other National Reports
BGP 05 REPORTS MENU	Reports
BGP 05 SEL ALL PATS	Selected Indicators w/All Communities
BGP 05 SEL PATIENT	
PANEL	Selected Indicators w/Patient Panel Population
BGP 05 SELECTED IND	
REPORT	Selected Indicators w/Community Specified
BGP 05 SITE	
PARAMETERS	Site Parameters
BGP 05 SYSTEM SETUP	System Setup
BGP 05 TAXONOMY	
CHECK	Taxonomy Check
BGP 05 TAXONOMY	
SETUP	Taxonomy Setup
BGP 05 UPLOAD FILES	Upload Report Files from Site
BGP 05 VIEW ANY	
TAXONOMY	View All CRS Taxonomies (Site-Defined&Hard Coded)
BGP 06 AREA	
CHILDHOOD WT/HT	Run AREA Childhood Height and Weight Data File
BGP 06 AREA ELDER	
REPORT	Run AREA Elder Care Report
BGP 06 AREA GPRA	Run AREA National GPRA Report
BGP 06 AREA GPU	Run AREA GPRA Performance Report
BGP 06 AREA HEDIS	Run Area HEDIS Report
BGP 06 AREA MENU	Area Options
BGP 06 CMS PERF RPT	CMS Performance Report
BGP 06 CMS TAXONOMY	
CHECK	Taxonomy Check-CMS Report
BGP 06 CMS TAXONOMY	
SETUP	Taxonomy Setup-CMS Report
BGP 06	
COMPREHENSIVE PAT	
LIST	Comprehensive National GPRA Patient List
BGP 06 CRS TAXONOMY	
SETUP	Taxonomy Setup-All CRS Reports
BGP 06 ELDER REPORT	Elder Care Report
BGP 06 ELDER	Taxonomy Check-Elder Care Report

Option Name	Description
TAXONOMY CHECK	
BGP 06 GPU GPRA	
PERFORMANCE	GPRA Performance Report
BGP 06 HEDIS AREA	
REPORT	Run AREA HEDIS Report
BGP 06 HEDIS REPORT	HEDIS Performance Report
BGP 06 HEDIS	
TAXONOMY CHECK	Taxonomy Check-HEDIS Report
BGP 06 HOSPITAL	
REPORTS	Hospital Reports
BGP 06 LAB TAX REPORT	Lab Taxonomy Report
BGP 06 LIST FILES	List files in a directory
BGP 06 LOCAL REPORTS	
MENU	Reports for Local Use: IHS Clinical Measures
BGP 06 LOCAL	
TAXONOMY CHECK	Taxonomy Check-Local CRS Reports
BGP 06 MAIN MENU	CRS 2006
BGP 06 MED TAX	
REPORT	Medication Taxonomy Report
BGP 06 MENU NATIONAL	National GPRA Reports
BGP 06 NATIONAL GPRA	
REPORT	National GPRA Report
BGP 06 NATIONAL PAT	
LISTS	National GPRA Report Patient List
BGP 06 NGR TAXONOMY	Taxonomy Check-National GPRA/GPRA Performance
CHECK	Rpts
BGP 06 NGR TAXONOMY	Taxonomy Setup-National GPRA/GPRA Performance
SETUP	Rpts
BGP 06 OTHER CMS	
REPORTS	Other CMS Reports
BGP 06 OTHER NTL RPTS	
MENU	Other National Reports
BGP 06 REPORTS MENU	Reports
BGP 06 SEL ALL PATS	Selected Measures w/All Communities
BGP 06 SEL PATIENT	
PANEL	Selected Measures w/Patient Panel Population
BGP 06 SELECTED IND	
REPORT	Selected Measures w/Community Specified
BGP 06 SITE	
PARAMETERS	Site Parameters
BGP 06 SYSTEM SETUP	System Setup
BGP 06 TAXONOMY	
CHECK	Taxonomy Check
BGP 06 TAXONOMY	
REPORTS	Taxonomy Reports
BGP 06 TAXONOMY	
SETUP	Taxonomy Setup
BGP 06 UPLOAD FILES	Upload Report Files from Site
BGP 06 VIEW ANY	
TAXONOMY	View All CRS Taxonomies (Site-Defined&Hard Coded)
BGP FAC TO AREA RPT	
PERIOD	Run Area Annual Report for Export to Area
BGP GPRA AREA	Run AREA GPRA Report (to be used at Area only)

Option Name	Description
REPORT	
BGP GPRA FAC TO AREA	
RPT	Run FY Report for Export to Area
BGP GPRA LOCAL	
REPORT	Run GPRA Report for Local Use
BGP GPRA UPLOAD FILE	Upload GPRA Data File from Site
BGP INDICATOR 32	
TESTS	Enter/Edit Indicator 32 Lab Tests
BGP TAXONOMY CHECK	Check for Taxonomies Required by the GRPA Report
BGPGRPC	BGPG Application Context
BGPMENU	IHS Clinical Reporting System (CRS) Main Menu

7.0 Internal Relations

All users should be given the access to the appropriate options and keys to them, as needed. All of the options in this system stand alone.

8.0 How to Generate Online Documentation

The file number range for this package is 90240 and 90370. The namespace is BGP. All templates, routines, screen forms, etc. begin with BGP.

This section describes some of the methods by which users can generate IHS RPMS CLINICAL Reporting system technical documentation. Online technical documentation pertaining to the IHS RPMS CLINICAL REPORTING Reporting software, in addition to that which is located in the help prompts and on the help screens throughout the IHS RPMS CLINICAL Reporting package, can be generated through the use of several Kernel options. These include, but are not limited to, the following:

- %INDEX
- Menu Management
- Inquire Option
- Print Option File
- VA FileMan
- Data Dictionary Utilities
- List File Attributes

Entering question marks at the “Select...Option” prompts can also provide users with valuable technical information. For example, a single question mark (?) lists all options that can be accessed from the current option. Entering two question marks (??) lists all options accessible from the current one, showing the formal name and lock for each. Three question marks (???) displays a brief description for each option in a menu, whereas an option name preceded by a question mark (?OPTION) shows extended help, if available, for that option.

For a more exhaustive option listing and further information about other utilities that supply online technical information, please consult the DHCP Kernel Reference manual.

8.1 %INDEX

This option analyzes the structure of a routine to determine in part if the routine adheres to RPMS Programming Standards. The %INDEX output can include the following components:

- Compiled list of errors and warnings
- Routine listing
- Local variables
- Global variables
- Naked globals

- Label references
- External references

By running %INDEX for a specified set of routines, you are afforded the opportunity to discover any deviations from RPMS Programming Standards that exist in the selected routines and to see how routines interact with one another (i.e., which routines call or are called by other routines).

To run %INDEX for the IHS RPMS CLINICAL Reporting package, specify the BGP namespace at the Routine(s)?> prompt.

8.2 Inquire Option

This menu management option provides the following information about a specified option:

- Option name
- Menu text
- Option description
- Type of option
- Lock (if any)

In addition, all items on the menu are listed for each menu option. To secure information about IHS RPMS CLINICAL REPORTING Reporting options, you must specify the BGP namespace.

8.3 Print Option File

This utility generates a listing of options from the Option file (#19). You can choose to print all of the entries in this file or you can specify a single option or range of options. For a list of IHS RPMS CLINICAL Reporting options, please refer to the Exported Options section of this manual.

8.4 List File Attributes

This VA FileMan option allows you to generate documentation pertaining to files and file structure. Using the Standard format of this option yields the following data dictionary information for a specified file:

- File name and description
- Identifiers
- Cross-references
- Files pointed to by the file specified
- Files that point to the file specified input, print, and sort templates

In addition, the following applicable data is supplied for each field in the file:

- Field name, number, title, and description
- Global location
- Help prompt
- Cross-references
- Input transform
- Date last edited
- Notes

Using the Global Map format of this option generates an output that lists the following information:

- All cross-references for the file selected
- Global location of each field in the file
- Input, print, and sort templates

For a comprehensive listing of CLINICAL package files, please refer to the files section of this manual.

9.0 SAC Exemptions

This program uses a standard global for export of data to the Area office. The global is called ^BGPDATA(and is killed at the root level in various routines.

The kill of these unsubscripted globals has been approved by the Standards and Conventions committee.

10.0 Glossary

Archiving	The storing of historical or little-used data off-line (often on tape).
Banner	A line of text with a user's name and domain.
Browser	An interactive application that displays ASCII text on a terminal that supports a scroll region. The text can be in the form of a word-processing field or sequential local or global array. The user is allowed to navigate freely within the document.
Callable Entry Points	Places in a routine that can be called from an application program.
Cross-reference	An indexing method whereby files can include pre-sorted lists of entries as part of the stored database. Cross-references (x-refs) facilitate look-up and reporting.
Entry Point	Entry point within a routine that is referenced by a "DO" or "GOTO" command from a routine internal to a package.
File	A set of related records or entries treated as a single unit.
FileMan	The database management system for RPMS.
Global	In MUMPS, global refers to a variable stored on disk (global variable) or the array to which the global variable may belong (global array).
ICD	International Classification of Diseases.
INDEX (%INDEX)	A Kernel utility used to verify routines and other MUMPS code associated with a package. Checking is done according to current ANSI MUMPS standards and RPMS programming standards. This tool can be invoked through an option or from direct mode (>D ^%INDEX).
IRM	Information Resource Management. The IHS personnel responsible for information systems management and security.

Kernel	The set of MUMPS software utilities that function as an intermediary between the host operating system and application packages, such as Laboratory and Pharmacy. The Kernel provides a standard and consistent user and programmer interface between application packages and the underlying MUMPS implementation. These utilities provide the foundation for RPMS.
Menu	A list of choices for computing activity. A menu is a type of option designed to identify a series of items (other options) for presentation to the user for selection. When displayed, menu-type options are preceded by the word “Select” and followed by the word “option” as in Select Menu Management option: (the menu’s select prompt).
Namespace	A unique set of 2 to 4 alpha characters that are assigned by the database administrator to a software application.
Option	An entry in the Option file. As an item on a menu, an option provides an opportunity for users to select it, thereby invoking the associated computing activity. Options may also be scheduled to run in the background, non-interactively, by TaskMan.
Patient Care Component (PCC)	The central repository for data in the Resource and Patient Management System (RPMS).
Queuing	Requesting that a job be processed at a later time rather than within the current session.
Routine	A program or sequence of instructions called by a program that may have some general or frequent use. MUMPS routines are groups of program lines that are saved, loaded, and called as a single unit via a specific name.
UCI	User Class Identification: a computing area.
Up-Hat (^)	A circumflex, also known as a “hat” or “caret,” that is used as a piece delimiter in a global. The up-hat is denoted as “^” and is typed by pressing Shift+6 on the keyboard.
Utility	A callable routine line tag or function. A universal routine usable by anyone.

Variable	A character or group of characters that refers to a value. MUMPS recognizes 3 types of variables: local variables, global variables, and special variables. Local variables exist in a partition of the main memory and disappear at sign-off. A global variable is stored on disk, potentially available to any user. Global variables usually exist as parts of global arrays.
----------	--

11.0 Contact Information

If you have any questions or comments regarding this distribution, please contact the OIT Help Desk by:

Phone: (505) 248-4371 or

(888) 830-7280

Fax: (505) 248-4363

Web: <http://www.rpms.ihs.gov/TechSupp.asp>

Email: ITSCHelp@mail.ihs.gov