iCare Tab in Electronic Health Record

In the Electronic Health Record (EHR), the Clinical Application Coordinators (CAC) are able to customize their view by adding a tab that will take them to a program outside the EHR. This function can be utilized to add a link to the iCare application. These instructions are intended for the CACs. If you'd like to see this tab in the EHR, please see your CAC.

Instructions:

- 1. Go to the Design Mode by right clicking in the blue tool bar at the top of the screen
- 2. Right click to the right of the tabs along the top and click on **Add Object.** Expand the options under Function, expand the options under **Tools** and go to **Program Launcher.**
- 3. Click the <Add> button



4. A plain blank square will appear on your screen, right click on it and go to Properties.

≱ IHS+EHR	~SCHILDT,HUNTER~	FORT HALL INDIAN HEA	ILTH SERVICE	_ 0 🗙
User Design	Patient Tools Help			
Patient Chart	eMaxx Portneuf Up-To-Date	Programs IRIS RPMS Clin	Pharm Communication FTH Intranet E-Mail	
2	RCIS/BH			
	Scheduling GUI	1-1		
	Narcotic Flow Sheet			
<u>ب</u>	Visual CPHAD			
SCHILDT, HU	INTER FORT-HALL-HC.PRT.IHS.	GOV FT HALL H CT	11Jun-2007 09:23	

5. Select the link that will launch you to the current version iCare by clicking on the 3 dotted lines (...) behind the EXENAME line

Properties for Program Launcher					
Property	Value				
ТОР	0				
LEFT	952				
HEIGHT	51				
WIDTH	89				
ALIGN	None				
ANCHORS	🔽 Top; Left 🗸 🗸				
CAPTION	Click here to launch iCare				
COLOR	×				
COMMANDLINE					
EXENAME	h Service\iCare Version 1.1.33\bgi-iCare.exe				
GLYPH					
LAYOUT	Image appears at the top.				
NUMGLYPHS	1				
EXENAME Name of executable file.					
ОК	Cancel Apply				

NOTE: These programs are normally located on your local PC. The exact path must be the same for every user. Like in the example the only users that will be able to use this shortcut are the ones with Microsoft Office 2003 installed since the path to the file is different from other versions. In such a case you may make another shortcut to get around that problem.

Open					? 🗙
Look jn:	🚞 iCare Version -	1.1.33	~ () 🦻 📂 🖽	•
My Recent Documents Desktop My Documents	 bqi-iCare.exe bqi-ultraGridToc care_Help.chm vcImage.ocx Infragistics2.Wi 	ıBar.dll n.v7.1.dll n.ultraWinListBar.v7.1.dll n.ultraWinGrid.ExcelExpor n.ultraWinTabbedMdi.v7.1 n.ultraWinTabbedMdi.v7.1 n.ultraWinTabcManager.v7 n.ultraWinToolbars.v7.1.cll n.ultraWinToolbars.v7.1.cl n.ultraWinTabControl.v7.1 n.ultraWinTabControl.v7.1.dl	t.v7.1.dll log.v7.1.dll dll '.1.dll ill dll L.dll I	 Infragistics2. Infragistics2. Infragistics2. Infragistics2. Infragistics2. Infragistics2. RichTextBoxI BMXNet20.dl bqi-iCare.exe MAIN.ico 	Win, UltraWinEditi Win, UltraWinData Win, UltraWinGrid Win, Misc. v7. 1. dll Shared. v7. 1. dll Excel. v7. 1. dll PrintControl. dll I a. config
	<]		>
	File <u>n</u> ame:	bqi-iCare.exe		~	<u>O</u> pen
My Network	Files of type:			~	Cancel

6. Enter the color and the caption (e.g. Click here to launch iCare) of your tab in the respective rows

Properties for Program Launcher					
Property	Value				
ТОР	0				
LEFT	952				
HEIGHT	51				
WIDTH	89				
ALIGN	None				
ANCHORS	🔽 Top; Left 💽				
CAPTION	Click here to launch iCare				
COLOR					
COMMANDLINE					
EXENAME	:h Service\iCare				
GLYPH					
LAYOUT	Image appears				
NUMGLYPHS					
	More Colors				
COLOR Background color of the control.					
ОК	Cancel Apply				

- 7. Click on the <Apply> button
- 8. Click on the <OK> button
- 9. Drag your tab to the desired location on your HER
- 10. Save your template by clicking on the **Design/Save as Template** menu option in the upper left corner of the screen.

Save As Template		×
Template:	Application 🔲 User	
GEBREMARIAM		*
%SPRTEST %TEST		^
GEBREMARIAM HEADER		

11. Exit Design Mode