



RESOURCE AND PATIENT MANAGEMENT SYSTEM

Clinical Reporting System (CRS) (BGP)

Technical Manual

Version 7.0
January 2007

Office of Information Technology (OIT)
Division of Information Resources
Albuquerque, New Mexico

PREFACE

This manual contains the technical manual for the CRS Clinical Reporting System version 7.0, which adds FY 2007 clinical performance measures to existing FY 2002 through FY 2006 measures.

The CRS Clinical Reporting System is an RPMS (Resource and Patient Management System) software application designed for national reporting as well as local and Area monitoring of clinical GPRA and developmental measures. CRS was first released for FY 2002 performance measures (as GPRA+) and is based on a design by the Aberdeen Area (GPRA2000).

The Government Performance and Results Act (GPRA) requires Federal agencies to report annually on how the agency measured up against the performance targets set in its annual Plan. IHS GPRA measures include measures for clinical prevention and treatment, quality of care, infrastructure, and administrative efficiency functions. The CRS Clinical Reporting System is the reporting tool used by the IHS Office of Planning and Evaluation to collect and report clinical performance results annually to the Department of Health and Human Services (DHHS) and to Congress.

Each year, an updated version of CRS software is released to reflect changes in the logic descriptions of the different denominators and numerators. Additional performance measures may also be added. Local facilities can run reports as often as they want to and can also use CRS to transmit data to their Area. The Area Office can use CRS to produce an aggregated Area report for either annual GPRA or Area Director Performance reports.

The CRS Clinical Reporting System will produce reports on demand from local RPMS databases for both GPRA and developmental clinical performance measures that are based on RPMS data. CRS is intended to eliminate the need for manual chart audits for evaluating and reporting clinical measures. Administrative and clinical users will be able to review individual or all measures at any time, and can:

- Identify potential data issues in their RPMS, i.e., missing or incorrect data.
- Monitor their site's performance against past national performance and upcoming agency goals.
- Identify specific areas where the facility is not meeting the measure in order to initiate business process or other changes.
- Quickly measure impact of process changes on performance measures.
- Identify IHS Areas meeting or exceeding measures to provide lessons learned.

To produce reports with comparable data across every facility, the GPRA measure definition was "translated" into programming code with the assistance of clinical subject matter experts. CRS uses pre-defined taxonomies to find data items in PCC to determine if a patient meets the performance measure criteria. Taxonomies contain

groups of codes (e.g., diagnoses or procedures) or site-specific terms. Each performance measure topic has one or more denominators and numerators defined.

CRS is intended for use by Area and site Quality Improvement staff, Compliance Officers, GPRA Coordinators, clinical staff such as physicians, nurses, nurse practitioners, and other providers, Area Directors, as well as any staff involved with quality assurance initiatives.

TABLE OF CONTENTS

1.0	Implementation and Maintenance	4
1.1	System Requirements.....	4
2.0	Security Keys	5
3.0	Routines	6
4.0	Files and Tables	56
4.1	File List	56
4.2	File Access	58
4.3	Cross References	61
4.4	Table File.....	79
5.0	Archiving and Purging	114
6.0	External Relations	115
6.1	External Calls.....	115
6.2	Callable Routines.....	198
6.3	Exported Options	198
7.0	Internal Relations	204
8.0	How to Generate Online Documentation	205
8.1	%INDEX.....	205
8.2	Inquire Option	206
8.3	Print Option File	206
8.4	List File Attributes	206
9.0	SAC Exemptions	208
10.0	Glossary	209
11.0	Contact Information	211

1.0 Implementation and Maintenance

The *Clinical Reporting System* occupies the BGP namespace. Options, security locks/keys, templates, routines, and globals are namespaced BGP.

1.1 System Requirements

- Kernel v8.0 or higher
- FileMan v21 or higher
- IHS Patient Dictionaries (AUPN) v99.1 through patch 17
- PCC Management Reports v3.0 through patch 18
- Taxonomy System v5.1 through patch 8
- XB/ZIB Utilities v3.0 through patch 11

2.0 Security Keys

The following table provides security key information.

Name	Descriptive Name	Description
BGPZAREA	BGP AREA REPORT	This key unlocks the area report options on the main BGPMENU. This key should only be assigned at the Area office level.
BGPZMENU	BGP MAIN MENU - CRS	This key unlocks the main CRS menu. It should be assigned to all users who need to run GPRA/CRS reports.
BGPZ PATIENT LISTS	BGP PATIENT LISTS	This key allows the holder to obtain patient lists. The patient lists contain sensitive patient/clinical information. This key should only be given to those users who have a need for such information.
BGPZ SITE PARAMETERS	BGP SITE PARAMETERS	This key should be given to those with authority to update the site parameters.
BGPZ TAXONOMY EDIT	BGP TAXONOMY EDIT	This key should be given only to users who will need to update taxonomies.

3.0 Routines

The following table provides information for routines and their descriptions.

Routine	Description
BGP3BAN	IHS/CMI/LAB - BANNER FOR GPRA 03
BGP3D1	IHS/CMI/LAB - calc indicators
BGP3D10	IHS/CMI/LAB - indicator J-4,J-5,J-6
BGP3D2	IHS/CMI/LAB - indicator 2
BGP3D3	IHS/CMI/LAB - indicator 3
BGP3D4	IHS/CMI/LAB - indicator 3
BGP3D5	IHS/CMI/LAB - indicator 31
BGP3D6	IHS/CMI/LAB - indicator 31
BGP3D7	IHS/CMI/LAB - indicator C
BGP3D8	IHS/CMI/LAB - indicator C
BGP3D9	IHS/CMI/LAB - indicator J
BGP3DA	IHS/CMI/LAB - IHS GPRA - report for area export
BGP3DAR	IHS/CMI/LAB - ihs area GPRA
BGP3DARA	IHS/CMI/LAB - ihs area AA
BGP3DARH	IHS/CMI/LAB - ihs area AA
BGP3DCI	IHS/CMI/LAB - IHS area GPRA
BGP3DG	IHS/CMI/LAB - IHS GPRA - report for area export
BGP3DH	IHS/CMI/LAB - cover page for gpra
BGP3DL	IHS/CMI/LAB - IHS GPRA - report for local use
BGP3DP	IHS/CMI/LAB - IHS gpra print
BGP3DP1	IHS/CMI/LAB - print ind 1
BGP3DP2	IHS/CMI/LAB - print ind 1
BGP3DP3	IHS/CMI/LAB - print ind 1
BGP3DP4	IHS/CMI/LAB - print ind 31
BGP3DP41	IHS/CMI/LAB - print ind 31
BGP3DP5	IHS/CMI/LAB - print ind H
BGP3DP50	IHS/CMI/LAB - print ind H
BGP3DP51	IHS/CMI/LAB - print ind H
BGP3DP6	IHS/CMI/LAB - print ind C-1
BGP3DP60	IHS/CMI/LAB - print ind C-1
BGP3DP7	IHS/CMI/LAB - print ind
BGP3DPH	IHS/CMI/LAB - AREA REPORT HEADER

Routine	Description
BGP3DS	IHS/CMI/LAB - IHS gpra print
BGP3DSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP3DSL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP3PARP	IHS/CMI/LAB - IHS gpra print
BGP3PDL	IHS/CMI/LAB - IHS gpra print
BGP3PDL1	IHS/CMI/LAB - print ind 1
BGP3PDL2	IHS/CMI/LAB - print ind 1
BGP3PDL3	IHS/CMI/LAB - print ind 1
BGP3PDL4	IHS/CMI/LAB - print ind 1
BGP3PDL5	IHS/CMI/LAB - print ind 1
BGP3PDL6	IHS/CMI/LAB - print ind 1
BGP3PDL7	IHS/CMI/LAB - print ind 1
BGP3PDL8	IHS/CMI/LAB - print ind 1
BGP3PDL9	IHS/CMI/LAB - print ind 1
BGP3PDL5	IHS/CMI/LAB - IHS gpra print
BGP3PDM1	IHS/CMI/LAB - print ind 1
BGP3PDM2	IHS/CMI/LAB - print ind 1
BGP3TXCH	IHS/CMI/LAB - TAXONOMY CHECK FOR FY03 GPRA REPORT
BGP3ULF	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP3UTL	IHS/CMI/LAB -
BGP3UTL1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP4BAN	IHS/CMI/LAB - BANNER FOR GPRA 04
BGP4D1	IHS/CMI/LAB - calc indicators
BGP4D10	IHS/CMI/LAB - calc indicators
BGP4D2	IHS/CMI/LAB - indicator 1,2,3,4
BGP4D21	IHS/CMI/LAB - indicator 6
BGP4D3	IHS/CMI/LAB - indicator 11
BGP4D31	IHS/CMI/LAB - indicator C
BGP4D4	IHS/CMI/LAB - indicator 3
BGP4D41	IHS/CMI/LAB - indicator 3
BGP4D5	IHS/CMI/LAB - indicator calc
BGP4D6	IHS/CMI/LAB - indicator 31
BGP4D61	IHS/CMI/LAB - indicator 31
BGP4D7	IHS/CMI/LAB - indicator 31
BGP4D71	IHS/CMI/LAB - indicator C
BGP4D8	IHS/CMI/LAB - indicator C

Routine	Description
BGP4D9	IHS/CMI/LAB - indicator J
BGP4DA	IHS/CMI/LAB - IHS AREA CLIN 04 REPORT DRIVER
BGP4DAP	IHS/CMI/LAB - IHS GPRA 04 SELECTED REPORT DRIVER ALL PATIENTS
BGP4DAR	IHS/CMI/LAB - ihs area GPRA
BGP4DARA	IHS/CMI/LAB - ihs area AA
BGP4DCHI	IHS/CMI/LAB - calculate HEDIS indicators
BGP4DCI	IHS/CMI/LAB - IHS area GPRA
BGP4DG	IHS/CMI/LAB - IHS GPRA 04 REPORT DRIVER
BGP4DH	IHS/CMI/LAB - cover page for gpra
BGP4DHE	IHS/CMI/LAB - IHS GPRA 04 REPORT DRIVER
BGP4DHE1	IHS/CMI/LAB - indicator HEDIS
BGP4DHEA	IHS/CMI/LAB - ihs area AA
BGP4DHEP	IHS/CMI/LAB - IHS HEDIS print
BGP4DHSL	IHS/CMI/LAB - FY 04 DISPLAY IND LISTS
BGP4DL	IHS/CMI/LAB - IHS GPRA 04 SELECTED REPORT DRIVER
BGP4DP	IHS/CMI/LAB - IHS gpra print
BGP4DP1	IHS/CMI/LAB - print ind 1
BGP4DP11	IHS/CMI/LAB - print ind 1
BGP4DP2	IHS/CMI/LAB - print ind 10
BGP4DP3	IHS/CMI/LAB - print ind 1
BGP4DP4	IHS/CMI/LAB - print ind 31
BGP4DP5	IHS/CMI/LAB - print ind H
BGP4DP50	IHS/CMI/LAB - print ind H
BGP4DP51	IHS/CMI/LAB - print ind H
BGP4DP6	IHS/CMI/LAB - print ind H
BGP4DP60	IHS/CMI/LAB - print ind H
BGP4DPH	IHS/CMI/LAB - AREA REPORT HEADER
BGP4DPP	IHS/CMI/LAB - IHS GPRA 04 SELECTED REPORT DRIVER ALL PATIENTS
BGP4DS	IHS/CMI/LAB - IHS gpra print
BGP4DSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP4DSL	IHS/CMI/LAB - FY 04 DISPLAY IND LISTS
BGP4DSP	IHS/CMI/LAB - IHS summary page
BGP4HEH	IHS/CMI/LAB - cover page for HEDIS
BGP4HEHH	IHS/CMI/LAB - cover page for HEDIS

Routine	Description
BGP4HEL	IHS/CMI/LAB - IHS gpra print
BGP4HEL1	IHS/CMI/LAB - print ind 1
BGP4HEP1	IHS/CMI/LAB - print ind 1
BGP4HES	IHS/CMI/LAB - IHS gpra print
BGP4HESL	IHS/CMI/LAB - IHS gpra print
BGP4HUTL	IHS/CMI/LAB -
BGP4LSTF	IHS/CMI/LAB - List 'BG04' files in pub
BGP4PARP	IHS/CMI/LAB - IHS gpra print
BGP4PDL	IHS/CMI/LAB - IHS gpra print
BGP4PDL0	IHS/CMI/LAB - print ind 1
BGP4PDL1	IHS/CMI/LAB - print ind 1
BGP4PDL2	IHS/CMI/LAB - print ind 1
BGP4PDL3	IHS/CMI/LAB - print ind 1
BGP4PDL4	IHS/CMI/LAB - print ind 1
BGP4PDL5	IHS/CMI/LAB - print ind 1
BGP4PDL7	IHS/CMI/LAB - print ind 1
BGP4PDL8	IHS/CMI/LAB - print ind 1
BGP4PDL9	IHS/CMI/LAB - print ind
BGP4PDLH	IHS/CMI/LAB - cover page for gpra del
BGP4PDLZ	IHS/CMI/LAB - IHS gpra print
BGP4PDLZ	IHS/CMI/LAB -
BGP4PHEP	IHS/CMI/LAB - IHS gpra print
BGP4TXCH	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 GPRA REPORT
BGP4TXHE	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 HEDIS REPORT
BGP4ULF	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP4UTL	IHS/CMI/LAB -
BGP4UTL1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP50P1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED 20 Dec 2004 9:23 AM
BGP5ASL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5BAN	IHS/CMI/LAB - BANNER FOR CRS 05 06 Jan 2005 9:02 AM
BGP5C1	IHS/CMI/LAB - calc CMS indicators 26 Sep 2004 11:28 AM 04 May 2005 2:38 PM
BGP5C11	IHS/CMI/LAB - calc CMS indicators
BGP5C12	IHS/CMI/LAB - calc CMS indicators 26 Sep 2004 11:28 AM
BGP5C13	IHS/CMI/LAB - calc CMS indicators 26 Sep 2004 11:28 AM
BGP5C2	IHS/CMI/LAB - calc CMS indicators 26 Sep 2004 11:28 AM

Routine	Description
BGP5C3	IHS/CMI/LAB - calc CMS indicators 26 Sep 2004 11:28 AM
BGP5CL	IHS/CMI/LAB - IHS GPRA 05 CMS REPORT DRIVER
BGP5CP	IHS/CMI/LAB - IHS gpra print
BGP5CSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5CTL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5CTS	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5CTXC	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP5D1	IHS/CMI/LAB - calc indicators 29 Apr 2005 7:38 PM
BGP5D10	IHS/CMI/LAB - calc indicators
BGP5D2	IHS/CMI/LAB - indicator 1,2,3,4
BGP5D21	IHS/CMI/LAB - indicator 6
BGP5D22	IHS/CMI/LAB - indicator I2
BGP5D3	IHS/CMI/LAB - indicator 11
BGP5D31	IHS/CMI/LAB - indicator C
BGP5D32	IHS/CMI/LAB - indicator C
BGP5D33	IHS/CMI/LAB - indicator C
BGP5D34	IHS/CMI/LAB - indicator C
BGP5D35	IHS/CMI/LAB - indicator C
BGP5D4	IHS/CMI/LAB - indicator 3
BGP5D41	IHS/CMI/LAB - indicator 3
BGP5D5	IHS/CMI/LAB - indicator calc
BGP5D6	IHS/CMI/LAB - indicator 31
BGP5D61	IHS/CMI/LAB - indicator 31
BGP5D7	IHS/CMI/LAB - indicator 31
BGP5D71	IHS/CMI/LAB - indicator C
BGP5D72	IHS/CMI/LAB - indicator 31
BGP5D73	IHS/CMI/LAB - indicator 31
BGP5D8	IHS/CMI/LAB - indicator C
BGP5D9	IHS/CMI/LAB - indicator J
BGP5DAP	IHS/CMI/LAB - IHS GPRA 04 SELECTED REPORT DRIVER ALL PATIENTS
BGP5DAR	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP5DCEI	IHS/CMI/LAB - calculate HEDIS indicators
BGP5DCHI	IHS/CMI/LAB - calculate HEDIS indicators
BGP5DCI	IHS/CMI/LAB - IHS area GPRA
BGP5DEL	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER

Routine	Description
BGP5DELA	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER
BGP5DELP	IHS/CMI/LAB - IHS HEDIS print
BGP5DESI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5DESL	IHS/CMI/LAB - FY 05 DISPLAY IND LISTS
BGP5DGPA	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP5DGPU	IHS/CMI/LAB - IHS AREA CLIN 05 REPORT DRIVER
BGP5DH	IHS/CMI/LAB - cover page for gpra
BGP5DHE	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER
BGP5DHE1	IHS/CMI/LAB - indicator HEDIS
BGP5DHEA	IHS/CMI/LAB - ihs area AA
BGP5DHEP	IHS/CMI/LAB - IHS HEDIS print
BGP5DHSL	IHS/CMI/LAB - FY 05 DISPLAY IND LISTS
BGP5DL	IHS/CMI/LAB - IHS GPRA 05 SELECTED REPORT DRIVER
BGP5DNG	IHS/CMI/LAB - IHS AREA CLIN 05 REPORT DRIVER
BGP5DP	IHS/CMI/LAB - IHS gpra print
BGP5DP1	IHS/CMI/LAB - print ind 1 21 Mar 2005 12:55 PM
BGP5DP11	IHS/CMI/LAB - print ind 1
BGP5DP2	IHS/CMI/LAB - print ind 10
BGP5DP3	IHS/CMI/LAB - print ind 1
BGP5DP4	IHS/CMI/LAB - print ind 31
BGP5DP5	IHS/CMI/LAB - print ind H
BGP5DP50	IHS/CMI/LAB - print ind H
BGP5DP51	IHS/CMI/LAB - print ind H
BGP5DP6	IHS/CMI/LAB - print ind H
BGP5DP60	IHS/CMI/LAB - print ind H
BGP5DP7	IHS/CMI/LAB - print ind 10
BGP5DPH	IHS/CMI/LAB - AREA REPORT HEADER
BGP5DPP	IHS/CMI/LAB - IHS GPRA 05 SELECTED REPORT DRIVER ALL PATIENTS
BGP5DS	IHS/CMI/LAB - IHS gpra print
BGP5DSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5DSL	IHS/CMI/LAB - FY 05 DISPLAY IND LISTS
BGP5DSP	IHS/CMI/LAB - IHS summary page
BGP5EL10	IHS/CMI/LAB - print ind 1
BGP5EL11	IHS/CMI/LAB - print ind 1
BGP5EL12	IHS/CMI/LAB - print ind 1

Routine	Description
BGP5EL13	IHS/CMI/LAB - print ind 1
BGP5EL14	IHS/CMI/LAB - print ind 1
BGP5EL15	IHS/CMI/LAB - print ind 1
BGP5EL16	IHS/CMI/LAB - print ind 1
BGP5EL17	IHS/CMI/LAB - print ind 1
BGP5EL18	IHS/CMI/LAB - print ind 1
BGP5EL19	IHS/CMI/LAB - print ind 1
BGP5EL1A	IHS/CMI/LAB - print ind 1
BGP5EL1B	IHS/CMI/LAB - print ind 1
BGP5EL1C	IHS/CMI/LAB - print ind 1
BGP5EL1D	IHS/CMI/LAB - print ind 1
BGP5EL1E	IHS/CMI/LAB - print ind 1
BGP5EL1F	IHS/CMI/LAB - print ind 1
BGP5EL1G	IHS/CMI/LAB - print ind 1
BGP5EL1H	IHS/CMI/LAB - print ind 1
BGP5EL1I	IHS/CMI/LAB - print ind 1
BGP5EL1J	IHS/CMI/LAB - print ind 1
BGP5EL1K	IHS/CMI/LAB - print ind 1
BGP5EL1L	IHS/CMI/LAB - print ind 1
BGP5EL2	IHS/CMI/LAB - indicator 1,2,3,4 17 Jan 2005 6:49 AM
BGP5EL3	IHS/CMI/LAB - indicator 1,2,3,4 05 Apr 2005 1:44 PM
BGP5EL4	IHS/CMI/LAB - indicator 1,2,3,4 05 Apr 2005 1:44 PM
BGP5ELH	IHS/CMI/LAB - cover page for ELDER
BGP5ELHH	IHS/CMI/LAB - cover page for ELDER
BGP5ELL	IHS/CMI/LAB - IHS gpra print
BGP5ELL1	IHS/CMI/LAB - print ind 1
BGP5ELL2	IHS/CMI/LAB - print ind
BGP5ELL3	IHS/CMI/LAB - print ind
BGP5ELL4	IHS/CMI/LAB - print ind
BGP5ELL5	IHS/CMI/LAB - print ind
BGP5ELL6	IHS/CMI/LAB - print ind
BGP5ELL7	IHS/CMI/LAB - print ind
BGP5ELL8	IHS/CMI/LAB - print ind
BGP5ELL9	IHS/CMI/LAB - print ind
BGP5ELLA	IHS/CMI/LAB - print ind
BGP5ELLB	IHS/CMI/LAB - print ind

Routine	Description
BGP5ELLC	IHS/CMI/LAB - print ind
BGP5ELLD	IHS/CMI/LAB - print ind
BGP5ELLE	IHS/CMI/LAB - print ind
BGP5ELLF	IHS/CMI/LAB - print ind
BGP5ELLG	IHS/CMI/LAB - print ind
BGP5ELLH	IHS/CMI/LAB - print ind
BGP5ELLI	IHS/CMI/LAB - print ind
BGP5ELLJ	IHS/CMI/LAB - print ind
BGP5ELLK	IHS/CMI/LAB - print ind
BGP5ELLL	IHS/CMI/LAB - print ind
BGP5ELLM	IHS/CMI/LAB - print ind
BGP5ELLN	IHS/CMI/LAB - print ind 1
BGP5ELLQ	IHS/CMI/LAB - print ind
BGP5ELP1	IHS/CMI/LAB - print ind 1
BGP5ELS	IHS/CMI/LAB - IHS gpra print
BGP5ELSL	IHS/CMI/LAB - IHS gpra print
BGP5ETL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5ETS	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5EUTL	IHS/CMI/LAB -
BGP5GAEL	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP5GAGP	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP5GAHE	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP5GALL	IHS/CMI/LAB - GUI COM REPORT
BGP5GAPU	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP5GCMS	IHS/CMI/LAB - GUI CMS REPORT
BGP5GCOM	IHS/CMI/LAB - GUI COM REPORT
BGP5GELD	IHS/CMI/LAB - GUI COM REPORT
BGP5GGPU	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP5GHED	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP5GLST	IHS/CMI/LAB - GUI CMS REPORT
BGP5GNPL	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP5GNTL	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP5GPP	IHS/CMI/LAB - GUI PP REPORT
BGP5GR	IHS/CMI/LAB - BGPG Visual CRS Reports 12/30/2004 12:29:35 PM
BGP5GRA	IHS/CMI/LAB - BGP Gui Area Reports 5/2/2005 8:38:59 PM
BGP5GT	IHS/CMI/LAB - BGPG Gui CRS Tables 2/2/2005 10:24:22 AM

Routine	Description
BGP5GU	IHS/CMI/LAB - BGP Gui Utilities 10/29/2004 3:28:39 PM
BGP5GUPL	IHS/CMI/LAB - GUI CMS REPORT
BGP5HEH	IHS/CMI/LAB - cover page for HEDIS
BGP5HEHH	IHS/CMI/LAB - cover page for HEDIS
BGP5HEL	IHS/CMI/LAB - IHS gpra print
BGP5HEL1	IHS/CMI/LAB - print ind 1
BGP5HEP1	IHS/CMI/LAB - print ind 1
BGP5HES	IHS/CMI/LAB - IHS gpra print
BGP5HESL	IHS/CMI/LAB - IHS gpra print
BGP5HUTL	IHS/CMI/LAB -
BGP5LSTF	IHS/CMI/LAB - List 'BG05' files in pub
BGP5NPL	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP5NPLD	IHS/CMI/LAB - IHS gpra print
BGP5NPLP	IHS/CMI/LAB - IHS gpra print
BGP5PARP	IHS/CMI/LAB - IHS gpra print
BGP5PDL	IHS/CMI/LAB - IHS gpra print
BGP5PDL0	IHS/CMI/LAB - print ind 1
BGP5PDL1	IHS/CMI/LAB - print ind 1
BGP5PDL2	IHS/CMI/LAB - print ind 1
BGP5PDL3	IHS/CMI/LAB - print ind 1
BGP5PDL4	IHS/CMI/LAB - print ind 1
BGP5PDL5	IHS/CMI/LAB - print ind 1
BGP5PDL7	IHS/CMI/LAB - print ind 1
BGP5PDL8	IHS/CMI/LAB - print ind 1
BGP5PDL9	IHS/CMI/LAB - print ind
BGP5PDLH	IHS/CMI/LAB - cover page for gpra del
BGP5PDL5	IHS/CMI/LAB - IHS gpra print
BGP5PHEL	IHS/CMI/LAB - IHS gpra print
BGP5PHEP	IHS/CMI/LAB - IHS gpra print
BGP5SDP	IHS/CMI/LAB - IHS summary page
BGP5TXCH	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP5TXCN	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP5TXEL	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP5TXHE	IHS/CMI/LAB - TAXONOMY CHECK FOR FY05 HEDIS REPORT
BGP5TXL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5TXS	IHS/CMI/LAB - DISPLAY IND LISTS

Routine	Description
BGP5TXV	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5TXV1	IHS/CMI/LAB - DISPLAY IND LISTS
BGP5ULF	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP5UTL	IHS/CMI/LAB -
BGP5UTL1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP6AHW	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP6ASL	IHS/CMI/LAB - DISPLAY IND LISTS 28 Apr 2006 2:10 PM
BGP6AX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXGB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXGC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6AXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005

Routine	Description
BGP6BAN	IHS/CMI/LAB - BANNER FOR CRS 05 06 Jan 2006 9:02 AM
BGP6BX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXBG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6BXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6C1	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM 04 May 2006 2:38 PM

Routine	Description
BGP6C11	IHS/CMI/LAB - calc CMS measures
BGP6C12	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP6C13	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP6C2	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP6C3	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP6C31	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP6CL	IHS/CMI/LAB - IHS GPRA 05 CMS REPORT DRIVER
BGP6CP	IHS/CMI/LAB - IHS gpra print
BGP6CSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6CTL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6CTS	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6CTXC	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT 18 Feb 2006 5:31 PM
BGP6D1	IHS/CMI/LAB - calc measures 29 Apr 2006 7:38 PM 01 May 2006 3:25 PM
BGP6D10	IHS/CMI/LAB - calc measures
BGP6D2	IHS/CMI/LAB - measure 1,2,3,4
BGP6D21	IHS/CMI/LAB - measure 6
BGP6D22	IHS/CMI/LAB - measure I2
BGP6D23	IHS/CMI/LAB - measure I2
BGP6D3	IHS/CMI/LAB - measure 11
BGP6D31	IHS/CMI/LAB - measure C
BGP6D32	IHS/CMI/LAB - measure C
BGP6D33	IHS/CMI/LAB - measure C
BGP6D34	IHS/CMI/LAB - measure C
BGP6D35	IHS/CMI/LAB - measure C
BGP6D4	IHS/CMI/LAB - measure 3
BGP6D41	IHS/CMI/LAB - measure 3
BGP6D42	IHS/CMI/LAB - measure 11
BGP6D5	IHS/CMI/LAB - measure calc
BGP6D51	IHS/CMI/LAB - measure I2
BGP6D52	IHS/CMI/LAB - measure 31
BGP6D6	IHS/CMI/LAB - measure 31
BGP6D61	IHS/CMI/LAB - measure 31
BGP6D7	IHS/CMI/LAB - measure 31
BGP6D71	IHS/CMI/LAB - measure C
BGP6D72	IHS/CMI/LAB - measure 31

Routine	Description
BGP6D73	IHS/CMI/LAB - measure 31
BGP6D8	IHS/CMI/LAB - measure C
BGP6D81	IHS/CMI/LAB - measure C
BGP6D82	IHS/CMI/LAB - measure C 14 Mar 2006 11:49 AM
BGP6D83	IHS/CMI/LAB - measure C
BGP6D84	IHS/CMI/LAB - measure C
BGP6D85	IHS/CMI/LAB - measure C
BGP6D86	IHS/CMI/LAB - measure C
BGP6D861	IHS/CMI/LAB - measure C
BGP6D862	IHS/CMI/LAB - measure C
BGP6D863	IHS/CMI/LAB - measure C
BGP6D9	IHS/CMI/LAB - measure J
BGP6DAP	IHS/CMI/LAB - IHS GPRA 04 SELECTED REPORT DRIVER ALL PATIENTS
BGP6DAR	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP6DCEI	IHS/CMI/LAB - calculate HEDIS measures
BGP6DCHA	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP6DCHI	IHS/CMI/LAB - calculate HEDIS measures
BGP6DCHW	IHS/CMI/LAB - calc measures 29 Apr 2006 7:38 PM
BGP6DCI	IHS/CMI/LAB - IHS area GPRA
BGP6DCL	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP6DCLD	IHS/CMI/LAB - IHS gpra print
BGP6DCLP	IHS/CMI/LAB - IHS gpra print
BGP6DEL	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER
BGP6DELA	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER
BGP6DELP	IHS/CMI/LAB - IHS HEDIS print
BGP6DESI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6DESL	IHS/CMI/LAB - FY 06 DISPLAY IND LISTS
BGP6DGPA	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP6DGPU	IHS/CMI/LAB - IHS AREA CLIN 05 REPORT DRIVER
BGP6DH	IHS/CMI/LAB - cover page for gpra 28 Apr 2006 11:30 AM
BGP6DHE	IHS/CMI/LAB - IHS GPRA 05 REPORT DRIVER
BGP6DHE1	IHS/CMI/LAB - measure HEDIS
BGP6DHEA	IHS/CMI/LAB - ihs area AA
BGP6DHEP	IHS/CMI/LAB - IHS HEDIS print
BGP6DHSL	IHS/CMI/LAB - FY 06 DISPLAY IND LISTS

Routine	Description
BGP6DL	IHS/CMI/LAB - IHS GPRA 05 SELECTED REPORT DRIVER
BGP6DLT	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP6DMT	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP6DNG	IHS/CMI/LAB - IHS AREA CLIN 05 REPORT DRIVER
BGP6DP	IHS/CMI/LAB - IHS gpra print
BGP6DP1	IHS/CMI/LAB - print ind 1 21 Mar 2006 12:55 PM
BGP6DP11	IHS/CMI/LAB - print ind 1
BGP6DP12	IHS/CMI/LAB - print ind 19 AGE DIST
BGP6DP1C	IHS/CMI/LAB - print ind 1 21 Mar 2006 12:55 PM
BGP6DP2	IHS/CMI/LAB - print ind 10
BGP6DP3	IHS/CMI/LAB - print ind 1
BGP6DP4	IHS/CMI/LAB - print ind 31
BGP6DP5	IHS/CMI/LAB - print ind H
BGP6DP50	IHS/CMI/LAB - print ind H
BGP6DP51	IHS/CMI/LAB - print ind H
BGP6DP6	IHS/CMI/LAB - print ind H
BGP6DP60	IHS/CMI/LAB - print ind H
BGP6DP7	IHS/CMI/LAB - print ind 10
BGP6DPH	IHS/CMI/LAB - AREA REPORT HEADER
BGP6DPP	IHS/CMI/LAB - IHS GPRA 05 SELECTED REPORT DRIVER ALL PATIENTS
BGP6DS	IHS/CMI/LAB - IHS gpra print
BGP6DSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6DSL	IHS/CMI/LAB - FY 06 DISPLAY IND LISTS
BGP6DSP	IHS/CMI/LAB - IHS summary page
BGP6DSTM	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP6DU	IHS/CMI/LAB - gpra utility calls
BGP6EL10	IHS/CMI/LAB - print ind 1
BGP6EL11	IHS/CMI/LAB - print ind 1
BGP6EL12	IHS/CMI/LAB - print ind 1
BGP6EL13	IHS/CMI/LAB - print ind 1
BGP6EL14	IHS/CMI/LAB - print ind 1
BGP6EL15	IHS/CMI/LAB - print ind 1
BGP6EL16	IHS/CMI/LAB - print ind 1
BGP6EL17	IHS/CMI/LAB - print ind 1
BGP6EL18	IHS/CMI/LAB - print ind 1

Routine	Description
BGP6EL19	IHS/CMI/LAB - print ind 1
BGP6EL1A	IHS/CMI/LAB - print ind 1
BGP6EL1B	IHS/CMI/LAB - print ind 1
BGP6EL1C	IHS/CMI/LAB - print ind 1
BGP6EL1D	IHS/CMI/LAB - print ind 1
BGP6EL1E	IHS/CMI/LAB - print ind 1
BGP6EL1F	IHS/CMI/LAB - print ind 1
BGP6EL1G	IHS/CMI/LAB - print ind 1
BGP6EL1H	IHS/CMI/LAB - print ind 1
BGP6EL1I	IHS/CMI/LAB - print ind 1
BGP6EL1J	IHS/CMI/LAB - print ind 1
BGP6EL1K	IHS/CMI/LAB - print ind 1
BGP6EL1L	IHS/CMI/LAB - print ind 1
BGP6EL1M	IHS/CMI/LAB - print ind 1
BGP6EL1N	IHS/CMI/LAB - print ind 1
BGP6EL2	IHS/CMI/LAB - measure 1,2,3,4 17 Jan 2006 6:49 AM
BGP6EL3	IHS/CMI/LAB - measure 1,2,3,4 05 Apr 2006 1:44 PM
BGP6EL4	IHS/CMI/LAB - measure 1,2,3,4 05 Apr 2006 1:44 PM
BGP6ELH	IHS/CMI/LAB - cover page for ELDER
BGP6ELHH	IHS/CMI/LAB - cover page for ELDER
BGP6ELL	IHS/CMI/LAB - IHS gpra print
BGP6ELL1	IHS/CMI/LAB - print ind 1
BGP6ELL2	IHS/CMI/LAB - print ind
BGP6ELL3	IHS/CMI/LAB - print ind
BGP6ELL4	IHS/CMI/LAB - print ind
BGP6ELL5	IHS/CMI/LAB - print ind
BGP6ELL6	IHS/CMI/LAB - print ind
BGP6ELL7	IHS/CMI/LAB - print ind
BGP6ELL8	IHS/CMI/LAB - print ind
BGP6ELL9	IHS/CMI/LAB - print ind
BGP6ELLA	IHS/CMI/LAB - print ind
BGP6ELLB	IHS/CMI/LAB - print ind
BGP6ELLC	IHS/CMI/LAB - print ind
BGP6ELLD	IHS/CMI/LAB - print ind
BGP6ELLE	IHS/CMI/LAB - print ind
BGP6ELLF	IHS/CMI/LAB - print ind

Routine	Description
BGP6ELLG	IHS/CMI/LAB - print ind
BGP6ELLH	IHS/CMI/LAB - print ind
BGP6ELLI	IHS/CMI/LAB - print ind
BGP6ELLJ	IHS/CMI/LAB - print ind
BGP6ELLK	IHS/CMI/LAB - print ind
BGP6ELLL	IHS/CMI/LAB - print ind
BGP6ELLM	IHS/CMI/LAB - print ind
BGP6ELLN	IHS/CMI/LAB - print ind 1
BGP6ELLO	IHS/CMI/LAB - print ind
BGP6ELLP	IHS/CMI/LAB - print ind
BGP6ELLQ	IHS/CMI/LAB - print ind
BGP6ELP1	IHS/CMI/LAB - print ind 1
BGP6ELS	IHS/CMI/LAB - IHS gpra print
BGP6ELSL	IHS/CMI/LAB - IHS gpra print
BGP6ETL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6EUTL	IHS/CMI/LAB -
BGP6EX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXBB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXHB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXHC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005

Routine	Description
BGP6EXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6EXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6FXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6FXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 07, 2005
BGP6FXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 08, 2005
BGP6FXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 08, 2005
BGP6FXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 08, 2005
BGP6FXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 08, 2005
BGP6FXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 08, 2005
BGP6FXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 08, 2005
BGP6FXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 08, 2005
BGP6FXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 08, 2005
BGP6FXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005
BGP6FXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005
BGP6FXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005
BGP6FXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005
BGP6FXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005
BGP6FXUB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005
BGP6FXUC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005
BGP6FXUD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005
BGP6FXUE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005

Routine	Description
BGP6FXUF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005
BGP6FXUG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005
BGP6FXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 12, 2005
BGP6FXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6FXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6FXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6FXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6GACW	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP6GAEL	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP6GAGP	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP6GAHE	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP6GALL	IHS/CMI/LAB - GUI COM REPORT
BGP6GAPU	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP6GCMP	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GCMS	IHS/CMI/LAB - GUI CMS REPORT 27 Apr 2006 10:30 PM
BGP6GCOM	IHS/CMI/LAB - GUI COM REPORT
BGP6GELD	IHS/CMI/LAB - GUI COM REPORT
BGP6GGPU	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GHED	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GLST	IHS/CMI/LAB - GUI CMS REPORT
BGP6GLTX	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GMTX	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GNPL	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GNST	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP6GNTL	IHS/CMI/LAB - IHS Diabetes Audit 2003 27 Apr 2006 10:56 PM
BGP6GPP	IHS/CMI/LAB - GUI PP REPORT
BGP6GR	IHS/CMI/LAB - BGPG Visual CRS Reports 12/30/2004 12:29:35 PM 27 Apr 2006 10:30 PM
BGP6GR1	IHS/CMI/LAB - BGPG Visual CRS Reports 12/30/2004 12:29:35 PM
BGP6GRA	IHS/CMI/LAB - BGP Gui Area Reports 5/2/2005 8:38:59 PM
BGP6GT	IHS/CMI/LAB - BGPG Gui CRS Tables 2/2/2005 10:24:22 AM
BGP6GU	IHS/CMI/LAB - BGP Gui Utilities 10/29/2004 3:28:39 PM 19 Sep 2005 5:28 PM 27 Apr 2006 10:28 PM
BGP6GUPL	IHS/CMI/LAB - GUI CMS REPORT
BGP6GX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 10, 2006
BGP6GXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005

Routine	Description
BGP6GXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6GXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6GXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6GXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6GXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6GXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6GXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 13, 2005
BGP6GXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 29, 2005
BGP6GXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 29, 2005
BGP6GXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 29, 2005
BGP6GXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 02, 2006
BGP6GXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 05, 2006
BGP6GXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXQB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXQC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXQD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXQE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXQF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXQG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXQH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXQI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 06, 2006
BGP6GXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 08, 2006
BGP6GXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 08, 2006
BGP6GXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 08, 2006
BGP6GXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 10, 2006
BGP6GXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 10, 2006
BGP6HEH	IHS/CMI/LAB - cover page for HEDIS
BGP6HEHH	IHS/CMI/LAB - cover page for HEDIS
BGP6HEL	IHS/CMI/LAB - IHS gpra print
BGP6HEL1	IHS/CMI/LAB - print ind 1
BGP6HEP1	IHS/CMI/LAB - print ind 1
BGP6HES	IHS/CMI/LAB - IHS gpra print 01 May 2006 4:31 PM
BGP6HESL	IHS/CMI/LAB - IHS gpra print 28 Apr 2006 3:42 PM

Routine	Description
BGP6HUTL	IHS/CMI/LAB -
BGP6HX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6HXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 21, 2005
BGP6JX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 09, 2006
BGP6JXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 10, 2006
BGP6JXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 10, 2006
BGP6JXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 10, 2006
BGP6JXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 10, 2006
BGP6JXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 10, 2006
BGP6JXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 30, 2006
BGP6JXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 30, 2006

Routine	Description
BGP6JXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JAN 30, 2006
BGP6JXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 16, 2006
BGP6JXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON FEB 02, 2006
BGP6JXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON FEB 02, 2006
BGP6JXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON FEB 02, 2006
BGP6JXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON FEB 06, 2006
BGP6JXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON FEB 09, 2006
BGP6JXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON FEB 20, 2006
BGP6JXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 22, 2006
BGP6JXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 22, 2006
BGP6JXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 23, 2006
BGP6JXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 27, 2006
BGP6JXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 27, 2006
BGP6JXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 27, 2006
BGP6JXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 28, 2006
BGP6JXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 09, 2006
BGP6LSTF	IHS/CMI/LAB - List 'BG06' files in pub
BGP6NPL	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP6NPLD	IHS/CMI/LAB - IHS gpra print
BGP6NPLP	IHS/CMI/LAB - IHS gpra print
BGP6PARP	IHS/CMI/LAB - IHS gpra print
BGP6PDL	IHS/CMI/LAB - IHS gpra print
BGP6PDL0	IHS/CMI/LAB - print ind 1
BGP6PDL1	IHS/CMI/LAB - print ind 1
BGP6PDL2	IHS/CMI/LAB - print ind 1
BGP6PDL3	IHS/CMI/LAB - print ind 1
BGP6PDL4	IHS/CMI/LAB - print ind 1
BGP6PDL5	IHS/CMI/LAB - print ind 1
BGP6PDL7	IHS/CMI/LAB - print ind 1
BGP6PDL8	IHS/CMI/LAB - print ind 1
BGP6PDL9	IHS/CMI/LAB - print ind
BGP6PDLA	IHS/CMI/LAB - print ind
BGP6PDLH	IHS/CMI/LAB - cover page for gpra del 02 May 2006 12:37 PM
BGP6PDL S	IHS/CMI/LAB - IHS gpra print
BGP6PHEL	IHS/CMI/LAB - IHS gpra print
BGP6PHEP	IHS/CMI/LAB - IHS gpra print

Routine	Description
BGP6SDP	IHS/CMI/LAB - IHS summary page
BGP6ULF	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP6UTL	IHS/CMI/LAB - 27 Apr 2006 11:01 PM
BGP6UTL1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP6XTCH	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP6XTCN	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP6XTEL	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP6XTHE	IHS/CMI/LAB - TAXONOMY CHECK FOR FY06 HEDIS REPORT
BGP6XTL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6XTS	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6XTV	IHS/CMI/LAB - DISPLAY IND LISTS
BGP6XTV1	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7ASL	IHS/CMI/LAB - DISPLAY IND LISTS 28 Apr 2006 2:10 PM
BGP7BAN	IHS/CMI/LAB - BANNER FOR CRS 07
BGP7C1	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM 04 May 2007 2:38 PM
BGP7C11	IHS/CMI/LAB - calc CMS measures
BGP7C12	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP7C13	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP7C2	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP7C3	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP7C31	IHS/CMI/LAB - calc CMS measures 26 Sep 2004 11:28 AM
BGP7CL	IHS/CMI/LAB - IHS GPRA 07 CMS REPORT DRIVER
BGP7CP	IHS/CMI/LAB - IHS gpra print
BGP7CSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7CTL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7CTS	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7CTXC	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT 18 Feb 2007 5:31 PM
BGP7D1	IHS/CMI/LAB - calc measures 29 Apr 2007 7:38 PM 14 Nov 2006 5:02 PM
BGP7D10	IHS/CMI/LAB - calc measures
BGP7D2	IHS/CMI/LAB - measure 1,2,3,4
BGP7D21	IHS/CMI/LAB - measure 6
BGP7D211	IHS/CMI/LAB - measure 6
BGP7D21A	IHS/CMI/LAB - measure 6
BGP7D22	IHS/CMI/LAB - measure I2

Routine	Description
BGP7D23	IHS/CMI/LAB - measure I2
BGP7D3	IHS/CMI/LAB - measure 11
BGP7D31	IHS/CMI/LAB - measure C
BGP7D32	IHS/CMI/LAB - measure C
BGP7D33	IHS/CMI/LAB - measure C
BGP7D34	IHS/CMI/LAB - measure C
BGP7D35	IHS/CMI/LAB - measure C
BGP7D4	IHS/CMI/LAB - measure 3
BGP7D41	IHS/CMI/LAB - measure 3
BGP7D42	IHS/CMI/LAB - measure 11
BGP7D5	IHS/CMI/LAB - measure calc
BGP7D51	IHS/CMI/LAB - measure I2
BGP7D52	IHS/CMI/LAB - measure 31
BGP7D6	IHS/CMI/LAB - measure 31
BGP7D61	IHS/CMI/LAB - measure 31
BGP7D7	IHS/CMI/LAB - measure 31
BGP7D71	IHS/CMI/LAB - measure C
BGP7D72	IHS/CMI/LAB - measure 31
BGP7D721	IHS/CMI/LAB - measure AHR.A
BGP7D722	IHS/CMI/LAB - measure AHR.A
BGP7D73	IHS/CMI/LAB - measure 31
BGP7D731	IHS/CMI/LAB - measure AHR.A
BGP7D732	IHS/CMI/LAB - measure AHR.A
BGP7D74	IHS/CMI/LAB - measure 1,2,3,4 05 Apr 2007 1:44 PM
BGP7D75	IHS/CMI/LAB - measure 31
BGP7D8	IHS/CMI/LAB - measure C
BGP7D81	IHS/CMI/LAB - measure C
BGP7D82	IHS/CMI/LAB - measure C 14 Mar 2007 11:49 AM
BGP7D83	IHS/CMI/LAB - measure C
BGP7D84	IHS/CMI/LAB - measure C
BGP7D85	IHS/CMI/LAB - measure C
BGP7D86	IHS/CMI/LAB - measure C
BGP7D861	IHS/CMI/LAB - measure C
BGP7D862	IHS/CMI/LAB - measure C
BGP7D863	IHS/CMI/LAB - measure C
BGP7D864	IHS/CMI/LAB - measure C

Routine	Description
BGP7D9	IHS/CMI/LAB - measure J
BGP7DAP	IHS/CMI/LAB - IHS GPRA 04 SELECTED REPORT DRIVER ALL PATIENTS
BGP7DAR	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP7DCEI	IHS/CMI/LAB - calculate HEDIS measures
BGP7DCHA	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP7DCHI	IHS/CMI/LAB - calculate HEDIS measures
BGP7DCHW	IHS/CMI/LAB - calc measures 29 Apr 2007 7:38 PM
BGP7DCI	IHS/CMI/LAB - IHS area GPRA
BGP7DCL	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP7DCLD	IHS/CMI/LAB - IHS gpra print
BGP7DCLP	IHS/CMI/LAB - IHS gpra print
BGP7DEL	IHS/CMI/LAB - IHS GPRA 07 REPORT DRIVER
BGP7DELA	IHS/CMI/LAB - IHS GPRA 07 REPORT DRIVER
BGP7DELP	IHS/CMI/LAB - IHS HEDIS print
BGP7DESI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7DESL	IHS/CMI/LAB - FY 07 DISPLAY IND LISTS
BGP7DGPA	IHS/CMI/LAB - ihs area GPRA 02 Sep 2004 1:11 PM
BGP7DGPU	IHS/CMI/LAB - IHS AREA CLIN 05 REPORT DRIVER
BGP7DH	IHS/CMI/LAB - cover page for gpra 28 Apr 2007 11:30 AM
BGP7DHE	IHS/CMI/LAB - IHS GPRA 07 REPORT DRIVER
BGP7DHE1	IHS/CMI/LAB - measure HEDIS
BGP7DHE2	IHS/CMI/LAB - measure HEDIS
BGP7DHEA	IHS/CMI/LAB - ihs area AA
BGP7DHEP	IHS/CMI/LAB - IHS HEDIS print
BGP7DHSL	IHS/CMI/LAB - FY 07 DISPLAY IND LISTS
BGP7DL	IHS/CMI/LAB - IHS GPRA 07 SELECTED REPORT DRIVER
BGP7DLT	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP7DMT	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP7DNE	IHS/CMI/LAB - NATL COMP EXPORT
BGP7DNG	IHS/CMI/LAB - NATL COMP EXPORT 13 Nov 2006 12:31 PM
BGP7DP	IHS/CMI/LAB - IHS gpra print
BGP7DP1	IHS/CMI/LAB - print ind 1 21 Mar 2007 12:55 PM
BGP7DP11	IHS/CMI/LAB - print ind 1
BGP7DP12	IHS/CMI/LAB - print ind 19 AGE DIST
BGP7DP19	IHS/CMI/LAB - print ind 1

Routine	Description
BGP7DP1C	IHS/CMI/LAB - print ind 1 21 Mar 2007 12:55 PM
BGP7DP2	IHS/CMI/LAB - print ind 10
BGP7DP3	IHS/CMI/LAB - print ind 1
BGP7DP4	IHS/CMI/LAB - print ind 31
BGP7DP5	IHS/CMI/LAB - print ind H
BGP7DP50	IHS/CMI/LAB - print ind H
BGP7DP51	IHS/CMI/LAB - print ind H
BGP7DP6	IHS/CMI/LAB - print ind H
BGP7DP60	IHS/CMI/LAB - print ind H
BGP7DP7	IHS/CMI/LAB - print ind 10
BGP7DPE	IHS/CMI/LAB - IHS GPRA 07 REPORT DRIVER
BGP7DPE1	IHS/CMI/LAB - calc measures 29 Apr 2007 7:38 PM 14 Nov 2006 5:02 PM
BGP7DPE2	IHS/CMI/LAB - calc measures 29 Apr 2007 7:38 PM 14 Nov 2006 5:02 PM
BGP7DPEA	IHS/CMI/LAB - IHS GPRA 07 REPORT DRIVER
BGP7DPED	IHS/CMI/LAB - IHS gpra print
BGP7DPEE	IHS/CMI/LAB - IHS gpra print
BGP7DPEL	IHS/CMI/LAB - IHS gpra print
BGP7DPEP	IHS/CMI/LAB - IHS gpra print
BGP7DPH	IHS/CMI/LAB - AREA REPORT HEADER
BGP7DPP	IHS/CMI/LAB - IHS GPRA 07 SELECTED REPORT DRIVER ALL PATIENTS
BGP7DPSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7DS	IHS/CMI/LAB - IHS gpra print
BGP7DSI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7DSL	IHS/CMI/LAB - FY 07 DISPLAY IND LISTS
BGP7DSP	IHS/CMI/LAB - IHS summary page
BGP7DSPN	IHS/CMI/LAB - IHS summary page
BGP7DSTM	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP7DU	IHS/CMI/LAB - gpra utility calls
BGP7EL10	IHS/CMI/LAB - print ind 1
BGP7EL11	IHS/CMI/LAB - print ind 1
BGP7EL12	IHS/CMI/LAB - print ind 1
BGP7EL13	IHS/CMI/LAB - print ind 1
BGP7EL14	IHS/CMI/LAB - print ind 1
BGP7EL15	IHS/CMI/LAB - print ind 1
BGP7EL16	IHS/CMI/LAB - print ind 1

Routine	Description
BGP7EL17	IHS/CMI/LAB - print ind 1
BGP7EL18	IHS/CMI/LAB - print ind 1
BGP7EL19	IHS/CMI/LAB - print ind 1
BGP7EL1A	IHS/CMI/LAB - print ind 1
BGP7EL1B	IHS/CMI/LAB - print ind 1
BGP7EL1C	IHS/CMI/LAB - print ind 1
BGP7EL1D	IHS/CMI/LAB - print ind 1
BGP7EL1E	IHS/CMI/LAB - print ind 1
BGP7EL1F	IHS/CMI/LAB - print ind 1
BGP7EL1G	IHS/CMI/LAB - print ind 1
BGP7EL1H	IHS/CMI/LAB - print ind 1
BGP7EL1I	IHS/CMI/LAB - print ind 1
BGP7EL1J	IHS/CMI/LAB - print ind 1
BGP7EL1K	IHS/CMI/LAB - print ind 1
BGP7EL1L	IHS/CMI/LAB - print ind 1
BGP7EL1M	IHS/CMI/LAB - print ind 1
BGP7EL1N	IHS/CMI/LAB - print ind 1
BGP7EL1O	IHS/CMI/LAB - print ind 1
BGP7EL1P	IHS/CMI/LAB - print ind
BGP7EL1Q	IHS/CMI/LAB - print ind
BGP7EL1R	IHS/CMI/LAB - print ind
BGP7EL2	IHS/CMI/LAB - measure 1,2,3,4 17 Jan 2007 6:49 AM
BGP7EL3	IHS/CMI/LAB - measure 1,2,3,4 05 Apr 2007 1:44 PM
BGP7EL31	IHS/CMI/LAB - measure 1,2,3,4 05 Apr 2007 1:44 PM
BGP7EL4	IHS/CMI/LAB - measure 1,2,3,4 05 Apr 2007 1:44 PM
BGP7ELH	IHS/CMI/LAB - cover page for ELDER
BGP7ELHH	IHS/CMI/LAB - cover page for ELDER
BGP7ELL	IHS/CMI/LAB - IHS gpra print
BGP7ELL1	IHS/CMI/LAB - print ind 1
BGP7ELL2	IHS/CMI/LAB - print ind
BGP7ELL3	IHS/CMI/LAB - print ind
BGP7ELL4	IHS/CMI/LAB - print ind
BGP7ELL5	IHS/CMI/LAB - print ind
BGP7ELL6	IHS/CMI/LAB - print ind
BGP7ELL7	IHS/CMI/LAB - print ind
BGP7ELL8	IHS/CMI/LAB - print ind

Routine	Description
BGP7ELL9	IHS/CMI/LAB - print ind
BGP7ELLA	IHS/CMI/LAB - print ind
BGP7ELLB	IHS/CMI/LAB - print ind
BGP7ELLC	IHS/CMI/LAB - print ind
BGP7ELLD	IHS/CMI/LAB - print ind
BGP7ELLE	IHS/CMI/LAB - print ind
BGP7ELLF	IHS/CMI/LAB - print ind
BGP7ELLG	IHS/CMI/LAB - print ind
BGP7EL LH	IHS/CMI/LAB - print ind
BGP7EL LI	IHS/CMI/LAB - print ind
BGP7EL LJ	IHS/CMI/LAB - print ind
BGP7EL LK	IHS/CMI/LAB - print ind
BGP7EL LL	IHS/CMI/LAB - print ind
BGP7EL LM	IHS/CMI/LAB - print ind
BGP7EL LN	IHS/CMI/LAB - print ind 1
BGP7EL LO	IHS/CMI/LAB - print ind
BGP7EL LP	IHS/CMI/LAB - print ind
BGP7EL LQ	IHS/CMI/LAB - print ind
BGP7EL LR	IHS/CMI/LAB - print ind
BGP7EL LS	IHS/CMI/LAB - print ind
BGP7EL LT	IHS/CMI/LAB - print ind
BGP7EL LU	IHS/CMI/LAB - print ind
BGP7EL P1	IHS/CMI/LAB - print ind 1
BGP7ELS	IHS/CMI/LAB - IHS gpra print
BGP7ELSL	IHS/CMI/LAB - IHS gpra print
BGP7ETL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7EUTL	IHS/CMI/LAB -
BGP7GACW	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP7GAEL	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP7GAGP	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP7GAHE	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP7GALL	IHS/CMI/LAB - GUI COM REPORT
BGP7GAPE	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP7GAPU	IHS/CMI/LAB - AREA NATIONAL GPRA REPORT
BGP7GCMP	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP7GCMS	IHS/CMI/LAB - GUI CMS REPORT 27 Apr 2007 10:30 PM

Routine	Description
BGP7GCOM	IHS/CMI/LAB - GUI COM REPORT
BGP7GELD	IHS/CMI/LAB - GUI COM REPORT
BGP7GGPU	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP7GHED	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP7GLST	IHS/CMI/LAB - GUI CMS REPORT
BGP7GLTX	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP7GMTX	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP7GNPL	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP7GNST	IHS/CMI/LAB - IHS Diabetes Audit 2003
BGP7GNTL	IHS/CMI/LAB - IHS Diabetes Audit 2003 27 Apr 2007 10:56 PM
BGP7GNXP	IHS/CMI/LAB - IHS Diabetes Audit 2003 27 Apr 2007 10:56 PM
BGP7GPED	IHS/CMI/LAB - GUI COM REPORT
BGP7GPP	IHS/CMI/LAB - GUI PP REPORT
BGP7GR	IHS/CMI/LAB - BGPG Visual CRS Reports 12/30/2004 12:29:35 PM 27 Apr 2007 10:30 PM
BGP7GR1	IHS/CMI/LAB - BGPG Visual CRS Reports 12/30/2004 12:29:35 PM
BGP7GR2	IHS/CMI/LAB - BGPG Visual CRS Reports 12/30/2004 12:29:35 PM 27 Apr 2007 10:30 PM
BGP7GT	IHS/CMI/LAB - BGPG Gui CRS Tables 2/2/2005 10:24:22 AM
BGP7GU	IHS/CMI/LAB - BGP Gui Utilities 10/29/2004 3:28:39 PM 19 Sep 2005 5:28 PM 27 Apr 2007 10:28 PM
BGP7GUPL	IHS/CMI/LAB - GUI CMS REPORT
BGP7GX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 27, 2006
BGP7GXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7GXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7GXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7GXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7GXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7GXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 20, 2006
BGP7GXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 20, 2006
BGP7GXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 20, 2006
BGP7GXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 21, 2006
BGP7GXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 24, 2006
BGP7G XK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 24, 2006
BGP7GXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXMB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006

Routine	Description
BGP7GXMC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXMD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXME	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXMF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXNB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXNC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXND	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXOB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXOC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXOD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXOE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXOF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXOG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXOH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXOI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXOJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 25, 2006
BGP7GXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 27, 2006
BGP7GXQB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 27, 2006
BGP7GXQC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 27, 2006
BGP7GXQD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 27, 2006
BGP7GXQE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 27, 2006
BGP7HEH	IHS/CMI/LAB - cover page for HEDIS
BGP7HEHH	IHS/CMI/LAB - cover page for HEDIS
BGP7HEL	IHS/CMI/LAB - IHS gpra print
BGP7HEL1	IHS/CMI/LAB - print ind 1
BGP7HEP1	IHS/CMI/LAB - print ind 1
BGP7HES	IHS/CMI/LAB - IHS gpra print 01 May 2007 4:31 PM
BGP7HESL	IHS/CMI/LAB - IHS gpra print 28 Apr 2007 3:42 PM
BGP7HUTL	IHS/CMI/LAB -
BGP7HX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006

Routine	Description
BGP7HXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7HXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7IX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7IXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXCB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006

Routine	Description
BGP7IXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXLB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXLC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXLD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXLE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXLF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXLG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXLH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7IXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 12, 2006
BGP7IXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 12, 2006
BGP7IXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 12, 2006
BGP7IXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 12, 2006
BGP7IXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXSS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXST	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXTB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006

Routine	Description
BGP7IXTC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXUB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXUC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXUD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXUE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 16, 2006
BGP7IXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 18, 2006
BGP7IXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 18, 2006
BGP7IXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 18, 2006
BGP7IXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7IXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 19, 2006
BGP7JX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7JXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXFB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXFC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXFD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXQB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXQC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXQD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXQE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006

Routine	Description
BGP7JXQF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXQG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7JXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7JXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7JXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7JXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 11, 2006
BGP7KX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXEB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXEC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXED	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXEE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXEF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXGB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXGC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXGD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXGE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXGF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXGG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXGH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXGI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXIB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXIC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006

Routine	Description
BGP7KXID	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7KXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7LSTF	IHS/CMI/LAB - List 'BG07' files in pub
BGP7LX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 13, 2006
BGP7LXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXAN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXBB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006

Routine	Description
BGP7LXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXCB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXDB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXDC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXEB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXEC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXED	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXEE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXEF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXFB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXFC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXFD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXFE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXFF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXFG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXFH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXFI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXFJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXGO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006

Routine	Description
BGP7LXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXHB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 28, 2006
BGP7LXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXIB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXIC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXID	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXJB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXJC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXJD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXKB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXKC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXKD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXKE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXKF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXKG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXKH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXKI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXKJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXKK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006

Routine	Description
BGP7LXLP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXLM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXNB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXNC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXOB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXOC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXOD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXOE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXOF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXOG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXOH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXOI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXOJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXOK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXPB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXPC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXPD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXPE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXPF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXPG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 29, 2006
BGP7LXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 13, 2006
BGP7LXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 30, 2006
BGP7LXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 30, 2006
BGP7LXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7LXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7LXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006

Routine	Description
BGP7LXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7LXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7LXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON SEP 10, 2006
BGP7MX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 20, 2006
BGP7MXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 20, 2006
BGP7MXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 20, 2006
BGP7NPL	IHS/CMI/LAB - national patient list 20 Dec 2004 9:24 AM
BGP7NPLD	IHS/CMI/LAB - IHS gpra print
BGP7NPLP	IHS/CMI/LAB - IHS gpra print
BGP7PARP	IHS/CMI/LAB - IHS gpra print
BGP7PD19	IHS/CMI/LAB - print ind
BGP7PDL	IHS/CMI/LAB - IHS gpra print
BGP7PDL0	IHS/CMI/LAB - print ind 1
BGP7PDL1	IHS/CMI/LAB - print ind 1
BGP7PDL2	IHS/CMI/LAB - print ind 1
BGP7PDL3	IHS/CMI/LAB - print ind 1
BGP7PDL4	IHS/CMI/LAB - print ind 1
BGP7PDL5	IHS/CMI/LAB - print ind 1
BGP7PDL7	IHS/CMI/LAB - print ind 1
BGP7PDL8	IHS/CMI/LAB - print ind 1
BGP7PDL9	IHS/CMI/LAB - print ind
BGP7PDLA	IHS/CMI/LAB - print ind
BGP7PDLH	IHS/CMI/LAB - cover page for gpra del 02 May 2007 12:37 PM
BGP7PDLN	IHS/CMI/LAB - IHS gpra print
BGP7PDL S	IHS/CMI/LAB - IHS gpra print
BGP7PESI	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7PESL	IHS/CMI/LAB - FY 06 DISPLAY IND LISTS
BGP7PHEL	IHS/CMI/LAB - IHS gpra print
BGP7PHEP	IHS/CMI/LAB - IHS gpra print
BGP7POS	IHS/CMI/LAB - NO DESCRIPTION PROVIDED 28 Jan 2005 1:34 PM
BGP7POS1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED 28 Jan 2005 1:34 PM
BGP7PEEL	IHS/CMI/LAB - IHS gpra print
BGP7PUTL	IHS/CMI/LAB -
BGP7SDP	IHS/CMI/LAB - IHS summary page
BGP7SDPN	IHS/CMI/LAB - IHS summary page
BGP7ULF	IHS/CMI/LAB - NO DESCRIPTION PROVIDED

Routine	Description
BGP7UTL	IHS/CMI/LAB - 27 Apr 2007 11:01 PM
BGP7UTL1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP7UTL2	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGP7XTCH	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP7XTCN	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP7XTEL	IHS/CMI/LAB - TAXONOMY CHECK FOR FY04 CRS REPORT
BGP7XTHE	IHS/CMI/LAB - TAXONOMY CHECK FOR FY07 HEDIS REPORT
BGP7XTL	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7XTS	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7XTV	IHS/CMI/LAB - DISPLAY IND LISTS
BGP7XTV1	IHS/CMI/LAB - DISPLAY IND LISTS
BGPAP	IHS/CMI/LAB - IHS gpra print
BGPAP1	IHS/CMI/LAB - print ind 1
BGPAP12	IHS/CMI/LAB - print ind 12
BGPAP13	IHS/CMI/LAB - print ind 13
BGPAP14	IHS/CMI/LAB - print ind 14
BGPAP1B	IHS/CMI/LAB - print ind 1B
BGPAP2	IHS/CMI/LAB - print ind 2
BGPAP22	IHS/CMI/LAB - print ind 22
BGPAP221	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGPAP22A	IHS/CMI/LAB - print ind 22
BGPAP22B	IHS/CMI/LAB - print ind 22
BGPAP22C	IHS/CMI/LAB - print ind 22
BGPAP23	IHS/CMI/LAB - print ind 23
BGPAP24	IHS/CMI/LAB - print ind 24
BGPAP3	IHS/CMI/LAB - print ind 3
BGPAP30	IHS/CMI/LAB - IHS gpra print
BGPAP30A	IHS/CMI/LAB - IHS gpra print
BGPAP30B	IHS/CMI/LAB - IHS gpra print
BGPAP30C	IHS/CMI/LAB - IHS gpra print
BGPAP4	IHS/CMI/LAB -print ind 4
BGPAP5	IHS/CMI/LAB - print ind 5
BGPAP6	IHS/CMI/LAB - print ind 6
BGPAP7	IHS/CMI/LAB - print ind 7
BGPAP8	IHS/CMI/LAB - print ind 8
BGPAPA	IHS/CMI/LAB - IHS gpra print

Routine	Description
BGPAPB	IHS/CMI/LAB - IHS gpra print
BGPAPC	IHS/CMI/LAB - IHS gpra print
BGPAPD	IHS/CMI/LAB - IHS gpra print
BGPAPH	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGPD	IHS/CMI/LAB - IHS GPRA - report for local use
BGPD0	IHS/CMI/LAB - DISPLAY IND LISTS
BGPD1	IHS/CMI/LAB - IHS area GPRA
BGPD12	IHS/CMI/LAB - indicator 12
BGPD13	IHS/CMI/LAB - indicator 13
BGPD14	IHS/CMI/LAB - indicator 14
BGPD1B	IHS/CMI/LAB - ind 1b
BGPD2	IHS/CMI/LAB - indicator 2
BGPD22	IHS/CMI/LAB - indicator 18
BGPD23	IHS/CMI/LAB - indicator 23
BGPD24	IHS/CMI/LAB - indicator 24
BGPD29	IHS/CMI/LAB - indicator 29
BGPD3	IHS/CMI/LAB - indicator 3
BGPD30	IHS/CMI/LAB - indicator 30
BGPD31	IHS/CMI/LAB - indicator 31
BGPD4	IHS/CMI/LAB - indicator 4
BGPD5	IHS/CMI/LAB - indicator 5
BGPD6	IHS/CMI/LAB - indicator 6
BGPD7	IHS/CMI/LAB - indicator 7
BGPD8	IHS/CMI/LAB - indicator 8
BGPDA	IHS/CMI/LAB - indicator A
BGPDAREA	IHS/CMI/LAB - ihs area GPRA
BGPDARP	IHS/CMI/LAB - IHS area GPRA
BGPDB	IHS/CMI/LAB - indicator B
BGPDC	IHS/CMI/LAB - indicator C
BGPDD	IHS/CMI/LAB - indicator D
BGPDFTA	IHS/CMI/LAB - IHS area GPRA
BGPDH	IHS/CMI/LAB - cover page for gpra
BGPDL	IHS/CMI/LAB - DISPLAY IND LISTS
BGPDP	IHS/CMI/LAB - IHS gpra print
BGPDP1	IHS/CMI/LAB - print ind 1
BGPDP12	IHS/CMI/LAB - print ind 12

Routine	Description
BGPDP13	IHS/CMI/LAB - print ind 13
BGPDP14	IHS/CMI/LAB - print ind 14
BGPDP1B	IHS/CMI/LAB - print ind 1B
BGPDP2	IHS/CMI/LAB - print ind 2
BGPDP22	IHS/CMI/LAB - print ind 22
BGPDP221	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGPDP22A	IHS/CMI/LAB - print ind 22
BGPDP22B	IHS/CMI/LAB - print ind 22
BGPDP22C	IHS/CMI/LAB - print ind 22
BGPDP23	IHS/CMI/LAB - print ind 23
BGPDP24	IHS/CMI/LAB - print ind 24
BGPDP29	IHS/CMI/LAB - IHS gpra print
BGPDP3	IHS/CMI/LAB - print ind 3
BGPDP30	IHS/CMI/LAB - IHS gpra print
BGPDP30A	IHS/CMI/LAB - IHS gpra print
BGPDP30B	IHS/CMI/LAB - IHS gpra print
BGPDP30C	IHS/CMI/LAB - IHS gpra print
BGPDP31	IHS/CMI/LAB - print ind 31
BGPDP32	IHS/CMI/LAB - print ind 32
BGPDP4	IHS/CMI/LAB -print ind 4
BGPDP5	IHS/CMI/LAB - print ind 5
BGPDP6	IHS/CMI/LAB - print ind 6
BGPDP7	IHS/CMI/LAB - print ind 7
BGPDP8	IHS/CMI/LAB - print ind 8
BGPDPA	IHS/CMI/LAB - IHS gpra print
BGPDPARP	IHS/CMI/LAB - IHS gpra print
BGPDPB	IHS/CMI/LAB - IHS gpra print
BGPDPC	IHS/CMI/LAB - IHS gpra print
BGPDPD	IHS/CMI/LAB - IHS gpra print
BGPDPH	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGPDS	IHS/CMI/LAB - IHS gpra print
BGPDT	IHS/CMI/LAB - 2001 DIABETES AUDIT
BGPDU	IHS/CMI/LAB - gpra utility calls
BGPDUP	IHS/CMI/LAB - NO DESCRIPTION PROVIDED
BGPGR	IHS/CMI/LAB - BGPG Visual CRS Reports 12/30/2004 12:29:35 PM
BGPGRA	IHS/CMI/LAB - BGP Gui Area Reports 5/2/2005 8:38:59 PM

Routine	Description
BGPGT	IHS/CMI/LAB - BGPG Gui CRS Tables 2/2/2005 10:24:22 AM
BGPGU	IHS/CMI/LAB - BGP Gui Utilities 10/29/2004 3:28:39 PM
BGPLTX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPLTXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPLTXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPLTXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPLTXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPLTXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPLTXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPLTXFB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPLTXFC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPLX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXBB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXCQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXEB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXEC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003

Routine	Description
BGPLXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXMB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPLXMC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 18, 2003
BGPMX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 20, 2006
BGPMXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 20, 2006
BGPMXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON NOV 20, 2006
BGPPOST	IHS/CMI/LAB - NO DESCRIPTION PROVIDED 28 Jan 2005 1:34 PM
BGPPOST1	IHS/CMI/LAB - NO DESCRIPTION PROVIDED 28 Jan 2005 1:34 PM
BGPTX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXEB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXEC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXED	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXEE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXEF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXEG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXGB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXKB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXKC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004

Routine	Description
BGPTXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXPB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXTB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXTC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXTD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXTE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXVB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 04, 2003
BGPTXVC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 04, 2003
BGPTXVD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 04, 2003
BGPTXVE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON DEC 04, 2003
BGPTXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXWB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004

Routine	Description
BGPTXXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON MAR 19, 2004
BGPTXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPTXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXAR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005

Routine	Description
BGPUTXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXZB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXZC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXZD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUTXZE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPUX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPUXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXAQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004

Routine	Description
BGPUXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXGB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXGC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 05, 2004
BGPUXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPUXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPUXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPUXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPUXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPUXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPVTX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005

Routine	Description
BGPVTXIB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXIC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVTXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPVX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPVXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPVXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPVXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPVXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPVXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPVXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 06, 2004
BGPVXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 08, 2004
BGPVXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 09, 2004
BGPVXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 09, 2004
BGPVXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 09, 2004
BGPVXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 09, 2004
BGPVXLB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 09, 2004
BGPVXLC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 09, 2004
BGPVXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 09, 2004
BGPVXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 09, 2004

Routine	Description
BGPVXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 13, 2004
BGPVXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 13, 2004
BGPVXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 13, 2004
BGPVXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON JUL 19, 2004
BGPVXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXUB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXUC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXUD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXUE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXV	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXW	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXY	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPVXZ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON AUG 03, 2004
BGPWTX	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 26, 2005
BGPWTXA	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBB	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXBG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXC	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXD	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXE	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXF	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXG	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXH	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXI	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXJ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXK	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXL	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXM	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005

Routine	Description
BGPWTXN	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXO	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXP	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXQ	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXR	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 21, 2005
BGPWTXS	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 24, 2005
BGPWTXT	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 24, 2005
BGPWTXU	IHS/CMI/LAB -CREATED BY ^ATXSTX ON APR 26, 2005

4.0 Files and Tables

4.1 File List

The following table provides information for files and tables.

File	Global	Filename
90240.01	^BGPD(BGP IHS GPRA DATA 02
90241.01	^BGPCTRL(BGP CONTROL FILE
90241.02	^BGPSITE(BGP SITE PARAMETERS
90242	^BGPINDC(BGP INDIVIDUAL INDICATOR DEFINITIONS
90242.01	^BGPIND(BGP INDICATORS
90243	^BGPDPDC(BGP 03 DATA CURRENT
90243.01	^BGPDPDP(BGP 03 DATA PREVIOUS
90243.02	^BGPDPDB(BGP 03 DATA BASELINE
90244.01	^BGPINDF(BGP 04 INDICATORS
90244.02	^BGPINDFC(BGP 04 INDIVIDUAL INDICATORS
90244.03	^BGPDPDCF(BGP 04 DATA CURRENT
90244.04	^BGPDPDPF(BGP 04 DATA PREVIOUS
90244.05	^BGPDPDBF(BGP 04 DATA BASELINE
90370.01	^BGPHEIF(BGP 04 HEDIS INDICATORS
90370.02	^BGPHEIIF(BGP 04 HEDIS INDIVIDUAL INDICATORS
90370.03	^BGPHEDCF(BGP 04 HEDIS DATA CURRENT
90370.04	^BGPHEDPF(BGP 04 HEDIS DATA PREVIOUS
90370.05	^BGPHEDBF(BGP 04 HEDIS DATA BASELINE
90371.01	^BGPINDV(BGP 05 INDICATORS
90371.02	^BGPINDVC(BGP 05 INDIVIDUAL INDICATORS
90371.03	^BGPDPDCV(BGP 05 DATA CURRENT
90371.04	^BGPDPDPV(BGP 05 DATA PREVIOUS
90371.05	^BGPDPDBV(BGP 05 DATA BASELINE
90371.06	^BGPVNPL(BGP 05 NATIONAL PATIENT LISTS
90371.08	^BGP TAXV(BGP 05 TAXONOMIES
90372.01	^BGPHEIV(BGP 05 HEDIS INDICATORS
90372.02	^BGPHEIIV(BGP 05 HEDIS INDIVIDUAL INDICATORS
90372.03	^BGPHEDCV(BGP 05 HEDIS DATA CURRENT
90372.04	^BGPHEDPV(BGP 05 HEDIS DATA PREVIOUS
90372.05	^BGPHEDBV(BGP 05 HEDIS DATA BASELINE

File	Global	Filename
90372.08	^BGPGUI(BGP 05 GUI REPORT OUTPUT
90373.01	^BGPCMSIF(BGP 05 CMS INDICATORS
90373.02	^BGPELIIV(BGP 05 ELDER INDIVIDUAL INDICATORS
90373.03	^BGPELDCV(BGP 05 ELDER DATA CURRENT
90373.04	^BGPELDPV(BGP 05 ELDER DATA PREVIOUS
90373.05	^BGPELDBV(BGP 05 ELDER DATA BASELINE
90374.01	^BGPINDS(BGP 06 INDICATORS
90374.02	^BGPINDSC(BGP 06 INDIVIDUAL INDICATORS
90374.03	^BGPDPDCS(BGP 06 DATA CURRENT
90374.04	^BGPDPDPS(BGP 06 DATA PREVIOUS
90374.05	^BGPDPDBS(BGP 06 DATA BASELINE
90374.06	^BGPSNPL(BGP 06 NATIONAL PATIENT LISTS
90374.08	^BGPTAXS(BGP 06 TAXONOMIES
90375.01	^BGPHEIS(BGP 06 HEDIS INDICATORS
90375.02	^BGPHEIIS(BGP 06 HEDIS INDIVIDUAL INDICATORS
90375.03	^BGPHEDCS(BGP 06 HEDIS DATA CURRENT
90375.04	^BGPHEDPS(BGP 06 HEDIS DATA PREVIOUS
90375.05	^BGPHEDBS(BGP 06 HEDIS DATA BASELINE
90375.08	^BGPGUI(BGP 06 GUI REPORT OUTPUT
90375.99	^BGPBMX(BGP BMX ADO SCHEMA
90376.01	^BGPCMSIS(BGP 06 CMS INDICATORS
90376.02	^BGPELIIS(BGP 06 ELDER INDIVIDUAL INDICATORS
90376.03	^BGPELDCS(BGP 06 ELDER DATA CURRENT
90376.04	^BGPELDPS(BGP 06 ELDER DATA PREVIOUS
90376.05	^BGPELDBS(BGP 06 ELDER DATA BASELINE
90376.06	^BGPELIS(BGP 06 ELDER INDICATORS
90530.01	^BGPINDA(BGP 07 INDICATORS
90530.02	^BGPINDAC(BGP 07 INDIVIDUAL INDICATORS
90530.03	^BGPDPDCA(BGP 07 DATA CURRENT
90530.04	^BGPDPDPA(BGP 07 DATA PREVIOUS
90530.05	^BGPDPDBA(BGP 07 DATA BASELINE
90530.06	^BGPANPL(BGP 07 NATIONAL PATIENT LISTS
90530.08	^BGPTAXA(BGP 07 TAXONOMIES
90530.09	^BGPTAXTA(BGP 07 TAXONOMY TEMP
90530.11	^BGPXPA(BGP 07 XP OUTPUT
90530.12	^BGPPEDECA(BGP 07 PATIENT ED DATA CURRENT

File	Global	Filename
90530.13	^BGPPEDPA(BGP 07 PATIENT ED DATA PREVIOUS
90530.14	^BGPPEDBA(BGP 07 PATIENT ED DATA BASELINE
90530.15	^BGPPEIA(BGP 07 PAT ED MEASURES
90531.01	^BGPHEIA(BGP 07 HEDIS INDICATORS
90531.02	^BGPHEIIA(BGP 07 HEDIS INDIVIDUAL INDICATORS
90531.03	^BGPHEIDCA(BGP 07 HEDIS DATA CURRENT
90531.04	^BGPHEIDPA(BGP 07 HEDIS DATA PREVIOUS
90531.05	^BGPHEIDBA(BGP 07 HEDIS DATA BASELINE
90531.08	^BGPGUIA(BGP 07 GUI REPORT OUTPUT
90532.01	^BGPCMSIA(BGP 07 CMS INDICATORS
90532.02	^BGPELIIA(BGP 07 ELDER INDIVIDUAL INDICATORS
90532.03	^BGPELDCA(BGP 07 ELDER DATA CURRENT
90532.04	^BGPELDPA(BGP 07 ELDER DATA PREVIOUS
90532.05	^BGPELDBA(BGP 07 ELDER DATA BASELINE
90532.06	^BGPELIA(BGP 07 ELDER INDICATORS

4.2 File Access

File (#)	Global	RD	WR	LYG	DD	DEL
90240.01	^BPGD	@	@	@	@	@
90241.01	^BGPCTRL	@	@	@	@	@
90241.02	^BGPSITE	M	M	M	@	@
90242	^BGPINDC	@	@	@	@	@
90242.01	^BGPIND	@	@	@	@	@
90243	^BGPGPDC	@	@	@	@	@
90243.01	^BGPGPDP	@	@	@	@	@
90243.02	^BGPGPDB	@	@	@	@	@
90244.01	^BGPINDF(0)	@	@	@	@	@
90244.02	^BGPINDFC(0)	@	@	@	@	@
90244.03	^BGPGPDCF(0)	@	@	@	@	@
90244.04	^BGPGPDPF(0)	@	@	@	@	@
90244.05	^BGPGPDBF(0)	@	@	@	@	@
90370.01	^BGPHEIF(0)	@	@	@	@	@
90370.02	^BGPHEIIF(0)	@	@	@	@	@
90370.03	^BGPHEIDCF(0)	@	@	@	@	@
90370.04	^BGPHEIDPF(0)	@	@	@	@	@

File (#)	Global	RD	WR	LYG	DD	DEL
90370.05	^BGPHEDBF(0)	@	@	@	@	@
90370.01	^BGPHEIF(0)	@	@	@	@	@
90370.02	^BGPHEIIF(0)	@	@	@	@	@
90370.03	^BGPHEDCF(0)	@	@	@	@	@
90370.04	^BGPHEDPF	@	@	@	@	@
90370.05	^BGPHEDBF	@	@	@	@	@
90371.01	^BGPINDV(0)	@	@	@	@	@
90371.02	^BGPINDVC(0)	@	@	@	@	@
90371.03	^BGPDPDCV(0)	@	@	@	@	@
90371.04	^BGPDPDPV(0)	@	@	@	@	@
90371.05	^BGPDPDBV(0)	@	@	@	@	@
90371.06	^BGPVNPL(0)	@	@	@	@	@
90371.08	^BGPTAXV(@	@	@	@	@
90372.01	^BGPHEIV(0)	@	@	@	@	@
90372.02	^BGPHEIIV(0)	@	@	@	@	@
90372.03	^BGPHEDCV(0)	@	@	@	@	@
90372.04	^BGPHEDPV	@	@	@	@	@
90372.05	^BGPHEDBV	@	@	@	@	@
90372.08	^BGPGUI(@	@	@	@	@
90373.01	^BGPCMSIF(@	@	@	@	@
90373.02	^BGPELIIV(@	@	@	@	@
90373.03	^BGPELDCV(@	@	@	@	@
90373.04	^BGPELDPV(@	@	@	@	@
90373.05	^BGPELDBV(@	@	@	@	@
90374.01	^BGPINDS(0)	@	@	@	@	@
90374.02	^BGPINDSC(0)	@	@	@	@	@
90374.03	^BGPDPDCS(0)	@	@	@	@	@
90374.04	^BGPDPDPS(0)	@	@	@	@	@
90374.05	^BGPDPDBS(0)	@	@	@	@	@
90374.06	^BGPSNPL(0)	@	@	@	@	@
90374.08	^BGPTAXS(@	@	@	@	@
90375.01	^BGPHEIS(0)	@	@	@	@	@
90375.02	^BGPHEIIS(0)	@	@	@	@	@
90375.03	^BGPHEDCS(0)	@	@	@	@	@
90375.04	^BGPHEDPS	@	@	@	@	@
90375.05	^BGPHEDBS	@	@	@	@	@

File (#)	Global	RD	WR	LYG	DD	DEL
90375.08	^BGPGUIS(@	@	@	@	@
90375.99	^BGPBMX(@	@	@	@	@
90376.01	^BGPCMSIS(@	@	@	@	@
90376.02	^BGPELIIS(@	@	@	@	@
90376.03	^BGPELDCS(@	@	@	@	@
90376.04	^BGPELDPS(@	@	@	@	@
90376.05	^BGPELDBS(@	@	@	@	@
90376.06	^BGPELIS(@	@	@	@	@
90530.01	^BGPINDA(@	@	@	@	@
90530.02	^BGPINDAC(@	@	@	@	@
90530.03	^BGPDPDCA(@	@	@	@	@
90530.04	^BGPDPDPA(@	@	@	@	@
90530.05	^BGPDPDBA(@	@	@	@	@
90530.06	^BGPANPL(@	@	@	@	@
90530.08	^BGPTAXA(@	@	@	@	@
90530.09	^BGPTAXTA(@	@	@	@	@
90530.11	^BGPXPA(@	@	@	@	@
90530.12	^BGPPEDCA(@	@	@	@	@
90530.13	^BGPPEDPA(@	@	@	@	@
90530.14	^BGPPEDBA(@	@	@	@	@
90530.15	^BGPPEIA(@	@	@	@	@
90531.01	^BGPHEIA(@	@	@	@	@
90531.02	^BGPHEIIA(@	@	@	@	@
90531.03	^BGPHEIDCA(M	M	M	@	M
90531.04	^BGPHEIDPA(@	@	@	@	@
90531.05	^BGPHEIDBA(@	@	@	@	@
90531.08	^BGPGUIA(@	@	@	@	@
90532.01	^BGPCMSIA(@	@	@	@	@
90532.02	^BGPELIIA(@	@	@	@	@
90532.03	^BGPELDCA(@	@	@	@	@
90532.04	^BGPELDPA(@	@	@	@	@
90532.05	^BGPELDDBA(@	@	@	@	@
90532.06	^BGPELIA(@	@	@	@	@

4.3 Cross References

90240.01 (BGP IHS GPRA DATA 02)

- .01 BEGINNING DATE
 - B Regular type cross reference
- 18000 PHN VISIT PRIMARY DX (90240.018)
 - .01 PHN VISIT PRIMARY DX
 - B Regular type cross reference
- 24000 CURRENT INJURIES (90240.24)
 - .01 CURRENT INJURIES
 - B Regular type cross reference
- 24000 E CODES (90240.2424)
 - .01 E CODES
 - B Regular type cross reference
- 25000 INJURIES 98 PREV (90240.25)
 - .01 INJURIES 98 PREV
 - B Regular type cross reference
- 25000 E CODE TALLY (90240.2525)
 - .01 E CODE TALLY
 - B Regular type cross reference
- 26000 INJURIES 98 (90240.26)
 - .01 INJURIES 98
 - B Regular type cross reference
- 26000 E CODES (90240.2626)
 - .01 E CODES
 - B Regular type cross reference
- 28000 COMMUNITIES (90240.28)
 - .01 COMMUNITIES
 - B Regular type cross reference
- 48000 PHN VISIT PRIMARY DX PR (90240.048)
 - .01 PHN VISIT PRIMARY DX PR
 - B Regular type cross reference
- 88000 PHN PRIMARY DIAGNOSIS BL (90240.088)
 - .01 PHN PRIMARY DIAGNOSIS BL
 - B Regular type cross reference

90241.01 (BGP CONTROL FILE)

- .01 FY
 - B Regular type cross reference
- 1101 CLINICS FOR CLINICAL USER (90241.0111)
 - .01 CLINIC CODE
 - B Regular type cross reference
- 1201 CLINIC CODES FOR 2ND VISIT (90241.0112)
 - .01 CLINIC CODES FOR 2ND VISIT
 - B Regular type cross reference

90241.02 (BGP SITE PARAMETERS)

- .01 LOCATION
B Regular type cross reference

90242 (BGP INDIVIDUAL INDICATOR DEFINITIONS)

- .01 INDICATOR
B Regular type cross reference
- .04 INDICATOR ID
C Regular type cross reference
AA MUMPS type cross reference
- .05 GPRA
AGPRA Regular type cross reference
- .06 AREA APP
AREA Regular type cross reference
- .07 EISS
AEISS Regular type cross reference
- .08 DENOM FIELD NUMBER
ADF Regular type cross reference
- .09 NUMERATOR FIELD NUMBER
ANF Regular type cross reference
- .12 DM INDICATOR?
ADM Regular type cross reference
- .13 CARDIO INDICATOR?
ACARD Regular type cross reference
- .14 EISS NUMERATOR PIECE
APO MUMPS type cross reference

90242.01 (BGP INDICATORS)

- .01 FY
B Regular type cross reference
AO MUMPS type cross reference
- .06 ORDER TO DISPLAY
AOTOO MUMPS type cross reference

90243 BGP 03 DATA CURRENT

- .01 BEGINNING DATE
B Regular type cross reference
- 9999 COMMUNITIES (90243.09999)
 - .01 COMMUNITIES
B Regular type cross reference
- 99999 PATIENTS EISS (90243.099999)
 - .01 PATIENTS EISS
B Regular type cross reference

90243.01 (BGP 03 Data Previous)

- .01 BEGINNING DATE
B Regular type cross reference
- 9999 COMMUNITIES (90243.10999)
 - .01 COMMUNITIES
B Regular type cross reference
- 99999 PATIENTS EISS (90243.109999)
 - .01 PATIENTS EISS
B Regular type cross reference

90243.02 (BGP 03 Data Baseline)

- .01 BEGINNING DATE
B Regular type cross reference
- 9999 COMMUNITIES (90243.11999)
 - .01 COMMUNITIES
B Regular type cross reference
- 99999 PATIENTS EISS (90243.119999)
 - .01 PATIENTS EISS
B Regular type cross reference

90244.01 BGP 04 INDICATORS

- .01 INDICATOR NUMBER
B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT
AO Regular type cross reference
- .03 TITLE
C Regular type cross reference
- .07 GPRA INDICATOR
GPRA Regular type cross reference
- .08 AREA INDICATOR
AREA Regular type cross reference
- .09 HEDIS INDICATOR
HEDIS Regular type cross reference
- .11 LOCAL REPORT
LOCAL Regular type cross reference

90244.02 BGP 04 INDIVIDUAL INDICATORS

- .01 INDICATOR NUMBER
B Regular type cross reference
- .04 INDICATOR ID
C Regular type cross reference
AA Mumps cross reference
- .05 GPRA INDICATOR
AGPRA Regular type cross reference
- .06 AREA APP
B Regular type cross reference

- .07 EISS
AEISS Regular type cross reference
- .08 DENOMINATOR FIELD
ADF Regular type cross reference
- .09 NUMERATOR FIELD
ANF Regular type cross reference
- .12 DM INDICATOR?
ADM Regular type cross reference
- .13 CARDIO INDICATOR?
ACARD Regular type cross reference
- .14 EISS PIECE
APO Mumps type cross reference
- 1202 WOMEN'S HEALTH
AWH Regular type cross reference
- 1203 PREVENTION RELATED?
APR Regular type cross reference
- 1204 ORDER
ORDER Mumps type cross reference
AB Mumps type cross reference
OR Regular type cross reference

90244.03 (BGP 04 Data Current)

- .01 BEGINNING DATE
B Regular type cross reference
- 9999 COMMUNITIES (90244.119999)
.01 COMMUNITIES
B Regular type cross reference
- 99999 PATIENTS EISS (90243.129999)
.01 PATIENTS EISS
B Regular type cross reference

90244.04 (BGP 04 Data Previous)

- .01 BEGINNING DATE
B Regular type cross reference
- 9999 COMMUNITIES (90244.13999)
.01 COMMUNITIES
B Regular type cross reference
- 99999 PATIENTS EISS (90243.139999)
.01 PATIENTS EISS
B Regular type cross reference

90244.05 (BGP 04 Data Baseline)

- .01 BEGINNING DATE
B Regular type cross reference
- 9999 COMMUNITIES (90244.14999)
.01 COMMUNITIES
B Regular type cross reference

99999 PATIENTS EISS (90244.1499999)
.01 PATIENTS EISS
B Regular type cross reference

90370.01 BGP 04 HEDIS INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.02 ORDER ON SELECTION/REPORT
AO Regular type cross reference
.03 TITLE
C Regular type cross reference
.07 GPRA INDICATOR
GPRA Regular type cross reference
.08 AREA INDICATOR
AREA Regular type cross reference
.09 HEDIS INDICATOR
HEDIS Regular type cross reference
.11 LOCAL REPORT
LOCAL Regular type cross reference

90370.02 BGP 04 HEDIS INDIVIDUAL INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.04 INDICATOR ID
C Regular type cross reference
AA Mumps cross reference
.05 GPRA INDICATOR
AGPRA Regular type cross reference
.06 AREA APP
B Regular type cross reference
.07 EISS
AEISS Regular type cross reference
.08 DENOMINATOR FIELD
ADF Regular type cross reference
.09 NUMERATOR FIELD
ANF Regular type cross reference
.12 DM INDICATOR?
ADM Regular type cross reference
.13 CARDIO INDICATOR?
ACARD Regular type cross reference
.14 EISS PIECE
APO Mumps type cross reference
1205 WOMEN'S HEALTH
AWH Regular type cross reference
1206 PREVENTION RELATED?
APR Regular type cross reference

1207 ORDER
ORDER Mumps type cross reference
AB Mumps type cross reference
OR Regular type cross reference

90370.03 (BGP 04 HEDIS Data Current)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90370.119999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90243.129999)
.01 PATIENTS EISS
B Regular type cross reference

90370.04 (BGP 04 HEDIS Data Previous)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90370.13999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90243.139999)
.01 PATIENTS EISS
B Regular type cross reference

90370.05 (BGP 04 HEDIS Data Baseline)

.01 BEGINNING DATE
B Regular type cross reference
9999 COMMUNITIES (90370.14999)
.01 COMMUNITIES
B Regular type cross reference
99999 PATIENTS EISS (90370.1499999)
.01 PATIENTS EISS
B Regular type cross reference

90371.01 BGP 05 INDICATORS

.01 INDICATOR NUMBER
B Regular type cross reference
.02 ORDER ON SELECTION/REPORT
AO Regular type cross reference
.03 TITLE
C Regular type cross reference
.07 GPRA INDICATOR
GPRA Regular type cross reference
.08 AREA INDICATOR
AREA Regular type cross reference
.09 HEDIS INDICATOR

- HEDIS Regular type cross reference
- .11 LOCAL REPORT
- LOCAL Regular type cross reference

90371.02 BGP 05 INDIVIDUAL INDICATORS

- .01 INDICATOR NUMBER
 - B Regular type cross reference
- .04 INDICATOR ID
 - C Regular type cross reference
 - AA Mumps cross reference
- .05 GPRA INDICATOR
 - AGPRA Regular type cross reference
- .06 AREA APP
 - B Regular type cross reference
- .07 EISS
 - AEISS Regular type cross reference
- .08 DENOMINATOR FIELD
 - ADF Regular type cross reference
- .09 NUMERATOR FIELD
 - ANF Regular type cross reference
- .12 DM INDICATOR?
 - ADM Regular type cross reference
- .13 CARDIO INDICATOR?
 - ACARD Regular type cross reference
- .14 EISS PIECE
 - APO Mumps type cross reference
- 1208 WOMEN'S HEALTH
 - AWH Regular type cross reference
- 1209 PREVENTION RELATED?
 - APR Regular type cross reference
- 1210 ORDER
 - ORDER Mumps type cross reference
 - AB Mumps type cross reference
 - OR Regular type cross reference

90371.03 (BGP 05 Data Current)

- .01 BEGINNING DATE
 - B Regular type cross reference
- 9999 COMMUNITIES (90371.119999)
 - .01 COMMUNITIES
 - B Regular type cross reference
- 99999 PATIENTS EISS (90243.129999)
 - .01 PATIENTS EISS
 - B Regular type cross reference

90371.04 (BGP 05 Data Previous)

- .01 BEGINNING DATE
B Regular type cross reference
- 9999 COMMUNITIES (90371.13999)
 - .01 COMMUNITIES
B Regular type cross reference
- 99999 PATIENTS EISS (90243.139999)
 - .01 PATIENTS EISS
B Regular type cross reference

90371.05 (BGP 05 Data Baseline)

- .01 BEGINNING DATE
B Regular type cross reference
- 9999 COMMUNITIES (90371.14999)
 - .01 COMMUNITIES
B Regular type cross reference
- 99999 PATIENTS EISS (90371.1499999)
 - .01 PATIENTS EISS
B Regular type cross reference

90371.08 BGP 05 TAXONOMIES

- .01 TAXONOMY NAME
B Regular type cross reference
- 1200 REPORTS
B Regular type cross reference

90372.01 BGP 05 HEDIS INDICATORS

- .01 INDICATOR NUMBER
B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT
AO Regular type cross reference
- .03 TITLE
C Regular type cross reference
- .07 GPRA INDICATOR
GPRA Regular type cross reference
- .08 AREA INDICATOR
AREA Regular type cross reference
- .09 HEDIS INDICATOR
HEDIS Regular type cross reference
- .11 LOCAL REPORT
LOCAL Regular type cross reference

90372.02 BGP 05 HEDIS INDIVIDUAL INDICATORS

- .01 INDICATOR NUMBER
B Regular type cross reference
- .04 INDICATOR ID

C Regular type cross reference
 AA Mumps cross reference
 .05 GPRA INDICATOR
 AGPRA Regular type cross reference
 .06 AREA APP
 B Regular type cross reference
 .07 EISS
 AEISS Regular type cross reference
 .08 DENOMINATOR FIELD
 ADF Regular type cross reference
 .09 NUMERATOR FIELD
 ANF Regular type cross reference
 .12 DM INDICATOR?
 ADM Regular type cross reference
 .13 CARDIO INDICATOR?
 ACARD Regular type cross reference
 .14 EISS PIECE
 APO Mumps type cross reference
 1211 WOMEN'S HEALTH
 AWH Regular type cross reference
 1212 PREVENTION RELATED?
 APR Regular type cross reference
 1213 ORDER
 ORDER Mumps type cross reference
 AB Mumps type cross reference
 OR Regular type cross reference

90372.03 (BGP 05 HEDIS Data Current)

.01 BEGINNING DATE
 B Regular type cross reference
 9999 COMMUNITIES (90372.119999)
 .01 COMMUNITIES
 B Regular type cross reference
 99999 PATIENTS EISS (90243.129999)
 .01 PATIENTS EISS
 B Regular type cross reference

90372.04 (BGP 05 HEDIS Data Previous)

.01 BEGINNING DATE
 B Regular type cross reference
 9999 COMMUNITIES (90372.13999)
 .01 COMMUNITIES
 B Regular type cross reference
 99999 PATIENTS EISS (90243.139999)
 .01 PATIENTS EISS
 B Regular type cross reference

90372.05 (BGP 05 HEDIS Data Baseline)

- .01 BEGINNING DATE
 - B Regular type cross reference
- 9999 COMMUNITIES (90372.14999)
 - .01 COMMUNITIES
 - B Regular type cross reference
- 99999 PATIENTS EISS (90372.1499999)
 - .01 PATIENTS EISS
 - B Regular type cross reference

90372.08 BGP 05 GUI REPORT OUTPUT

- .01 NAME
 - B Regular type cross reference
- .02 USER WHO REQUESTED REPORT
 - AUSR Regular type cross reference
- .03 START TIME
 - AC Regular type cross reference

90373.01 BGP 05 CMS INDICATORS

- .01 INDICATOR NUMBER
 - B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT
 - AO Regular type cross reference
- .03 TITLE
 - C Regular type cross reference

90373.02 BGP 05 ELDER INDIVIDUAL INDICATORS

- .01 INDICATOR
 - B Regular type cross reference
- .04 INDICATOR ID
 - C Regular type cross reference
- 1204 ORDER
 - ABC MUMPS
 - AB MUMPS
 - OR Regular type cross reference
 - AP MUMPS

90373.03 BGP 05 ELDER DATA CURRENT

- .01 BEGINNING DATE
 - B Regular type cross reference

90373.04 BGP 05 ELDER DATA PREVIOUS

- .01 BEGINNING DATE
 - B Regular type cross reference

90373.05 BGP 05 ELDER DATA BASELINE

- .01 BEGINNING DATE
 - B Regular type cross reference

90373.06 BGP 05 ELDER INDICATORS

- .01 INDICATOR NUMBER
 - B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT
 - AO Regular type cross reference
- .03 TITLE
 - C Regular type cross reference
- .07 NATIONAL REPORT
 - GPRA Regular type cross reference
- .11 LOCAL REPORT
 - LOCAL Regular type cross reference
- .01 DENOMINATOR DEFINITIONS
 - B Regular type cross reference
- .01 NUMERATOR DEFINITIONS
 - B Regular type cross reference

90374.01 BGP 06 INDICATORS

- .01 INDICATOR NUMBER
 - B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT
 - AO Regular type cross reference
- .03 TITLE
 - C Regular type cross reference
- .07 NATIONAL REPORT
 - GPRA Regular type cross reference
- .11 LOCAL REPORT
 - LOCAL Regular type cross reference
- .01 DENOMINATOR DEFINITIONS
 - B Regular type cross reference
- .01 NUMERATOR DEFINITIONS
 - B Regular type cross reference

90374.02 BGP 06 INDIVIDUAL INDICATORS

- .01 INDICATOR
 - B Regular type cross reference
- .04 INDICATOR ID
 - C Regular type cross reference
- .05 NATIONAL
 - AGPRA Regular type cross reference
- .12 DM INDICATOR?
 - ADM Regular type cross reference
- .13 CARDIO INDICATOR?
 - ACARD Regular type cross reference

1202 WOMEN'S HEALTH
 AWH Regular type cross reference
1203 ELDER RELATED
 AEL Regular type cross reference
1204 ORDER
 ABC MUMPS
 AB MUMPS
 OR Regular type cross reference

90374.03 BGP 06 DATA CURRENT

.01 BEGINNING DATE
 B Regular type cross reference

90374.04 BGP 06 DATA PREVIOUS

.01 BEGINNING DATE
 B Regular type cross reference

90374.05 BGP 06 DATA BASELINE

.01 BEGINNING DATE
 B Regular type cross reference

90374.06 BGP 06 NATIONAL PATIENT LISTS

.01 POINTER
 B Regular type cross reference

90374.08 BGP 06 TAXONOMIES

.01 TAXONOMY NAME
 B Regular type cross reference
.01 REPORTS
 B Regular type cross reference
.01 COMMUNITIES
 B Regular type cross reference
.01 PATIENTS EISS
 B Regular type cross reference
.01 COMMUNITIES
 B Regular type cross reference
.01 PATIENTS EISS
 B Regular type cross reference
.01 COMMUNITIES
 B Regular type cross reference
.01 PATIENTS EISS
 B Regular type cross reference

90375.01 BGP 06 HEDIS INDICATORS

.01 INDICATOR NUMBER
 B Regular type cross reference

- .02 ORDER ON SELECTION/REPORT
 - AO Regular type cross reference
- .03 TITLE
 - C Regular type cross reference
- .09 HEDIS INDICATOR
 - HEDIS Regular type cross reference
- .01 DENOMINATOR DEFINITIONS
 - B Regular type cross reference
- .01 NUMERATOR DEFINITIONS
 - B Regular type cross reference

90375.02 BGP 06 HEDIS INDIVIDUAL INDICATORS

- .01 INDICATOR
 - B Regular type cross reference
- .04 INDICATOR ID
 - C Regular type cross reference
 - AA MUMPS
- .08 DENOM FIELD NUMBER
 - ADF Regular type cross reference
- .09 NUMERATOR FIELD NUMBER
 - ANF Regular type cross reference
- 1204 ORDER
 - AB MUMPS
 - OR Regular type cross reference

90375.03 BGP 06 HEDIS DATA CURRENT

- .01 BEGINNING DATE
 - B Regular type cross reference

90375.04 BGP 06 HEDIS DATA PREVIOUS

- .01 BEGINNING DATE
 - B Regular type cross reference

90375.05 BGP 06 HEDIS DATA BASELINE

- .01 BEGINNING DATE
 - B Regular type cross reference

90375.08 BGP 06 GUI REPORT OUTPUT

- .01 NAME
 - B Regular type cross reference
- .03 START TIME
 - AC Regular type cross reference
 - AUSR MUMPS
- .01 COMMUNITIES
 - B Regular type cross reference
- .01 COMMUNITIES

- B Regular type cross reference
- .01 COMMUNITIES
- B Regular type cross reference

90375.99 BGP BMX ADO SCHEMA

- .01 NAME
- B Regular type cross reference

90376.01 BGP 06 CMS INDICATORS

- .01 INDICATOR NUMBER
- B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT
- AO Regular type cross reference
- .03 TITLE
- C Regular type cross reference

90376.02 BGP 06 ELDER INDIVIDUAL INDICATORS

- .01 INDICATOR
- B Regular type cross reference
- .04 INDICATOR ID
- C Regular type cross reference
- 1204 ORDER
- ABC MUMPS
- AB MUMPS
- OR Regular type cross reference
- AP MUMPS

90376.03 BGP 06 ELDER DATA CURRENT

- .01 BEGINNING DATE
- B Regular type cross reference

90376.04 BGP 06 ELDER DATA PREVIOUS

- .01 BEGINNING DATE
- B Regular type cross reference

90376.05 BGP 06 ELDER DATA BASELINE

- .01 BEGINNING DATE
- B Regular type cross reference

90376.06 BGP 06 ELDER INDICATORS

- .01 INDICATOR NUMBER
- B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT
- AO Regular type cross reference
- .03 TITLE
- C Regular type cross reference

- .07 NATIONAL REPORT
 - GPRA Regular type cross reference
- .11 LOCAL REPORT
 - LOCAL Regular type cross reference
- .01 DENOMINATOR DEFINITIONS
 - B Regular type cross reference
- .01 NUMERATOR DEFINITIONS
 - B Regular type cross reference
- .01 COMMUNITIES
 - B Regular type cross reference
- .01 PATIENTS EISS
 - B Regular type cross reference
- .01 COMMUNITIES
 - B Regular type cross reference
- .01 PATIENTS EISS
 - B Regular type cross reference
- .01 COMMUNITIES
 - B Regular type cross reference
- .01 PATIENTS EISS
 - B Regular type cross reference

90530.01 BGP 07 INDICATORS

- .01 INDICATOR NUMBER
 - B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT
 - AO Regular type cross reference
- .03 TITLE
 - C Regular type cross reference
- .07 NATIONAL REPORT
 - GPRA Regular type cross reference
- .11 LOCAL REPORT
 - LOCAL Regular type cross reference
- .01 DENOMINATOR DEFINITIONS
 - B Regular type cross reference
- .01 NUMERATOR DEFINITIONS
 - B Regular type cross reference

90530.02 BGP 07 INDIVIDUAL INDICATORS

- .01 INDICATOR
 - B Regular type cross reference
- .04 INDICATOR ID
 - C Regular type cross reference
- .05 NATIONAL
 - AGPRA Regular type cross reference
- .12 DM INDICATOR?
 - ADM Regular type cross reference

- .13 CARDIO INDICATOR?
 - ACARD Regular type cross reference
- 1202 WOMEN'S HEALTH
 - AWH Regular type cross reference
- 1203 ELDER RELATED
 - AEL Regular type cross reference
- 1204 ORDER
 - ABC MUMPS
 - AB MUMPS
 - OR Regular type cross reference

90530.03 BGP 07 DATA CURRENT

- .01 BEGINNING DATE
 - B Regular type cross reference

90530.04 BGP 07 DATA PREVIOUS

- .01 BEGINNING DATE
 - B Regular type cross reference

90530.05 BGP 07 DATA BASELINE

- .01 BEGINNING DATE
 - B Regular type cross reference

90530.06 BGP 07 NATIONAL PATIENT LISTS

- .01 POINTER
 - B Regular type cross reference

90530.08 BGP 07 TAXONOMIES

- .01 TAXONOMY NAME
 - B Regular type cross reference
- .01 REPORTS
 - B Regular type cross reference
- .01 COMMUNITIES
 - B Regular type cross reference
- .01 PATIENTS EISS
 - B Regular type cross reference
- .01 COMMUNITIES
 - B Regular type cross reference
- .01 PATIENTS EISS
 - B Regular type cross reference
- .01 COMMUNITIES
 - B Regular type cross reference
- .01 PATIENTS EISS
 - B Regular type cross reference

90531.01 BGP 07 HEDIS INDICATORS

- .01 INDICATOR NUMBER
 - B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT
 - AO Regular type cross reference
- .03 TITLE
 - C Regular type cross reference
- .09 HEDIS INDICATOR
 - HEDIS Regular type cross reference
- .01 DENOMINATOR DEFINITIONS
 - B Regular type cross reference
- .01 NUMERATOR DEFINITIONS
 - B Regular type cross reference

90531.02 BGP 07 HEDIS INDIVIDUAL INDICATORS

- .01 INDICATOR
 - B Regular type cross reference
- .04 INDICATOR ID
 - C Regular type cross reference
 - AA MUMPS
- .08 DENOM FIELD NUMBER
 - ADF Regular type cross reference
- .09 NUMERATOR FIELD NUMBER
 - ANF Regular type cross reference
- 1204 ORDER
 - AB MUMPS
 - OR Regular type cross reference

90531.03 BGP 07 HEDIS DATA CURRENT

- .01 BEGINNING DATE
 - B Regular type cross reference

90531.04 BGP 07 HEDIS DATA PREVIOUS

- .01 BEGINNING DATE
 - B Regular type cross reference

90531.05 BGP 07 HEDIS DATA BASELINE

- .01 BEGINNING DATE
 - B Regular type cross reference

90531.08 BGP 07 GUI REPORT OUTPUT

- .01 NAME
 - B Regular type cross reference
- .03 START TIME
 - AC Regular type cross reference
 - AUSR MUMPS

- .01 COMMUNITIES
 - B Regular type cross reference
- .01 COMMUNITIES
 - B Regular type cross reference
- .01 COMMUNITIES
 - B Regular type cross reference

90531.99 BGP BMX ADO SCHEMA

- .01 NAME
 - B Regular type cross reference

90532.01 BGP 07 CMS INDICATORS

- .01 INDICATOR NUMBER
 - B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT
 - AO Regular type cross reference
- .03 TITLE
 - C Regular type cross reference

90532.02 BGP 07 ELDER INDIVIDUAL INDICATORS

- .01 INDICATOR
 - B Regular type cross reference
- .04 INDICATOR ID
 - C Regular type cross reference
- 1204 ORDER
 - ABC MUMPS
 - AB MUMPS
 - OR Regular type cross reference
 - AP MUMPS

90532.03 BGP 07 ELDER DATA CURRENT

- .01 BEGINNING DATE
 - B Regular type cross reference

90532.04 BGP 07 ELDER DATA PREVIOUS

- .01 BEGINNING DATE
 - B Regular type cross reference

90532.05 BGP 07 ELDER DATA BASELINE

- .01 BEGINNING DATE
 - B Regular type cross reference

90532.06 BGP 07 ELDER INDICATORS

- .01 INDICATOR NUMBER
 - B Regular type cross reference
- .02 ORDER ON SELECTION/REPORT

AO Regular type cross reference
 .03 TITLE
 C Regular type cross reference
 .07 NATIONAL REPORT
 GPRA Regular type cross reference
 .11 LOCAL REPORT
 LOCAL Regular type cross reference
 .01 DENOMINATOR DEFINITIONS
 B Regular type cross reference
 .01 NUMERATOR DEFINITIONS
 B Regular type cross reference
 .01 COMMUNITIES
 B Regular type cross reference
 .01 PATIENTS EISS
 B Regular type cross reference
 .01 COMMUNITIES
 B Regular type cross reference
 .01 PATIENTS EISS
 B Regular type cross reference
 .01 COMMUNITIES
 B Regular type cross reference
 .01 PATIENTS EISS
 B Regular type cross reference

4.4 Table File

Note: The following 22 files have over 1,000 fields per file. These files are used to hold calculated data values for the GPRA reports. The table files are not listed here to conserve space in this document. If you want a listing, then use the ^XBFLD or FM data dictionary listing utility to list the data dictionary in detail.

File: 90240.01 BGP IHS GPRA DATA 02

Global: ^BGPD(

FILE: 90243 BGP 03 DATA CURRENT

GLOBAL: ^BGPDPDC(

FILE: 90243.01 BGP 03 DATA PREVIOUS

GLOBAL: ^BGPDPDP(

FILE: 90243.02 BGP 03 DATA BASELINE

GLOBAL: ^BGPDPDB(

FILE: 90244.03 BGP 04 DATA CURRENT

GLOBAL: ^BGPDPDCF(

FILE: 90244.04 BGP 04 DATA PREVIOUS
GLOBAL: ^BGPDPDF(

FILE: 90244.05 BGP 04 DATA BASELINE
GLOBAL: ^BGPHEDBF(

FILE: 90370.03 BGP 04 HEDIS DATA CURRENT
GLOBAL: ^BGPHEDCF(

FILE: 90370.04 BGP 04 HEDIS DATA PREVIOUS
GLOBAL: ^BGPHEDPF(

FILE: 90370.05 BGP 04 HEDIS DATA BASELINE
GLOBAL: ^BGPHEDBV(

FILE: 90371.03 BGP 05 DATA CURRENT
GLOBAL: ^BGPDPDCV(

FILE: 90371.04 BGP 05 DATA PREVIOUS
GLOBAL: ^BGPDPDPV(

FILE: 90371.05 BGP 05 DATA BASELINE
GLOBAL: ^BGPHEDBV(

FILE: 90372.03 BGP 05 HEDIS DATA CURRENT
GLOBAL: ^BGPHEDCV(

FILE: 90372.04 BGP 05 HEDIS DATA PREVIOUS
GLOBAL: ^BGPHEDPV(

FILE: 90372.05 BGP 05 HEDIS DATA BASELINE
GLOBAL: ^BGPHEDBV(

FILE: 90373.03 BGP 05 ELDER DATA CURRENT
GLOBAL: ^BGPPELDCV(

FILE: 90373.04 BGP 05 ELDER DATA PREVIOUS
GLOBAL: ^BGPPELDPV(

FILE: 90373.05 BGP 05 ELDER DATA BASELINE
GLOBAL: ^BGPPELDBV(

FILE: 90374.03 BGP 06 DATA CURRENT
GLOBAL: ^BGPDPDCV(

FILE: 90374.04 BGP 06 DATA PREVIOUS
GLOBAL: ^BGPDPDPV(

FILE: 90374.05 BGP 06 DATA BASELINE

GLOBAL: ^BGPHEDBV(

**FILE: 90375.03 BGP 06 HEDIS DATA CURRENT
GLOBAL: ^BGPHEDCV(**

**FILE: 90375.04 BGP 06 HEDIS DATA PREVIOUS
GLOBAL: ^BGPHEDPV(**

**FILE: 90375.05 BGP 06 HEDIS DATA BASELINE
GLOBAL: ^BGPHEDBV(**

**FILE: 90376.03 BGP 06 ELDER DATA CURRENT
GLOBAL: ^BGPELDCV(**

**FILE: 90376.04 BGP 06 ELDER DATA PREVIOUS
GLOBAL: ^BGPELDPV(**

**FILE: 90376.05 BGP 06 ELDER DATA BASELINE
GLOBAL: ^BGPELDBV(**

**FILE: 90530.03 BGP 07 DATA CURRENT
GLOBAL: ^BGPGPDCV(**

**FILE: 90530.04 BGP 07 DATA PREVIOUS
GLOBAL: ^BGPGPDPV(**

**FILE: 90530.05 BGP 07 DATA BASELINE
GLOBAL: ^BGPHEDBV(**

**FILE: 90531.03 BGP 07 HEDIS DATA CURRENT
GLOBAL: ^BGPHEDCV(**

**FILE: 90531.04 BGP 07 HEDIS DATA PREVIOUS
GLOBAL: ^BGPHEDPV(**

**FILE: 90531.05 BGP 07 HEDIS DATA BASELINE
GLOBAL: ^BGPHEDBV(**

**FILE: 90532.03 BGP 07 ELDER DATA CURRENT
GLOBAL: ^BGPELDCV(**

**FILE: 90532.04 BGP 07 ELDER DATA PREVIOUS
GLOBAL: ^BGPELDPV(**

**FILE: 90532.05 BGP 07 ELDER DATA BASELINE
GLOBAL: ^BGPELDBV(**

File: 90241.01 BGP CONTROL FILE**Global: ^BGPCTRL(**

Field #	Field Name	Subscript	Piece	Type
.01	FY	D0,0	1	F
1101	Clinics for Clinical User (90241.0111)			
.01	Clinic Code	D0,11,D1,0	1	F
.019	Clinic Name	Computed		
1201	Clinic Codes for 2ND Visit (90241.0112)			
.01	Clinic Codes for 2ND Visit	D0,12,D1,0	1	F

FILE: 90241.02 BGP SITE PARAMETERS**GLOBAL: ^BGPSITE(**

Field #	Field Name	Subscript	Piece	Type
.01	LOCATION	D0,0	1	P
.02	HOME LOCATION	"	2	P
.03	AREA EXPORT?	"	3	S
.04	EISS EXPORT?	"	4	S

FILE: 90242 BGP INDIVIDUAL INDICATOR DEFINITIONS**GLOBAL: ^BGPINDC(**

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.02	INDICATOR SECOND LEVEL	D0,0	2	F
.03	INDICATOR THIRD LEVEL	"	3	F
.04	INDICATOR ID	"	4	F
.05	GPRA	"	5	S
.06	AREA APP	"	6	S
.07	EISS	"	7	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.11	EISS DENOM PIECE	"	11	N
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS NUMERATOR PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F

Field #	Field Name	Subscript	Piece	Type
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1100	DENOMINATOR (90242.11)			
.01	DENOMINATOR	D0,11,D1,0	1	W
1101	GPRA DENOM LINE TITLE	D0,11	1	F
1102	GPRA DENOM TITLE 2	"	2	F
1103	GPRA DENOM TITLE 3	"	3	F
1104	AREA DENOM TITLE 1	"	4	F
1105	AREA DENOM TITLE 2	"	5	F
1106	AREA DENOM TITLE 3	"	6	F
1200	NUMERATRO (90242.12)			
.01	NUMERATRO	D0,12,D1,0	1	W
1300	LOGIC (90242.13)			
.01	LOGIC	D0,13,D1,0	1	W

FILE: 90242.01 BGP INDICATORS
GLOBAL: ^BGPIND(

Field #	Field Name	Subscript	Piece	Type
.01	FY	D0,0	1	F
.02	INDICATOR NUMBER	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	GPRA INDICATOR	"	7	S
.08	AREA INDICATOR	"	8	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	LOGIC USED (90242.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	DESCRIPTION (90242.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	LOCAL REPORT PRINT HEADER (90242.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90242.0153)			

Field #	Field Name	Subscript	Piece	Type
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
5400	GPRA REPORT PRINT HEADER (90242.0154)			
.01	GPRA REPORT PRINT HEADER	D0,54,D1,0	1	W
5500	AREA DIRECTORS REPORT HEADER (90242.0155)			
.01	AREA DIRECTORS REPORT HEADER	D0,55,D1,0	1	W
7100	REPORT LIST 1ST PAGE (90242.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90242.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90244.01 BGP 04 INDICATORS
GLOBAL: ^BGPINDF(

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	GPRA INDICATOR	"	7	S
.08	AREA INDICATOR	"	8	S
.09	HEDIS INDICATOR	"	9	S
.11	LOCAL REPORT	"	11	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90244.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90244.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90244.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90244.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90244.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N

Field #	Field Name	Subscript	Piece	Type
.02	REPORTS	"	2	F
1	DENOM DESCRIPTION (90244.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D 2,0	1	W
6200	NUMERATOR DEFINITIONS (90244.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
1	NUMERATOR DESCRIPTION (90244.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D 2,0	1	W
7100	REPORT LIST 1ST PAGE (90244.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90244.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90244.02 BGP 04 INDIVIDUAL INDICATORS
GLOBAL: ^BGPINDFC(

Field #	Field Name	Subscript	Piece	Type
.001	NUMBER			N
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.02	INDICATOR SECOND LEVEL	D0,0	2	F
.03	INDICATOR THIRD LEVEL	"	3	F
.04	INDICATOR ID	"	4	F
.05	GPRA	"	5	S
.06	AREA APP	"	6	S
.07	EISS	"	7	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.11	EISS DENOM PIECE	"	11	N
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS NUMERATOR PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F

Field #	Field Name	Subscript	Piece	Type
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1101	GPRA DENOM LINE TITLE	D0,11	1	F
1102	GPRA DENOM TITLE 2	"	2	F
1103	GPRA DENOM TITLE 3	"	3	F
1104	AREA DENOM TITLE 1	"	4	F
1105	AREA DENOM TITLE 2	"	5	F
1106	AREA DENOM TITLE 3	"	6	F
1201	HEDIS?	D0,12	1	S
1202	WOMEN'S HEALTH	"	2	S
1203	PREVENTION RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	NATIONAL PREVIOUS	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N

FILE: 90370.01 BGP 04 HEDIS INDICATORS
GLOBAL: ^BGPHEIF(

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.09	HEDIS INDICATOR	"	9	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90244.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90244.0151)			

Field #	Field Name	Subscript	Piece	Type
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90244.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90244.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90244.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
1	DENOM DESCRIPTION (90244.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D 2,0	1	W
6200	NUMERATOR DEFINITIONS (90244.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
1	NUMERATOR DESCRIPTION (90244.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D 2,0	1	W
7100	REPORT LIST 1ST PAGE (90244.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90244.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90370.02 BGP 04 HEDIS INDIVIDUAL INDICATORS
GLOBAL: ^BGPHEIIF(

Field #	Field Name	Subscript	Piece	Type
.001	NUMBER			N
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.02	INDICATOR SECOND LEVEL	D0,0	2	F
.03	INDICATOR THIRD LEVEL	"	3	F
.04	INDICATOR ID	"	4	F
.07	EISS	"	7	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.11	EISS DENOM PIECE	"	11	N
.13	CARDIO INDICATOR?	"	13	S
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F

Field #	Field Name	Subscript	Piece	Type
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1101	GPRA DENOM LINE TITLE	D0,11	1	F
1102	GPRA DENOM TITLE 2	"	2	F
1103	GPRA DENOM TITLE 3	"	3	F
1104	AREA DENOM TITLE 1	"	4	F
1105	AREA DENOM TITLE 2	"	5	F
1106	AREA DENOM TITLE 3	"	6	F
1201	HEDIS?	D0,12	1	S
1202	WOMEN'S HEALTH	"	2	S
1203	PREVENTION RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	NATIONAL PREVIOUS	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N

FILE: 90371.01 BGP 05 INDICATORS
GLOBAL: ^BGPINDV(

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	NATIONAL REPORT	"	7	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K

Field #	Field Name	Subscript	Piece	Type
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90371.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90371.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90371.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90371.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90371.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	“	2	F
.03	SEL	“	3	F
1	DENOM DESCRIPTION (90371.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D 2,0	1	W
6200	NUMERATOR DEFINITIONS (90371.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	NUMERATOR DESCRIPTION (90371.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D 2,0	1	W
7100	REPORT LIST 1ST PAGE (90371.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90371.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1, 0	1	W

FILE: 90371.02 BGP 05 INDIVIDUAL INDICATORS
GLOBAL: ^BGPINDVC(

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.04	INDICATOR ID	D0,0	4	F
.05	NATIONAL	"	5	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
.22	ADD BANG	"	22	S
.23	NATIONAL LIST	"	23	S
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1202	WOMEN'S HEALTH	D0,12	2	S
1203	ELDER RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	GOAL 2010	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N
1407	TITLE 2	"	7	F

**FILE: 90371.06 BGP 05 NATIONAL PATIENT LISTS
GLOBAL: ^BGPVNPL(**

Field #	Field Name	Subscript	Piece	Type
.01	POINTER	D0,0	1	P
.02	INDIVIDUAL INDICATOR ID	"	2	F
.03	LIST TITLE	"	3	F
1101	LIST TITLE FOR REPORT (90371.61101)			
.01	LIST TITLE FOR REPORT	D0,11,D1,0	1	W
1201	SET LIST VALUE	D0,12		K

**FILE:90371.08 BGP 05 TAXONOMIES
GLOBAL: ^BGPTAXV(**

Field #	Field Name	Subscript	Piece	Type
.01	TAXONOMY NAME	D0,0	1	F
.02	TYPE	"	2	S
.03	LAST CHECKED/EXPORTED	"	3	D
.04	EDITABLE	"	4	S
1100	DESCRIPTION (90371.0811)	"	5	W
1200	REPORTS (90371.0812)	"	1	
.01	REPORTS	"	1	S

**FILE: 90372.01 BGP 05 HEDIS INDICATORS
GLOBAL: ^BGPHEIV(**

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.09	HEDIS INDICATOR	"	9	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90372.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W

Field #	Field Name	Subscript	Piece	Type
5100	INDICATOR DESCRIPTION (90372.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90372.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90372.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90372.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
1	DENOM DESCRIPTION (90372.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D 2,0	1	W
6200	NUMERATOR DEFINITIONS (90372.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
1	NUMERATOR DESCRIPTION (90372.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D 2,0	1	W
7100	REPORT LIST 1ST PAGE (90372.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90372.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

**FILE: 90372.02 BGP 05 HEDIS INDIVIDUAL INDICATORS
GLOBAL: ^BGPHEIIV(**

Field #	Field Name	Subscript	Piece	Type
.001	NUMBER			N
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.02	CORRESPONDING GPRA	D0,0	2	F
.03	CORRESPONDING GPRA ORDER	"	3	F
.04	INDICATOR ID	"	4	F
.07	EISS	"	7	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.11	EISS DENOM PIECE	"	11	N

Field #	Field Name	Subscript	Piece	Type
.13	CARDIO INDICATOR?	"	13	S
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1101	GPRA DENOM LINE TITLE	D0,11	1	F
1102	GPRA DENOM TITLE 2	"	2	F
1103	GPRA DENOM TITLE 3	"	3	F
1104	AREA DENOM TITLE 1	"	4	F
1105	AREA DENOM TITLE 2	"	5	F
1106	AREA DENOM TITLE 3	"	6	F
1201	HEDIS?	D0,12	1	S
1202	WOMEN'S HEALTH	"	2	S
1203	PREVENTION RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	NATIONAL PREVIOUS	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N

**FILE: 90373.01 BGP 05 CMS INDICATORS
GLOBAL: ^BGPCMSV(**

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.06	ORDER TO DISPLAY	"	6	N
1	M CODE TO EXECUTE	D0,1		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
5100	INDICATOR DESCRIPTION (90372.0151)			

Field #	Field Name	Subscript	Piece	Type
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	LIST 2 HEADER (90372.0152)			
.01	LIST 2 HEADER	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90372.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
5400	LIST 1 HEADER (90372.0154)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
5500	LIST 1 TITLE (90372.0155)			
.01	LIST 1 TITLE	D0,72,D1,0	1	W

**FILE: 90372.08 BGP 05 GUI REPORT OUTPUT
GLOBAL: ^BGPGUI(**

Field #	Field Name	Subscript	Piece	Type
.01	NAME	D0,0	1	F
.02	USER WHO REQUESTED REPORT	"	2	P
.03	START TIME	"	3	D
.04	END TIME	"	4	D
.05	TYPE OF REPORT	"	5	F
.06	STATUS	"	6	S
.07	TYPE OF OUTPUT	"	7	S
1100	PRINT FORMATTED OUTPUT			
.01	OUTPUT	D0,11,D1,0	1	W
1200	DELIMITED FORMATTED OUTPUT	D0,12,d1,0	1	W

**FILE: 90373.06 BGP 05 ELDER INDICATORS
GLOBAL: ^BGPELIV(**

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	NATIONAL REPORT	"	7	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K

Field #	Field Name	Subscript	Piece	Type
1100	FY LOGIC (90371.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90371.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90371.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90371.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90371.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	DENOM DESCRIPTION (90371.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D 2,0	1	W
6200	NUMERATOR DEFINITIONS (90371.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	NUMERATOR DESCRIPTION (90371.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D 2,0	1	W
7100	REPORT LIST 1ST PAGE (90371.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90371.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90373.02 BGP 05 ELDER INDIVIDUAL INDICATORS
GLOBAL: ^BGPELIIC(

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.04	INDICATOR ID	D0,0	4	F
.05	NATIONAL	"	5	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N

Field #	Field Name	Subscript	Piece	Type
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
.22	ADD BANG	"	22	S
.23	NATIONAL LIST	"	23	S
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1202	WOMEN'S HEALTH	D0,12	2	S
1203	ELDER RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	GOAL 2010	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N
1407	TITLE 2	"	7	F

**FILE: 90374.01 BGP 06 INDICATORS
GLOBAL: ^BGPINDS(**

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	NATIONAL REPORT	"	7	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K

Field #	Field Name	Subscript	Piece	Type
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90371.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90371.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90371.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90371.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90371.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	“	2	F
.03	SEL	“	3	F
1	DENOM DESCRIPTION (90371.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D 2,0	1	W
6200	NUMERATOR DEFINITIONS (90371.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	NUMERATOR DESCRIPTION (90371.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D 2,0	1	W
7100	REPORT LIST 1ST PAGE (90371.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90371.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

FILE: 90374.02 BGP 06 INDIVIDUAL INDICATORS
GLOBAL: ^BGPINDSC(

FIELD #	Field Name	Subscript	Piece	Type
.01	INDICATOR	D0,0	1	P
.04	INDICATOR ID	D0,0	4	F
.05	NATIONAL	"	5	S
.06	GRANT1 PIECE	D0,0	6	N
.07	GRANT2 PIECE	D0,0	7	N
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
.22	ADD BANG	"	22	S
.23	NATIONAL LIST	"	23	S
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
3	M CODE FOR COMP PL	D0,3		K
1201	COLUMN TEXT FOR CMP PL	D0,12	1	F
1202	WOMEN'S HEALTH	D0,12	2	S
1203	ELDER RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	GOAL 2010	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N
1407	TITLE 2	"	7	F
1408	GOAL 06	"	8	F
1409	SECOND LINE GOAL	"	9	F
1410	NATIONAL CURRENT 2	"	10	F
1411	SECOND LINE 2010 GOAL	"	11	F

**FILE: 90375.06 BGP 06 NATIONAL PATIENT LISTS
GLOBAL: ^BGPSNPL(**

Field #	Field Name	Subscript	Piece	Type
.01	POINTER	D0,0	1	P
.02	INDIVIDUAL INDICATOR ID	"	2	F
.03	LIST TITLE	"	3	F
1101	LIST TITLE FOR REPORT (90371.61101)			
.01	LIST TITLE FOR REPORT	D0,11,D1,0	1	W
1201	SET LIST VALUE	D0,12		K

**FILE:90375.08 BGP 06 TAXONOMIES
GLOBAL: ^BGPTAXS(**

Field #	Field Name	Subscript	Piece	Type
.01	TAXONOMY NAME	D0,0	1	F
.02	TYPE	"	2	S
.03	LAST CHECKED/EXPORTED	"	3	D
.04	EDITABLE	"	4	S
1100	DESCRIPTION (90371.0811)	"	5	W
1200	REPORTS (90371.0812)	"	1	
.01	REPORTS	"	1	S

**FILE: 90375.01 BGP 05 HEDIS INDICATORS
GLOBAL: ^BGPHEIS(**

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.09	HEDIS INDICATOR	"	9	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90372.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION			

Field #	Field Name	Subscript	Piece	Type
	(90372.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90372.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90372.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90372.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
1	DENOM DESCRIPTION (90372.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D 2,0	1	W
6200	NUMERATOR DEFINITIONS (90372.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
1	NUMERATOR DESCRIPTION (90372.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D 2,0	1	W
7100	REPORT LIST 1ST PAGE (90372.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90372.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

**FILE: 90375.02 BGP 06 HEDIS INDIVIDUAL INDICATORS
GLOBAL: ^BGPHEIIS(**

Field #	Field Name	Subscript	Piece	Type
.001	NUMBER			N
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.02	CORRESPONDING GPRA	D0,0	2	F
.03	CORRESPONDING GPRA ORDER	"	3	F
.04	INDICATOR ID	"	4	F
.07	EISS	"	7	S
.08	DENOM FIELD NUMBER	"	8	N
.09	NUMERATOR FIELD NUMBER	"	9	N
.11	EISS DENOM PIECE	"	11	N

Field #	Field Name	Subscript	Piece	Type
.13	CARDIO INDICATOR?	"	13	S
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1101	GPRA DENOM LINE TITLE	D0,11	1	F
1102	GPRA DENOM TITLE 2	"	2	F
1103	GPRA DENOM TITLE 3	"	3	F
1104	AREA DENOM TITLE 1	"	4	F
1105	AREA DENOM TITLE 2	"	5	F
1106	AREA DENOM TITLE 3	"	6	F
1201	HEDIS?	D0,12	1	S
1202	WOMEN'S HEALTH	"	2	S
1203	PREVENTION RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	NATIONAL PREVIOUS	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N

**FILE: 90376.01 BGP 06 CMS INDICATORS
GLOBAL: ^BGPCMSS(**

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.06	ORDER TO DISPLAY	"	6	N
1	M CODE TO EXECUTE	D0,1		K
3	PRINT M CODE	D0,3		K
4	DELIMITED FILE PRINT	D0,4		K
5100	INDICATOR DESCRIPTION (90372.0151)			

Field #	Field Name	Subscript	Piece	Type
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	LIST 2 HEADER (90372.0152)			
.01	LIST 2 HEADER	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90372.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
5400	LIST 1 HEADER (90372.0154)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
5500	LIST 1 TITLE (90372.0155)			
.01	LIST 1 TITLE	D0,72,D1,0	1	W

FILE: 90375.08 BGP 06 GUI REPORT OUTPUT
GLOBAL: ^BPGUIS(

Field #	Field Name	Subscript	Piece	Type
.01	NAME	D0,0	1	F
.02	USER WHO REQUESTED REPORT	"	2	P
.03	START TIME	"	3	D
.04	END TIME	"	4	D
.05	TYPE OF REPORT	"	5	F
.06	STATUS	"	6	S
.07	TYPE OF OUTPUT	"	7	S
1100	PRINT FORMATTED OUTPUT			
.01	OUTPUT	D0,11,D1,0	1	W
1200	DELIMITED FORMATTED OUTPUT	D0,12,d1,0	1	W

FILE: 90376.06 BGP 06 ELDER INDICATORS
GLOBAL: ^BGPELIS(

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR NUMBER	D0,0	1	N
.02	ORDER ON SELECTION/REPORT	"	2	F
.03	TITLE	"	3	F
.04	SELECTION TITLE	"	4	F
.05	LIST SELECTION TITLE	"	5	F
.06	ORDER TO DISPLAY	"	6	N
.07	NATIONAL REPORT	"	7	S
.11	LOCAL REPORT	"	11	S
.12	NO LIST	"	12	S
1	M CODE TO EXECUTE	D0,1		K
2	SET LIST VALUE	D0,2		K
3	PRINT M CODE	D0,3		K

Field #	Field Name	Subscript	Piece	Type
4	DELIMITED FILE PRINT	D0,4		K
1100	FY LOGIC (90371.0111)			
.01	LOGIC USED	D0,11,D1,0	1	W
5100	INDICATOR DESCRIPTION (90371.0151)			
.01	DESCRIPTION	D0,51,D1,0	1	W
5200	PAST PERFORMANCE (90371.0152)			
.01	REPORT PRINT SECTION	D0,52,D1,0	1	W
5300	REPORT PRINT 2ND PAGE (90371.0153)			
.01	REPORT PRINT 2ND PAGE	D0,53,D1,0	1	W
6100	DENOMINATOR DEFINITIONS (90371.0161)			
.01	DENOMINATOR DEFINITIONS	D0,61,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	DENOM DESCRIPTION (90371.01611)			
.01	DENOM DESCRIPTION	D0,61,D1,1,D 2,0	1	W
6200	NUMERATOR DEFINITIONS (90371.0162)			
.01	NUMERATOR DEFINITIONS	D0,62,D1,0	1	N
.02	REPORTS	"	2	F
.03	SEL	"	3	F
1	NUMERATOR DESCRIPTION (90371.01621)			
.01	NUMERATOR DESCRIPTION	D0,62,D1,1,D 2,0	1	W
7100	REPORT LIST 1ST PAGE (90371.0171)			
.01	REPORT LIST 1ST PAGE	D0,71,D1,0	1	W
7200	LIST PRINT SECOND PAGE (90371.0172)			
.01	LIST PRINT SECOND PAGE	D0,72,D1,0	1	W

**FILE: 90376.02 BGP 06 ELDER INDIVIDUAL INDICATORS
GLOBAL: ^BGPELIIS(**

Field #	Field Name	Subscript	Piece	Type
.01	INDICATOR	D0,0	1	P
.019	INDICATOR	COMPUTED		
.04	INDICATOR ID	D0,0	4	F
.05	NATIONAL	"	5	S
.08	DENOM FIELD NUMBER	"	8	N

Field #	Field Name	Subscript	Piece	Type
.09	NUMERATOR FIELD NUMBER	"	9	N
.12	DM INDICATOR?	"	12	S
.13	CARDIO INDICATOR?	"	13	S
.14	EISS PIECE	"	14	N
.15	LINE TITLE 1	"	15	F
.16	LINE TITLE 2	"	16	F
.17	LOCAL DENOM TITLE 1	"	17	F
.18	LOCAL DENOM TITLE 2	"	18	F
.19	LINE TITLE 3	"	19	F
.21	LOCAL DENOM TITLE 3	"	21	F
.22	ADD BANG	"	22	S
.23	NATIONAL LIST	"	23	S
1	M CODE NUMERATOR	D0,1		K
2	M CODE DENOMINATOR	D0,2		K
1202	WOMEN'S HEALTH	D0,12	2	S
1203	ELDER RELATED	"	3	S
1204	ORDER	"	4	F
1205	ADD GPRA TO 1ST LINE	"	5	S
1206	ADD GPRA TO SECOND LINE	"	6	S
1207	ADD GPRA TO 3RD LINE	"	7	S
1401	SUMMARY?	D0,14	1	S
1402	NATIONAL CURRENT	"	2	F
1403	GOAL 2010	"	3	F
1404	TITLE	"	4	F
1405	CATEGORY	"	5	S
1406	ORDER ON SUMMARY SHEET	"	6	N
1407	TITLE 2	"	7	F

FILE: 90530.01 BGP 07 INDICATORS
GLOBAL: ^BGPINDA(

File #	Field Name	Subscript	Piece	Type
.01	INDICATOR NUMBER	0	1	RNJ3,0X
.02	ORDER ON SELECTION/REPORT	0	2	F
.03	TITLE	0	3	F
.04	SELECTION TITLE	0	4	F
.05	LIST SELECTION TITLE	0	5	F
.06	ORDER TO DISPLAY	0	6	NJ4,0
.07	NATIONAL REPORT	0	7	S
.11	LOCAL REPORT	0	11	S
.12	NO LIST	0	12	S
.13	# DENOMINATORS	0	13	NJ2,0
1	M CODE TO EXECUTE	1	E1,245	K
2	SET LIST VALUE	2	E1,245	K

File #	Field Name	Subscript	Piece	Type
3	PRINT M CODE	3	E1,245	K
4	DELIMITED FILE PRINT	4	E1,245	K
1100	FY LOGIC	11	0	90530.0111
5100	INDICATOR DESCRIPTION	51	0	90530.0151
5200	PAST PERFORMANCE	52	0	90530.0152
5300	REPORT PRINT 2ND PAGE	53	0	90530.0153
5400	NATIONAL ONLY	54	0	90530.54
6100	DENOMINATOR DEFINITIONS	61	0	90530.0161
6200	NUMERATOR DEFINITIONS	62	0	90530.0162
7100	REPORT LIST 1ST PAGE	71	0	90530.0171

FILE: 90530.02 BGP 07 INDIVIDUAL INDICATORS
GLOBAL: ^BGPINDAC(

Field #	Field Name	Subscript	Piece	Type
01	INDICATOR	0	1	RP90530.0 1'
.04	INDICATOR ID	0	4	F
.05	NATIONAL	0	5	S
.06	CRSNT1 PIECE	0	6	NJ3,0
.07	CRSNT2 PIECE	0	7	NJ3,0
.08	DENOM FIELD NUMBER	0	8	NJ10,3
.09	NUMERATOR FIELD NUMBER	0	9	NJ10,3
.11	PAGE BREAK DENOM NGR	0	11	S
.12	DM INDICATOR?	0	12	S
.13	CARDIO INDICATOR?	0	13	S
.14	GPRANT1 PIECE	0	14	NJ5,0
.15	LINE TITLE 1	0	15	F
.16	LINE TITLE 2	0	16	F
.17	LOCAL DENOM TITLE 1	0	17	F
.18	LOCAL DENOM TITLE 2	0	18	F
.19	LINE TITLE 3	0	19	F
.21	LOCAL DENOM TITLE 3	0	21	F
.22	ADD BANG	0	22	S
.23	NATIONAL LIST	0	23	S
.24	PAGE BREAK DENOM	0	24	S
1	M CODE NUMERATOR	1	E1,245	K
2	M CODE DENOMINATOR	2	E1,245	K
3	M CODE FOR COMP PL	3	E1,245	K
1201	COLUMN TEXT FOR CMP PL	12	1	F
1202	WOMEN'S HEALTH	12	2	S
1203	ELDER RELATED	12	3	S
1204	ORDER	12	4	F

Field #	Field Name	Subscript	Piece	Type
1205	ADD GPRA TO 1ST LINE	12	5	S
1206	ADD GPRA TO SECOND LINE	12	6	S
1207	ADD GPRA TO 3RD LINE	12	7	S
1208	PAGE BREAK NUMER COM	12	8	S
1401	SUMMARY?	14	1	S
1402	NATIONAL CURRENT	14	2	F
1403	GOAL 2010	14	3	F
1404	TITLE	14	4	F
1405	CATEGORY	14	5	P90241.03'
1406	ORDER ON SUMMARY SHEET	14	6	NJ4,0
1407	TITLE 2	14	7	F
1408	GOAL 07	14	8	F
1409	SECOND LINE GOAL	14	9	F
1410	NATIONAL CURRENT 2	14	10	F
1411	SECOND LINE 2010 GOAL	14	11	F
1412	TITLE 3	14	12	F
1501	NON GPRA SUMMARY?	15	1	S
1502	NG NATIONAL CURRENT	15	2	F
1503	NG GOAL 2010	15	3	F
1504	NG TITLE	15	4	F
1505	NG CATEGORY	15	5	P90241.03'
1506	NG ORDER ON SUMMARY SHEET	15	6	NJ4,0
1507	NG TITLE 2	15	7	F
1508	NG GOAL 07	15	8	F
1509	NG SECOND LINE GOAL	15	9	F
1510	NG NATIONAL CURRENT 2	15	10	F

**FILE: 90530.06 BGP 07 NATIONAL PATIENT LISTS
GLOBAL: ^BGPANPL(**

Field #	Field Name	Subscript	Piece	Type
.01	POINTER	0	1	RP90530.0 1'
.02	INDIVIDUAL INDICATOR ID	0	2	F
.03	LIST TITLE	0	3	F
1101	LIST TITLE FOR REPORT	11	0	90530.6110 1

**FILE: 90530.08 BGP 07 TAXONOMIES
GLOBAL: ^BGPTAXA(**

Field #	Field Name	Subscript	Piece	Type
.01	TAXONOMY NAME	0	1	RF
.02	TYPE	0	2	S
.03	LAST CHECKED/EXPORTED	0	3	D
.04	EDITABLE?	0	4	S
1100	DESCRIPTION	11	0	90530.0811
1200	REPORTS	12	0	90530.0812 S

**FILE: 90530.09 BGP 07 TAXONOMY TEMP
GLOBAL: ^BGPTAXTA(**

Field #	Field Name	Subscript	Piece	Type
.01	NAME	0	1	RFX
.02	BRIEF DESCRIPTION	0	2	F
.03	LOCATION	0	3	P9999999.06'
.04	NODELETE	0	4	S
.05	CREATOR	0	5	P200'
.06	STARTING DATE FOR PATIENT LIST	0	6	DX
.07	BULLETIN	0	7	P3.6'
.08	CREATED BY TAXONOMY SYSTEM	0	8	S
.09	DATE CREATED	0	9	D
.11	BULLETIN SERVICE CATEGORY	0	11	S
.12	QMAN LINK	0	12	P9009071'
.13	NON-CANONIC	0	13	S
.14	X-REF	0	14	F
.15	FILE	0	15	P1'
.16	RANGE	0	16	S
.17	SEND BULLETIN IF RECENTLY SENT	0	17	S

Field #	Field Name	Subscript	Piece	Type
.18	LOCATION OF ENCOUNTER	0	18	P9999999.06'
.19	SEND BULLETIN-1ST VISIT ONLY	0	19	S
.21	VISIT TYPE	0	21	F
.22	READ ONLY	0	22	S
1101	EXTENDED DESCRIPTION	11	0	90530.10101
2101	VALUE	21	0	90530.11101A
3101	PRE BULLETIN CODE	31	E1,245	K

FILE: 90530.11 BGP 07 XP OUTPUT
GLOBAL: ^BGPXPA(

.01	USER_DATE/TIME RUN	0	1	RF
.02	BEGIN DATE	0	2	D
.03	END DATE	0	3	D
.04	COMMUNITY TAXONOMY	0	4	F
.05	MFI TAXONOMY	0	5	F
1101	DATA RECORDS	11	0	90530.111101A
1201	EXPORT FILE NAMES	12	0	90530.111201A
1301	COMMUNITIES	13	0	90530.111301A

FILE: 90530.15 BGP 07 PAT ED MEASURES
GLOBAL: ^BGPPEIA(

.01	MEASURE NUMBER	0	1	RNJ2,0
.02	TITLE	0	2	F
.03	ORDER TO DISPLAY	0	3	NJ2,0
.04	LIST SELECTION TITLE	0	4	F
1101	FY LOGIC	11	0	90530.151101
6100	DENOMINATOR DEFINITIONS	61	0	90530.1561A
6200	NUMERATOR	62	0	90530.1562A

7100	LIST PRINT 1ST PAGE	71	0	90530.1571
------	------------------------	----	---	------------

FILE: 90531.01 BGP 07 HEDIS INDICATORS
GLOBAL: ^BGPHEIA(

.01	INDICATOR NUMBER	0	1	RNJ3,0X
.02	ORDER ON SELECTION/REPORT	0	2	F
.03	TITLE	0	3	F
.04	SELECTION TITLE	0	4	F
.05	LIST SELECTION TITLE	0	5	F
.06	ORDER TO DISPLAY	0	6	NJ4,0
.09	HEDIS INDICATOR	0	9	S
.12	NO LIST	0	12	S
1	M CODE TO EXECUTE	1	E1,245	K
2	SET LIST VALUE	2	E1,245	K
3	PRINT M CODE	3	E1,245	K
4	DELIMITED FILE PRINT	4	E1,245	K
1100	FY LOGIC	11	0	90531.0111
5100	INDICATOR DESCRIPTION	51	0	90531.0151
5200	PAST PERFORMANCE	52	0	90531.0152
5300	REPORT PRINT 2ND PAGE	53	0	90531.0153
6100	DENOMINATOR DEFINITIONS	61	0	90531.0161
6200	NUMERATOR DEFINITIONS	62	0	90531.0162
7100	REPORT LIST 1ST PAGE	71	0	90531.0171

**FILE: 90531.02 BGP 07 HEDIS INDIVIDUAL INDICATORS
GLOBAL: ^BGPHEIIA(**

.001	NUMBER			NJ7,0
.01	INDICATOR	0	1	RP90531.01'
.019	INDICATOR			CJ8
.02	CORRESPONDING GPRA	0	2	F
.03	CORRESPONDING GPRA ORDER	0	3	F
.04	INDICATOR ID	0	4	F
.07	EISS	0	7	S
.08	DENOM FIELD NUMBER	0	8	NJ10,3
.09	NUMERATOR FIELD NUMBER	0	9	NJ10,3
.11	EISS DENOM PIECE	0	11	NJ4,0
.13	CARDIO INDICATOR?	0	13	S
.15	LINE TITLE 1	0	15	F
.16	LINE TITLE 2	0	16	F
.17	LOCAL DENOM TITLE 1	0	17	F
.18	LOCAL DENOM TITLE 2	0	18	F
.19	LINE TITLE 3	0	19	F
.21	LOCAL DENOM TITLE 3	0	21	F
1	M CODE NUMERATOR	1	E1,245	K
2	M CODE DENOMINATOR	2	E1,245	K
1101	GPRA DENOM LINE TITLE	11	1	F
1102	GPRA DENOM TITLE 2	11	2	F
1103	GPRA DENOM TITLE 3	11	3	F
1104	AREA DENOM TITLE 1	11	4	F
1105	AREA DENOM TITLE 2	11	5	F
1106	AREA DENOM TITLE 3	11	6	F

1201	HEDIS?	12	1	S
1202	WOMEN'S HEALTH	12	2	S
1203	PREVENTION RELATED	12	3	S
1204	ORDER	12	4	F
1205	ADD GPRA TO 1ST LINE	12	5	S
1206	ADD GPRA TO SECOND LINE	12	6	S
1207	ADD GPRA TO 3RD LINE	12	7	S
1401	SUMMARY?	14	1	S
1402	NATIONAL CURRENT	14	2	F
1403	NATIONAL PREVIOUS	14	3	F
1404	TITLE	14	4	F
1405	CATEGORY	14	5	S

FILE: 90531.08 BGP 07 GUI REPORT OUTPUT
GLOBAL: ^BGPGUIA(

.01	NAME	0	1	RF
.02	USER WHO REQUESTED REPORT	0	2	P200'
.03	START TIME	0	3	D
.04	END TIME	0	4	D
.05	TYPE OF REPORT	0	5	F
.06	STATUS	0	6	S
.07	TYPE OF OUTPUT	0	7	S
1100	PRINT FORMATTED OUTPUT	11	0	90531.0811

**FILE: 90532.01 BGP 07 CMS INDICATORS
GLOBAL: ^BGPCMSIA(**

.01	INDICATOR NUMBER	0	1	RNJ3,0X
.02	ORDER ON SELECTION/REPORT	0	2	F
.03	TITLE	0	3	F
.04	SELECTION TITLE	0	4	F
.06	ORDER TO DISPLAY	0	6	NJ4,0
1	M CODE TO EXECUTE	1	E1,245	K
3	PRINT M CODE	3	E1,245	K
4	DELIMITED FILE PRINT	4	E1,245	K
5100	INDICATOR DESCRIPTION	51	0	90532.0151
5200	LIST 2 HEADER	52	0	90532.52
5300	REPORT PRINT 2ND PAGE	53	0	90532.0153
5400	LIST 1 HEADER	54	0	90532.54

**FILE: 90532.02 BGP 07 ELDER INDIVIDUAL INDICATORS
GLOBAL: ^BGPELIIA(**

.01	INDICATOR	0	1	RP90532.06'
.019	INDICATOR			CJ8
.04	INDICATOR ID	0	4	F
.08	DENOM FIELD NUMBER	0	8	NJ10,3
.09	NUMERATOR FIELD NUMBER	0	9	NJ10,3
.15	LINE TITLE 1	0	15	F
.16	LINE TITLE 2	0	16	F
.17	LOCAL DENOM TITLE 1	0	17	F
.18	LOCAL DENOM TITLE 2	0	18	F
.19	LINE TITLE 3	0	19	F
.21	LOCAL DENOM TITLE 3	0	21	F
.22	ADD BANG	0	22	S
1	M CODE	1	E1,245	K

	NUMERATOR			
2	M CODE DENOMINATOR	2	E1,245	K

FILE: 90532.06 BGP 07 ELDER INDICATORS
GLOBAL: ^BGPELIA(

.01	INDICATOR NUMBER	0	1	RNJ3,0X
.02	ORDER ON SELECTION/REPORT	0	2	F
.03	TITLE	0	3	F
.04	SELECTION TITLE	0	4	F
.05	LIST SELECTION TITLE	0	5	F
.06	ORDER TO DISPLAY	0	6	NJ4,0
.07	NATIONAL REPORT	0	7	S
.11	LOCAL REPORT	0	11	S
.12	NO LIST	0	12	S
1	M CODE TO EXECUTE	1	E1,245	K
2	SET LIST VALUE	2	E1,245	K
3	PRINT M CODE	3	E1,245	K
4	DELIMITED FILE PRINT	4	E1,245	K
1100	FY LOGIC	11	0	90532.0611
5100	INDICATOR DESCRIPTION	51	0	90532.0651
5200	PAST PERFORMANCE	52	0	90532.0652
5300	REPORT PRINT 2ND PAGE	53	0	90532.0653
6100	DENOMINATOR DEFINITIONS	61	0	90532.0661
6200	NUMERATOR DEFINITIONS	62	0	90532.0662
7100	REPORT LIST 1ST PAGE	71	0	90532.0671
7200	LIST PRINT SECOND PAGE	72	0	90532.0672

5.0 Archiving and Purging

There is no archiving and purging in this package.

6.0 External Relations

6.1 External Calls

This package calls the following documented entry points:

```
Routine          is Invoked by:
  ^%DT
BGP3UTL1 ,BGP4UTL1 ,BGP5UTL1 ,BGP6UTL1 ,BGP7UTL1 , |dd90240.01
              |dd90240.0197, |dd90243, |dd90243.01, |dd90243.02
              |dd90244.03, |dd90244.04, |dd90244.05, |dd90370.03
              |dd90370.04, |dd90370.05, |dd90371.03, |dd90371.04
              |dd90371.05, |dd90371.08, |dd90372.03
              |dd90372.04, |dd90372.05, |dd90372.08, |dd90373.03
              |dd90373.04, |dd90373.05, |dd90374.03, |dd90374.04
              |dd90374.05, |dd90374.08, |dd90374.09, |dd90375.03
              |dd90375.04, |dd90375.05, |dd90375.08, |dd90376.03
              |dd90376.04, |dd90376.05, |dd90530.03, |dd90530.04
              |dd90530.05, |dd90530.08, |dd90530.09, |dd90530.11
              |dd90531.03, |dd90531.04, |dd90531.05, |dd90531.08
              |dd90532.03, |dd90532.04, |dd90532.05
  COMMA^%DTC
BGP3DP ,BGP3DP1 ,BGP3DP2 ,BGP3DP3 ,BGP3DP4 ,BGP3DP41 ,BGP3DP5
BGP3DP50 ,BGP3DP51 ,BGP3DP6 ,BGP3DP60 ,BGP3DP7 ,BGP3DPH
BGP3PDL ,BGP3PDL1 ,BGP3PDL2 ,BGP3PDL3 ,BGP3PDL4 ,BGP3PDL5
BGP3PDL6 ,BGP3PDL7 ,BGP3PDL8 ,BGP3PDL9 ,BGP3PDM1 ,BGP3PDM2
BGP4DHEP ,BGP4DP ,BGP4DP1
BGP4DP11 ,BGP4DP2 ,BGP4DP3 ,BGP4DP4 ,BGP4DP5 ,BGP4DP50
BGP4DP51 ,BGP4DP6 ,BGP4DP60 ,BGP4DPH ,BGP4HEL ,BGP4HEL1
BGP4HEP1 ,BGP4PDL ,BGP4PDL0 ,BGP4PDL1 ,BGP4PDL2 ,BGP4PDL3
BGP4PDL4 ,BGP4PDL5 ,BGP4PDL7 ,BGP4PDL8 ,BGP4PDL9 ,BGP5CP
BGP5DELP ,BGP5DHEP ,BGP5DP
BGP5DP1 ,BGP5DP11 ,BGP5DP2 ,BGP5DP3 ,BGP5DP4 ,BGP5DP5
BGP5DP50 ,BGP5DP51 ,BGP5DP6 ,BGP5DP60 ,BGP5DP7 ,BGP5DPH
BGP5EL10 ,BGP5EL11 ,BGP5EL12 ,BGP5EL13 ,BGP5EL14 ,BGP5EL15
BGP5EL16 ,BGP5EL17 ,BGP5EL18 ,BGP5EL19 ,BGP5EL1A ,BGP5EL1B
BGP5EL1C ,BGP5EL1D ,BGP5EL1E
BGP5EL1F ,BGP5EL1G ,BGP5EL1H ,BGP5EL1I ,BGP5EL1J ,BGP5EL1K
BGP5EL1L ,BGP5ELL ,BGP5ELL1 ,BGP5ELL2 ,BGP5ELL3 ,BGP5ELL4
BGP5ELL5 ,BGP5ELL6 ,BGP5ELL7 ,BGP5ELL8 ,BGP5ELL9 ,BGP5ELLA
BGP5ELLB ,BGP5ELLC ,BGP5ELLD ,BGP5ELLE ,BGP5ELLF ,BGP5ELLG
BGP5EL LH ,BGP5EL LI
BGP5EL LJ ,BGP5EL LK ,BGP5EL LL ,BGP5EL LM ,BGP5EL LN ,BGP5EL LQ
BGP5EL P1 ,BGP5HEL ,BGP5HEL1 ,BGP5HEP1 ,BGP5PDL ,BGP5PDL0
BGP5PDL1 ,BGP5PDL2 ,BGP5PDL3 ,BGP5PDL4 ,BGP5PDL5 ,BGP5PDL7
BGP5PDL8 ,BGP5PDL9 ,BGP6CP ,BGP6DELP ,BGP6DHEP ,BGP6DP
BGP6DP11 ,BGP6DP12 ,BGP6DP1C
```

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 416

Routine	is Invoked by:
COMMA^%DTC	BGP6DP2, BGP6DP3, BGP6DP4, BGP6DP5, BGP6DP50, BGP6DP51 BGP6DP6, BGP6DP60, BGP6DP7, BGP6DPH, BGP6EL10, BGP6EL11 BGP6EL12, BGP6EL13, BGP6EL14, BGP6EL15, BGP6EL16, BGP6EL17 BGP6EL18, BGP6EL19, BGP6EL1A, BGP6EL1B, BGP6EL1C, BGP6EL1D BGP6EL1E, BGP6EL1F, BGP6EL1G BGP6EL1H, BGP6EL1I, BGP6EL1J, BGP6EL1K, BGP6EL1L, BGP6EL1M BGP6EL1N, BGP6ELL, BGP6ELL1, BGP6ELL2, BGP6ELL3, BGP6ELL4 BGP6ELL5, BGP6ELL6, BGP6ELL7, BGP6ELL8, BGP6ELL9, BGP6ELLA BGP6ELLB, BGP6ELLC, BGP6ELLD, BGP6ELLE, BGP6ELLF, BGP6ELLG BGP6EL LH, BGP6EL LI BGP6ELLJ, BGP6ELLK, BGP6ELLL, BGP6ELLM, BGP6ELLN, BGP6ELLO BGP6ELLP, BGP6ELLQ, BGP6ELP1, BGP6HEL, BGP6HEL1, BGP6HEP1 BGP6PDL, BGP6PDL0, BGP6PDL1, BGP6PDL2, BGP6PDL3, BGP6PDL4 BGP6PDL5, BGP6PDL7, BGP6PDL8, BGP6PDL9, BGP6PDLA, BGP7CP BGP7DELP, BGP7DHEP, BGP7DP BGP7DP11, BGP7DP12, BGP7DP19, BGP7DP1C, BGP7DP2, BGP7DP3 BGP7DP4, BGP7DP5, BGP7DP50, BGP7DP51, BGP7DP6, BGP7DP60 BGP7DP7, BGP7DPED, BGP7DPEP, BGP7DPH, BGP7EL10, BGP7EL11 BGP7EL12, BGP7EL13, BGP7EL14, BGP7EL15, BGP7EL16, BGP7EL17 BGP7EL18, BGP7EL19, BGP7EL1A BGP7EL1B, BGP7EL1C, BGP7EL1D, BGP7EL1E, BGP7EL1F, BGP7EL1G BGP7EL1H, BGP7EL1I, BGP7EL1J, BGP7EL1K, BGP7EL1L, BGP7EL1M BGP7EL1N, BGP7EL1O, BGP7EL1P, BGP7EL1Q, BGP7EL1R, BGP7ELL BGP7ELL1, BGP7ELL2, BGP7ELL3, BGP7ELL4, BGP7ELL5, BGP7ELL6 BGP7ELL7, BGP7ELL8 BGP7ELL9, BGP7ELLA, BGP7ELLB, BGP7ELLC, BGP7ELLD, BGP7ELLE BGP7ELLF, BGP7ELLG, BGP7EL LH, BGP7EL LI, BGP7ELLJ, BGP7ELLK BGP7ELLL, BGP7ELLM, BGP7ELLN, BGP7ELLO, BGP7ELLP, BGP7ELLQ BGP7ELLR, BGP7ELLS, BGP7ELLT, BGP7ELLU, BGP7ELP1, BGP7HEL BGP7HEL1, BGP7HEP1 BGP7PD19, BGP7PDL, BGP7PDL0, BGP7PDL1, BGP7PDL2, BGP7PDL3 BGP7PDL4, BGP7PDL5, BGP7PDL7, BGP7PDL8, BGP7PDL9, BGP7PDLA BGP7PDLN, BGPAP, BGPAP1, BGPAP12, BGPAP13, BGPAP14, BGPAP1B BGPAP2, BGPAP22, BGPAP22A, BGPAP22B, BGPAP22C, BGPAP23 BGPAP24, BGPAP3, BGPAP30 BGPAP30A, BGPAP30B, BGPAP30C, BGPAP4, BGPAP5, BGPAP6, BGPAP7 BGPAP8, BGPAPA, BGPAPB, BGPAPC, BGPAPD, BGPAPH, BGPDP, BGPDP1 BGPDP12, BGPDP13, BGPDP14, BGPDP1B, BGPDP2, BGPDP22, BGPDP22A BGPDP22B, BGPDP22C, BGPDP23, BGPDP24, BGPDP29, BGPDP3, BGPDP30 BGPDP30A BGPDP30B, BGPDP30C, BGPDP31, BGPDP32, BGPDP4, BGPDP5, BGPDP6 BGPDP7, BGPDP8, BGPDPA, BGPDPARP, BGPDPB, BGPDPD, BGPDPD BGPDPH BGP6GR, BGP6GRA, BGP6GU, BGP7GR, BGP7GR2, BGP7GRA, BGP7GU BGP7GU BGP3DA, BGP3DAR, BGP3DARA, BGP3DARH, BGP3DG, BGP3DL, BGP4DA BGP4DAP, BGP4DAR, BGP4DARA, BGP4DG, BGP4DHE, BGP4DHEA, BGP4DL BGP4DPP, BGP5CL, BGP5DAP, BGP5DAR, BGP5DEL, BGP5DELA, BGP5DGPA BGP5DGPU, BGP5DHE, BGP5DHEA, BGP5DL, BGP5DNG, BGP5DPP, BGP5NPL BGP6AHW
DEBUG^%Serenji	BGP6GR, BGP6GRA, BGP6GU, BGP7GR, BGP7GR2, BGP7GRA, BGP7GU BGP7GU
^%ZIS	BGP3DA, BGP3DAR, BGP3DARA, BGP3DARH, BGP3DG, BGP3DL, BGP4DA BGP4DAP, BGP4DAR, BGP4DARA, BGP4DG, BGP4DHE, BGP4DHEA, BGP4DL

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 417

Routine is Invoked by:

```

^%ZIS
BGP6CL, BGP6DAP, BGP6DAR, BGP6DCL, BGP6DEL, BGP6DELA, BGP6DGPA
BGP6DGPU, BGP6DHE, BGP6DHEA, BGP6DL, BGP6DLT, BGP6DMT, BGP6DNG
BGP6DPP, BGP6DSTM, BGP6NPL, BGP7CL, BGP7DAP, BGP7DAR, BGP7DCL
BGP7DEL, BGP7DELA, BGP7DGPA, BGP7DGPU, BGP7DHE, BGP7DHEA
BGP7DL, BGP7DLT
BGP7DMT, BGP7DNE, BGP7DNG, BGP7DPE, BGP7DPEA, BGP7DPP
BGP7DSTM, BGP7NPL, BGP7D, BGP7DAREA, BGP7DARP, BGP7DFTA
HOME^%ZIS
BGP3TXCH, BGP3ULF, BGP4LSTF, BGP4TXCH, BGP4TXHE, BGP4ULF
BGP50P1, BGP5CTXC, BGP5GUPL, BGP5LSTF, BGP5TXCH, BGP5TXCN
BGP5TXEL, BGP5TXHE, BGP5ULF, BGP6CTXC, BGP6DLT, BGP6DMT
BGP6GUPL, BGP6LSTF, BGP6ULF, BGP6XTCH, BGP6XTCN, BGP6XTEL
BGP6XTHE, BGP7CTXC, BGP7DAR
BGP7DLT, BGP7DMT, BGP7GUPL, BGP7LSTF, BGP7ULF, BGP7XTCH
BGP7XTCN, BGP7XTEL, BGP7XTHE, BGP7DARP, BGP7DFTA, BGP7DUP
^%ZISC
BGP3DA, BGP3DAR, BGP3DARA, BGP3DARH, BGP3DG, BGP3DL, BGP3DP
BGP3PARP, BGP3ULF, BGP4DA, BGP4DAP, BGP4DAR, BGP4DARA, BGP4DG
BGP4DHE, BGP4DHEA, BGP4DHEP, BGP4DL, BGP4DP, BGP4DPP, BGP4PARP
BGP4PHEP, BGP4ULF, BGP4UTL, BGP5CL, BGP5DAP, BGP5DAR, BGP5DEL
BGP5DELA
BGP5DELP, BGP5DGPA, BGP5DGPU, BGP5DHE, BGP5DHEA, BGP5DHEP
BGP5DL, BGP5DNG, BGP5DP, BGP5DPP, BGP5GUPL, BGP5NPL, BGP5PARP
BGP5PHEL, BGP5PHEP, BGP5ULF, BGP5UTL, BGP6AHW, BGP6CL, BGP6DAP
BGP6DAR, BGP6DCHA, BGP6DCL, BGP6DEL, BGP6DELA, BGP6DELP
BGP6DGPA, BGP6DGPU
BGP6DHE, BGP6DHEA, BGP6DHEP, BGP6DL, BGP6DLT, BGP6DMT, BGP6DNG
BGP6DP, BGP6DPP, BGP6DSTM, BGP6GUPL, BGP6NPL, BGP6PARP
BGP6PHEL, BGP6PHEP, BGP6ULF, BGP6UTL, BGP7CL, BGP7DAP, BGP7DAR
BGP7DCHA, BGP7DCL, BGP7DEL, BGP7DELA, BGP7DELP, BGP7DGPA
BGP7DGPU, BGP7DHE
BGP7DHEA, BGP7DHEP, BGP7DL, BGP7DLT, BGP7DMT, BGP7DNE, BGP7DNG
BGP7DP, BGP7DPE, BGP7DPEA, BGP7DPEP, BGP7DPP, BGP7DSTM
BGP7GUPL, BGP7NPL, BGP7PARP, BGP7PHEL, BGP7PHEP, BGP7PPEL
BGP7ULF, BGP7UTL, BGP7D, BGP7DAREA, BGP7DARP, BGP7DFTA, BGP7DUP
$$LIST^%ZISH
BGP4LSTF, BGP5GLST, BGP5GU, BGP5LSTF, BGP6GLST, BGP6GU
BGP6LSTF, BGP7GLST, BGP7GU, BGP7LSTF, BGP7GU
$$OPEN^%ZISH
BGP3ULF, BGP4ULF, BGP4UTL, BGP5GUPL, BGP5ULF, BGP5UTL, BGP6DNG
BGP6GUPL, BGP6ULF, BGP6UTL, BGP7DCHA, BGP7DNE, BGP7DNG
BGP7GUPL, BGP7ULF, BGP7UTL, BGP7DUP
$$STATUS^%ZISH
BGP3ULF, BGP4ULF, BGP5GUPL, BGP5ULF, BGP6GUPL, BGP6ULF
BGP7GUPL, BGP7ULF, BGP7DUP

```

^%ZTLOAD BGP3DA , BGP3DAR , BGP3DARA , BGP3DARH , BGP3DG , BGP3DL , BGP4DA
BGP4DAP , BGP4DAR , BGP4DARA , BGP4DG , BGP4DHE , BGP4DHEA , BGP4DL
BGP4DPP , BGP5CL , BGP5DAP , BGP5DAR , BGP5DEL , BGP5DELA , BGP5DGPA
BGP5DGPU , BGP5DHE , BGP5DHEA , BGP5DL , BGP5DNG , BGP5DPP
BGP5GAEL , BGP5GAGP
BGP5GAHE , BGP5GALL , BGP5GAPU , BGP5GCMS , BGP5GCOM , BGP5GELD
BGP5GGPU , BGP5GHED , BGP5GLST , BGP5GNPL , BGP5GNTL , BGP5GPP
BGP5NPL , BGP6AHW , BGP6CL , BGP6DAP , BGP6DAR , BGP6DCL , BGP6DEL
BGP6DELA , BGP6DGPA , BGP6DGPU , BGP6DHE , BGP6DHEA , BGP6DL
BGP6DLT , BGP6DMT , BGP6DNG

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 418

```

Routine          is Invoked by:
  ^%ZTLOAD      BGP6DPP, BGP6DSTM, BGP6GACW, BGP6GAEL, BGP6GAGP, BGP6GAHE
                BGP6GALL, BGP6GAPU, BGP6GCMP, BGP6GCMS, BGP6GCOM, BGP6GELD
                BGP6GGPU, BGP6GHED, BGP6GLST, BGP6GLTX, BGP6GMTX, BGP6GNPL

BGP6GNST, BGP6GNTL, BGP6GPP, BGP6NPL, BGP7CL, BGP7DAP, BGP7DAR
                BGP7DCL, BGP7DEL
                BGP7DELA, BGP7DGPA, BGP7DGPU, BGP7DHE, BGP7DHEA, BGP7DL

BGP7DLT, BGP7DMT, BGP7DNE, BGP7DNG, BGP7DPE, BGP7DPEA, BGP7DPP
                BGP7DSTM, BGP7GACW, BGP7GAEL, BGP7GAGP, BGP7GAHE, BGP7GALL
                BGP7GAPE, BGP7GAPU, BGP7GCMP, BGP7GCMS, BGP7GCOM, BGP7GELD
                BGP7GGPU, BGP7GHED
                BGP7GLST, BGP7GLTX, BGP7GMTX, BGP7GNPL, BGP7GNST, BGP7GNTL

BGP7GNXP, BGP7GPED, BGP7GPP, BGP7NPL, BGP7D, BGP7DAREA, BGP7DARP
                BGP7DFTA
  ^APCDCHKJ     BGP5C1, BGP6C1, BGP7C1
  $$START1^APCLDF
BGP3D1, BGP3D10, BGP3D2, BGP3D3, BGP3D4, BGP3D5, BGP3D6, BGP3D7
                BGP3D8, BGP3D9, BGP4D1, BGP4D10, BGP4D2, BGP4D21, BGP4D3
                BGP4D31, BGP4D4, BGP4D41, BGP4D5, BGP4D6, BGP4D61, BGP4D7
                BGP4D71, BGP4D8, BGP4D9, BGP5C1, BGP5C11, BGP5C13, BGP5C2
                BGP5C3, BGP5D1, BGP5D10
                BGP5D2, BGP5D21, BGP5D22, BGP5D3, BGP5D31, BGP5D33, BGP5D34
                BGP5D35, BGP5D4, BGP5D41, BGP5D5, BGP5D6, BGP5D61, BGP5D7
                BGP5D71, BGP5D72, BGP5D8, BGP5D9, BGP5DHE1, BGP5EL2, BGP5EL3
                BGP5EL4, BGP6C1, BGP6C11, BGP6C12, BGP6C13, BGP6C2, BGP6C3
                BGP6C31, BGP6D1, BGP6D10
                BGP6D2, BGP6D21, BGP6D22, BGP6D23, BGP6D3, BGP6D31, BGP6D33
                BGP6D34, BGP6D35, BGP6D4, BGP6D41, BGP6D5, BGP6D51, BGP6D52
                BGP6D6, BGP6D61, BGP6D7, BGP6D71, BGP6D72, BGP6D8, BGP6D81
                BGP6D82, BGP6D83, BGP6D84, BGP6D85, BGP6D861, BGP6D862
                BGP6D863, BGP6D9, BGP6DCHW

BGP6DHE1, BGP6EL2, BGP6EL3, BGP6EL4, BGP7C1, BGP7C11, BGP7C12
                BGP7C13, BGP7C2, BGP7C3, BGP7C31, BGP7D1, BGP7D10, BGP7D2

BGP7D21, BGP7D211, BGP7D21A, BGP7D22, BGP7D23, BGP7D3, BGP7D31
                BGP7D33, BGP7D34, BGP7D35, BGP7D4, BGP7D41, BGP7D5, BGP7D51
                BGP7D52, BGP7D6
                BGP7D61, BGP7D7, BGP7D71, BGP7D72, BGP7D721, BGP7D722

BGP7D731, BGP7D732, BGP7D74, BGP7D75, BGP7D8, BGP7D81, BGP7D82
                BGP7D83, BGP7D84, BGP7D85, BGP7D862, BGP7D9, BGP7DCHW
                BGP7DCLP, BGP7DHE1, BGP7DHE2, BGP7DPE1, BGP7DPE2, BGP7EL2
                BGP7EL3, BGP7EL4, BGP7D1, BGP7D12
                BGP7D13, BGP7D14, BGP7D1B, BGP7D2, BGP7D22, BGP7D24, BGP7D29, BGP7D3
                BGP7D31, BGP7D4, BGP7D5, BGP7D6, BGP7D7, BGP7DA, BGP7DB, BGP7DC, BGP7DD

  ^APCLTAX0     |opt
  $$CLINIC^APCLV
BGP3D1, BGP3D2, BGP3D3, BGP3D4, BGP3D6, BGP3D9, BGP4D1, BGP4D10
                BGP4D2, BGP4D21, BGP4D41, BGP4D7, BGP4D71, BGP5D1, BGP5D10

```

BGP5D2 ,BGP5D21 ,BGP5D22 ,BGP5D41 ,BGP5D7 ,BGP5D71 ,BGP5DHE1
BGP5EL2 ,BGP5EL4 ,BGP6C12 ,BGP6C2 ,BGP6C3 ,BGP6D1 ,BGP6D10
BGP6D2 ,BGP6D21
BGP6D22 ,BGP6D23 ,BGP6D41 ,BGP6D51 ,BGP6D7 ,BGP6D71 ,BGP6D83
BGP6D84 ,BGP6DHE1 ,BGP6EL2 ,BGP6EL4 ,BGP7C12 ,BGP7C2 ,BGP7C3
BGP7D1 ,BGP7D10 ,BGP7D2 ,BGP7D21 ,BGP7D22 ,BGP7D23 ,BGP7D41

BGP7D51 ,BGP7D7 ,BGP7D71 ,BGP7D83 ,BGP7D84 ,BGP7DHE1 ,BGP7EL2
BGP7EL4 ,BGPD14

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 419

```

Routine          is Invoked by:
  $$CLINIC^APCLV    BGPD2, BGPD22, BGPD3, BGPD8, BGPDD
  $$DSCHDATE^APCLV
BGP5D22, BGP5D72, BGP5EL3, BGP6D22, BGP6D51, BGP6EL3, BGP7D22
  BGP7D51, BGP7EL3
  $$PRIMPOV^APCLV   BGP3D2, BGP3D3, BGP3D8, BGP4D21, BGP4D31, BGP4D71, BGP5C1
  BGP5C2, BGP5C3, BGP5D21, BGP5D22, BGP5D31, BGP5D71, BGP5DHE1
  BGP6C1, BGP6C2, BGP6C3, BGP6D21, BGP6D22, BGP6D31, BGP6D51
  BGP6D71, BGP6DHE1, BGP7C1, BGP7C2, BGP7C3, BGP7D21, BGP7D211
  BGP7D22, BGP7D31
  BGP7D51, BGP7D71, BGP7DHE1, BGPD14, BGPD8, BGPDD
  $$PRIMPROV^APCLV  BGP3D1, BGP3D2, BGP4D1, BGP4D10, BGP4D21, BGP4D71, BGP5D1

BGP5D10, BGP5D21, BGP5D71, BGP5DHE1, BGP6D1, BGP6D10, BGP6D21

BGP6D22, BGP6D51, BGP6D71, BGP6DHE1, BGP7D1, BGP7D10, BGP7D21
  BGP7D22, BGP7D51, BGP7D71, BGP7DCLP, BGP7DHE1, BGPD14, BGPD2
  BGPDD
  $$ICD^ATXCHK      BGP3D1, BGP3D6, BGP3D7, BGP3D8, BGP3D9, BGP4D1, BGP4D10
  BGP4D31, BGP4D5, BGP4D7, BGP4D8, BGP4D9, BGP5C1, BGP5C11
  BGP5C13, BGP5C2, BGP5C3, BGP5D1, BGP5D10, BGP5D22, BGP5D31
  BGP5D4, BGP5D5, BGP5D61, BGP5D7, BGP5D72, BGP5D8, BGP5D9
  BGP5DHE1, BGP5EL3, BGP5UTL1
  BGP6C1, BGP6C11, BGP6C13, BGP6C2, BGP6C3, BGP6C31, BGP6D1
  BGP6D10, BGP6D22, BGP6D31, BGP6D4, BGP6D5, BGP6D51, BGP6D61
  BGP6D7, BGP6D72, BGP6D8, BGP6D81, BGP6D82, BGP6D83, BGP6D84

BGP6D862, BGP6D9, BGP6DHE1, BGP6DU, BGP6EL3, BGP6UTL1, BGP7C1
  BGP7C11, BGP7C13
  BGP7C2, BGP7C3, BGP7C31, BGP7D1, BGP7D10, BGP7D22, BGP7D31
  BGP7D4, BGP7D5, BGP7D51, BGP7D61, BGP7D7, BGP7D72, BGP7D721
  BGP7D722, BGP7D731, BGP7D732, BGP7D74, BGP7D75, BGP7D8

BGP7D81, BGP7D82, BGP7D83, BGP7D84, BGP7D862, BGP7D9, BGP7DHE1
  BGP7DU, BGP7EL3, BGP7UTL1
  BGPD2, BGPDU

  BULL^ATXSTX2
BGP6AXA, BGP6AXB, BGP6AXC, BGP6AXD, BGP6AXE, BGP6AXF, BGP6AXG

BGP6AXH, BGP6AXI, BGP6AXJ, BGP6AXK, BGP6AXL, BGP6AXM, BGP6AXN

BGP6AXO, BGP6AXP, BGP6AXQ, BGP6AXR, BGP6AXS, BGP6AXT, BGP6AXU

BGP6AXV, BGP6AXW, BGP6AXX, BGP6AXY, BGP6AXZ, BGP6BXA, BGP6BXB
  BGP6BXC

BGP6BXD, BGP6BXE, BGP6BXF, BGP6BXG, BGP6BXH, BGP6BXI, BGP6BXJ

BGP6BXK, BGP6BXL, BGP6BXM, BGP6BXN, BGP6BXO, BGP6BXP, BGP6BXQ

BGP6BXR, BGP6BXS, BGP6BXT, BGP6BXU, BGP6BXV, BGP6BXW, BGP6BXX

BGP6BXY, BGP6BXZ, BGP6EXA, BGP6EXB, BGP6EXC, BGP6EXD, BGP6EXE

```

BGP6EXF

BGP6EXG , BGP6EXH , BGP6EXI , BGP6EXJ , BGP6EXK , BGP6EXL , BGP6EXM

BGP6EXN , BGP6EXO , BGP6EXP , BGP6EXQ , BGP6EXR , BGP6EXS , BGP6EXT

BGP6EXU , BGP6EXV , BGP6EXW , BGP6EXX , BGP6EXY , BGP6EXZ , BGP6FXA

BGP6FXB , BGP6FXC , BGP6FXD , BGP6FXE , BGP6FXF , BGP6FXG , BGP6FXH

BGP6FXI

BGP6FXJ , BGP6FXK , BGP6FXL , BGP6FXM , BGP6FXN , BGP6FXO , BGP6FXP

BGP6FXQ , BGP6FXR , BGP6FXS , BGP6FXT , BGP6FXU , BGP6FXV , BGP6FXW

BGP6FXX , BGP6FXY , BGP6FXZ , BGP6GXA , BGP6GXB , BGP6GXC , BGP6GXD

BGP6GXE , BGP6GXG , BGP6GXH , BGP6GXI , BGP6GXJ , BGP6GXK , BGP6GXL

BGP6GXM

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
420

Routine is Invoked by:
BULL^ATXSTX2
BGP6GXN, BGP6GXO, BGP6GXP, BGP6GXQ, BGP6GXR, BGP6GXS, BGP6GXT
BGP6GXU, BGP6GXX, BGP6GXY, BGP6HXA, BGP6HXB, BGP6HXC, BGP6HXD
BGP6HXE, BGP6HXF, BGP6HXG, BGP6HXH, BGP6HXI, BGP6HXJ, BGP6HXK
BGP6HXL, BGP6HXM, BGP6HYN, BGP6HXO, BGP6HXP, BGP6HXQ, BGP6HXR
BGP6HXS
BGP6HXT, BGP6HXU, BGP6HXV, BGP6HXW, BGP6HXX, BGP6HXY, BGP6HXZ
BGP6JXA, BGP6JXB, BGP6JXC, BGP6JXD, BGP6JXE, BGP6JXF, BGP6JXG
BGP6JXH, BGP6JXI, BGP6JXJ, BGP6JXK, BGP6JXL, BGP6JXM, BGP6JXN
BGP6JXO, BGP6JXP, BGP6JXQ, BGP6JXR, BGP6JXS, BGP6JXT, BGP6JXU
BGP6JXV
BGP6JXW, BGP7GXA, BGP7GXB, BGP7GXC, BGP7GXD, BGP7GXE, BGP7GXF
BGP7GXG, BGP7GXH, BGP7GXI, BGP7GXJ, BGP7GXK, BGP7GXL, BGP7GXM
BGP7GXN, BGP7GXO, BGP7GXP, BGP7GXQ, BGP7HXA, BGP7HXB, BGP7HXC
BGP7HXD, BGP7HXE, BGP7HXF, BGP7HXG, BGP7HXH, BGP7HXI, BGP7HXJ
BGP7HXK
BGP7HXL, BGP7HXM, BGP7HYN, BGP7HXO, BGP7HXP, BGP7HXQ, BGP7HXR
BGP7HXS, BGP7HXT, BGP7HXU, BGP7HXV, BGP7HXW, BGP7HXX, BGP7HXY
BGP7HXZ, BGP7IXA, BGP7IXB, BGP7IXC, BGP7IXD, BGP7IXE, BGP7IXF
BGP7IXG, BGP7IXH, BGP7IXI, BGP7IXJ, BGP7IXK, BGP7IXL, BGP7IXM
BGP7IXN
BGP7IXO, BGP7IXP, BGP7IXQ, BGP7IXR, BGP7IXS, BGP7IXT, BGP7IXU
BGP7IXV, BGP7IXW, BGP7IXX, BGP7IXY, BGP7IXZ, BGP7JXA, BGP7JXB
BGP7JXC, BGP7JXD, BGP7JXE, BGP7JXF, BGP7JXG, BGP7JXH, BGP7JXI
BGP7JXJ, BGP7JXK, BGP7JXL, BGP7JXM, BGP7JXN, BGP7JXO, BGP7JXP
BGP7JXQ
BGP7JXR, BGP7JXS, BGP7JXT, BGP7JXU, BGP7JXV, BGP7JXW, BGP7JXX
BGP7JXY, BGP7JXZ, BGP7KXA, BGP7KXB, BGP7KXC, BGP7KXD, BGP7KXE
BGP7KXF, BGP7KXG, BGP7KXH, BGP7KXI, BGP7KXJ, BGP7KXK, BGP7KXL

BGP7KXM, BGP7KXN, BGP7KXO, BGP7KXP, BGP7KXQ, BGP7KXR, BGP7KXS
BGP7KXT

BGP7KXU, BGP7KXV, BGP7KXW, BGP7KXX, BGP7KXY, BGP7KXZ, BGP7LXA

BGP7LXB, BGP7LXC, BGP7LXD, BGP7LXE, BGP7LXF, BGP7LXG, BGP7LXH

BGP7LXI, BGP7LXJ, BGP7LXK, BGP7LXL, BGP7LXM, BGP7LXN, BGP7LXO

BGP7LXP, BGP7LXR, BGP7LXS, BGP7LXT, BGP7LXU, BGP7LXV, BGP7LXW
BGP7LXX

BGP7LXY, BGP7LXZ, BGP7MXA, BGP7MXB, BGPLTXA, BGPLTXB, BGPLTXC
BGPLTXD, BGPLTXE, BGPLTXF, BGPLXA, BGPLXB, BGPLXC, BGPLXD

BGPLXE, BGPLXF, BGPLXG, BGPLXH, BGPLXI, BGPLXJ, BGPLXK, BGPLXL

BGPLXM, BGPMXA, BGPMXB, BGPTXA, BGPTXB, BGPTXC, BGPTXD, BGPTXE
BGPTXF, BGPTXG

BGPTXH, BGPTXI, BGPTXJ, BGPTXK, BGPTXL, BGPTXM, BGPTXN, BGPTXO

BGPTXP, BGPTXQ, BGPTXR, BGPTXS, BGPTXT, BGPTXU, BGPTXV, BGPTXW
BGPTXX, BGPTXY, BGPTXZ, BGPUTA, BGPUTB, BGPUTC, BGPUTD

BGPUTE, BGPUTF, BPUTXG, BPUTXH, BPUTXI, BPUTXJ, BPUTXK
BPUTXL, BPUTXM

BGPUN, BPUTO, BPUTP, BPUTQ, BPUTR, BPUTS, BPUTT
BPUTU, BPUTV, BPUTW, BPUTX, BPUTY, BPUTZ, BPUTA

BPUTB, BPUTC, BPUTD, BPUTE, BPUXF, BPUTG, BPUXH, BPUTI

BPUTJ, BPUTK, BPUTL, BPUTM, BPUTN, BPUXO, BPUTP, BPUTQ
BPUXR, BPUTS

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 421

Routine is Invoked by:

BULL^ATXSTX2
BGPUXT, BGPUXU, BGPUXV, BGPUXW, BGPUXX, BGPUXY, BGPUXZ, BGPVTXA
BGPVTXB, BGPVTXC, BGPVTXD, BGPVTXE, BGPVTXF, BGPVTXG, BGPVTXH
BGPVTXI, BGPVTXJ, BGPVTXK, BGPVTXL, BGPVTXM, BGPVTXN, BGPVTXO
BGPVTXP, BGPVTXQ, BGPVTXR, BGPVTXS, BGPVTXT, BGPVTXU, BGPVTXV
BGPVTXW
BGPVTXX, BGPVTTY, BGPVTXZ, BGPVXA, BGPVXB, BGPVXC, BGPVXD
BGPVXE, BGPVXF, BGPVXG, BGPVXH, BGPVXI, BGPVXJ, BGPVXK, BGPVXL
BGPVXM, BGPVXN, BGPVXO, BGPVXP, BGPVXQ, BGPVXR, BGPVXS, BGPVXT
BGPVXU, BGPVXV, BGPVXW, BGPVXX, BGPVXY, BGPVXZ, BGPWTXA
BGPWTXB, BGPWTXC, BGPWTXD
BGPWTXE, BGPWTXF, BGPWTXG, BGPWTXH, BGPWTXI, BGPWTXJ, BGPWTXK
BGPWTXL, BGPWTXM, BGPWTXN, BGPWTXO, BGPWTXP, BGPWTXQ, BGPWTXR
BGPWTXS, BGPWTXT, BGPWTXU
KILL^ATXSTX2
BGP6AXA, BGP6AXB, BGP6AXC, BGP6AXD, BGP6AXE, BGP6AXF, BGP6AXG
BGP6AXH, BGP6AXI, BGP6AXJ, BGP6AXK, BGP6AXL, BGP6AXM, BGP6AXN
BGP6AXO, BGP6AXP, BGP6AXQ, BGP6AXR, BGP6AXS, BGP6AXT, BGP6AXU
BGP6AXV, BGP6AXW, BGP6AXX, BGP6AXY, BGP6AXZ, BGP6BXA, BGP6BXB
BGP6BXC
BGP6BXD, BGP6BXE, BGP6BXF, BGP6BXG, BGP6BXH, BGP6BXI, BGP6BXJ
BGP6BXK, BGP6BXL, BGP6BXM, BGP6BXN, BGP6BXO, BGP6BXP, BGP6BXQ
BGP6BXR, BGP6BXS, BGP6BXT, BGP6BXU, BGP6BXV, BGP6BXW, BGP6BXX
BGP6BXY, BGP6BXZ, BGP6EXA, BGP6EXB, BGP6EXC, BGP6EXD, BGP6EXE
BGP6EXF
BGP6EXG, BGP6EXH, BGP6EXI, BGP6EXJ, BGP6EXK, BGP6EXL, BGP6EXM
BGP6EXN, BGP6EXO, BGP6EXP, BGP6EXQ, BGP6EXR, BGP6EXS, BGP6EXT
BGP6EXU, BGP6EXV, BGP6EXW, BGP6EXX, BGP6EXY, BGP6EXZ, BGP6FXA
BGP6FXB, BGP6FXC, BGP6FXD, BGP6FXE, BGP6FXF, BGP6FXG, BGP6FXH
BGP6FXI
BGP6FXJ, BGP6FXK, BGP6FXL, BGP6FXM, BGP6FXN, BGP6FXO, BGP6FXP
BGP6FXQ, BGP6FXR, BGP6FXS, BGP6FXT, BGP6FXU, BGP6FXV, BGP6FXW

BGP6FXX, BGP6FXY, BGP6FXZ, BGP6GXA, BGP6GXB, BGP6GXC, BGP6GXD
BGP6GXE, BGP6GXG, BGP6GXH, BGP6GXI, BGP6GXJ, BGP6GXK, BGP6GXL
BGP6GXM
BGP6GXN, BGP6GXO, BGP6GXP, BGP6GXQ, BGP6GXR, BGP6GXS, BGP6GXT
BGP6GXU, BGP6GXX, BGP6GXY, BGP6HXA, BGP6HXB, BGP6HXC, BGP6HXD
BGP6HXE, BGP6HXF, BGP6HXG, BGP6HXH, BGP6HXI, BGP6HXJ, BGP6HXK
BGP6HXL, BGP6HXM, BGP6HYN, BGP6HXO, BGP6HXP, BGP6HXQ, BGP6HXR
BGP6HXS
BGP6HXT, BGP6HXU, BGP6HXV, BGP6HXW, BGP6HXX, BGP6HXY, BGP6HXZ
BGP6JXA, BGP6JXB, BGP6JXC, BGP6JXD, BGP6JXE, BGP6JXF, BGP6JXG
BGP6JXH, BGP6JXI, BGP6JXJ, BGP6JXK, BGP6JXL, BGP6JXM, BGP6JXN
BGP6JXO, BGP6JXP, BGP6JXQ, BGP6JXR, BGP6JXS, BGP6JXT, BGP6JXU
BGP6JXV
BGP6JXW, BGP7GXA, BGP7GXB, BGP7GXC, BGP7GXD, BGP7GXE, BGP7GXF
BGP7GXG, BGP7GXH, BGP7GXI, BGP7GXJ, BGP7GXK, BGP7GXL, BGP7GXM
BGP7GXN, BGP7GXO, BGP7GXP, BGP7GXQ, BGP7HXA, BGP7HXB, BGP7HXC
BGP7HXD, BGP7HXE, BGP7HXF, BGP7HXG, BGP7HXH, BGP7HXI, BGP7HXJ
BGP7HXK
BGP7HXL, BGP7HXM, BGP7HYN, BGP7HXO, BGP7HXP, BGP7HXQ, BGP7HXR
BGP7HXS, BGP7HXT, BGP7HXU, BGP7HXV, BGP7HXW, BGP7HXX, BGP7HXY
BGP7HXZ, BGP7IXA, BGP7IXB, BGP7IXC, BGP7IXD, BGP7IXE, BGP7IXF
BGP7IXG, BGP7IXH, BGP7IXI, BGP7IXJ, BGP7IXK, BGP7IXL, BGP7IXM
BGP7IXN

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
422

Routine is Invoked by:

KILL^ATXSTX2
BGP7IXO, BGP7IXP, BGP7IXQ, BGP7IXR, BGP7IXS, BGP7IXT, BGP7IXU
BGP7IXV, BGP7IXW, BGP7IXX, BGP7IXY, BGP7IXZ, BGP7JXA, BGP7JXB
BGP7JXC, BGP7JXD, BGP7JXE, BGP7JXF, BGP7JXG, BGP7JXH, BGP7JXI
BGP7JXJ, BGP7JXK, BGP7JXL, BGP7JXM, BGP7JXN, BGP7JXO, BGP7JXP
BGP7JXQ
BGP7JXR, BGP7JXS, BGP7JXT, BGP7JXU, BGP7JXV, BGP7JXW, BGP7JXX
BGP7JXY, BGP7JXZ, BGP7KXA, BGP7KXB, BGP7KXC, BGP7KXD, BGP7KXE
BGP7KXF, BGP7KXG, BGP7KXH, BGP7KXI, BGP7KXJ, BGP7KXK, BGP7KXL
BGP7KXM, BGP7KXN, BGP7KXO, BGP7KXP, BGP7KXQ, BGP7KXR, BGP7KXS
BGP7KXT
BGP7KXU, BGP7KXV, BGP7KXW, BGP7KXX, BGP7KXY, BGP7KXZ, BGP7LXA
BGP7LXB, BGP7LXC, BGP7LXD, BGP7LXE, BGP7LXF, BGP7LXG, BGP7LXH
BGP7LXI, BGP7LXJ, BGP7LXK, BGP7LXL, BGP7LXM, BGP7LXN, BGP7LXO
BGP7LXP, BGP7LXR, BGP7LXS, BGP7LXT, BGP7LXU, BGP7LXV, BGP7LXW
BGP7LXX
BGP7LXY, BGP7LXZ, BGP7MXA, BGP7MXB, BGPLTXA, BGPLTXB, BGPLTXC
BGPLTXD, BGPLTXE, BGPLTXF, BGPLXA, BGPLXB, BGPLXC, BGPLXD
BGPLXE, BGPLXF, BGPLXG, BGPLXH, BGPLXI, BGPLXJ, BGPLXK, BGPLXL
BGPLXM, BGPMXA, BGPMXB, BGPTXA, BGPTXB, BGPTXC, BGPTXD, BGPTXE
BGPTXF, BGPTXG
BGPTXH, BGPTXI, BGPTXJ, BGPTXK, BGPTXL, BGPTXM, BGPTXN, BGPTXO
BGPTXP, BGPTXQ, BGPTXR, BGPTXS, BGPTXT, BGPTXU, BGPTXV, BGPTXW
BGPTXX, BGPTXY, BGPTXZ, BGPUTXA, BGPUTXB, BGPUTXC, BGPUTXD
BGPUTXE, BGPUTXF, BGPUTXG, BGPUTXH, BGPUTXI, BGPUTXJ, BGPUTXK
BGPUTXL, BGPUTXM
BGPUTXN, BGPUTXO, BGPUTXP, BGPUTXQ, BGPUTXR, BGPUTXS, BGPUTXT
BGPUTXU, BGPUTXV, BGPUTXW, BGPUTXX, BGPUTXY, BGPUTXZ, BGPUXA
BGPUXB, BGPUXC, BGPUXD, BGPUXE, BGPUXF, BGPUXG, BGPUXH, BGPUXI
BGPUXJ, BGPUXK, BGPUXL, BGPUXM, BGPUXN, BGPUXO, BGPUXP, BGPUXQ
BGPUXR, BGPUXS

BGPUXT , BGPUXU , BGPUXV , BGPUXW , BGPUXX , BGPUXY , BGPUXZ , BGPVTXA
BGPVTXB , BGPVTXC , BGPVTXD , BGPVTXE , BGPVTXF , BGPVTXG , BGPVTXH
BGPVTXI , BGPVTXJ , BGPVTXK , BGPVTXL , BGPVTXM , BGPVTXN , BGPVTXO
BGPVTXP , BGPVTXQ , BGPVTXR , BGPVTXS , BGPVTXT , BGPVTXU , BGPVTXV
BGPVTXW
BGPVTXX , BGPVTXY , BGPVTXZ , BGPVXA , BGPVXB , BGPVXC , BGPVXD
BGPVXE , BGPVXF , BGPVXG , BGPVXH , BGPVXI , BGPVXJ , BGPVXK , BGPVXL
BGPVXM , BGPVXN , BGPVXO , BGPVXP , BGPVXQ , BGPVXR , BGPVXS , BGPVXT
BGPVXU , BGPVXV , BGPVXW , BGPVXX , BGPVXY , BGPVXZ , BGPWTXA
BGPWTXB , BGPWTXC , BGPWTXD
BGPWTXE , BGPWTXF , BGPWTXG , BGPWTXH , BGPWTXI , BGPWTXJ , BGPWTXK
BGPWTXL , BGPWTXM , BGPWTXN , BGPWTXO , BGPWTXP , BGPWTXQ , BGPWTXR
BGPWTXS , BGPWTXT , BGPWTXU
TAX^ATXSTX2
BGP6AXA , BGP6AXB , BGP6AXC , BGP6AXD , BGP6AXE , BGP6AXF , BGP6AXG
BGP6AXH , BGP6AXI , BGP6AXJ , BGP6AXK , BGP6AXL , BGP6AXM , BGP6AXN
BGP6AXO , BGP6AXP , BGP6AXQ , BGP6AXR , BGP6AXS , BGP6AXT , BGP6AXU
BGP6AXV , BGP6AXW , BGP6AXX , BGP6AXY , BGP6AXZ , BGP6BXA , BGP6BXB
BGP6BXC
BGP6BXD , BGP6BXE , BGP6BXF , BGP6BXG , BGP6BXH , BGP6BXI , BGP6BXJ
BGP6BXK , BGP6BXL , BGP6BXM , BGP6BXN , BGP6BXO , BGP6BXP , BGP6BXQ
BGP6BXR , BGP6BXS , BGP6BXT , BGP6BXU , BGP6BXV , BGP6BXW , BGP6BXX
BGP6BXY , BGP6BXZ , BGP6EXA , BGP6EXB , BGP6EXC , BGP6EXD , BGP6EXE
BGP6EXF

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
423

Routine is Invoked by:
TAX^ATXSTX2
BGP6EXG, BGP6EXH, BGP6EXI, BGP6EXJ, BGP6EXK, BGP6EXL, BGP6EXM
BGP6EXN, BGP6EXO, BGP6EXP, BGP6EXQ, BGP6EXR, BGP6EXS, BGP6EXT
BGP6EXU, BGP6EXV, BGP6EXW, BGP6EXX, BGP6EXY, BGP6EXZ, BGP6FXA
BGP6FXB, BGP6FXC, BGP6FXD, BGP6FXE, BGP6FXF, BGP6FXG, BGP6FXH
BGP6FXI
BGP6FXJ, BGP6FXK, BGP6FXL, BGP6FXM, BGP6FXN, BGP6FXO, BGP6FXP
BGP6FXQ, BGP6FXR, BGP6FXS, BGP6FXT, BGP6FXU, BGP6FXV, BGP6FXW
BGP6FXX, BGP6FXY, BGP6FXZ, BGP6GXA, BGP6GXB, BGP6GXC, BGP6GXD
BGP6GXE, BGP6GXG, BGP6GXH, BGP6GXI, BGP6GXJ, BGP6GXK, BGP6GXL
BGP6GXM
BGP6GXN, BGP6GXO, BGP6GXP, BGP6GXQ, BGP6GXR, BGP6GXS, BGP6GXT
BGP6GXU, BGP6GXX, BGP6GXY, BGP6HXA, BGP6HXB, BGP6HXC, BGP6HXD
BGP6HXE, BGP6HXF, BGP6HXG, BGP6HXH, BGP6HXI, BGP6HXJ, BGP6HXK
BGP6HXL, BGP6HXM, BGP6HYN, BGP6HXO, BGP6HXP, BGP6HXQ, BGP6HXR
BGP6HXS
BGP6HXT, BGP6HXU, BGP6HXV, BGP6HXW, BGP6HXX, BGP6HXY, BGP6HXZ
BGP6JXA, BGP6JXB, BGP6JXC, BGP6JXD, BGP6JXE, BGP6JXF, BGP6JXG
BGP6JXH, BGP6JXI, BGP6JXJ, BGP6JXK, BGP6JXL, BGP6JXM, BGP6JXN
BGP6JXO, BGP6JXP, BGP6JXQ, BGP6JXR, BGP6JXS, BGP6JXT, BGP6JXU
BGP6JXV
BGP6JXW, BGP7GXA, BGP7GXB, BGP7GXC, BGP7GXD, BGP7GXE, BGP7GXF
BGP7GXG, BGP7GXH, BGP7GXI, BGP7GXJ, BGP7GXK, BGP7GXL, BGP7GXM
BGP7GXN, BGP7GXO, BGP7GXP, BGP7GXQ, BGP7HXA, BGP7HXB, BGP7HXC
BGP7HXD, BGP7HXE, BGP7HXF, BGP7HXG, BGP7HXH, BGP7HXI, BGP7HXJ
BGP7HXK
BGP7HXL, BGP7HXM, BGP7HYN, BGP7HXO, BGP7HXP, BGP7HXQ, BGP7HXR
BGP7HXS, BGP7HXT, BGP7HXU, BGP7HXV, BGP7HXW, BGP7HXX, BGP7HXY
BGP7HXZ, BGP7IXA, BGP7IXB, BGP7IXC, BGP7IXD, BGP7IXE, BGP7IXF

BGP7IXG, BGP7IXH, BGP7IXI, BGP7IXJ, BGP7IXK, BGP7IXL, BGP7IXM
BGP7IXN

BGP7IXO, BGP7IXP, BGP7IXQ, BGP7IXR, BGP7IXS, BGP7IXT, BGP7IXU

BGP7IXV, BGP7IXW, BGP7IXX, BGP7IXY, BGP7IXZ, BGP7JXA, BGP7JXB

BGP7JXC, BGP7JXD, BGP7JXE, BGP7JXF, BGP7JXG, BGP7JXH, BGP7JXI

BGP7JXJ, BGP7JXK, BGP7JXL, BGP7JXM, BGP7JXN, BGP7JXO, BGP7JXP
BGP7JXQ

BGP7JXR, BGP7JXS, BGP7JXT, BGP7JXU, BGP7JXV, BGP7JXW, BGP7JXX

BGP7JXY, BGP7JXZ, BGP7KXA, BGP7KXB, BGP7KXC, BGP7KXD, BGP7KXE

BGP7KXF, BGP7KXG, BGP7KXH, BGP7KXI, BGP7KXJ, BGP7KXK, BGP7KXL

BGP7KXM, BGP7KXN, BGP7KXO, BGP7KXP, BGP7KXQ, BGP7KXR, BGP7KXS
BGP7KXT

BGP7KXU, BGP7KXV, BGP7KXW, BGP7KXX, BGP7KXY, BGP7KXZ, BGP7LXA

BGP7LXB, BGP7LXC, BGP7LXD, BGP7LXE, BGP7LXF, BGP7LXG, BGP7LXH

BGP7LXI, BGP7LXJ, BGP7LXK, BGP7LXL, BGP7LXM, BGP7LXN, BGP7LXO

BGP7LXP, BGP7LXR, BGP7LXS, BGP7LXT, BGP7LXU, BGP7LXV, BGP7LXW
BGP7LXX

BGP7LXY, BGP7LXZ, BGP7MXA, BGP7MXB, BGPLTXA, BGPLTXB, BGPLTXC
BGPLTXD, BGPLTXE, BGPLTXF, BGPLXA, BGPLXB, BGPLXC, BGPLXD

BGPLXE, BGPLXF, BGPLXG, BGPLXH, BGPLXI, BGPLXJ, BGPLXK, BGPLXL

BGPLXM, BGPMXA, BGPMXB, BGPTXA, BGPTXB, BGPTXC, BGPTXD, BGPTXE
BGPTXF, BGPTXG

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 424

Routine is Invoked by:
 TAX^ATXSTX2
 BGPTXH, BGPTXI, BGPTXJ, BGPTXK, BGPTXL, BGPTXM, BGPTXN, BGPTXO
 BGPTXP, BGPTXQ, BGPTXR, BGPTXS, BGPTXT, BGPTXU, BGPTXV, BGPTXW
 BGPTXX, BGPTXY, BGPTXZ, BGPputXA, BGPputXB, BGPputXC, BGPputXD
 BGPputXE, BGPputXF, BGPputXG, BGPputXH, BGPputXI, BGPputXJ, BGPputXK
 BGPputXL, BGPputXM
 BGPputXN, BGPputXO, BGPputXP, BGPputXQ, BGPputXR, BGPputXS, BGPputXT
 BGPputXU, BGPputXV, BGPputXW, BGPputXX, BGPputXY, BGPputXZ, BGPputXA
 BGPputXB, BGPputXC, BGPputXD, BGPputXE, BGPputXF, BGPputXG, BGPputXH, BGPputXI
 BGPputXJ, BGPputXK, BGPputXL, BGPputXM, BGPputXN, BGPputXO, BGPputXP, BGPputXQ
 BGPputXR, BGPputXS
 BGPputXT, BGPputXU, BGPputXV, BGPputXW, BGPputXX, BGPputXY, BGPputXZ, BGPputXA
 BGPputXB, BGPputXC, BGPputXD, BGPputXE, BGPputXF, BGPputXG, BGPputXH
 BGPputXI, BGPputXJ, BGPputXK, BGPputXL, BGPputXM, BGPputXN, BGPputXO
 BGPputXP, BGPputXQ, BGPputXR, BGPputXS, BGPputXT, BGPputXU, BGPputXV
 BGPputXW
 BGPputXX, BGPputXY, BGPputXZ, BGPputXA, BGPputXB, BGPputXC, BGPputXD
 BGPputXE, BGPputXF, BGPputXG, BGPputXH, BGPputXI, BGPputXJ, BGPputXK, BGPputXL
 BGPputXM, BGPputXN, BGPputXO, BGPputXP, BGPputXQ, BGPputXR, BGPputXS, BGPputXT
 BGPputXU
 BGPputXX, BGPputXY, BGPputXZ, BGPputXA, BGPputXB, BGPputXC, BGPputXD
 BGPputXE, BGPputXF, BGPputXG, BGPputXH, BGPputXI, BGPputXJ, BGPputXK, BGPputXL
 BGPputXM, BGPputXN, BGPputXO, BGPputXP, BGPputXQ, BGPputXR, BGPputXS, BGPputXT
 BGPputXU
 BGPputXX, BGPputXY, BGPputXZ, BGPputXA, BGPputXB, BGPputXC, BGPputXD
 BGPputXE, BGPputXF, BGPputXG, BGPputXH, BGPputXI, BGPputXJ, BGPputXK
 BGPputXL, BGPputXM, BGPputXN, BGPputXO, BGPputXP, BGPputXQ, BGPputXR
 BGPputXS, BGPputXT, BGPputXU
 \$\$AGE^AUPNPAT
 BGP3D1, BGP3D3, BGP4D1, BGP4D10, BGP4D7, BGP5C1, BGP5C2, BGP5C3
 BGP5D1, BGP5D10, BGP5D6, BGP5D7, BGP6C1, BGP6C2, BGP6C3, BGP6D1
 BGP6D10, BGP6D23, BGP6D6, BGP6D7, BGP6D83, BGP6D84, BGP6D85
 BGP6DCHW, BGP7C1, BGP7C2, BGP7C3, BGP7D1, BGP7D10, BGP7D21
 BGP7D23, BGP7D6
 BGP7D7, BGP7D83, BGP7D84, BGP7D85, BGP7DCHW, BGP7DPE1, BGP7D1
 BGP7D2, BGP7D22
 \$\$BEN^AUPNPAT
 BGP3D1, BGP4D1, BGP4D10, BGP5D1, BGP5D10, BGP6D1, BGP6D10
 BGP7D1, BGP7D10, BGP7DPE1, BGP7D1
 \$\$COMMRES^AUPNPAT
 BGP6DCHW, BGP7DCHW, BGP7DNE
 \$\$DOB^AUPNPAT
 BGP3D1, BGP3D2, BGP3D3, BGP3D4, BGP3DCI, BGP4D1, BGP4D10
 BGP4D2, BGP4D3, BGP4D4, BGP5C1, BGP5C11, BGP5C13, BGP5C2
 BGP5D1, BGP5D10, BGP5D2, BGP5D22, BGP5D3, BGP5D32, BGP5D33

BGP5D34 ,BGP5D35 ,BGP5D4 ,BGP5D61 ,BGP5D71 ,BGP5D72 ,BGP5D9
BGP5EL2 ,BGP5EL3 ,BGP6C1
BGP6C11 ,BGP6C13 ,BGP6C2 ,BGP6D1 ,BGP6D10 ,BGP6D2 ,BGP6D21
BGP6D22 ,BGP6D3 ,BGP6D32 ,BGP6D33 ,BGP6D34 ,BGP6D35 ,BGP6D4

BGP6D61 ,BGP6D71 ,BGP6D72 ,BGP6D82 ,BGP6D85 ,BGP6D862 ,BGP6D9
BGP6DHE1 ,BGP6EL2 ,BGP6EL3 ,BGP7C1 ,BGP7C11 ,BGP7C13 ,BGP7C2
BGP7D1 ,BGP7D10

BGP7D2 ,BGP7D21 ,BGP7D211 ,BGP7D21A ,BGP7D22 ,BGP7D3 ,BGP7D32
BGP7D33 ,BGP7D34 ,BGP7D35 ,BGP7D4 ,BGP7D61 ,BGP7D71 ,BGP7D72
BGP7D721 ,BGP7D731 ,BGP7D732 ,BGP7D75 ,BGP7D82 ,BGP7D85
BGP7D862 ,BGP7D9 ,BGP7DHE1 ,BGP7EL2 ,BGP7EL3 ,BGPD1 ,BGPD24
BGP3D1 ,BGP4D1 ,BGP4D10 ,BGP5D1 ,BGP5D10 ,BGP6CP ,BGP6D1
BGP6D10 ,BGP7CP ,BGP7D1 ,BGP7D10 ,BGP7DPE1 ,BGPD1

\$\$DOD^AUPNPAT

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 425

```

Routine          is Invoked by:
  $$HRN^AUPNPAT
BGP3DS, BGP3PDL, BGP4DS, BGP4HES, BGP4HESL, BGP4PDL, BGP5C1
      BGP5C2, BGP5C3, BGP5DS, BGP5ELS, BGP5ELSL, BGP5HES, BGP5HESL

BGP5NPLD, BGP5NPLP, BGP5PDL, BGP6C1, BGP6C2, BGP6C3, BGP6DCLD
      BGP6DCLP, BGP6DS, BGP6ELS, BGP6ELSL, BGP6HES, BGP6HESL
      BGP6NPLD, BGP6NPLP

BGP6PDL, BGP7C1, BGP7C2, BGP7C3, BGP7DCLD, BGP7DCLP, BGP7DPEE
      BGP7DPEL, BGP7DS, BGP7ELS, BGP7ELSL, BGP7HES, BGP7HESL
      BGP7NPLD, BGP7NPLP, BGP7PDL, BGP7PDS
  $$TRIBE^AUPNPAT  BGP6DCHW, BGP7DCHW
  KILL^AUPNPAT     BGP3DA, BGP3DAR, BGP3DARA, BGP3DARH, BGP3DG, BGP3DL, BGP4DA

BGP4DAP, BGP4DAR, BGP4DARA, BGP4DG, BGP4DHE, BGP4DHEA, BGP4DL

BGP4DPP, BGP5CL, BGP5DAP, BGP5DAR, BGP5DEL, BGP5DELA, BGP5DGPA
      BGP5DGPU, BGP5DHE, BGP5DHEA, BGP5DL, BGP5DNG, BGP5DPP
      BGP5GAEL, BGP5GAGP
      BGP5GAHE, BGP5GALL, BGP5GAPU, BGP5GCMS, BGP5GCOM, BGP5GELD
      BGP5GGPU, BGP5GHED, BGP5GLST, BGP5GNPL, BGP5GNTL, BGP5GPP

BGP5NPL, BGP6AHW, BGP6CL, BGP6DAP, BGP6DAR, BGP6DCHA, BGP6DCL
      BGP6DEL, BGP6DELA, BGP6DGPA, BGP6DGPU, BGP6DHE, BGP6DHEA
      BGP6DL, BGP6DNG, BGP6DPP
      BGP6DSTM, BGP6GACW, BGP6GAEL, BGP6GAGP, BGP6GAHE, BGP6GALL
      BGP6GAPU, BGP6GCMP, BGP6GCMS, BGP6GCOM, BGP6GELD, BGP6GGPU
      BGP6GHED, BGP6GLST, BGP6GLTX, BGP6GMTX, BGP6GNPL, BGP6GNST

BGP6GNTL, BGP6GPP, BGP6NPL, BGP7CL, BGP7DAP, BGP7DAR, BGP7DCHA
      BGP7DCL, BGP7DEL
      BGP7DELA, BGP7DGPA, BGP7DGPU, BGP7DHE, BGP7DHEA, BGP7DL
      BGP7DNE, BGP7DNG, BGP7DPE, BGP7DPEA, BGP7DPP, BGP7DSTM
      BGP7GACW, BGP7GAEL, BGP7GAGP, BGP7GAHE, BGP7GALL, BGP7GAPE
      BGP7GAPU, BGP7GCMP, BGP7GCMS, BGP7GCOM, BGP7GELD, BGP7GGPU
      BGP7GHED, BGP7GLST, BGP7GLTX
      BGP7GMTX, BGP7GNPL, BGP7GNST, BGP7GNTL, BGP7GNXP, BGP7GPED
      BGP7GPP, BGP7NPL, BGP7D, BGP7DAREA, BGP7DARP, BGP7DFTA

  ^BGP3BAN          | opt
  ^BGP3D1           BGP3DA, BGP3DG, BGP3DL
  $$LASTHF^BGP3D6  BGP3D8
  $$PED^BGP3D6     BGP3D9
  $$FIRSTIHD^BGP3D9 BGP3D1
  $$V2IHD^BGP3D9  BGP3D1
  ^BGP3DA          | opt
  ^BGP3DAR         | opt
  ^BGP3DARA        | opt
  CALCIND^BGP3DCI  BGP3D1
  ^BGP3DG          | opt
  ^BGP3DH          BGP3DP
  AREACP^BGP3DH    BGP3PARP
  ^BGP3DL          | opt

```

^BGP3DP	BGP3DA , BGP3DG , BGP3DL
PRINT1^BGP3DP	BGP3PARP
PI^BGP3DP1	BGP3DP3 , BGP3DP4 , BGP3DP5 , BGP3DP6
I1AGE^BGP3DP41	BGP3DP4

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 426

Routine	is Invoked by:
I1AGEP^BGP3DP50	BGP3DP5, BGP3DP51
I1AGE^BGP3DP51	BGP3DP5
I1AGE^BGP3DP60	BGP3DP6
H1^BGP3DPH	BGP3DP1, BGP3DP4, BGP3DP5, BGP3DP6, BGP3DP7
H2^BGP3DPH	BGP3DP1, BGP3DP6, BGP3DP7
H3^BGP3DPH	BGP3DP1
HEADER^BGP3DPH	BGP3DH, BGP3DP, BGP3DP1, BGP3DP3, BGP3DP4, BGP3DP41, BGP3DP5 BGP3DP50, BGP3DP51, BGP3DP6, BGP3DP60, BGP3DP7
^BGP3DS	BGP3DP
EN^BGP3DSI	BGP3DL
EN^BGP3DSL	BGP3DL
PT^BGP3DSL	BGP3DA, BGP3DAR, BGP3DARA, BGP3DARH, BGP3DG, BGP3DL
RT^BGP3DSL	BGP3DL
PRINT^BGP3PARP	BGP3DAR, BGP3DARA, BGP3DARH
^BGP3PDL	BGP3DP, BGP3PARP
H1^BGP3PDL1	BGP3PDL3, BGP3PDL4, BGP3PDL6, BGP3PDL7
PI^BGP3PDL1	BGP3PDL3, BGP3PDL4, BGP3PDL6, BGP3PDL7
I1AGE^BGP3PDL5	BGP3PDL4
I1AGEP^BGP3PDL8	BGP3PDL7, BGP3PDL9
I1AGE^BGP3PDL9	BGP3PDL7
^BGP3PDL5	BGP3PDL
I1AGE^BGP3PDM1	BGP3PDL6
^BGP3TXCH	opt
TAXCHK^BGP3TXCH	BGP3DA, BGP3DG, BGP3DL
^BGP3ULF	opt
\$\$DATE^BGP3UTL	
BGP3D10, BGP3D2, BGP3D3, BGP3D4, BGP3D5, BGP3D6, BGP3D7, BGP3D8	BGP3D9, BGP4D31
EISSGS^BGP3UTL	BGP3DA
GETIEN^BGP3UTL	BGP3ULF
GS^BGP3UTL	BGP3DA, BGP3DG
REPORT^BGP3UTL	BGP3DA, BGP3DG, BGP3DL
STMP^BGP3UTL	BGP3DCI
\$\$REFUSAL^BGP3UTL1	BGP3D2, BGP3D3, BGP3D4, BGP3D6, BGP3D8
^BGP4BAN	opt
^BGP4D1	BGP4DA, BGP4DAP, BGP4DG, BGP4DHE, BGP4DL
^BGP4D10	BGP4DPP
\$\$HDL^BGP4D2	BGP4D41, BGP4D9
\$\$HGBA1C^BGP4D2	BGP4DHE1
\$\$LDL^BGP4D2	BGP4D41, BGP4D9, BGP4DHE1
\$\$LIPID^BGP4D2	BGP4D41, BGP4D9
\$\$TRIG^BGP4D2	BGP4D41, BGP4D9
\$\$EYE^BGP4D21	BGP4DHE1
\$\$POSUR^BGP4D21	BGP4DHE1
\$\$PNEU^BGP4D31	BGP4D3
\$\$HYSTER^BGP4D4	BGP4D3
\$\$MEANBP^BGP4D41	BGP4D9
\$\$CHL^BGP4D5	BGP4DHE1
\$\$BMI^BGP4D6	BGP4D71
\$\$OB^BGP4D6	BGP4D71

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 427

Routine	is Invoked by:
\$\$OW^BGP4D6	BGP4D71
\$\$FOB^BGP4D61	BGP4DHE1
\$\$IBH^BGP4D61	BGP4DHE1
\$\$RECT^BGP4D61	BGP4DHE1
\$\$DENT^BGP4D7	BGP4D71
\$\$DX^BGP4D7	BGP4D71
\$\$TOBACCO^BGP4D7	BGP4D71
^BGP4DA	opt
^BGP4DAP	opt
^BGP4DAR	opt
^BGP4DARA	opt
CALCIND^BGP4DCHI	BGP4D1
CALCIND^BGP4DCI	BGP4D1 , BGP4D10
^BGP4DG	opt
^BGP4DH	BGP4DP
AREACP^BGP4DH	BGP4PARP
^BGP4DHE	opt
^BGP4DHEA	opt
^BGP4DHEP	BGP4DHE
PRINT1^BGP4DHEP	BGP4PHEP
EN^BGP4DHSL	BGP4DHE
PT^BGP4DHSL	BGP4DHE , BGP4DHEA
RT^BGP4DHSL	BGP4DHE
^BGP4DL	opt
^BGP4DP	BGP4DA , BGP4DAP , BGP4DG , BGP4DL , BGP4DPP
PRINT1^BGP4DP	BGP4PARP
SETEXCEL^BGP4DP	BGP4PARP
PI^BGP4DP1	BGP4DP3 , BGP4DP4 , BGP4DP5 , BGP4DP6
SETEXA^BGP4DP1	BGP4DP2 , BGP4DP3 , BGP4PDL1 , BGP4PDL2 , BGP4PDL3
SETN^BGP4DP1	BGP4DP2 , BGP4DP3 , BGP4PDL0 , BGP4PDL1 , BGP4PDL2 , BGP4PDL3
I1AGE^BGP4DP11	BGP4DP1
I1AGEP^BGP4DP50	BGP4DP5 , BGP4DP51
I1AGEP^BGP4DP60	BGP4DP6
H1^BGP4DPH	BGP4DP1 , BGP4DP4 , BGP4DP5 , BGP4DP6 , BGP4HEP1
H2^BGP4DPH	BGP4DP1 , BGP4HEP1
H3^BGP4DPH	BGP4DP11
HEADER^BGP4DPH	
BGP4DHEP , BGP4DP , BGP4DP1 , BGP4DP11 , BGP4DP3 , BGP4DP4 , BGP4DP5	
	BGP4DP50 , BGP4DP51 , BGP4DP6 , BGP4DP60 , BGP4DS , BGP4DSP
	BGP4HEP1 , BGP4HES
^BGP4DPP	opt
^BGP4DS	BGP4DP
EN^BGP4DSI	BGP4DAP , BGP4DL , BGP4DPP
EN^BGP4DSL	BGP4DAP , BGP4DL , BGP4DPP
PT^BGP4DSL	BGP4DA , BGP4DAP , BGP4DAR , BGP4DARA , BGP4DG , BGP4DL , BGP4DPP
RT^BGP4DSL	BGP4DAP , BGP4DL , BGP4DPP
^BGP4DSP	BGP4DP
^BGP4HEH	BGP4DHEP
AREACP^BGP4HEH	BGP4PHEP
^BGP4HEHH	BGP4HEL

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 428

Routine	is Invoked by:
^BGP4HEL	BGP4DHEP ,BGP4PHEP
^BGP4HES	BGP4DHEP
^BGP4HESL	BGP4HEL
GETIEN^BGP4HUTL	BGP4ULF
GS^BGP4HUTL	BGP4DHE
REPORT^BGP4HUTL	BGP4DHE
STMP^BGP4HUTL	BGP4DCHI
^BGP4LSTF	opt
PRINT^BGP4PARP	BGP4DAR ,BGP4DARA
^BGP4PDL	BGP4DP ,BGP4PARP
HEADER^BGP4PDL	BGP4PDLZ
H1^BGP4PDL1	BGP4PDL0 ,BGP4PDL3 ,BGP4PDL4 ,BGP4PDL7
PI^BGP4PDL1	BGP4PDL0 ,BGP4PDL3 ,BGP4PDL4 ,BGP4PDL7
I1AGEP^BGP4PDL5	BGP4PDL0
I1AGEP^BGP4PDL8	BGP4PDL7
I1AGE^BGP4PDL9	BGP4PDL1
^BGP4PDLH	BGP4PDL
^BGP4PDL5	BGP4PDL
PRINT^BGP4PHEP	BGP4DHEA
^BGP4TXCH	opt
TAXCHK^BGP4TXCH	BGP4DA ,BGP4DAP ,BGP4DG ,BGP4DL ,BGP4DPP
TAXCHK^BGP4TXHE	BGP4DHE
^BGP4ULF	opt
\$\$D^BGP4UTL	BGP4DA ,BGP4DAR ,BGP4DARA ,BGP4DG
\$\$DATE^BGP4UTL	BGP4D2 ,BGP4D21 ,BGP4D3 ,BGP4D31 ,BGP4D4 ,BGP4D41 ,BGP4D5 BGP4D61 ,BGP4D7 ,BGP4D71 ,BGP4D9 ,BGP4DHE1 ,BGP4DP
\$\$LZERO^BGP4UTL	BGP4DAR ,BGP4DARA

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>

\$\$RZERO^BGP4UTL	BGP4DA ,BGP4DAR ,BGP4DARA ,BGP4DG
EXCELGS^BGP4UTL	BGP4DA ,BGP4DAR ,BGP4DARA ,BGP4DG
GETIEN^BGP4UTL	BGP4ULF
GS^BGP4UTL	BGP4DA ,BGP4DG
JRNL^BGP4UTL	BGP4D1 ,BGP4D10
REPORT^BGP4UTL	BGP4DA ,BGP4DAP ,BGP4DG ,BGP4DL ,BGP4DPP
STMP^BGP4UTL	BGP4DCI
\$\$REFUSAL^BGP4UTL1	BGP4D21 ,BGP4D3 ,BGP4D31 ,BGP4D4 ,BGP4D5 ,BGP4D61 ,BGP4D8
EN^BGP5ASL	BGP5DAR ,BGP5DELA ,BGP5DGPA ,BGP5DHEA
GET^BGP5ASL	BGP5GU ,BGP5GU
^BGP5BAN	opt
^BGP5C1	BGP5CL ,BGP5GCMS
GETMEDS^BGP5C1	BGP5C11 ,BGP5C13 ,BGP5C3
ACEIALG1^BGP5C11	BGP5C2
ACEICON1^BGP5C11	BGP5C2
ACEIRX1^BGP5C11	BGP5C2
EN^BGP5C11	BGP5C1
^BGP5C12	BGP5C11
EN^BGP5C12	BGP5C2 ,BGP5C3
ARBALG1^BGP5C13	BGP5C11

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 429

Routine	is Invoked by:
ARBCON1^BGP5C13	BGP5C11
ARBRX1^BGP5C13	BGP5C11
^BGP5CL	opt
^BGP5CP	BGP5CL
EP^BGP5CSI	BGP5CL
EP^BGP5CTL	BGP5CTS, BGP5ETS
^BGP5CTS	opt
^BGP5CTXC	opt
TAXCHK^BGP5CTXC	BGP5CL
^BGP5D1	
BGP5DAP, BGP5DEL, BGP5DELA, BGP5DGPU, BGP5DHE, BGP5DL, BGP5DNG	BGP5GALL, BGP5GCOM, BGP5GELD, BGP5GGPU, BGP5GHED, BGP5GNPL
	BGP5GNTL, BGP5NPL
^BGP5D10	BGP5DPP, BGP5GPP
\$\$HGBA1C^BGP5D2	BGP5D22, BGP5DHE1
\$\$LDL^BGP5D2	BGP5D22, BGP5D41, BGP5D72, BGP5D73, BGP5DHE1, BGP5EL3
\$\$MEANBP^BGP5D2	BGP5D22
\$\$DENTSRV^BGP5D21	BGP5DHE1, BGP5EL2
\$\$DEP^BGP5D21	BGP5D41, BGP5EL3
\$\$DEPSCR^BGP5D21	BGP5D41, BGP5EL3
\$\$DNKA^BGP5D21	BGP5DHE1
\$\$EYE^BGP5D21	BGP5D22, BGP5EL2
\$\$GFR^BGP5D21	BGP5D22, BGP5EL2
\$\$POSUR^BGP5D21	BGP5D22, BGP5EL2
\$\$FLU^BGP5D3	BGP5EL3
\$\$HOSP^BGP5D31	BGP5EL3
\$\$LAST^BGP5D31	BGP5EL3
\$\$PNEU^BGP5D31	BGP5D3, BGP5EL3
\$\$V2ASTH^BGP5D31	BGP5EL3
\$\$IMMREF^BGP5D32	BGP5D33, BGP5D34, BGP5D35
GETIMMS^BGP5D32	BGP5D33, BGP5D34, BGP5D35
\$\$MMR^BGP5D33	BGP5D32
\$\$OPV^BGP5D33	BGP5D32
\$\$DTAP^BGP5D34	BGP5D32
\$\$HEP^BGP5D35	BGP5D32
\$\$HIB^BGP5D35	BGP5D32
\$\$VAR^BGP5D35	BGP5D32
\$\$MAM^BGP5D4	BGP5EL3
\$\$MAS^BGP5D4	BGP5EL3
\$\$MEANBP^BGP5D41	BGP5D9, BGP5EL3
\$\$CHL^BGP5D5	BGP5DHE1
\$\$DVDX^BGP5D5	BGP5EL3
\$\$DVEX^BGP5D5	BGP5EL3
\$\$DVPED^BGP5D5	BGP5EL3
\$\$REFDV^BGP5D5	BGP5EL3
\$\$BMI^BGP5D6	BGP5D41, BGP5D71, BGP5EL3
\$\$OB^BGP5D6	BGP5D71, BGP5EL3
\$\$OW^BGP5D6	BGP5D71, BGP5EL3
\$\$CRC^BGP5D61	BGP5DHE1, BGP5EL3
\$\$FOB^BGP5D61	BGP5DHE1, BGP5EL3

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 430

Routine	is Invoked by:
\$\$IBH^BGP5D61	BGP5DHE1, BGP5EL3
\$\$RECT^BGP5D61	BGP5DHE1, BGP5EL3
\$\$REF^BGP5D61	BGP5DHE1, BGP5EL3
\$\$DENT^BGP5D7	BGP5D41, BGP5D71, BGP5EL3
\$\$DX^BGP5D7	BGP5D71, BGP5EL3
\$\$PED^BGP5D7	BGP5D41, BGP5EL3
\$\$PREG^BGP5D7	BGP5D8
\$\$TOBACCO^BGP5D7	BGP5D41, BGP5D71, BGP5EL3
\$\$MEDNUTR^BGP5D71	BGP5D41
\$\$OTHREL^BGP5D71	BGP5D41
\$\$SPECEX^BGP5D71	BGP5D41
\$\$SPECNUTR^BGP5D71	BGP5D41
\$\$CHOL^BGP5D73	BGP5EL3
\$\$FIRSTIHD^BGP5D9	BGP5D21, BGP5D41, BGP5D73, BGP5DHE1
\$\$V2IHD^BGP5D9	BGP5D21, BGP5D41, BGP5D73, BGP5DHE1
^BGP5DAP	opt
^BGP5DAR	opt
CALCIND^BGP5DCEI	BGP5D1
CALCIND^BGP5DCHI	BGP5D1
CALCIND^BGP5DCI	BGP5D1, BGP5D10
^BGP5DEL	opt
^BGP5DELA	opt
^BGP5DELP	BGP5DEL, BGP5DELA
PRINT1^BGP5DELP	BGP5PHEL
EN^BGP5DESI	BGP5DEL
EN^BGP5DESL	BGP5DEL, BGP5DELA
PT^BGP5DESL	BGP5DEL, BGP5DELA
RT^BGP5DESL	BGP5DEL, BGP5DELA
^BGP5DGPA	opt
^BGP5DGPU	opt
^BGP5DH	BGP5DP
AREACP^BGP5DH	BGP5PARP
^BGP5DHE	opt
^BGP5DHEA	opt
^BGP5DHEP	BGP5DHE
PRINT1^BGP5DHEP	BGP5PHEP
EN^BGP5DHSL	BGP5DHE
PT^BGP5DHSL	BGP5DHE, BGP5DHEA
RT^BGP5DHSL	BGP5DHE
^BGP5DL	opt
^BGP5DNG	opt
^BGP5DP	BGP5DAP, BGP5DGPU, BGP5DL, BGP5DNG, BGP5DPP, BGP5NPL
PRINT1^BGP5DP	BGP5PARP
SETEXCEL^BGP5DP	BGP5PARP
PI^BGP5DP1	BGP5DP3, BGP5DP4, BGP5DP5, BGP5DP6
PI1^BGP5DP1	BGP5DP4, BGP5DP5
SETEXA^BGP5DP1	BGP5DP2, BGP5DP3, BGP5DP7, BGP5PDL1, BGP5PDL2, BGP5PDL3
SETN^BGP5DP1	BGP5DP2, BGP5DP3, BGP5DP7, BGP5PDL0, BGP5PDL1, BGP5PDL2, BGP5PDL3

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 431

Routine	is Invoked by:
I1AGE^BGP5DP11	BGP5DP1
I1AGEP^BGP5DP50	BGP5DP5, BGP5DP51
I1AGEP^BGP5DP60	BGP5DP6
H1^BGP5DPH	BGP5DP1, BGP5DP4, BGP5DP5, BGP5DP6, BGP5ELP1, BGP5HEP1
H2^BGP5DPH	BGP5DP1, BGP5ELP1, BGP5HEP1
H3^BGP5DPH	BGP5DP11
HEADER^BGP5DPH	BGP5DELP, BGP5DHEP, BGP5DP, BGP5DP1, BGP5DP11, BGP5DP3 BGP5DP4, BGP5DP5, BGP5DP50, BGP5DP51, BGP5DP6, BGP5DP60 BGP5DS, BGP5DSP, BGP5EL10, BGP5EL11, BGP5EL12, BGP5EL13 BGP5EL14, BGP5EL15, BGP5EL16, BGP5EL17, BGP5EL18, BGP5EL19 BGP5EL1A, BGP5EL1B, BGP5EL1C BGP5EL1D, BGP5EL1E, BGP5EL1F, BGP5EL1G, BGP5EL1H, BGP5EL1I BGP5EL1J, BGP5EL1K, BGP5EL1L, BGP5ELP1, BGP5ELS, BGP5HEP1 BGP5HES, BGP5SDP
^BGP5DPP	opt
^BGP5DS	BGP5DP
EN^BGP5DSI	BGP5DAP, BGP5DL, BGP5DPP, BGP5NPL
EN^BGP5DSL	BGP5DAP, BGP5DL, BGP5DPP
PT^BGP5DSL	
BGP5DAP, BGP5DAR, BGP5DGPA, BGP5DGPU, BGP5DL, BGP5DNG, BGP5DPP	BGP5NPL
RT^BGP5DSL	BGP5DAP, BGP5DL, BGP5DPP, BGP5NPL
^BGP5DSP	BGP5DP
I1AGE^BGP5EL10	BGP5ELP1
I1AGE^BGP5EL11	BGP5ELP1
I1AGE^BGP5EL12	BGP5ELP1
I1AGE^BGP5EL13	BGP5ELP1
I1AGE^BGP5EL14	BGP5ELP1
I1AGE^BGP5EL15	BGP5ELP1
I1AGE^BGP5EL16	BGP5ELP1
I1AGE^BGP5EL17	BGP5ELP1
I1AGE^BGP5EL18	BGP5ELP1
I1AGE^BGP5EL19	BGP5ELP1
I1AGE^BGP5EL1A	BGP5ELP1
I1AGE^BGP5EL1B	BGP5ELP1
I1AGE^BGP5EL1C	BGP5ELP1
I1AGE^BGP5EL1D	BGP5ELP1
I1AGE^BGP5EL1E	BGP5ELP1
I1AGE^BGP5EL1F	BGP5ELP1
I1AGE^BGP5EL1G	BGP5ELP1
I1AGE^BGP5EL1H	BGP5ELP1
I1AGE^BGP5EL1I	BGP5ELP1
I1AGE^BGP5EL1J	BGP5ELP1
I1AGE^BGP5EL1K	BGP5ELP1
I22^BGP5EL1L	BGP5ELP1
\$\$FRACTURE^BGP5EL3	BGP5D22
\$\$FUNCTION^BGP5EL4	BGP5EL3
\$\$PHNV^BGP5EL4	BGP5EL3
\$\$TXBMD^BGP5EL4	BGP5D22, BGP5EL3
^BGP5ELH	BGP5DELP

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
432

Routine	is Invoked by:
AREACP^BGP5ELH	BGP5PHEL
^BGP5ELHH	BGP5ELL
^BGP5ELL	BGP5DELP ,BGP5PHEL
SETN^BGP5ELL1	BGP5ELLN
I1AGE^BGP5ELL2	BGP5ELL1
I1AGE^BGP5ELL3	BGP5ELL1
I1AGE^BGP5ELL4	BGP5ELL1
I1AGE^BGP5ELL5	BGP5ELL1
I1AGE^BGP5ELL6	BGP5ELL1
I1AGE^BGP5ELL7	BGP5ELL1
I1AGE^BGP5ELL8	BGP5ELL1
I1AGE^BGP5ELL9	BGP5ELL1
I1AGE^BGP5ELLA	BGP5ELL1
I1AGE^BGP5ELLB	BGP5ELL1
I1AGE^BGP5ELLC	BGP5ELL1
I1AGE^BGP5ELLD	BGP5ELL1
I1AGE^BGP5ELLE	BGP5ELL1
I1AGE^BGP5ELLF	BGP5ELL1
I1AGE^BGP5ELLG	BGP5ELL1
I1AGE^BGP5ELLH	BGP5ELL1
I1AGE^BGP5ELLI	BGP5ELL1
I1AGE^BGP5ELLJ	BGP5ELL1
I1AGE^BGP5ELLK	BGP5ELL1
I1AGE^BGP5ELLL	BGP5ELL1
I1AGE^BGP5ELLM	BGP5ELL1
I22^BGP5ELLN	BGP5ELL1
I1AGEP^BGP5ELLQ	BGP5ELL3
SETN^BGP5ELP1	BGP5EL1L
^BGP5ELS	BGP5DELP
^BGP5ELSL	BGP5ELL
GS^BGP5EUTL	BGP5DEL ,BGP5DELA ,BGP5GELD
REPORT^BGP5EUTL	BGP5DEL ,BGP5GELD
STMP^BGP5EUTL	BGP5DCEI
EP^BGP5GAEL	BGP5GRA ,BGP5GRA
EP^BGP5GAGP	BGP5GRA ,BGP5GRA
EP^BGP5GAHE	BGP5GRA ,BGP5GRA
EP^BGP5GALL	BGP5GR ,BGP5GR
EP^BGP5GAPU	BGP5GRA ,BGP5GRA
EP^BGP5GCMS	BGP5GR ,BGP5GR
EP^BGP5GCOM	BGP5GR ,BGP5GR
EP^BGP5GELD	BGP5GR ,BGP5GR
EP^BGP5GGPU	BGP5GR ,BGP5GR
EP^BGP5GHED	BGP5GR ,BGP5GR
EP^BGP5GLST	BGP5GU ,BGP5GU
EP^BGP5GNPL	BGP5GR ,BGP5GR
EP^BGP5GNTL	BGP5GR ,BGP5GR
EP^BGP5GPP	BGP5GR ,BGP5GR
EP^BGP5GUPL	BGP5GU ,BGP5GU
^BGP5HEH	BGP5DHEP

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 433

Routine	is Invoked by:
AREACP^BGP5HEH	BGP5PHEP
^BGP5HEHH	BGP5HEL
^BGP5HEL	BGP5DHEP, BGP5PHEP
^BGP5HES	BGP5DHEP
^BGP5HESL	BGP5HEL
GETIEN^BGP5HUTL	BGP5GUPL, BGP5ULF
GS^BGP5HUTL	BGP5DHE, BGP5GHED
REPORT^BGP5HUTL	BGP5DHE, BGP5GHED
STMP^BGP5HUTL	BGP5DCHI
^BGP5LSTF	opt
^BGP5NPL	opt
NPL1^BGP5NPLD	BGP5PDL
NPL1^BGP5NPLP	BGP5DP
PRINT^BGP5PARP	BGP5DAR, BGP5DGPA
^BGP5PDL	BGP5DP, BGP5PARP
\$\$\$SB^BGP5PDL1	BGP5D6, BGP5EL3
H1^BGP5PDL1	BGP5PDL0, BGP5PDL3, BGP5PDL4, BGP5PDL7
PI^BGP5PDL1	BGP5PDL0, BGP5PDL3, BGP5PDL4, BGP5PDL7
PI1^BGP5PDL1	BGP5PDL4, BGP5PDL7
I1AGEP^BGP5PDL5	BGP5PDL0
I1AGEP^BGP5PDL8	BGP5PDL7
I1AGE^BGP5PDL9	BGP5PDL1
^BGP5PDLH	BGP5PDL
^BGP5PDL5	BGP5PDL
PRINT^BGP5PHEL	BGP5DELA
PRINT^BGP5PHEP	BGP5DHEA
^BGP5SDP	BGP5DSP
^BGP5TXCH	opt
TAXCHK^BGP5TXCH	BGP5DAP, BGP5DL, BGP5DPP
^BGP5TXCN	opt
TAXCHK^BGP5TXCN	BGP5DGPU, BGP5DNG, BGP5NPL
^BGP5TXEL	opt
TAXCHK^BGP5TXEL	BGP5DEL
^BGP5TXHE	opt
TAXCHK^BGP5TXHE	BGP5DHE
EP^BGP5TXL	BGP5TXS
^BGP5TXV	opt
EP^BGP5TXV1	BGP5CTS, BGP5TXV
^BGP5ULF	opt
\$\$\$D^BGP5UTL	BGP5DAR, BGP5DNG, BGP5GAGP
\$\$\$DATE^BGP5UTL	BGP5ASL, BGP5C1, BGP5C11, BGP5C12, BGP5C13, BGP5C2, BGP5C3 BGP5D2, BGP5D21, BGP5D22, BGP5D3, BGP5D31, BGP5D4, BGP5D41 BGP5D5, BGP5D6, BGP5D61, BGP5D7, BGP5D71, BGP5D72, BGP5D73 BGP5D8, BGP5D9, BGP5DGPA, BGP5DHE1, BGP5DP, BGP5EL2, BGP5EL3 BGP5EL4, BGP6DHE1 BGP7DHE1
\$\$\$LZERO^BGP5UTL	BGP5DAR, BGP5GAGP
<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>	
\$\$\$RZERO^BGP5UTL	BGP5DAR, BGP5DNG, BGP5GAGP

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 434

Routine	is Invoked by:
EXCELGS^BGP5UTL	BGP5DAR, BGP5DNG, BGP5GAGP
GETIEN^BGP5UTL	BGP5GUPL, BGP5ULF
GS^BGP5UTL	BGP5DGPU, BGP5DNG, BGP5GGPU, BGP5GNTL
JRNL^BGP5UTL	BGP5D1, BGP5D10
REPORT^BGP5UTL	BGP5DAP, BGP5DGPU, BGP5DL, BGP5DNG, BGP5DPP, BGP5GALL BGP5GCOM, BGP5GGPU, BGP5GNPL, BGP5GNTL, BGP5GPP, BGP5NPL
STMP^BGP5UTL	BGP5DCI
\$\$FIRSTPRC^BGP5UTL1	BGP5EL3
\$\$LASTDX^BGP5UTL1	BGP5D72
\$\$LASTPRC^BGP5UTL1	BGP5D72, BGP5DHE1, BGP5EL4
\$\$LASTPRCI^BGP5UTL1	BGP5D3, BGP5D61
\$\$RADREF^BGP5UTL1	BGP5D61
\$\$REFUSAL^BGP5UTL1	BGP5D21, BGP5D3, BGP5D31, BGP5D4, BGP5D5, BGP5D61, BGP5D8 BGP5DHE1
EN^BGP6ASL	BGP6DAR, BGP6DCHA, BGP6DELA, BGP6DGPA, BGP6DHEA
GET^BGP6ASL	BGP6GU
^BGP6AX	BGPPOST
^BGP6AXA	BGP6AX
^BGP6AXB	BGP6AX
^BGP6AXC	BGP6AX
^BGP6AXD	BGP6AX
^BGP6AXF	BGP6AX
^BGP6AXG	BGP6AX
^BGP6AXGB	BGP6AXG
^BGP6AXGC	BGP6AXG
^BGP6AXH	BGP6AX
^BGP6AXI	BGP6AX
^BGP6AXJ	BGP6AX
^BGP6AXM	BGP6AX
^BGP6AXN	BGP6AX
^BGP6AXO	BGP6AX
^BGP6AXP	BGP6AX
^BGP6AXS	BGP6AX
^BGP6AXT	BGP6AX
^BGP6AXU	BGP6AX
^BGP6AXV	BGP6AX
^BGP6AXW	BGP6AX
^BGP6AXX	BGP6AX
^BGP6AXY	BGP6AX
^BGP6AXZ	BGP6AX
^BGP6BAN	opt
^BGP6BX	BGPPOST
^BGP6BXA	BGP6BX
^BGP6BXB	BGP6BX
^BGP6BXBB	BGP6BXB
^BGP6BXBC	BGP6BXB
^BGP6BXBD	BGP6BXB
^BGP6BXBE	BGP6BXB
^BGP6BXBF	BGP6BXB

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
435

Routine	is Invoked by:
^BGP6BXBG	BGP6BXB
^BGP6BXC	BGP6BX
^BGP6BXD	BGP6BX
^BGP6BXE	BGP6BX
^BGP6BXF	BGP6BX
^BGP6BXG	BGP6BX
^BGP6BXH	BGP6BX
^BGP6BXJ	BGP6BX
^BGP6BXL	BGP6BX
^BGP6BXM	BGP6BX
^BGP6BXN	BGP6BX
^BGP6BXO	BGP6BX
^BGP6BXQ	BGP6BX
^BGP6BXR	BGP6BX
^BGP6BXS	BGP6BX
^BGP6BXT	BGP6BX
^BGP6BXU	BGP6BX
^BGP6BXV	BGP6BX
^BGP6BXW	BGP6BX
^BGP6BXX	BGP6BX
^BGP6BXY	BGP6BX
^BGP6C1	BGP6CL, BGP6GCMS
GETMEDS^BGP6C1	BGP6C11, BGP6C12, BGP6C13, BGP6C3
ACEIALG1^BGP6C11	BGP6C2
ACEICON1^BGP6C11	BGP6C2
ACEIRX1^BGP6C11	BGP6C2
EN^BGP6C11	BGP6C1
^BGP6C12	BGP6C11
CESS1^BGP6C12	BGP6C3
EN^BGP6C12	BGP6C2, BGP6C3
SMOKER^BGP6C12	BGP6C1
SMOKER1^BGP6C12	BGP6C3
ARBALG1^BGP6C13	BGP6C11
ARBCON1^BGP6C13	BGP6C11
ARBRX1^BGP6C13	BGP6C11
ABGP01^BGP6C31	BGP6C3
^BGP6CL	opt
^BGP6CP	BGP6CL
EP^BGP6CSI	BGP6CL
EP^BGP6CTL	BGP6CTS
^BGP6CTS	opt
INIT^BGP6CTS	BGP6GU
^BGP6CTXC	opt
TAXCHK^BGP6CTXC	BGP6CL
^BGP6D1	
BGP6AHW, BGP6DAP, BGP6DCL, BGP6DEL, BGP6DGPU, BGP6DHE, BGP6DL	
	BGP6DNG, BGP6DSTM, BGP6GALL, BGP6GCMP, BGP6GCOM, BGP6GELD
	BGP6GGPU, BGP6GHED, BGP6GNPL, BGP6GNST, BGP6GNTL, BGP6NPL
^BGP6D10	BGP6DPP, BGP6GPP
\$\$HGBA1C^BGP6D2	BGP6D22, BGP6DHE1

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 436

Routine	is Invoked by:
\$\$LDL^BGP6D2	
BGP6D22, BGP6D23, BGP6D41, BGP6D72, BGP6D73, BGP6DHE1, BGP6EL4	
\$\$MEANBP^BGP6D2	BGP6D22
\$\$DENTSRV^BGP6D21	BGP6DHE1, BGP6EL2
\$\$DEP^BGP6D21	BGP6D23, BGP6D41, BGP6EL3
\$\$DEPEDU^BGP6D21	BGP6EL3
\$\$DEPREF^BGP6D21	BGP6D23, BGP6EL3
\$\$DEPSCR^BGP6D21	BGP6D23, BGP6D41, BGP6EL3
\$\$DNKA^BGP6D21	BGP6D22, BGP6DHE1
\$\$EYE^BGP6D21	BGP6D22, BGP6EL2
\$\$GFR^BGP6D21	BGP6D22, BGP6EL2
\$\$LOINC^BGP6D21	BGP6C3, BGP6C31, BGP6D8, BGP6D82, BGP6D83, BGP6D84, BGP6D861 BGP6D862
\$\$POSUR^BGP6D21	BGP6D22, BGP6D23, BGP6EL2
\$\$FLU^BGP6D3	BGP6EL3
\$\$HOSP^BGP6D31	BGP6EL3
\$\$LAST^BGP6D31	BGP6EL3
\$\$PNEU^BGP6D31	BGP6D3, BGP6EL3
\$\$V2ASTH^BGP6D31	BGP6EL3
\$\$IMMREF^BGP6D32	BGP6D33, BGP6D34, BGP6D35, BGP6D85
GETIMMS^BGP6D32	BGP6D33, BGP6D34, BGP6D35, BGP6D85
\$\$MMR^BGP6D33	BGP6D32
\$\$OPV^BGP6D33	BGP6D32
\$\$DTAP^BGP6D34	BGP6D32
\$\$HEP^BGP6D35	BGP6D32, BGP6D85
\$\$HIB^BGP6D35	BGP6D32
\$\$PNEUMO^BGP6D35	BGP6D32
\$\$VAR^BGP6D35	BGP6D32, BGP6D85
\$\$MAM^BGP6D4	BGP6EL3
\$\$MAS^BGP6D4	BGP6EL3
\$\$LIFE^BGP6D41	BGP6D23
\$\$MEANBP^BGP6D41	BGP6D9, BGP6EL3
\$\$OSTEOSCR^BGP6D42	BGP6EL3
\$\$CHL^BGP6D5	BGP6DHE1
\$\$DVDX^BGP6D5	BGP6EL3
\$\$DVEX^BGP6D5	BGP6EL3
\$\$DVPED^BGP6D5	BGP6EL3
\$\$REFDV^BGP6D5	BGP6EL3
\$\$EAPT^BGP6D52	BGP6D51
\$\$BMI^BGP6D6	BGP6D23, BGP6D41, BGP6D71, BGP6EL3
\$\$OB^BGP6D6	BGP6D71, BGP6EL3
\$\$OW^BGP6D6	BGP6D71, BGP6EL3
\$\$REF^BGP6D6	BGP6D23, BGP6EL3
\$\$CRC^BGP6D61	BGP6DHE1, BGP6EL3
\$\$FOB^BGP6D61	BGP6DHE1, BGP6EL3
\$\$IBH^BGP6D61	BGP6DHE1, BGP6EL3
\$\$RECT^BGP6D61	BGP6DHE1, BGP6EL3
\$\$REF^BGP6D61	BGP6DHE1, BGP6EL3
\$\$DENT^BGP6D7	
BGP6C12, BGP6C2, BGP6D23, BGP6D41, BGP6D71, BGP6DHE1, BGP6EL3	
\$\$DX^BGP6D7	BGP6C12, BGP6C2, BGP6D23, BGP6D71, BGP6DHE1, BGP6EL3

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 437

Routine	is Invoked by:
\$\$PED^BGP6D7	BGP6D23, BGP6D41, BGP6EL3
\$\$PREG^BGP6D7	BGP6D8
\$\$TOBACCO^BGP6D7	
BGP6C12, BGP6C2, BGP6D23, BGP6D41, BGP6D71, BGP6DHE1, BGP6EL3	
\$\$MEDNUTR^BGP6D71	BGP6D41
\$\$OTHREL^BGP6D71	BGP6D41
\$\$SPECEX^BGP6D71	BGP6D41
\$\$SPECNUTR^BGP6D71	BGP6D41
\$\$CHOL^BGP6D73	BGP6EL4
\$\$CLASS^BGP6D82	BGP7D863
IRAA^BGP6D82	BGP6EL3
\$\$CBC^BGP6D861	BGP6D86
\$\$CBC4^BGP6D861	BGP6D86
\$\$CBC6^BGP6D861	BGP6D86
\$\$GOLDLAB^BGP6D861	BGP6D86
\$\$LFT^BGP6D861	BGP6D86
\$\$LFT6^BGP6D861	BGP6D86
\$\$POT^BGP6D861	BGP6D86
\$\$SERUM12^BGP6D861	BGP6D86
\$\$SERUM6^BGP6D861	BGP6D86
\$\$UP4^BGP6D861	BGP6D86
\$\$RHEUAR^BGP6D862	BGP6D86
\$\$MEDSPRE^BGP6D863	BGP6D86
\$\$FIRSTIHD^BGP6D9	BGP6D21, BGP6D41, BGP6D73, BGP6DHE1
\$\$V2IHD^BGP6D9	BGP6D21, BGP6D41, BGP6D73, BGP6DHE1
^BGP6DAP	opt
^BGP6DAR	opt
CALCIND^BGP6DCEI	BGP6D1
^BGP6DCHA	opt
\$\$D^BGP6DCHA	BGP6GACW
GS^BGP6DCHA	BGP6GACW
CALCIND^BGP6DCHI	BGP6D1
GEN^BGP6DCHW	BGP6D1
CALCIND^BGP6DCI	BGP6D1, BGP6D10
^BGP6DCL	opt
CPPL1^BGP6DCLD	BGP6PDL
CPPL1^BGP6DCLP	BGP6DP
^BGP6DEL	opt
^BGP6DELA	opt
^BGP6DELP	BGP6DEL
PRINT1^BGP6DELP	BGP6PHEL
EN^BGP6DESI	BGP6DEL
EN^BGP6DESL	BGP6DEL, BGP6DELA
PT^BGP6DESL	BGP6DEL, BGP6DELA
RT^BGP6DESL	BGP6DEL, BGP6DELA
^BGP6DGPA	opt
^BGP6DGPU	opt
^BGP6DH	BGP6DP
AREACP^BGP6DH	BGP6PARP
EOP^BGP6DH	BGP6GUPL, BGP6LSTF, BGP6ULF

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 438

Routine	is Invoked by:
^BGP6DHE	opt
^BGP6DHEA	opt
^BGP6DHEP	BGP6DHE
PRINT1^BGP6DHEP	BGP6PHEP
EN^BGP6DHSL	BGP6DHE
PT^BGP6DHSL	BGP6DHE, BGP6DHEA
RT^BGP6DHSL	BGP6DHE
^BGP6DL	opt
^BGP6DLT	opt
^BGP6DMT	opt
^BGP6DNG	opt
HWSF^BGP6DNG	BGP6GNTL
^BGP6DP	
BGP6AHW, BGP6DAP, BGP6DCL, BGP6DGPU, BGP6DL, BGP6DNG, BGP6DPP	BGP6DSTM, BGP6NPL
PRINT1^BGP6DP	BGP6PARP
SETEXCEL^BGP6DP	BGP6PARP
I1AGE^BGP6DP11	BGP6DP1C
I1AGE^BGP6DP12	BGP6DP1
I1AGE^BGP6DP1C	BGP6DP1
PI^BGP6DP1C	BGP6DP1, BGP6DP3, BGP6DP4, BGP6DP5, BGP6DP6
PI1^BGP6DP1C	BGP6DP4, BGP6DP5
SETEXA^BGP6DP1C	BGP6DP2, BGP6DP3, BGP6DP7, BGP6PDL1, BGP6PDL2, BGP6PDL3
SETN^BGP6DP1C	BGP6DP2, BGP6DP3, BGP6DP7, BGP6PDL0, BGP6PDL1, BGP6PDL2 BGP6PDL3
I1AGEP^BGP6DP50	BGP6DP5, BGP6DP51
I1AGEP^BGP6DP60	BGP6DP6
H1^BGP6DPH	BGP6DP1, BGP6DP1C, BGP6DP4, BGP6DP5, BGP6DP6, BGP6ELP1 BGP6HEP1
H2^BGP6DPH	BGP6DP1C, BGP6ELP1, BGP6HEP1
H3^BGP6DPH	BGP6DP11
H4^BGP6DPH	BGP6DP12
HEADER^BGP6DPH	BGP6DELP, BGP6DHEP, BGP6DP, BGP6DP1, BGP6DP11, BGP6DP12 BGP6DP1C, BGP6DP3, BGP6DP4, BGP6DP5, BGP6DP50, BGP6DP51 BGP6DP6, BGP6DP60, BGP6DSP, BGP6DSTM, BGP6EL10, BGP6EL11 BGP6EL12, BGP6EL13, BGP6EL14, BGP6EL15, BGP6EL16, BGP6EL17 BGP6EL18, BGP6EL19, BGP6EL1A BGP6EL1B, BGP6EL1C, BGP6EL1D, BGP6EL1E, BGP6EL1F, BGP6EL1G BGP6EL1H, BGP6EL1I, BGP6EL1J, BGP6EL1K, BGP6EL1L, BGP6EL1M BGP6EL1N, BGP6ELP1, BGP6HEP1, BGP6SDP
^BGP6DPP	opt
^BGP6DS	BGP6DP
EN^BGP6DSI	BGP6AHW, BGP6DAP, BGP6DL, BGP6DPP, BGP6DSTM, BGP6NPL
EN^BGP6DSL	BGP6DAP, BGP6DL, BGP6DPP
PT^BGP6DSL	
BGP6AHW, BGP6DAP, BGP6DAR, BGP6DCL, BGP6DGPA, BGP6DGPU, BGP6DL	BGP6DNG, BGP6DPP, BGP6DSTM, BGP6NPL
RT^BGP6DSL	
BGP6AHW, BGP6DAP, BGP6DCL, BGP6DL, BGP6DPP, BGP6DSTM, BGP6NPL	
^BGP6DSP	BGP6DP
^BGP6DSTM	opt
CT^BGP6DSTM	BGP6DP, BGP6PDL

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 439

Routine	is Invoked by:
\$\$CPT^BGP6DU	BGP6C12, BGP6D21, BGP6D22, BGP6D3, BGP6D4, BGP6D42, BGP6D5
BGP6D61, BGP6D72, BGP6D8, BGP6D83, BGP6D84, BGP6D861, BGP6DHE1	
\$\$CPTI^BGP6DU	BGP6EL4
\$\$PLCODE^BGP6DU	BGP6C1, BGP6C11, BGP6C13, BGP6C3, BGP6D2, BGP6D21, BGP6D23
\$\$PLTAX^BGP6DU	BGP6D31, BGP6D61, BGP6D71, BGP6D72, BGP6D73, BGP6D8, BGP6D82
\$\$RAD^BGP6DU	BGP6D861, BGP6D9, BGP6DHE1, BGP6EL2
\$\$TRAN^BGP6DU	BGP6D34
\$\$TRANI^BGP6DU	BGP6D33, BGP6D34, BGP6D35, BGP6D85
\$\$WH^BGP6DU	BGP6D42, BGP6D61, BGP6EL4
I1AGE^BGP6EL10	BGP6C12
I1AGE^BGP6EL11	BGP6C1, BGP6C11, BGP6C3
I1AGE^BGP6EL12	BGP6D3, BGP6D4
I1AGE^BGP6EL13	BGP6ELP1
I1AGE^BGP6EL14	BGP6ELP1
I1AGE^BGP6EL15	BGP6ELP1
I1AGE^BGP6EL16	BGP6ELP1
I1AGE^BGP6EL17	BGP6ELP1
I1AGE^BGP6EL18	BGP6ELP1
I1AGE^BGP6EL19	BGP6ELP1
I1AGE^BGP6EL1A	BGP6ELP1
I1AGE^BGP6EL1B	BGP6ELP1
I1AGE^BGP6EL1C	BGP6ELP1
I1AGE^BGP6EL1D	BGP6ELP1
I1AGE^BGP6EL1E	BGP6ELP1
I1AGE^BGP6EL1F	BGP6ELP1
I1AGE^BGP6EL1G	BGP6ELP1
I1AGE^BGP6EL1H	BGP6ELP1
I1AGE^BGP6EL1I	BGP6ELP1
I1AGE^BGP6EL1J	BGP6ELP1
I1AGE^BGP6EL1K	BGP6ELP1
IELDPHA^BGP6EL1L	BGP6ELP1
I1AGE^BGP6EL1M	BGP6ELP1
I1AGE^BGP6EL1N	BGP6ELP1
\$\$FRACTURE^BGP6EL3	BGP6D22
\$\$FUNCTION^BGP6EL4	BGP6EL3
\$\$PHNV^BGP6EL4	BGP6EL3
\$\$TXBMD^BGP6EL4	BGP6D22, BGP6EL3
^BGP6ELH	BGP6DELP
AREACP^BGP6ELH	BGP6PHEL
^BGP6ELHH	BGP6ELL
^BGP6ELL	BGP6DELP, BGP6PHEL
SETN^BGP6ELL1	BGP6ELLN
I1AGE^BGP6ELL2	BGP6ELLN
I1AGE^BGP6ELL3	BGP6ELL1
I1AGE^BGP6ELL4	BGP6ELL1
I1AGE^BGP6ELL5	BGP6ELL1

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 440

Routine	is Invoked by:
I1AGE^BGP6ELL6	BGP6ELL1
I1AGE^BGP6ELL7	BGP6ELL1
I1AGE^BGP6ELL8	BGP6ELL1
I1AGE^BGP6ELL9	BGP6ELL1
I1AGE^BGP6ELLA	BGP6ELL1
I1AGE^BGP6ELLB	BGP6ELL1
I1AGE^BGP6ELLC	BGP6ELL1
I1AGE^BGP6ELLD	BGP6ELL1
I1AGE^BGP6ELLE	BGP6ELL1
I1AGE^BGP6ELLF	BGP6ELL1
I1AGE^BGP6ELLG	BGP6ELL1
I1AGE^BGP6EL LH	BGP6ELL1
I1AGE^BGP6EL LI	BGP6ELL1
I1AGE^BGP6EL LJ	BGP6ELL1
I1AGE^BGP6EL LK	BGP6ELL1
I1AGE^BGP6EL LL	BGP6ELL1
I1AGE^BGP6EL LM	BGP6ELL1
IELDPHA^BGP6EL LN	BGP6ELL1
I1AGE^BGP6EL LO	BGP6ELL1
I1AGE^BGP6EL LP	BGP6ELL1
I1AGEP^BGP6EL LQ	BGP6ELL3
SETN^BGP6EL P1	BGP6EL1L
^BGP6ELS	BGP6DELP
^BGP6ELSL	BGP6ELL
GS^BGP6EUTL	BGP6DEL , BGP6DELA , BGP6GELD
REPORT^BGP6EUTL	BGP6DEL , BGP6GELD
STMP^BGP6EUTL	BGP6DCEI
^BGP6EX	BGP6POST
^BGP6EXA	BGP6EX
^BGP6EXBB	BGP6EXB
^BGP6EXD	BGP6EX
^BGP6EXE	BGP6EX
^BGP6EXG	BGP6EX
^BGP6EXHB	BGP6EXH
^BGP6EXHC	BGP6EXH
^BGP6EXI	BGP6EX
^BGP6EXJ	BGP6EX
^BGP6EXK	BGP6EX
^BGP6EXL	BGP6EX
^BGP6EXM	BGP6EX
^BGP6EXN	BGP6EX
^BGP6EXO	BGP6EX
^BGP6EXP	BGP6EX
^BGP6EXQ	BGP6EX
^BGP6EXS	BGP6EX
^BGP6EXT	BGP6EX
^BGP6EXU	BGP6EX
^BGP6EXV	BGP6EX
^BGP6EXW	BGP6EX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 441

Routine	is Invoked by:
^BGP6EXX	BGP6EX
^BGP6EXY	BGP6EX
^BGP6EXZ	BGP6EX
^BGP6FX	BGPPOST
^BGP6FXA	BGP6FX
^BGP6FXB	BGP6FX
^BGP6FXC	BGP6FX
^BGP6FXD	BGP6FX
^BGP6FXE	BGP6FX
^BGP6FXF	BGP6FX
^BGP6FXG	BGP6FX
^BGP6FXH	BGP6FX
^BGP6FXI	BGP6FX
^BGP6FXJ	BGP6FX
^BGP6FXK	BGP6FX
^BGP6FXL	BGP6FX
^BGP6FXM	BGP6FX
^BGP6FXN	BGP6FX
^BGP6FXO	BGP6FX
^BGP6FXP	BGP6FX
^BGP6FXQ	BGP6FX
^BGP6FXR	BGP6FX
^BGP6FXS	BGP6FX
^BGP6FXT	BGP6FX
^BGP6FXU	BGP6FX
^BGP6FXUB	BGP6FXU
^BGP6FXUC	BGP6FXU
^BGP6FXUD	BGP6FXU
^BGP6FXUE	BGP6FXU
^BGP6FXUF	BGP6FXU
^BGP6FXUG	BGP6FXU
^BGP6FXV	BGP6FX
^BGP6FXW	BGP6FX
^BGP6FXX	BGP6FX
^BGP6FXY	BGP6FX
^BGP6FXZ	BGP6FX
EP^BGP6GACW	BGP6GRA
EP^BGP6GAEL	BGP6GRA
EP^BGP6GAGP	BGP6GRA
EP^BGP6GAHE	BGP6GRA
EP^BGP6GALL	BGP6GR
EP^BGP6GAPU	BGP6GRA
EP^BGP6GCMP	BGP6GR
EP^BGP6GCMS	BGP6GR
EP^BGP6GCOM	BGP6GR
EP^BGP6GELD	BGP6GR1
EP^BGP6GGPU	BGP6GR
EP^BGP6GHED	BGP6GR
EP^BGP6GLST	BGP6GU

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
442

Routine	is Invoked by:
EP^BGP6GLTX	BGP6GU
EP^BGP6GMTX	BGP6GU
EP^BGP6GNPL	BGP6GR
EP^BGP6GNST	BGP6GR
EP^BGP6GNTL	BGP6GR
EP^BGP6GPP	BGP6GR
CATSTR^BGP6GR	BGP6GR1
ELD^BGP6GR1	BGP6GR
\$\$SEVEN^BGP6GU	BGPGU
EP^BGP6GUPL	BGP6GU
^BGP6GX	BGPPOST
^BGP6GXA	BGP6GX
^BGP6GXB	BGP6GX
^BGP6GXC	BGP6GX
^BGP6GXD	BGP6GX
^BGP6GXE	BGP6GX
^BGP6GXG	BGP6GX
^BGP6GXH	BGP6GX
^BGP6GXI	BGP6GX
^BGP6GXJ	BGP6GX
^BGP6GXK	BGP6GX
^BGP6GXL	BGP6GX
^BGP6GXM	BGP6GX
^BGP6GXN	BGP6GX
^BGP6GXO	BGP6GX
^BGP6GXP	BGP6GX
^BGP6GXQ	BGP6GX
^BGP6GXQB	BGP6GXQ
^BGP6GXQC	BGP6GXQ
^BGP6GXQD	BGP6GXQ
^BGP6GXQE	BGP6GXQ
^BGP6GXQF	BGP6GXQ
^BGP6GXQG	BGP6GXQ
^BGP6GXQH	BGP6GXQ
^BGP6GXQI	BGP6GXQ
^BGP6GXR	BGP6GX
^BGP6GXS	BGP6GX
^BGP6GXT	BGP6GX
^BGP6GXU	BGP6GX
^BGP6GXX	BGP6GX
^BGP6GXY	BGP6GX
^BGP6HEH	BGP6DHEP
AREACP^BGP6HEH	BGP6PHEP
^BGP6HEHH	BGP6HEL
^BGP6HEL	BGP6DHEP , BGP6PHEP
^BGP6HES	BGP6DHEP
^BGP6HESL	BGP6HEL
GETIEN^BGP6HUTL	BGP6GUPL , BGP6ULF
GS^BGP6HUTL	BGP6DHE , BGP6GHED

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
443

Routine	is Invoked by:
REPORT^BGP6HUTL	BGP6DHE,BGP6GHED
STMP^BGP6HUTL	BGP6DCHI
^BGP6HX	BGPPOST
^BGP6HXA	BGP6HX
^BGP6HXB	BGP6HX
^BGP6HXC	BGP6HX
^BGP6HXD	BGP6HX
^BGP6HXE	BGP6HX
^BGP6HXF	BGP6HX
^BGP6HXG	BGP6HX
^BGP6HXI	BGP6HX
^BGP6HXJ	BGP6HX
^BGP6H XK	BGP6HX
^BGP6HXL	BGP6HX
^BGP6HXM	BGP6HX
^BGP6H XN	BGP6HX
^BGP6H XO	BGP6HX
^BGP6H XP	BGP6HX
^BGP6H XQ	BGP6HX
^BGP6H XR	BGP6HX
^BGP6H XS	BGP6HX
^BGP6H XT	BGP6HX
^BGP6H XU	BGP6HX
^BGP6H XV	BGP6HX
^BGP6H XW	BGP6HX
^BGP6H XX	BGP6HX
^BGP6H XY	BGP6HX
^BGP6H XZ	BGP6HX
^BGP6JX	BGPPOST
^BGP6JXA	BGP6JX
^BGP6JXB	BGP6JX
^BGP6JXC	BGP6JX
^BGP6JXD	BGP6JX
^BGP6JXE	BGP6JX
^BGP6JXF	BGP6JX
^BGP6JXG	BGP6JX
^BGP6JXH	BGP6JX
^BGP6JXI	BGP6JX
^BGP6JXJ	BGP6JX
^BGP6JXK	BGP6JX
^BGP6JXL	BGP6JX
^BGP6JXM	BGP6JX
^BGP6JXN	BGP6JX
^BGP6JXO	BGP6JX
^BGP6JXP	BGP6JX
^BGP6JXQ	BGP6JX
^BGP6JXR	BGP6JX
^BGP6JXS	BGP6JX
^BGP6JXT	BGP6JX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 444

Routine	is Invoked by:
^BGP6JXU	BGP6JX
^BGP6JXV	BGP6JX
^BGP6JXW	BGP6JX
^BGP6LSTF	opt
^BGP6NPL	opt
NPL1^BGP6NPLD	BGP6PDL
S^BGP6NPLD	BGP6DSTM
NPL1^BGP6NPLP	BGP6DP
PRINT^BGP6PARP	BGP6DAR, BGP6DGPA
^BGP6PDL	BGP6DP, BGP6PARP
\$\$SB^BGP6PDL1	BGP6D23, BGP6D6, BGP6EL3
H1^BGP6PDL1	BGP6PDL0, BGP6PDL3, BGP6PDL4, BGP6PDL7
PI^BGP6PDL1	BGP6PDL0, BGP6PDL3, BGP6PDL4, BGP6PDL7
PI1^BGP6PDL1	BGP6PDL4, BGP6PDL7
I1AGEP^BGP6PDL5	BGP6PDL0
I1AGEP^BGP6PDL8	BGP6PDL7
I1AGE^BGP6PDL9	BGP6PDL1
I1AGE^BGP6PDLA	BGP6PDL1
^BGP6PDLH	BGP6PDL
^BGP6PDLs	BGP6PDL
PRINT^BGP6PHEL	BGP6DELA
PRINT^BGP6PHEP	BGP6DHEA
^BGP6SDP	BGP6DSP
^BGP6ULF	opt
\$\$D^BGP6UTL	BGP6DAR, BGP6DNG, BGP6GAGP, BGP6GNTL
\$\$DATE^BGP6UTL	BGP6ASL, BGP6C1, BGP6C11, BGP6C12, BGP6C13, BGP6C2, BGP6C3 BGP6C31, BGP6D2, BGP6D21, BGP6D22, BGP6D23, BGP6D3, BGP6D31 BGP6D4, BGP6D41, BGP6D42, BGP6D5, BGP6D51, BGP6D52, BGP6D6 BGP6D61, BGP6D7, BGP6D71, BGP6D72, BGP6D73, BGP6D8, BGP6D81 BGP6D82, BGP6D83, BGP6D84
BGP6D86, BGP6D862, BGP6D9, BGP6DGPA, BGP6DHE1, BGP6DP, BGP6EL2	BGP6EL3, BGP6EL4, BGP7D61
\$\$LZERO^BGP6UTL	BGP6DAR, BGP6GAGP
<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>	
\$\$RZERO^BGP6UTL	BGP6DAR, BGP6DNG, BGP6GACW, BGP6GAGP, BGP6GNTL
EXCELGS^BGP6UTL	BGP6DAR, BGP6DNG, BGP6GAGP
GETIEN^BGP6UTL	BGP6GUPL, BGP6ULF
GS^BGP6UTL	BGP6DGPU, BGP6DNG, BGP6GGPU, BGP6GNTL
JRNL^BGP6UTL	BGP6D1, BGP6D10
REPORT^BGP6UTL	
BGP6AHW, BGP6DAP, BGP6DCL, BGP6DGPU, BGP6DL, BGP6DNG, BGP6DPP	BGP6DSTM, BGP6GALL, BGP6GCMF, BGP6GCOM, BGP6GGPU, BGP6GNPL BGP6GNST, BGP6GNTL, BGP6GPP, BGP6NPL
STMP^BGP6UTL	BGP6DCI
XTMP^BGP6UTL	BGP6C1, BGP6D1
\$\$FIRSTPRC^BGP6UTL1	BGP6EL3
\$\$LASTDX^BGP6UTL1	BGP6D23, BGP6D42, BGP6D72, BGP6EL3, BGP6EL4
\$\$LASTDXI^BGP6UTL1	BGP6D42
\$\$LASTPRC^BGP6UTL1	BGP6C12, BGP6C2, BGP6D3, BGP6D61, BGP6D72, BGP6DHE1, BGP6EL4

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 445

Routine	is	Invoked by:
\$LASTPRCI^BGP6UTL1		BGP6D3, BGP6D42, BGP6D61, BGP6D83, BGP6D84
\$RADREF^BGP6UTL1		BGP6D42, BGP6D61
\$REFUSAL^BGP6UTL1		BGP6D21, BGP6D22, BGP6D3, BGP6D31, BGP6D4, BGP6D5, BGP6D6 BGP6D61, BGP6D8, BGP6DHE1, BGP7D61
^BGP6XTCH		opt
TAXCHK^BGP6XTCH		BGP6DAP, BGP6DL, BGP6DPP
^BGP6XTCN		opt
TAXCHK^BGP6XTCN		BGP6AHW, BGP6DCL, BGP6DGPU, BGP6DNG, BGP6DSTM, BGP6NPL
^BGP6XTEL		opt
TAXCHK^BGP6XTEL		BGP6DEL
^BGP6XTHE		opt
TAXCHK^BGP6XTHE		BGP6DHE
EP^BGP6XTL		BGP6XTS
^BGP6XTV		opt
INIT^BGP6XTV		BGP6GU
EP^BGP6XTV1		BGP6CTS, BGP6XTV
EN^BGP7ASL		BGP7DAR, BGP7DCHA, BGP7DELA, BGP7DGPA, BGP7DHEA, BGP7DPEA
GET^BGP7ASL		BGP7GU
^BGP7BAN		opt
^BGP7C1		BGP7CL, BGP7GCMS
ASAALLEG^BGP7C1		BGP7D721
GETMEDS^BGP7C1		BGP7C11, BGP7C12, BGP7C13, BGP7C3
ACEIALG1^BGP7C11		BGP7C2, BGP7D722, BGP7D732
ACEICON1^BGP7C11		BGP7C2
ACEIRX1^BGP7C11		BGP7C2
EN^BGP7C11		BGP7C1
^BGP7C12		BGP7C11
CESS1^BGP7C12		BGP7C2, BGP7C3
EN^BGP7C12		BGP7C2, BGP7C3
SMOKER^BGP7C12		BGP7C1
SMOKER1^BGP7C12		BGP7C2, BGP7C3
ARBALG1^BGP7C13		BGP7C11
ARBCON1^BGP7C13		BGP7C11
ARBRX1^BGP7C13		BGP7C11
ABGPO1^BGP7C31		BGP7C3
^BGP7CL		opt
^BGP7CP		BGP7CL
EP^BGP7CSI		BGP7CL
EP^BGP7CTL		BGP7CTS
^BGP7CTS		opt
INIT^BGP7CTS		BGP7GU
^BGP7CTXC		opt
TAXCHK^BGP7CTXC		BGP7CL
^BGP7D1		
BGP7DAP, BGP7DCL, BGP7DEL, BGP7DGPU, BGP7DHE, BGP7DL, BGP7DNG		BGP7DSTM, BGP7GALL, BGP7GCMP, BGP7GCOM, BGP7GELD, BGP7GGPU BGP7GHED, BGP7GNPL, BGP7GNST, BGP7GNTL, BGP7NPL
\$\$ACTUP^BGP7D1		BGP7DNE

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>

\$\$V2^BGP7D1		BGP7D721
---------------	--	----------

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 446

Routine is Invoked by:

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>

PROCCY^BGP7D1 BGP7DNE
 ^BGP7D10 BGP7DPP, BGP7GPP
 \$\$HGBA1C^BGP7D2 BGP7D22, BGP7DHE1
 \$\$LDL^BGP7D2
 BGP7D22, BGP7D23, BGP7D41, BGP7D72, BGP7D73, BGP7DHE1, BGP7EL4

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>

\$\$LOINC^BGP7D2 BGP7D722
 \$\$MEANBP^BGP7D2 BGP7D22
 \$\$DENTSRV^BGP7D21 BGP7DHE1, BGP7EL2
 \$\$DEP^BGP7D21 BGP7D23, BGP7D41, BGP7EL3
 \$\$DEPEDU^BGP7D21 BGP7EL3
 \$\$DEPREF^BGP7D21 BGP7D23, BGP7EL3
 \$\$DEPSCR^BGP7D21 BGP7D23, BGP7D41, BGP7EL3
 \$\$DNKA^BGP7D21 BGP7D22, BGP7DHE1
 \$\$EYE^BGP7D21 BGP7D22, BGP7EL2
 \$\$LOINC^BGP7D21 BGP7C3, BGP7C31, BGP7D8, BGP7D82, BGP7D83, BGP7D84, BGP7D861
 BGP7D862, BGP7D864
 \$\$ESRD^BGP7D211 BGP7D21, BGP7D22, BGP7D23, BGP7D9, BGP7EL2
 \$\$GFR^BGP7D211 BGP7D21, BGP7D22, BGP7D23, BGP7EL2
 \$\$QUANTUP^BGP7D211 BGP7D21, BGP7D22, BGP7D23, BGP7EL2
 \$\$CREAT^BGP7D22 BGP7D82, BGP7D86
 \$\$FLU^BGP7D3 BGP7EL3
 \$\$HOSP^BGP7D31 BGP7EL31
 \$\$LAST^BGP7D31 BGP7EL31
 \$\$PNEU^BGP7D31 BGP7D3, BGP7EL3
 \$\$V2ASTH^BGP7D31 BGP7EL31
 \$\$IMMREF^BGP7D32 BGP7D33, BGP7D34, BGP7D35, BGP7D85
 GETIMMS^BGP7D32 BGP7D33, BGP7D34, BGP7D35, BGP7D85
 \$\$MMR^BGP7D33 BGP7D32
 \$\$OPV^BGP7D33 BGP7D32
 \$\$DTAP^BGP7D34 BGP7D32
 \$\$HEP^BGP7D35 BGP7D32, BGP7D85
 \$\$HIB^BGP7D35 BGP7D32
 \$\$PNEUMO^BGP7D35 BGP7D32
 \$\$VAR^BGP7D35 BGP7D32, BGP7D85
 \$\$MAM^BGP7D4 BGP7EL3
 \$\$MAS^BGP7D4 BGP7EL3
 \$\$LIFE^BGP7D41 BGP7D23
 \$\$MEANBP^BGP7D41 BGP7D9, BGP7EL3
 \$\$OSTEOSCR^BGP7D42 BGP7EL3
 \$\$CHL^BGP7D5 BGP7DHE1
 \$\$DVDX^BGP7D5 BGP7EL3
 \$\$DVEX^BGP7D5 BGP7EL3
 \$\$DVPED^BGP7D5 BGP7EL3
 \$\$REFDV^BGP7D5 BGP7EL3
 \$\$EAPT^BGP7D52 BGP7D51
 \$\$BMI^BGP7D6 BGP7D23, BGP7D41, BGP7D71, BGP7EL3

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 447

Routine	is Invoked by:
\$\$OB^BGP7D6	BGP7D71, BGP7EL3
\$\$OW^BGP7D6	BGP7D71, BGP7EL3
\$\$REF^BGP7D6	BGP7D23, BGP7EL3
\$\$CRC^BGP7D61	BGP7DHE1, BGP7EL3
\$\$FOB^BGP7D61	BGP7DHE1, BGP7EL3
\$\$IBH^BGP7D61	BGP7DHE1, BGP7EL3
\$\$RECT^BGP7D61	BGP7DHE1
\$\$REF^BGP7D61	BGP7DHE1, BGP7EL3
\$\$DENT^BGP7D7	
BGP7C12, BGP7C2, BGP7D23, BGP7D41, BGP7D71, BGP7DHE1, BGP7EL3	
\$\$DX^BGP7D7	BGP7C12, BGP7C2, BGP7D23, BGP7D71, BGP7DHE1, BGP7EL3
\$\$PED^BGP7D7	BGP7D23, BGP7D41, BGP7EL3
\$\$PREG^BGP7D7	BGP7D722, BGP7D732, BGP7D8
\$\$TOBACCO^BGP7D7	
BGP7C12, BGP7C2, BGP7D23, BGP7D41, BGP7D71, BGP7DHE1, BGP7EL3	
\$\$MEDNUTR^BGP7D71	BGP7D41
\$\$OTHREL^BGP7D71	BGP7D41
\$\$SPECEX^BGP7D71	BGP7D41
\$\$SPECNUTR^BGP7D71	BGP7D41
\$\$AMA^BGP7D72	BGP7D731
\$\$BETAALG1^BGP7D72	BGP7D721, BGP7D731, BGP7D75
<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>	
\$\$EXPIRED^BGP7D72	BGP7D731
<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>	
\$\$TRANS^BGP7D72	BGP7D731
<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>	
\$\$ASA^BGP7D721	BGP7D731
\$\$ASAALLEG^BGP7D721	BGP7D731, BGP7D75
<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>	
\$\$ASAREF^BGP7D721	BGP7D731, BGP7D75
\$\$BETA^BGP7D721	BGP7D731
\$\$BETACONT^BGP7D721	BGP7D731, BGP7D75
\$\$BETAREF^BGP7D721	BGP7D731, BGP7D75
\$\$IHD^BGP7D721	BGP7D21, BGP7D41, BGP7D73
\$\$ACEALG^BGP7D722	BGP7D721, BGP7D731, BGP7D75
\$\$ACECONT^BGP7D722	BGP7D721, BGP7D731, BGP7D75
<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>	
\$\$ACEREF^BGP7D722	BGP7D721, BGP7D731, BGP7D75
\$\$ACERX^BGP7D722	BGP7D721, BGP7D731
\$\$STATALG^BGP7D722	BGP7D721, BGP7D731, BGP7D75
\$\$STATCON^BGP7D722	BGP7D721, BGP7D731, BGP7D75
<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>	
\$\$STATREF^BGP7D722	BGP7D721, BGP7D731, BGP7D75
\$\$STATRX^BGP7D722	BGP7D721, BGP7D731
\$\$CHOL^BGP7D73	BGP7EL4

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 448

Routine	is	Invoked by:
\$\$DAYS^BGP7D82		BGP7D52, BGP7D721, BGP7D722, BGP7D732, BGP7D863
IRAA^BGP7D82		BGP7EL3
\$\$CBC^BGP7D861		BGP7D86
\$\$CBC4^BGP7D861		BGP7D86
\$\$CBC6^BGP7D861		BGP7D86
\$\$GLUCOSE^BGP7D861		BGP7D86
\$\$LFT^BGP7D861		BGP7D86
\$\$LFT6^BGP7D861		BGP7D86
\$\$POT^BGP7D861		BGP7D86
\$\$SERUM12^BGP7D861		BGP7D86
\$\$SERUM6^BGP7D861		BGP7D86
\$\$UP4^BGP7D861		BGP7D86
\$\$RHEUAR^BGP7D862		BGP7D86
\$\$MEDSPRE^BGP7D863		BGP7D86
\$\$GOLDLAB^BGP7D864		BGP7D86
\$\$FIRSTIHD^BGP7D9		BGP7DHE1
\$\$V2IHD^BGP7D9		BGP7DHE1
^BGP7DAP		opt
^BGP7DAR		opt
CALCIND^BGP7DCEI		BGP7D1
^BGP7DCHA		opt
\$\$D^BGP7DCHA		BGP7GACW
GS^BGP7DCHA		BGP7GACW
CALCIND^BGP7DCHI		BGP7D1
\$\$STATE^BGP7DCHW		BGP7DNE
\$\$UID^BGP7DCHW		BGP7DNE
GEN^BGP7DCHW		BGP7D1
CALCIND^BGP7DCI		BGP7D1, BGP7D10
^BGP7DCL		opt
CPPL1^BGP7DCLD		BGP7PDL
\$\$LAST^BGP7DCLP		BGP7DCLD
CPPL1^BGP7DCLP		BGP7DP
^BGP7DEL		opt
^BGP7DELA		opt
^BGP7DELP		BGP7DEL
PRINT1^BGP7DELP		BGP7PHEL
EN^BGP7DESI		BGP7DEL
EN^BGP7DESL		BGP7DEL
PT^BGP7DESL		BGP7DEL, BGP7DELA
RT^BGP7DESL		BGP7DEL
^BGP7DGPA		opt
^BGP7DGPU		opt

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>

ENDDATE^BGP7DGPU	BGP7DAP, BGP7DL, BGP7DPP
^BGP7DH	BGP7DP, BGP7DPEP
AREACP^BGP7DH	BGP7PARP, BGP7PPEL
EOP^BGP7DH	BGP7GUPL, BGP7LSTF, BGP7ULF
^BGP7DHE	opt

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 449

Routine	is Invoked by:
\$\$POSUR^BGP7DHE2	BGP7DHE1
^BGP7DHEA	opt
^BGP7DHEP	BGP7DHE
PRINT1^BGP7DHEP	BGP7PHEP
EN^BGP7DHSL	BGP7DHE
PT^BGP7DHSL	BGP7DHE, BGP7DHEA
RT^BGP7DHSL	BGP7DHE
^BGP7DL	opt
^BGP7DLT	opt
^BGP7DMT	opt
^BGP7DNE	opt
PRINT^BGP7DNE	BGP7DNE
PROC^BGP7DNE	BGP7DNE, BGP7GNXP
^BGP7DNG	opt
HWSF^BGP7DNG	BGP7GNTL
^BGP7DP	
BGP7DAP, BGP7DCL, BGP7DGPU, BGP7DL, BGP7DNG, BGP7DPP, BGP7DSTM	
	BGP7NPL
PRINT1^BGP7DP	BGP7PARP
SETEXCEL^BGP7DP	BGP7PARP
I1AGE^BGP7DP11	BGP7DP1C
I1AGE^BGP7DP12	BGP7DP1
I1AGE^BGP7DP19	BGP7DP1
\$\$V^BGP7DP1C	BGP7DP11, BGP7DP12, BGP7DP19, BGP7DP1C, BGP7DP2, BGP7DP3 BGP7DP4, BGP7DP5, BGP7DP50, BGP7DP51, BGP7DP6, BGP7DP60 BGP7DP7, BGP7PD19, BGP7PDL0, BGP7PDL1, BGP7PDL2, BGP7PDL3 BGP7PDL4, BGP7PDL5, BGP7PDL7, BGP7PDL8, BGP7PDL9, BGP7PDLA
I1AGE^BGP7DP1C	BGP7DP1
PI^BGP7DP1C	BGP7DP1, BGP7DP3, BGP7DP4, BGP7DP5, BGP7DP6
PI1^BGP7DP1C	BGP7DP4, BGP7DP5
SETEXA^BGP7DP1C	BGP7DP2, BGP7DP3, BGP7DP7, BGP7PDL1, BGP7PDL2, BGP7PDL3
SETN^BGP7DP1C	BGP7DP2, BGP7DP3, BGP7DP7, BGP7PDL1, BGP7PDL2, BGP7PDL3
I1AGEP^BGP7DP50	BGP7DP5, BGP7DP51
I1AGEP^BGP7DP60	BGP7DP6
^BGP7DPE	opt
^BGP7DPE1	BGP7DPE, BGP7GPED
5^BGP7DPE2	BGP7DPE1
6^BGP7DPE2	BGP7DPE1
^BGP7DPEA	opt
^BGP7DPED	BGP7PPEL
<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>	
DEL^BGP7DPED	BGP7DPEP
^BGP7DPEE	BGP7DPED
^BGP7DPEL	BGP7DPEP
^BGP7DPEP	BGP7DPE
PRINT1^BGP7DPEP	BGP7PPEL
H1^BGP7DPH	BGP7DP1, BGP7DP1C, BGP7DP4, BGP7DP5, BGP7DP6, BGP7DPEP BGP7ELP1, BGP7HEP1
H2^BGP7DPH	BGP7DP1C, BGP7DPEP, BGP7ELP1, BGP7HEP1

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 450

Routine	is Invoked by:
H3^BGP7DPH	BGP7DP11
H4^BGP7DPH	BGP7DP12
H9^BGP7DPH	BGP7DP19
HEADER^BGP7DPH	BGP7DELP, BGP7DHEP, BGP7DP, BGP7DP1, BGP7DP11, BGP7DP12 BGP7DP19, BGP7DP1C, BGP7DP3, BGP7DP4, BGP7DP5, BGP7DP50 BGP7DP51, BGP7DP6, BGP7DP60, BGP7DSP, BGP7DSPN, BGP7DSTM BGP7EL10, BGP7EL11, BGP7EL12, BGP7EL13, BGP7EL14, BGP7EL15 BGP7EL16, BGP7EL17, BGP7EL18 BGP7EL19, BGP7EL1A, BGP7EL1B, BGP7EL1C, BGP7EL1D, BGP7EL1E BGP7EL1F, BGP7EL1G, BGP7EL1H, BGP7EL1I, BGP7EL1J, BGP7EL1K BGP7EL1L, BGP7EL1M, BGP7EL1N, BGP7EL1O, BGP7EL1P, BGP7EL1Q BGP7EL1R, BGP7ELP1, BGP7HEP1, BGP7SDP, BGP7SDPN
^BGP7DPP	opt
^BGP7DS	BGP7DP
EN^BGP7DSI	BGP7DAP, BGP7DL, BGP7DPP, BGP7DSTM, BGP7NPL
EN^BGP7DSL	BGP7DAP, BGP7DL, BGP7DPP
PT^BGP7DSL	
BGP7DAP, BGP7DAR, BGP7DCL, BGP7DGPA, BGP7DGPU, BGP7DL, BGP7DNG	BGP7DPP, BGP7DSTM, BGP7NPL
RT^BGP7DSL	BGP7DAP, BGP7DCL, BGP7DL, BGP7DPP, BGP7DSTM, BGP7NPL
^BGP7DSP	BGP7DP
^BGP7DSPN	BGP7DP
^BGP7DSTM	opt
CT^BGP7DSTM	BGP7DP, BGP7PDL
\$\$CPT^BGP7DU	
BGP7C12, BGP7D21, BGP7D211, BGP7D21A, BGP7D22, BGP7D3, BGP7D4	BGP7D42, BGP7D5, BGP7D61, BGP7D7, BGP7D72, BGP7D8, BGP7D83 BGP7D84, BGP7D861, BGP7DHE1, BGP7EL4
\$\$CPTI^BGP7DU	BGP7C1, BGP7C11, BGP7C12, BGP7C13, BGP7C3, BGP7D2, BGP7D211 BGP7D21A, BGP7D23, BGP7D31, BGP7D61, BGP7D71, BGP7D72 BGP7D721, BGP7D73, BGP7D731, BGP7D732, BGP7D75, BGP7D8 BGP7D82, BGP7D861, BGP7D864, BGP7D9, BGP7DHE2, BGP7EL2
\$\$PLCODE^BGP7DU	BGP7D34
\$\$PLTAX^BGP7DU	BGP7D33, BGP7D34, BGP7D35, BGP7D85
\$\$RAD^BGP7DU	BGP7D42, BGP7D61, BGP7EL4
\$\$TRAN^BGP7DU	
BGP7C12, BGP7D21, BGP7D211, BGP7D21A, BGP7D22, BGP7D3, BGP7D4	BGP7D42, BGP7D5, BGP7D61, BGP7D7, BGP7D72, BGP7D8, BGP7D83 BGP7D84, BGP7D861, BGP7DHE1, BGP7EL4
\$\$TRANI^BGP7DU	BGP7C1, BGP7C11, BGP7C12, BGP7C13, BGP7C3, BGP7D2, BGP7D21 BGP7D211, BGP7D21A, BGP7D23, BGP7D31, BGP7D61, BGP7D71 BGP7D72, BGP7D721, BGP7D73, BGP7D731, BGP7D732, BGP7D75
BGP7D8, BGP7D82, BGP7D861, BGP7D864, BGP7D9, BGP7DHE2, BGP7EL2	
\$\$WH^BGP7DU	BGP7D3, BGP7D4
I1AGE^BGP7EL10	BGP7ELP1
I1AGE^BGP7EL11	BGP7ELP1
I1AGE^BGP7EL12	BGP7ELP1
I1AGE^BGP7EL13	BGP7ELP1
I1AGE^BGP7EL14	BGP7ELP1
I1AGE^BGP7EL15	BGP7ELP1
I1AGE^BGP7EL16	BGP7ELP1

I1AGE^BGP7EL17

BGP7ELP1

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 451

Routine	is Invoked by:
I1AGE^BGP7EL18	BGP7ELP1
I1AGE^BGP7EL19	BGP7ELP1
I1AGE^BGP7EL1A	BGP7ELP1
I1AGE^BGP7EL1B	BGP7ELP1
I1AGE^BGP7EL1C	BGP7ELP1
I1AGE^BGP7EL1D	BGP7ELP1
I1AGE^BGP7EL1E	BGP7ELP1
I1AGE^BGP7EL1F	BGP7ELP1
I1AGE^BGP7EL1G	BGP7ELP1
I1AGE^BGP7EL1H	BGP7ELP1
I1AGE^BGP7EL1I	BGP7ELP1
I1AGE^BGP7EL1J	BGP7ELP1
I1AGE^BGP7EL1K	BGP7ELP1
IELDPHA^BGP7EL1L	BGP7ELP1
I1AGE^BGP7EL1M	BGP7ELP1
I1AGE^BGP7EL1N	BGP7ELP1
I1AGE^BGP7EL1O	BGP7ELP1
H3^BGP7EL1P	BGP7EL1Q, BGP7EL1R
I1AGE^BGP7EL1P	BGP7ELP1
^BGP7EL1Q	BGP7EL1P
^BGP7EL1R	BGP7EL1Q
\$\$FRACTURE^BGP7EL3	BGP7D22
IELDASA^BGP7EL31	BGP7EL3
IELDFSA^BGP7EL31	BGP7EL3
IELDPHA^BGP7EL31	BGP7EL3
\$\$FUNCTION^BGP7EL4	BGP7EL31
\$\$PHNV^BGP7EL4	BGP7EL31
\$\$TXBMD^BGP7EL4	BGP7D22, BGP7EL3
^BGP7ELH	BGP7DELP
AREACP^BGP7ELH	BGP7PHEL
^BGP7ELHH	BGP7ELL
^BGP7ELL	BGP7DELP, BGP7PHEL
SETN^BGP7ELL1	BGP7ELLN
I1AGE^BGP7ELL2	BGP7ELL1
I1AGE^BGP7ELL3	BGP7ELL1
I1AGE^BGP7ELL4	BGP7ELL1
I1AGE^BGP7ELL5	BGP7ELL1
I1AGE^BGP7ELL6	BGP7ELL1
I1AGE^BGP7ELL7	BGP7ELL1
I1AGE^BGP7ELL8	BGP7ELL1
I1AGE^BGP7ELL9	BGP7ELL1
I1AGE^BGP7ELLA	BGP7ELL1
I1AGE^BGP7ELLB	BGP7ELL1
I1AGE^BGP7ELLC	BGP7ELL1
I1AGE^BGP7ELLD	BGP7ELL1
I1AGE^BGP7ELLE	BGP7ELL1
I1AGE^BGP7ELLF	BGP7ELL1
I1AGE^BGP7ELLG	BGP7ELL1
I1AGE^BGP7ELLH	BGP7ELL1

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
452

Routine	is Invoked by:
I1AGE^BGP7ELLI	BGP7ELL1
I1AGE^BGP7ELLJ	BGP7ELL1
I1AGE^BGP7ELLK	BGP7ELL1
I1AGE^BGP7ELLL	BGP7ELL1
I1AGE^BGP7ELLM	BGP7ELL1
IELDPHA^BGP7ELLN	BGP7ELL1
I1AGE^BGP7ELLO	BGP7ELL1
I1AGE^BGP7ELLP	BGP7ELL1
I1AGEP^BGP7ELLQ	BGP7ELL3
I1AGE^BGP7ELLR	BGP7ELL1
I1AGE^BGP7ELLS	BGP7ELL1
I1AGEP^BGP7ELLT	BGP7ELLS
^BGP7ELLU	BGP7ELLT
BL^BGP7ELLV	BGP7ELLJ
SETN^BGP7ELP1	BGP7EL1L
^BGP7ELS	BGP7DELP
^BGP7ELSL	BGP7ELL
GS^BGP7EUTL	BGP7DEL , BGP7GELD
REPORT^BGP7EUTL	BGP7DEL , BGP7GELD
STMP^BGP7EUTL	BGP7DCEI
EP^BGP7GACW	BGP7GRA
EP^BGP7GAEL	BGP7GRA
EP^BGP7GAGP	BGP7GRA
EP^BGP7GAHE	BGP7GRA
EP^BGP7GALL	BGP7GR
EP^BGP7GAPE	BGP7GRA
EP^BGP7GAPU	BGP7GRA
EP^BGP7GCMP	BGP7GR
EP^BGP7GCMS	BGP7GR
EP^BGP7GCOM	BGP7GR
EP^BGP7GELD	BGP7GR1
EP^BGP7GGPU	BGP7GR
EP^BGP7GHED	BGP7GR
EP^BGP7GLST	BGP7GU
EP^BGP7GLTX	BGP7GU
EP^BGP7GMTX	BGP7GU
EP^BGP7GNPL	BGP7GR
EP^BGP7GNST	BGP7GR
EP^BGP7GNTL	BGP7GR
EP^BGP7GNXP	BGP7GR2
EP^BGP7GPED	BGP7GR2
EP^BGP7GPP	BGP7GR
CATSTR^BGP7GR	BGP7GR1
ELD^BGP7GR1	BGP7GR
EP^BGP7GUPL	BGP7GU
^BGP7GX	BGP7POS
^BGP7GXA	BGP7GX
^BGP7GXB	BGP7GX
^BGP7GXC	BGP7GX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
453

Routine	is Invoked by:
^BGP7GXD	BGP7GX
^BGP7GXE	BGP7GX
^BGP7GXF	BGP7GX
^BGP7GXG	BGP7GX
^BGP7GXH	BGP7GX
^BGP7GXI	BGP7GX
^BGP7GXJ	BGP7GX
^BGP7GXK	BGP7GX
^BGP7GXL	BGP7GX
^BGP7GXM	BGP7GX
^BGP7GXMB	BGP7GXM
^BGP7GXMC	BGP7GXM
^BGP7GXMD	BGP7GXM
^BGP7GXME	BGP7GXM
^BGP7GXMF	BGP7GXM
^BGP7GXN	BGP7GX
^BGP7GXNB	BGP7GXN
^BGP7GXNC	BGP7GXN
^BGP7GXND	BGP7GXN
^BGP7GXO	BGP7GX
^BGP7GXOB	BGP7GXO
^BGP7GXOC	BGP7GXO
^BGP7GXOD	BGP7GXO
^BGP7GXOE	BGP7GXO
^BGP7GXOF	BGP7GXO
^BGP7GXOG	BGP7GXO
^BGP7GXOH	BGP7GXO
^BGP7GXOI	BGP7GXO
^BGP7GXOJ	BGP7GXO
^BGP7GXP	BGP7GX
^BGP7GXQ	BGP7GX
^BGP7GXQB	BGP7GXQ
^BGP7GXQC	BGP7GXQ
^BGP7GXQD	BGP7GXQ
^BGP7GXQE	BGP7GXQ
^BGP7HEH	BGP7DHEP
AREACP^BGP7HEH	BGP7PHEP
^BGP7HEHH	BGP7HEL
^BGP7HEL	BGP7DHEP , BGP7PHEP
^BGP7HES	BGP7DHEP
^BGP7HESL	BGP7HEL
GETIEN^BGP7HUTL	BGP7GUPL , BGP7ULF
GS^BGP7HUTL	BGP7DHE , BGP7GHED
REPORT^BGP7HUTL	BGP7DHE , BGP7GHED
STMP^BGP7HUTL	BGP7DCHI
^BGP7HX	BGP7POS
^BGP7HXA	BGP7HX
^BGP7HXB	BGP7HX
^BGP7HXC	BGP7HX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
454

Routine	is Invoked by:
^BGP7HXD	BGP7HX
^BGP7HXE	BGP7HX
^BGP7HXF	BGP7HX
^BGP7HXG	BGP7HX
^BGP7HXH	BGP7HX
^BGP7HXI	BGP7HX
^BGP7HXJ	BGP7HX
^BGP7HXX	BGP7HX
^BGP7HXL	BGP7HX
^BGP7HXM	BGP7HX
^BGP7HXN	BGP7HX
^BGP7HXO	BGP7HX
^BGP7HXP	BGP7HX
^BGP7HXQ	BGP7HX
^BGP7HXR	BGP7HX
^BGP7HXS	BGP7HX
^BGP7HXT	BGP7HX
^BGP7HXU	BGP7HX
^BGP7HXV	BGP7HX
^BGP7HXW	BGP7HX
^BGP7HXX	BGP7HX
^BGP7HXY	BGP7HX
^BGP7HXZ	BGP7HX
^BGP7IX	BGP7POS
^BGP7IXA	BGP7IX
^BGP7IXB	BGP7IX
^BGP7IXC	BGP7IX
^BGP7IXCB	BGP7IXC
^BGP7IXD	BGP7IX
^BGP7IXE	BGP7IX
^BGP7IXF	BGP7IX
^BGP7IXG	BGP7IX
^BGP7IXH	BGP7IX
^BGP7IXI	BGP7IX
^BGP7IXJ	BGP7IX
^BGP7IXK	BGP7IX
^BGP7IXLB	BGP7IXL
^BGP7IXLC	BGP7IXL
^BGP7IXLD	BGP7IXL
^BGP7IXLE	BGP7IXL
^BGP7IXLF	BGP7IXL
^BGP7IXLG	BGP7IXL
^BGP7IXLH	BGP7IXL
^BGP7IXM	BGP7IX
^BGP7IXN	BGP7IX
^BGP7IXO	BGP7IX
^BGP7IXP	BGP7IX
^BGP7IXQ	BGP7IX
^BGP7IXR	BGP7IX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
455

Routine	is Invoked by:
^BGP7IXS	BGP7IX
^BGP7IXSB	BGP7IXS
^BGP7IXSC	BGP7IXS
^BGP7IXSD	BGP7IXS
^BGP7IXSE	BGP7IXS
^BGP7IXSF	BGP7IXS
^BGP7IXSG	BGP7IXS
^BGP7IXSH	BGP7IXS
^BGP7IXSI	BGP7IXS
^BGP7IXSJ	BGP7IXS
^BGP7IXSK	BGP7IXS
^BGP7IXSL	BGP7IXS
^BGP7IXSM	BGP7IXS
^BGP7IXSN	BGP7IXS
^BGP7IXSO	BGP7IXS
^BGP7IXSP	BGP7IXS
^BGP7IXSQ	BGP7IXS
^BGP7IXSR	BGP7IXS
^BGP7IXSS	BGP7IXS
^BGP7IXST	BGP7IXS
^BGP7IXT	BGP7IX
^BGP7IXTB	BGP7IXT
^BGP7IXTC	BGP7IXT
^BGP7IXU	BGP7IX
^BGP7IXUB	BGP7IXU
^BGP7IXUC	BGP7IXU
^BGP7IXUD	BGP7IXU
^BGP7IXUE	BGP7IXU
^BGP7IXV	BGP7IX
^BGP7IXW	BGP7IX
^BGP7IXX	BGP7IX
^BGP7IXY	BGP7IX
^BGP7IXZ	BGP7IX
^BGP7JX	BGP7POS
^BGP7JXA	BGP7JX
^BGP7JXB	BGP7JX
^BGP7JXC	BGP7JX
^BGP7JXD	BGP7JX
^BGP7JXE	BGP7JX
^BGP7JXF	BGP7JX
^BGP7JXFB	BGP7JXF
^BGP7JXFC	BGP7JXF
^BGP7JXFD	BGP7JXF
^BGP7JXG	BGP7JX
^BGP7JXH	BGP7JX
^BGP7JXI	BGP7JX
^BGP7JXJ	BGP7JX
^BGP7JXK	BGP7JX
^BGP7JXL	BGP7JX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
456

Routine	is Invoked by:
^BGP7JXM	BGP7JX
^BGP7JXN	BGP7JX
^BGP7JXO	BGP7JX
^BGP7JXP	BGP7JX
^BGP7JXQ	BGP7JX
^BGP7JXQB	BGP7JXQ
^BGP7JXQC	BGP7JXQ
^BGP7JXQD	BGP7JXQ
^BGP7JXQE	BGP7JXQ
^BGP7JXQF	BGP7JXQ
^BGP7JXQG	BGP7JXQ
^BGP7JXR	BGP7JX
^BGP7JXS	BGP7JX
^BGP7JXT	BGP7JX
^BGP7JXU	BGP7JX
^BGP7JXV	BGP7JX
^BGP7JXW	BGP7JX
^BGP7JXX	BGP7JX
^BGP7JXY	BGP7JX
^BGP7JXZ	BGP7JX
^BGP7KX	BGP7POS
^BGP7KXA	BGP7KX
^BGP7KXB	BGP7KX
^BGP7KXC	BGP7KX
^BGP7KXD	BGP7KX
^BGP7KXEB	BGP7KXE
^BGP7KXEC	BGP7KXE
^BGP7KXED	BGP7KXE
^BGP7KXEE	BGP7KXE
^BGP7KXEF	BGP7KXE
^BGP7KXF	BGP7KX
^BGP7KXG	BGP7KX
^BGP7KXGB	BGP7KXG
^BGP7KXGC	BGP7KXG
^BGP7KXGD	BGP7KXG
^BGP7KXGE	BGP7KXG
^BGP7KXGF	BGP7KXG
^BGP7KXGG	BGP7KXG
^BGP7KXGH	BGP7KXG
^BGP7KXGI	BGP7KXG
^BGP7KXH	BGP7KX
^BGP7KXIB	BGP7KXI
^BGP7KXIC	BGP7KXI
^BGP7KXID	BGP7KXI
^BGP7KXJ	BGP7KX
^BGP7KXK	BGP7KX
^BGP7KXL	BGP7KX
^BGP7KXN	BGP7KX
^BGP7KXO	BGP7KX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
457

Routine	is Invoked by:
^BGP7KXP	BGP7KX
^BGP7KXQ	BGP7KX
^BGP7KXR	BGP7KX
^BGP7KXS	BGP7KX
^BGP7KXT	BGP7KX
^BGP7KXU	BGP7KX
^BGP7KXV	BGP7KX
^BGP7KXW	BGP7KX
^BGP7KXX	BGP7KX
^BGP7KXY	BGP7KX
^BGP7KXZ	BGP7KX
^BGP7LSTF	opt
^BGP7LX	BGP7POS
^BGP7LXA	BGP7LX
^BGP7LXAB	BGP7LXA
^BGP7LXAC	BGP7LXA
^BGP7LXAD	BGP7LXA
^BGP7LXAE	BGP7LXA
^BGP7LXAF	BGP7LXA
^BGP7LXAG	BGP7LXA
^BGP7LXAH	BGP7LXA
^BGP7LXAI	BGP7LXA
^BGP7LXAJ	BGP7LXA
^BGP7LXAK	BGP7LXA
^BGP7LXAL	BGP7LXA
^BGP7LXAM	BGP7LXA
^BGP7LXAN	BGP7LXA
^BGP7LXB	BGP7LX
^BGP7LXBB	BGP7LXB
^BGP7LXC	BGP7LX
^BGP7LXCB	BGP7LXC
^BGP7LXD	BGP7LX
^BGP7LXDB	BGP7LXD
^BGP7LXDC	BGP7LXD
^BGP7LXE	BGP7LX
^BGP7LXEB	BGP7LXE
^BGP7LXEC	BGP7LXE
^BGP7LXED	BGP7LXE
^BGP7LXEE	BGP7LXE
^BGP7LXEF	BGP7LXE
^BGP7LXF	BGP7LX
^BGP7LXFB	BGP7LXF
^BGP7LXFC	BGP7LXF
^BGP7LXFD	BGP7LXF
^BGP7LXFE	BGP7LXF
^BGP7LXFF	BGP7LXF
^BGP7LXFG	BGP7LXF
^BGP7LXFH	BGP7LXF
^BGP7LXFI	BGP7LXF

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
458

Routine	is Invoked by:
^BGP7LXFJ	BGP7LXF
^BGP7LXG	BGP7LX
^BGP7LXGB	BGP7LXG
^BGP7LXGC	BGP7LXG
^BGP7LXGD	BGP7LXG
^BGP7LXGE	BGP7LXG
^BGP7LXGF	BGP7LXG
^BGP7LXGG	BGP7LXG
^BGP7LXGH	BGP7LXG
^BGP7LXGI	BGP7LXG
^BGP7LXGJ	BGP7LXG
^BGP7LXGK	BGP7LXG
^BGP7LXGL	BGP7LXG
^BGP7LXGM	BGP7LXG
^BGP7LXGN	BGP7LXG
^BGP7LXGO	BGP7LXG
^BGP7LXH	BGP7LX
^BGP7LXHB	BGP7LXH
^BGP7LXI	BGP7LX
^BGP7LXIB	BGP7LXI
^BGP7LXIC	BGP7LXI
^BGP7LXID	BGP7LXI
^BGP7LXJ	BGP7LX
^BGP7LXJB	BGP7LXJ
^BGP7LXJC	BGP7LXJ
^BGP7LXJD	BGP7LXJ
^BGP7LXK	BGP7LX
^BGP7LXKB	BGP7LXK
^BGP7LXKC	BGP7LXK
^BGP7LXKD	BGP7LXK
^BGP7LXKE	BGP7LXK
^BGP7LXKF	BGP7LXK
^BGP7LXKG	BGP7LXK
^BGP7LXKH	BGP7LXK
^BGP7LXKI	BGP7LXK
^BGP7LXKJ	BGP7LXK
^BGP7LXKK	BGP7LXK
^BGP7LXL	BGP7LX
^BGP7LXLB	BGP7LXL
^BGP7LXLC	BGP7LXL
^BGP7LXLD	BGP7LXL
^BGP7LXLE	BGP7LXL
^BGP7LXLF	BGP7LXL
^BGP7LXLG	BGP7LXL
^BGP7LXLH	BGP7LXL
^BGP7LXLI	BGP7LXL
^BGP7LXLJ	BGP7LXL
^BGP7LXLK	BGP7LXL
^BGP7LXLL	BGP7LXL

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
459

Routine	is Invoked by:
^BGP7LXLM	BGP7LXL
^BGP7LXLN	BGP7LXL
^BGP7LXLLO	BGP7LXL
^BGP7LXLQP	BGP7LXL
^BGP7LXLQ	BGP7LXL
^BGP7LXLRL	BGP7LXL
^BGP7LXLRS	BGP7LXL
^BGP7LXLRT	BGP7LXL
^BGP7LXLRLU	BGP7LXL
^BGP7LXLRLV	BGP7LXL
^BGP7LXLRLW	BGP7LXL
^BGP7LXLM	BGP7LX
^BGP7LXLN	BGP7LX
^BGP7LXNLB	BGP7LXN
^BGP7LXNLC	BGP7LXN
^BGP7LXLO	BGP7LX
^BGP7LXLOB	BGP7LXO
^BGP7LXLOC	BGP7LXO
^BGP7LXLOD	BGP7LXO
^BGP7LXLOE	BGP7LXO
^BGP7LXLOF	BGP7LXO
^BGP7LXLOG	BGP7LXO
^BGP7LXLOH	BGP7LXO
^BGP7LXLOI	BGP7LXO
^BGP7LXLOJ	BGP7LXO
^BGP7LXLOK	BGP7LXO
^BGP7LXP	BGP7LX
^BGP7LXPB	BGP7LXP
^BGP7LXPC	BGP7LXP
^BGP7LXPD	BGP7LXP
^BGP7LXPE	BGP7LXP
^BGP7LXPF	BGP7LXP
^BGP7LXPG	BGP7LXP
^BGP7LXR	BGP7LX
^BGP7LXS	BGP7LX
^BGP7LXT	BGP7LX
^BGP7LXU	BGP7LX
^BGP7LXV	BGP7LX
^BGP7LXW	BGP7LX
^BGP7LXX	BGP7LX
^BGP7LXY	BGP7LX
^BGP7LXZ	BGP7LX
^BGP7MX	BGP7POS
^BGP7MXA	BGP7MX
^BGP7MXB	BGP7MX
^BGP7NPL	opt
NPL1^BGP7NPLD	BGP7PDL
S^BGP7NPLD	BGP7DSTM
NPL1^BGP7NPLP	BGP7DP

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
460

Routine	is Invoked by:
PRINT^BGP7PARP	BGP7DAR, BGP7DGPA
^BGP7PD19	BGP7PDL1
I1AGE^BGP7PD19	BGP7PDL1
^BGP7PDL	BGP7DP, BGP7PARP
HEADER^BGP7PDL	BGP7PDLN
\$\$SB^BGP7PDL1	BGP7D23, BGP7D6, BGP7EL3
H1^BGP7PDL1	BGP7DPED, BGP7PDL0, BGP7PDL3, BGP7PDL4, BGP7PDL7

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>

H2^BGP7PDL1	BGP7DPED
PI^BGP7PDL1	BGP7PDL0, BGP7PDL3, BGP7PDL4, BGP7PDL7
PI1^BGP7PDL1	BGP7PDL4, BGP7PDL7

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>

SETN^BGP7PDL1	BGP7PDL0
I1AGEP^BGP7PDL5	BGP7PDL0
I1AGEP^BGP7PDL8	BGP7PDL7
I1AGE^BGP7PDL9	BGP7PDL1
I1AGE^BGP7PDLA	BGP7PDL1
^BGP7PDLH	BGP7DPED, BGP7PDL

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>

NONSUM^BGP7PDLN	BGP7PDL
^BGP7PDL5	BGP7PDL
EN^BGP7PESI	BGP7DPE
EN^BGP7PESL	BGP7DPE
PT^BGP7PESL	BGP7DPE, BGP7DPEA
RT^BGP7PESL	BGP7DPE
PRINT^BGP7PHEL	BGP7DELA
PRINT^BGP7PHEP	BGP7DHEA
CLTAX^BGP7POS1	BGP7POS
DRUGS^BGP7POS1	BGP7POS
LAB^BGP7POS1	BGP7POS
MHTAX^BGP7POS1	BGP7POS
PRE^BGP7POS1	BGP7POS
PRVTAX^BGP7POS1	BGP7POS
^BGP7PPEL	BGP7DPEA
GS^BGP7PUTL	BGP7DPE, BGP7GPED
REPORT^BGP7PUTL	BGP7DPE, BGP7GPED
^BGP7SDP	BGP7DP
^BGP7SDPN	BGP7DP
^BGP7ULF	opt
\$\$D^BGP7UTL	BGP7DAR, BGP7DNE, BGP7DNG, BGP7GAGP, BGP7GNTL, BGP7GNXP
\$\$DATE^BGP7UTL	BGP7ASL, BGP7C1, BGP7C11, BGP7C12, BGP7C13, BGP7C2, BGP7C3 BGP7C31, BGP7D2, BGP7D21, BGP7D22, BGP7D23, BGP7D3, BGP7D31 BGP7D4, BGP7D41, BGP7D42, BGP7D5, BGP7D51, BGP7D52, BGP7D6

BGP7D61, BGP7D7, BGP7D71, BGP7D72, BGP7D721, BGP7D722, BGP7D73
BGP7D731, BGP7D732
BGP7D74, BGP7D75, BGP7D8, BGP7D81, BGP7D82, BGP7D83, BGP7D84
BGP7D86, BGP7D862, BGP7D9, BGP7DCLP, BGP7DGPA, BGP7DHE1
BGP7DP, BGP7EL2, BGP7EL3, BGP7EL31, BGP7EL4

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 461

Routine	is Invoked by:
\$\$EDT^BGP7UTL	BGP7DNE
\$\$LZERO^BGP7UTL	BGP7DAR, BGP7DCHA, BGP7DNE, BGP7DNG, BGP7GAGP
<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>	
\$\$RZERO^BGP7UTL	BGP7DAR, BGP7DNG, BGP7GACW, BGP7GAGP, BGP7GNTL
EXCELGS^BGP7UTL	BGP7DAR, BGP7DNG, BGP7GAGP
GETIEN^BGP7UTL	BGP7GUPL, BGP7ULF
GS^BGP7UTL	BGP7DGPU, BGP7DNG, BGP7GGPU, BGP7GNTL
JRNL^BGP7UTL	BGP7D1, BGP7D10, BGP7DNE
REPORT^BGP7UTL	
BGP7DAP, BGP7DCL, BGP7DGPU, BGP7DL, BGP7DNG, BGP7DPP, BGP7DSTM	BGP7GALL, BGP7GCMP, BGP7GCOM, BGP7GGPU, BGP7GNPL, BGP7GNST
	BGP7GNTL, BGP7GPP, BGP7NPL
STMP^BGP7UTL	BGP7DCI
XTMP^BGP7UTL	BGP7C1, BGP7D1, BGP7DPE1
\$\$CPTREFT^BGP7UTL1	BGP7D61
\$\$FIRSTPRC^BGP7UTL1	BGP7EL3
\$\$LASTDX^BGP7UTL1	
BGP7D23, BGP7D42, BGP7D7, BGP7D72, BGP7D721, BGP7D74, BGP7EL3	BGP7EL4
\$\$LASTDXI^BGP7UTL1	BGP7D21, BGP7D42, BGP7D721, BGP7D731, BGP7D75
\$\$LASTPRC^BGP7UTL1	BGP7C12, BGP7C2, BGP7D3, BGP7D5, BGP7D61, BGP7D72, BGP7DHE1
	BGP7EL4
\$\$LASTPRCI^BGP7UTL1	BGP7D3, BGP7D42, BGP7D61, BGP7D83, BGP7D84
\$\$RADREF^BGP7UTL1	BGP7D42, BGP7D61
\$\$REFUSAL^BGP7UTL1	BGP7D21, BGP7D22, BGP7D3, BGP7D31, BGP7D4, BGP7D5, BGP7D6
	BGP7D61, BGP7D74, BGP7D8, BGP7DHE1
GETMEDS^BGP7UTL2	BGP7D22, BGP7D721, BGP7D722, BGP7D731, BGP7D732, BGP7D74
	BGP7D75, BGP7D81, BGP7D82, BGP7D83, BGP7D84, BGP7D863
	BGP7D864
^BGP7XTCH	opt
TAXCHK^BGP7XTCH	BGP7DAP, BGP7DL, BGP7DPP
^BGP7XTCN	opt
TAXCHK^BGP7XTCN	BGP7DCL, BGP7DGPU, BGP7DNE, BGP7DNG, BGP7DSTM, BGP7NPL
^BGP7XTEL	opt
TAXCHK^BGP7XTEL	BGP7DEL
^BGP7XTHE	opt
TAXCHK^BGP7XTHE	BGP7DHE
EP^BGP7XTL	BGP7XTS
^BGP7XTV	opt
INIT^BGP7XTV	BGP7GU
EP^BGP7XTV1	BGP7CTS, BGP7XTV
^BGPAP	BGPDPARP
I1^BGPAP1	BGPAP
I12^BGPAP12	BGPAP
I13^BGPAP13	BGPAP
I14^BGPAP14	BGPAP
I1B^BGPAP1B	BGPAP
I2A^BGPAP2	BGPAP
I2B^BGPAP2	BGPAP
I2C^BGPAP2	BGPAP

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
462

Routine	is Invoked by:
I22^BGPAP22	BGPAP
I23^BGPAP23	BGPAP
I24^BGPAP24	BGPAP
I3A^BGPAP3	BGPAP
I3B^BGPAP3	BGPAP
I3C^BGPAP3	BGPAP
I4A^BGPAP4	BGPAP
I4B^BGPAP4	BGPAP
I4C^BGPAP4	BGPAP
I5A^BGPAP5	BGPAP
I5B^BGPAP5	BGPAP
I5C^BGPAP5	BGPAP
I6^BGPAP6	BGPAP
I6A^BGPAP6	BGPAP
I7^BGPAP7	BGPAP
I8^BGPAP8	BGPAP
IA^BGPAPA	BGPAP
IB^BGPAPB	BGPAP
IC^BGPAPC	BGPAP
ID^BGPAPD	BGPAP
H2^BGPAPH	BGPAP13, BGPAP22B, BGPAP22C, BGPAP30, BGPAP30A, BGPAP30B BGPAP30C
HEADER^BGPAPH	BGPAP1, BGPAP12, BGPAP13, BGPAP14, BGPAP1B, BGPAP2, BGPAP22 BGPAP22A, BGPAP22B, BGPAP22C, BGPAP23, BGPAP24, BGPAP3
BGPAP30, BGPAP30A, BGPAP30B, BGPAP30C, BGPAP4, BGPAP5, BGPAP6 BGPAP7, BGPAP8, BGPAPA, BGPAPC, BGPAPD, BGPAPH	
^BGPD	opt
EN^BGPD0	BGPD
^BGPD1	BGPD, BGPDARP, BGPDFTA
\$\$DM^BGPD1	BGPD24, BGPDA, BGPDD
I12^BGPD12	BGPD1
I13^BGPD13	BGPD1
I14^BGPD14	BGPD1
I1B^BGPD1B	BGPD1
I2A^BGPD2	BGPD1
I2B^BGPD2	BGPD1
I2C^BGPD2	BGPD1
I22^BGPD22	BGPD1
I24^BGPD24	BGPD1
I29^BGPD29	BGPD1
I3A^BGPD3	BGPD1
I3B^BGPD3	BGPD1
I3C^BGPD3	BGPD1
I30^BGPD30	BGPD1
I4A^BGPD4	BGPD1
I4B^BGPD4	BGPD1
I4C^BGPD4	BGPD1
I5A^BGPD5	BGPD1
I5B^BGPD5	BGPD1

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
463

Routine	is Invoked by:
I5C^BGPD5	BGPD1
I6^BGPD6	BGPD1
I6A^BGPD6	BGPD1
I7^BGPD7	BGPD1
I8^BGPD8	BGPD1
IA^BGPDA	BGPD1
^BGPDAREA	opt
^BGPDARP	opt
IB^BGPDB	BGPD1
IC^BGPDC	BGPD1
ID^BGPDD	BGPD1
^BGPDFTA	opt
^BGPDPH	BGPDP
EOP^BGPDPH	BGP3ULF , BGP4LSTF , BGP4ULF , BGP5GUPL , BGP5LSTF , BGP5ULF BGPDUP
EN^BGPDL	BGPD , BGPDARP , BGPDFTA
^BGPDP	BGPD , BGPDARP , BGPDFTA
PRINT1^BGPDP	BGPDARP
I1^BGPDP1	BGPDP
I12^BGPDP12	BGPDP
I13^BGPDP13	BGPDP
I14^BGPDP14	BGPDP
I1B^BGPDP1B	BGPDP
I2A^BGPDP2	BGPDP
I2B^BGPDP2	BGPDP
I2C^BGPDP2	BGPDP
I22^BGPDP22	BGPDP
POV^BGPDP221	BGPDP22
AAGE^BGPDP22A	BGPDP22
MAGE^BGPDP22B	BGPDP22
FAGE^BGPDP22C	BGPDP22
I23^BGPDP23	BGPDP
I24^BGPDP24	BGPDP
I29^BGPDP29	BGPDP
I3A^BGPDP3	BGPDP
I3B^BGPDP3	BGPDP
I3C^BGPDP3	BGPDP
I30^BGPDP30	BGPDP
^BGPDP30A	BGPDP30
^BGPDP30B	BGPDP30
^BGPDP30C	BGPDP30
I4A^BGPDP4	BGPDP
I4B^BGPDP4	BGPDP
I4C^BGPDP4	BGPDP
I5A^BGPDP5	BGPDP
I5B^BGPDP5	BGPDP
I5C^BGPDP5	BGPDP
I6^BGPDP6	BGPDP
I6A^BGPDP6	BGPDP

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 464

Routine	is Invoked by:
I7^BGPDP7	BGPDP
I8^BGPDP8	BGPDP
IA^BGPDPA	BGPDP
PRINT^BGPDPARP	BGPDAREA
IB^BGPDPB	BGPDP
IC^BGPDPD	BGPDP
ID^BGPDPD	BGPDP
AREACP^BGPDPH	BGPDPARP
H1^BGPDPH	BGPDP1, BGPDP12, BGPDP13, BGPDP14, BGPDP1B, BGPDP2, BGPDP22
	BGPDP22B, BGPDP22C, BGPDP23, BGPDP24, BGPDP29, BGPDP3, BGPDP30
	BGPDP30A, BGPDP30B, BGPDP30C, BGPDP31, BGPDP32, BGPDP4, BGPDP5
	BGPDP6, BGPDP7, BGPDP8, BGPDPA, BGPDPB, BGPDPC, BGPDPD
H2^BGPDPH	BGPDP1, BGPDP12, BGPDP13, BGPDP14, BGPDP1B, BGPDP2, BGPDP22
	BGPDP22B, BGPDP22C, BGPDP23, BGPDP24, BGPDP29, BGPDP3, BGPDP30
	BGPDP30A, BGPDP30B, BGPDP30C, BGPDP31, BGPDP32, BGPDP4, BGPDP5
	BGPDP6, BGPDP7, BGPDP8, BGPDPA, BGPDPB, BGPDPC, BGPDPD
H3^BGPDPH	BGPDP1
H4^BGPDPH	BGPDP1B
H6^BGPDPH	BGPDPC
HEADER^BGPDPH	BGP3DPH, BGPAP1, BGPAP12, BGPAP13, BGPAP14, BGPAP1B, BGPAP2
	BGPAP22, BGPAP23, BGPAP24, BGPAP3, BGPAP4, BGPAP5, BGPAP6
	BGPAP7, BGPAP8, BGPAPB, BGPAPD, BGPDP1, BGPDP12, BGPDP13
	BGPDP14, BGPDP1B, BGPDP2, BGPDP22, BGPDP22A, BGPDP22B
	BGPDP22C, BGPDP23, BGPDP24
	BGPDP29, BGPDP3, BGPDP30, BGPDP30A, BGPDP30B, BGPDP30C
	BGPDP31, BGPDP32, BGPDP4, BGPDP5, BGPDP6, BGPDP7, BGPDP8
	BGPDPA, BGPDPARP, BGPDPB, BGPDPC, BGPDPD, BGPDPH
^BGPDS	BGPDP
TAXCHK^BGPDT	BGPD, BGPDARP, BGPDFTA, opt
^BGPDU	BGP3D6, BGP4D61, BGPDB
\$\$CPT^BGPDU	
BGP3D3, BGP3D4, BGP3D8, BGP4D3, BGP4D4, BGP4D5, BGP4D61, BGP4D8	
	BGP5D22, BGP5D3, BGP5D4, BGP5D5, BGP5D61, BGP5D72, BGP5D8
	BGP5DHE1, BGP5EL4, BGPD24, BGPD6, BGPD7
\$\$CPTI^BGPDU	BGP3D2, BGP3D4, BGP3D6, BGP3D7, BGP3D9, BGP4D2, BGP4D21
	BGP4D31, BGP4D61, BGP4D7, BGP4D71, BGP4D8, BGP5D2, BGP5D21
	BGP5D31, BGP5D61, BGP5D7, BGP5D71, BGP5D72, BGP5D73, BGP5D8
	BGP5D9, BGP5EL2, BGPD4, BGPDB
\$\$PLCODE^BGPDU	BGP5D34
\$\$PLTAX^BGPDU	BGP5D33, BGP5D34, BGP5D35
\$\$RAD^BGPDU	BGP5D61, BGP5EL4
\$\$WH^BGPDU	BGP3D3, BGP4D3, BGP4D4, BGP5D3, BGP5D4, BGPD6, BGPD7
BANNER^BGPDU	opt
^BGPDUP	opt
CATSTR^BGPGR	BGP5GRA, BGP6GRA, BGP7GRA, BGP6GRA
^BGPLTXA	BGPLTX
^BGPLTXB	BGPLTX
^BGPLTXC	BGPLTX

^BGPLTXD
^BGPLTXE

BGPLTX
BGPLTX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
465

Routine	is Invoked by:
^BGPLTXF	BGPLTX
^BGPLTXFB	BGPLTXF
^BGPLTXFC	BGPLTXF
^BGPLXA	BGPLX
^BGPLXB	BGPLX
^BGPLXBB	BGPLXB
^BGPLXC	BGPLX
^BGPLXCB	BGPLXC
^BGPLXCC	BGPLXC
^BGPLXCD	BGPLXC
^BGPLXCE	BGPLXC
^BGPLXCF	BGPLXC
^BGPLXCG	BGPLXC
^BGPLXCH	BGPLXC
^BGPLXCI	BGPLXC
^BGPLXCJ	BGPLXC
^BGPLXCK	BGPLXC
^BGPLXCL	BGPLXC
^BGPLXCM	BGPLXC
^BGPLXCN	BGPLXC
^BGPLXCO	BGPLXC
^BGPLXCP	BGPLXC
^BGPLXCQ	BGPLXC
^BGPLXD	BGPLX
^BGPLXE	BGPLX
^BGPLXEB	BGPLXE
^BGPLXEC	BGPLXE
^BGPLXF	BGPLX
^BGPLXG	BGPLX
^BGPLXH	BGPLX
^BGPLXI	BGPLX
^BGPLXJ	BGPLX
^BGPLXK	BGPLX
^BGPLXL	BGPLX
^BGPLXM	BGPLX
^BGPLXMB	BGPLXM
^BGPLXMC	BGPLXM
^BGPMXA	BGPMX
^BGPMXB	BGPMX
CLTAX^BGPPPOST1	BGPPPOST
DRUGS^BGPPPOST1	BGPPPOST
LAB^BGPPPOST1	BGPPPOST
MHTAX^BGPPPOST1	BGPPPOST
PRE^BGPPPOST1	BGPPPOST
PRVTAX^BGPPPOST1	BGPPPOST
^BGPTXA	BGPTX
^BGPTXB	BGPTX
^BGPTXC	BGPTX
^BGPTXD	BGPTX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
466

Routine	is Invoked by:
^BGPTXE	BGPTX
^BGPTXEB	BGPTXE
^BGPTXEC	BGPTXE
^BGPTXED	BGPTXE
^BGPTXEE	BGPTXE
^BGPTXEF	BGPTXE
^BGPTXEG	BGPTXE
^BGPTXF	BGPTX
^BGPTXG	BGPTX
^BGPTXGB	BGPTXG
^BGPTXH	BGP50P1, BGPTX
^BGPTXI	BGPTX
^BGPTXJ	BGPTX
^BGPTXK	BGPTX
^BGPTXL	BGPTX
^BGPTXM	BGPTX
^BGPTXN	BGPTX
^BGPTXO	BGPTX
^BGPTXP	BGPTX
^BGPTXPB	BGPTXP
^BGPTXQ	BGPTX
^BGPTXR	BGPTX
^BGPTXS	BGPTX
^BGPTXT	BGPTX
^BGPTXU	BGPTX
^BGPTXV	BGPTX
^BGPTXW	BGPTX
^BGPTXX	BGPTX
^BGPTXY	BGPTX
^BGPTXZ	BGPTX
^BGPUTXA	BGPUTX
^BGPUTXAB	BGPUTXA
^BGPUTXAC	BGPUTXA
^BGPUTXAD	BGPUTXA
^BGPUTXAE	BGPUTXA
^BGPUTXAF	BGPUTXA
^BGPUTXAG	BGPUTXA
^BGPUTXAH	BGPUTXA
^BGPUTXAI	BGPUTXA
^BGPUTXAJ	BGPUTXA
^BGPUTXAK	BGPUTXA
^BGPUTXAL	BGPUTXA
^BGPUTXAM	BGPUTXA
^BGPUTXAN	BGPUTXA
^BGPUTXAO	BGPUTXA
^BGPUTXAP	BGPUTXA
^BGPUTXAQ	BGPUTXA
^BGPUTXAR	BGPUTXA
^BGPUTXB	BGPUTX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
467

Routine	is Invoked by:
^BGPUTCX	BGPUTCX
^BGPUTCXD	BGPUTCX
^BGPUTCXE	BGPUTCX
^BGPUTCXF	BGPUTCX
^BGPUTCXG	BGPUTCX
^BGPUTCXH	BGPUTCX
^BGPUTCXI	BGPUTCX
^BGPUTCXJ	BGPUTCX
^BGPUTCXK	BGPUTCX
^BGPUTCXL	BGPUTCX
^BGPUTCXM	BGPUTCX
^BGPUTCXN	BGPUTCX
^BGPUTCXO	BGPUTCX
^BGPUTCXP	BGPUTCX
^BGPUTCXQ	BGPUTCX
^BGPUTCXR	BGPUTCX
^BGPUTCXS	BGPUTCX
^BGPUTCXT	BGPUTCX
^BGPUTCXU	BGPUTCX
^BGPUTCXV	BGPUTCX
^BGPUTCXW	BGPUTCX
^BGPUTCXX	BGPUTCX
^BGPUTCXY	BGPUTCX
^BGPUTCXZ	BGPUTCX
^BGPUTCXZB	BGPUTCXZ
^BGPUTCXZC	BGPUTCXZ
^BGPUTCXZD	BGPUTCXZ
^BGPUTCXZE	BGPUTCXZ
^BGPUXA	BGPUX
^BGPUXAB	BGPUXA
^BGPUXAC	BGPUXA
^BGPUXAD	BGPUXA
^BGPUXAE	BGPUXA
^BGPUXAF	BGPUXA
^BGPUXAG	BGPUXA
^BGPUXAH	BGPUXA
^BGPUXAI	BGPUXA
^BGPUXAJ	BGPUXA
^BGPUXAK	BGPUXA
^BGPUXAL	BGPUXA
^BGPUXAM	BGPUXA
^BGPUXAN	BGPUXA
^BGPUXAO	BGPUXA
^BGPUXAP	BGPUXA
^BGPUXAQ	BGPUXA
^BGPUXB	BGPUX
^BGPUXC	BGPUX
^BGPUXD	BGPUX
^BGPUXE	BGPUX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
468

Routine	is Invoked by:
^BGPUXF	BGPUX
^BGPUXG	BGPUX
^BGPUXGB	BGPUXG
^BGPUXGC	BGPUXG
^BGPUXH	BGPUX
^BGPUXI	BGPUX
^BGPUXJ	BGPUX
^BGPUXK	BGPUX
^BGPUXL	BGPUX
^BGPUXM	BGPUX
^BGPUXN	BGPUX
^BGPUXO	BGPUX
^BGPUXP	BGPUX
^BGPUXQ	BGPUX
^BGPUXR	BGPUX
^BGPUXS	BGPUX
^BGPUXT	BGPUX
^BGPUXU	BGPUX
^BGPUXV	BGPUX
^BGPUXW	BGPUX
^BGPUXY	BGPUX
^BGPUXZ	BGPUX
^BGPVTXA	BGPVTX
^BGPVTXB	BGPVTX
^BGPVTXC	BGPVTX
^BGPVTXD	BGPVTX
^BGPVTXE	BGPVTX
^BGPVTXF	BGPVTX
^BGPVTXG	BGPVTX
^BGPVTXH	BGPVTX
^BGPVTXI	BGPVTX
^BGPVTXIB	BGPVTXI
^BGPVTXIC	BGPVTXI
^BGPVTXJ	BGPVTX
^BGPVTXK	BGPVTX
^BGPVTXL	BGPVTX
^BGPVTXM	BGPVTX
^BGPVTXN	BGPVTX
^BGPVTXO	BGPVTX
^BGPVTXP	BGPVTX
^BGPVTXQ	BGPVTX
^BGPVTXR	BGPVTX
^BGPVTXS	BGPVTX
^BGPVTXT	BGPVTX
^BGPVTXU	BGPVTX
^BGPVTXV	BGPVTX
^BGPVTXW	BGPVTX
^BGPVXXX	BGPVTX
^BGPVXTY	BGPVTX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page
469

Routine	is Invoked by:
^BGPVTXZ	BGPVTX
^BGPVXA	BGPVX
^BGPVXB	BGPVX
^BGPVXC	BGPVX
^BGPVXD	BGPVX
^BGPVXE	BGPVX
^BGPVXF	BGPVX
^BGPVXH	BGPVX
^BGPVXI	BGPVX
^BGPVXJ	BGPVX
^BGPVXK	BGPVX
^BGPVXL	BGPVX
^BGPVXLB	BGPVXL
^BGPVXLC	BGPVXL
^BGPVXM	BGPVX
^BGPVXN	BGPVX
^BGPVXO	BGPVX
^BGPVXP	BGPVX
^BGPVXQ	BGPVX
^BGPVXR	BGPVX
^BGPVXS	BGPVX
^BGPVXU	BGPVX
^BGPVXUB	BGPVXU
^BGPVXUC	BGPVXU
^BGPVXUD	BGPVXU
^BGPVXUE	BGPVXU
^BGPVXV	BGPVX
^BGPVXW	BGPVX
^BGPVXX	BGPVX
^BGPVXY	BGPVX
^BGPVXZ	BGPVX
^BGPWTXA	BGPWTX
^BGPWTXB	BGPWTX
^BGPWTXBB	BGPWTXB
^BGPWTXBC	BGPWTXB
^BGPWTXBD	BGPWTXB
^BGPWTXBE	BGPWTXB
^BGPWTXBF	BGPWTXB
^BGPWTXBG	BGPWTXB
^BGPWTXC	BGPWTX
^BGPWTXD	BGPWTX
^BGPWTXE	BGPWTX
^BGPWTXF	BGPWTX
^BGPWTXG	BGPWTX
^BGPWTXH	BGPWTX
^BGPWTXI	BGPWTX
^BGPWTXJ	BGPWTX
^BGPWTXK	BGPWTX
^BGPWTXL	BGPWTX

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 470

Routine	is Invoked by:
^BGPWTXM	BGPWTX
^BGPWTXN	BGPWTX
^BGPWTXO	BGPWTX
^BGPWTXP	BGPWTX
^BGPWTXQ	BGPWTX
^BGPWTXR	BGPWTX
^BGPWTXS	BGPWTX
^BGPWTXT	BGPWTX
^BGPWTXU	BGPWTX
IMMFORC^BIRPC	BGPD23
^DIC	BGP3DA, BGP3DG, BGP3DL, BGP3UTL, BGP4DA, BGP4DG, BGP4DHE BGP4DL, BGP4DPP, BGP4HUTL, BGP4UTL, BGP5CL, BGP5CTL, BGP5DEL
	BGP5DGPU, BGP5DHE, BGP5DL, BGP5DNG, BGP5DPP, BGP5ETL, BGP5EUTL
	BGP5HUTL, BGP5NPL, BGP5TXL, BGP5UTL, BGP6AHW, BGP6CL, BGP6CTL BGP6DAP, BGP6DCL
	BGP6DEL, BGP6DGPU, BGP6DHE, BGP6DL, BGP6DNG, BGP6DPP, BGP6DSTM
	BGP6ETL, BGP6EUTL, BGP6HUTL, BGP6NPL, BGP6UTL, BGP6XTL, BGP7CL
	BGP7CTL, BGP7DAP, BGP7DCL, BGP7DEL, BGP7DGPU, BGP7DHE, BGP7DL BGP7DNE, BGP7DNG, BGP7DPE, BGP7DPP, BGP7DSTM, BGP7ETL BGP7EUTL, BGP7GNXP BGP7HUTL, BGP7NPL, BGP7POS, BGP7POS1, BGP7PUTL, BGP7UTL BGP7XTL, BGPD, BGPDARP, BGPDFTA, BGPPOST, BGPPOST1 dd90241.0111, dd90241.02
IX^DIC	BGP7POS, BGPPOST
FILE^DICN	BGP3ULF, BGP4ULF, BGP5CTL, BGP5ETL, BGP5GAEL, BGP5GAGP BGP5GAHE, BGP5GALL, BGP5GAPU, BGP5GCMS, BGP5GCOM, BGP5GELD BGP5GGPU, BGP5GHED, BGP5GLST, BGP5GNPL, BGP5GNTL, BGP5GPP BGP5GUPL, BGP5TXL, BGP5ULF, BGP6CTL, BGP6ETL, BGP6GAEL BGP6GAGP, BGP6GAHE, BGP6GALL BGP6GAPU, BGP6GCMP, BGP6GCMS, BGP6GCOM, BGP6GELD, BGP6GGPU BGP6GHED, BGP6GLST, BGP6GLTX, BGP6GMTX, BGP6GNPL, BGP6GNST BGP6GNTL, BGP6GPP, BGP6GUPL, BGP6ULF, BGP6XTL, BGP7CTL BGP7ETL, BGP7GAEL, BGP7GAGP, BGP7GAHE, BGP7GALL, BGP7GAPE BGP7GAPU, BGP7GCMP, BGP7GCMS BGP7GCOM, BGP7GELD, BGP7GGPU, BGP7GHED, BGP7GLST, BGP7GLTX BGP7GMTX, BGP7GNPL, BGP7GNST, BGP7GNTL, BGP7GNXP, BGP7GPED BGP7GPP, BGP7GUPL, BGP7ULF, BGP7XTL, BGPDUP, BGPPOST
DT^DICRW	BGP50P1
^DIE	BGP5CTL, BGP5ETL, BGP5GAEL, BGP5GAGP, BGP5GAHE, BGP5GALL BGP5GAPU, BGP5GCMS, BGP5GCOM, BGP5GELD, BGP5GGPU, BGP5GHED BGP5GLST, BGP5GNPL, BGP5GNTL, BGP5GPP, BGP5TXL, BGP6CTL BGP6DSTM, BGP6ETL, BGP6GACW, BGP6GAEL, BGP6GAGP, BGP6GAHE BGP6GALL, BGP6GAPU, BGP6GCMP BGP6GCMS, BGP6GCOM, BGP6GELD, BGP6GGPU, BGP6GHED, BGP6GLST BGP6GLTX, BGP6GMTX, BGP6GNPL, BGP6GNST, BGP6GNTL, BGP6GPP BGP6XTL, BGP7CTL, BGP7DSTM, BGP7ETL, BGP7GACW, BGP7GAEL BGP7GAGP, BGP7GAHE, BGP7GALL, BGP7GAPE, BGP7GAPU, BGP7GCMP

BGP7GCMS , BGP7GCOM , BGP7GELD
BGP7GGPU , BGP7GHED , BGP7GLST , BGP7GLTX , BGP7GMTX , BGP7GNPL
BGP7GNST , BGP7GNTL , BGP7GNXP , BGP7GPED , BGP7GPP , BGP7POS1
BGP7XTL , BGPPOST1

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 471

Routine	is Invoked by:
FILE^DIE	BGP5GU, BGP6GU, BGP7GU, BGP7GU
UPDATE^DIE	BGP5GU, BGP6GU, BGP7GU, BGP7GU
^DIK	BGP3DA, BGP3DG, BGP3DL, BGP3ULF, BGP4DA, BGP4DAP, BGP4DG BGP4DHE, BGP4DL, BGP4DPP, BGP4ULF, BGP50P1, BGP5DAP, BGP5DEL
	BGP5DGPU, BGP5DHE, BGP5DL, BGP5DNG, BGP5DPP, BGP5GU, BGP5GUPL
	BGP5NPL, BGP5ULF, BGP6AHW, BGP6DAP, BGP6DCL, BGP6DEL, BGP6DGPU BGP6DHE, BGP6DL
	BGP6DNG, BGP6DPP, BGP6DSTM, BGP6GU, BGP6GUPL, BGP6NPL, BGP6ULF
	BGP7DAP, BGP7DCL, BGP7DEL, BGP7DGPU, BGP7DHE, BGP7DL, BGP7DNE
	BGP7DNG, BGP7DPE, BGP7DPP, BGP7DSTM, BGP7GU, BGP7GUPL, BGP7NPL BGP7POS, BGP7POS1, BGP7ULF, BGP7D, BGP7DARP, BGP7DFTA, BGP7DUP BGP7GU, BGP7POST BGP7POST1
IX1^DIK	BGP3ULF, BGP4ULF, BGP5GUPL, BGP5ULF, BGP6DSTM, BGP6GUPL BGP6ULF, BGP7DSTM, BGP7GUPL, BGP7POS, BGP7POS1, BGP7ULF BGP7DUP, BGP7POST, BGP7POST1
IXALL^DIK	BGP7D
^DIM	dd90241.551101, dd90242, dd90242.01, dd90244.01 dd90244.02, dd90370.01, dd90370.02, dd90371.01 dd90371.02, dd90371.06, dd90372.01, dd90372.02 dd90373.01, dd90373.02, dd90373.06, dd90374.01 dd90374.02, dd90374.06, dd90374.09, dd90375.01 dd90375.02, dd90376.01, dd90376.02, dd90376.06 dd90530.01, dd90530.02, dd90530.06, dd90530.09 dd90531.01, dd90531.02, dd90532.01, dd90532.02 dd90532.06
\$\$GET1^DIQ	BGP5GU, BGP6DCHW, BGP6GU, BGP7DCHW, BGP7GU, BGP7GU
^DIR	BGP3DA, BGP3DAR, BGP3DARA, BGP3DARH, BGP3DG, BGP3DH, BGP3DL
	BGP3DP, BGP3DPH, BGP3DS, BGP3DSI, BGP3DSL, BGP3PDL, BGP3TXCH BGP3ULF, BGP4DA, BGP4DAP, BGP4DAR, BGP4DARA, BGP4DG, BGP4DH
	BGP4DHE, BGP4DHEA, BGP4DHEP, BGP4DHSL, BGP4DL, BGP4DP, BGP4DPH BGP4DPP, BGP4DS BGP4DSI, BGP4DSL, BGP4HEH, BGP4HEHH, BGP4HES, BGP4HESL BGP4LSTF, BGP4PDL, BGP4TXCH, BGP4TXHE, BGP4ULF, BGP5ASL BGP5CL, BGP5CP, BGP5CSI, BGP5CTL, BGP5CTS, BGP5CTX, BGP5DAP BGP5DAR, BGP5DEL, BGP5DELA, BGP5DELP, BGP5DESI, BGP5DESL BGP5DGPA, BGP5DGPU, BGP5DH
	BGP5DHE, BGP5DHEA, BGP5DHEP, BGP5DHSL, BGP5DL, BGP5DNG, BGP5DP
	BGP5DPH, BGP5DPP, BGP5DS, BGP5DSI, BGP5DSL, BGP5ELH, BGP5ELHH BGP5ELS, BGP5ELSL, BGP5ETL, BGP5ETS, BGP5HEH, BGP5HEHH BGP5HES, BGP5HESL, BGP5LSTF, BGP5NPL, BGP5NPLP, BGP5PDL BGP5TXCH, BGP5TXCN BGP5TXEL, BGP5TXHE, BGP5TXL, BGP5TXS, BGP5TXV, BGP5ULF

BGP6AHW , BGP6ASL , BGP6CL , BGP6CP , BGP6CSI , BGP6CTL , BGP6CTS
BGP6CTXC , BGP6DAP , BGP6DAR , BGP6DCHA , BGP6DCL , BGP6DCLP
BGP6DEL , BGP6DELA , BGP6DELP , BGP6DESI , BGP6DESL , BGP6DGPA
BGP6DGPU , BGP6DH , BGP6DHE
BGP6DHEA , BGP6DHEP , BGP6DHSL , BGP6DL , BGP6DLT , BGP6DMT
BGP6DNG , BGP6DP , BGP6DPH , BGP6DPP , BGP6DS , BGP6DSI , BGP6DSL
BGP6DSTM , BGP6ELH , BGP6ELHH , BGP6ELS , BGP6ELSL , BGP6ETL
BGP6HEH , BGP6HEHH , BGP6HES , BGP6HESL , BGP6LSTF , BGP6NPL
BGP6NPLP , BGP6PDLs , BGP6ULF

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 472

```

Routine          is Invoked by:
  ^DIR           BGP6XTCH, BGP6XTCN, BGP6XTEL, BGP6XTHE, BGP6XTL, BGP6XTS
                BGP6XTV, BGP7ASL, BGP7CL, BGP7CP, BGP7CSI, BGP7CTL, BGP7CTS
                BGP7CTXC, BGP7DAP, BGP7DAR, BGP7DCHA, BGP7DCL, BGP7DCLP
                BGP7DEL, BGP7DELA, BGP7DELP, BGP7DESI, BGP7DESL, BGP7DGPA
                BGP7DGPU, BGP7DH, BGP7DHE
                BGP7DHEA, BGP7DHEP, BGP7DHSL, BGP7DL, BGP7DLT, BGP7DMT
                BGP7DNE, BGP7DNG, BGP7DP, BGP7DPE, BGP7DPEA, BGP7DPEL

BGP7DPEP, BGP7DPH, BGP7DPP, BGP7DPSI, BGP7DS, BGP7DSI, BGP7DSL
                BGP7DSTM, BGP7ELH, BGP7ELHH, BGP7ELS, BGP7ELSL, BGP7ETL
                BGP7HEH, BGP7HEHH, BGP7HES
                BGP7HESL, BGP7LSTF, BGP7NPL, BGP7NPLP, BGP7PDLS, BGP7PESI
                BGP7PESL, BGP7ULF, BGP7XTCH, BGP7XTCN, BGP7XTEL, BGP7XTHE

BGP7XTL, BGP7XTS, BGP7XTV, BGPAP, BGPAPH, BGPD, BGPD0, BGPDAREA
                BGPDARP, BGPDFTA, BGPDPH, BGPDL, BGPDP, BGPDPH, BGPDS, BGPDT
                BGPDUP
  EN^DIU2       BGP50P1

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>
  $$CODEN^ICDCODE
BGP5UTL1, BGP6D51, BGP6DU, BGP6UTL1, BGP7D51, BGP7D61, BGP7DU
                BGP7UTL1, BGPDP22, BGPDU

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>
  $$ICDDX^ICDCODE   BGP3D10, BGP3D5, BGP3D9, BGP4D6, BGP4D7, BGP5C1, BGP5C11
                    BGP5C12, BGP5C13, BGP5C2, BGP5C3, BGP5D3, BGP5D41, BGP5D5
                    BGP5D6, BGP5D7, BGP5D72, BGP5UTL1, BGP6C1, BGP6C11, BGP6C12
                    BGP6C13, BGP6C2, BGP6C3, BGP6D23, BGP6D3, BGP6D41, BGP6D5
                    BGP6D6, BGP6D7, BGP6D72

BGP6D82, BGP6D83, BGP6D84, BGP6D862, BGP6UTL1, BGP7C1, BGP7C11
                BGP7C12, BGP7C13, BGP7C2, BGP7C3, BGP7D23, BGP7D3, BGP7D41
                BGP7D5, BGP7D6, BGP7D7, BGP7D72, BGP7D722, BGP7D732, BGP7D74

BGP7D82, BGP7D83, BGP7D84, BGP7D862, BGP7UTL1, BGPDP22, BGPDP23
                BGPDP29, BGPDP8
                BGPDP22

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>
  $$ICDOP^ICDCODE   BGP3D3, BGP4D4, BGP5C11, BGP5C2, BGP5C3, BGP5D3, BGP5D61

BGP5DHE1, BGP5UTL1, BGP6C11, BGP6C2, BGP6C3, BGP6D61, BGP6DHE1

BGP6UTL1, BGP7C11, BGP7C2, BGP7C3, BGP7D61, BGP7DHE1, BGP7UTL1
                BGPDP6

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>
  $$CODEN^ICPTCOD   BGP3D2, BGP3D4, BGP3D6, BGP3D7, BGP3D9, BGP4D2, BGP4D21
                    BGP4D31, BGP4D61, BGP4D7, BGP4D71, BGP4D8, BGP5D2, BGP5D21
                    BGP5D31, BGP5D4, BGP5D61, BGP5D7, BGP5D71, BGP5D72, BGP5D73
                    BGP5D8, BGP5D9, BGP5EL2, BGP6C1, BGP6C11, BGP6C13, BGP6C3

```

BGP6D2 ,BGP6D21 ,BGP6D23
BGP6D31 ,BGP6D4 ,BGP6D61 ,BGP6D71 ,BGP6D72 ,BGP6D73 ,BGP6D8
BGP6D82 ,BGP6D861 ,BGP6D9 ,BGP6DHE1 ,BGP6DU ,BGP6EL2 ,BGP7C1
BGP7C11 ,BGP7C12 ,BGP7C13 ,BGP7C3 ,BGP7D2 ,BGP7D21 ,BGP7D211

BGP7D21A ,BGP7D23 ,BGP7D31 ,BGP7D4 ,BGP7D61 ,BGP7D71 ,BGP7D72
BGP7D721 ,BGP7D73

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 473

Routine is Invoked by:

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>

\$\$CODEN^ICPTCOD BGP7D731, BGP7D732, BGP7D75, BGP7D8, BGP7D82, BGP7D861
BGP7D864, BGP7D9, BGP7DHE2, BGP7DU, BGP7EL2, BGPD4, BGPDB
BGPDD, BGPDU

<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>

\$\$CPT^ICPTCOD BGP4D4, BGP4D61, BGP5C11, BGP5C2, BGP5C3, BGP5D33, BGP5D34
BGP5D35, BGP5D4, BGP5EL3, BGP6C11, BGP6C2, BGP6C3, BGP6C31
BGP6D33, BGP6D34, BGP6D35, BGP6D4, BGP6D85, BGP6DU, BGP6EL3
BGP7C11, BGP7C2, BGP7C3, BGP7C31, BGP7D33, BGP7D34, BGP7D35
BGP7D4, BGP7D85, BGP7DU
BGP7EL3, BGPDU

EN^VALM

BGP3DSI, BGP3DSL, BGP4DHSL, BGP4DSI, BGP4DSL, BGP5ASL, BGP5CSI
BGP5CTL, BGP5CTS, BGP5DESI, BGP5DESL, BGP5DHSL, BGP5DSI

BGP5DSL, BGP5ETL, BGP5ETS, BGP5TXL, BGP5TXS, BGP5TXV, BGP5TXV1
BGP6ASL, BGP6CSI, BGP6CTL, BGP6CTS, BGP6DESI, BGP6DESL
BGP6DHSL, BGP6DSI

BGP6DSL, BGP6ETL, BGP6XTL, BGP6XTS, BGP6XTV, BGP6XTV1, BGP7ASL
BGP7CSI, BGP7CTL, BGP7CTS, BGP7DESI, BGP7DESL, BGP7DHSL
BGP7DPSI, BGP7DSI, BGP7DSL, BGP7ETL, BGP7PESI, BGP7PESL
BGP7XTL, BGP7XTS, BGP7XTV, BGP7XTV1, BGPD0, BGPLD

TERM^VALM0

BGP3DSI, BGP3DSL, BGP4DHSL, BGP4DSI, BGP4DSL, BGP5ASL, BGP5CSI
BGP5CTL, BGP5CTS, BGP5DESI, BGP5DESL, BGP5DHSL, BGP5DSI

BGP5DSL, BGP5ETL, BGP5ETS, BGP5TXL, BGP5TXS, BGP5TXV, BGP5TXV1
BGP6ASL, BGP6CSI, BGP6CTL, BGP6CTS, BGP6DESI, BGP6DESL
BGP6DHSL, BGP6DSI

BGP6DSL, BGP6ETL, BGP6XTL, BGP6XTS, BGP6XTV, BGP6XTV1, BGP7ASL
BGP7CSI, BGP7CTL, BGP7CTS, BGP7DESI, BGP7DESL, BGP7DHSL
BGP7DPSI, BGP7DSI, BGP7DSL, BGP7ETL, BGP7NPL, BGP7PESI
BGP7PESL, BGP7XTL, BGP7XTS, BGP7XTV, BGP7XTV1, BGPD0, BGPLD

CLEAR^VALM1

BGP3DSI, BGP3DSL, BGP4DHSL, BGP4DSI, BGP4DSL, BGP5ASL, BGP5CSI
BGP5CTL, BGP5CTS, BGP5DESI, BGP5DESL, BGP5DHSL, BGP5DSI

BGP5DSL, BGP5ETL, BGP5ETS, BGP5TXL, BGP5TXS, BGP5TXV, BGP5TXV1
BGP6ASL, BGP6CSI, BGP6CTL, BGP6CTS, BGP6DESI, BGP6DESL
BGP6DHSL, BGP6DSI

BGP6DSL, BGP6ETL, BGP6XTL, BGP6XTS, BGP6XTV, BGP6XTV1, BGP7ASL
BGP7CSI, BGP7CTL, BGP7CTS, BGP7DESI, BGP7DESL, BGP7DHSL
BGP7DPSI, BGP7DSI, BGP7DSL, BGP7ETL, BGP7PESI, BGP7PESL
BGP7XTL, BGP7XTS, BGP7XTV, BGP7XTV1, BGPD0, BGPLD

FULL^VALM1

BGP3DSI, BGP3DSL, BGP4DHSL, BGP4DSI, BGP4DSL, BGP5ASL, BGP5CSI
BGP5CTL, BGP5CTS, BGP5DESI, BGP5DESL, BGP5DHSL, BGP5DSI

BGP5DSL, BGP5ETL, BGP5ETS, BGP5TXL, BGP5TXS, BGP5TXV, BGP5TXV1
BGP6ASL, BGP6CSI, BGP6CTL, BGP6CTS, BGP6DESI, BGP6DESL
BGP6DHSL, BGP6DSI

BGP6DSL, BGP6ETL, BGP6XTL, BGP6XTS, BGP6XTV, BGP6XTV1, BGP7ASL
BGP7CSI, BGP7CTL, BGP7CTS, BGP7DESI, BGP7DESL, BGP7DHSL
BGP7DPSI, BGP7DSI, BGP7DSL, BGP7ETL, BGP7PESI, BGP7PESL
BGP7XTL, BGP7XTS, BGP7XTV, BGP7XTV1, BGP7D0, BGP7DL

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 474

Routine	is Invoked by:
^XBDBQUE	BGP3DA, BGP3DG, BGP3DL, BGP3TXCH, BGP4DA, BGP4DAP, BGP4DG BGP4DHE, BGP4DL, BGP4DPP, BGP4TXCH, BGP4TXHE, BGP5CTXC
BGP5DAP, BGP5DEL, BGP5DELA, BGP5DGPU, BGP5DHE, BGP5DL, BGP5DNG	BGP5DPP, BGP5NPL, BGP5TXCH, BGP5TXCN, BGP5TXEL, BGP5TXHE BGP6AHW, BGP6CTXC, BGP6DAP
BGP6DCL, BGP6DEL, BGP6DELA, BGP6DGPU, BGP6DHE, BGP6DL, BGP6DNG	BGP6DPP, BGP6DSTM, BGP6NPL, BGP6XTCH, BGP6XTCN, BGP6XTEL BGP6XTHE, BGP7CTXC, BGP7DAP, BGP7DCL, BGP7DEL, BGP7DELA
BGP7DGPU, BGP7DHE, BGP7DL, BGP7DNG, BGP7DPE, BGP7DPEA, BGP7DPP	BGP7DSTM, BGP7NPL
BGP7XTCH, BGP7XTCN, BGP7XTEL, BGP7XTHE, BGPD, BGPDAREA, BGPDT	
\$\$VAL^XBIDIQ1	BGP3D6, BGP3D9, BGP3DAR, BGP3DARA, BGP3DARH, BGP4D5, BGP4D7 BGP4DAR, BGP4DARA, BGP4DH, BGP4DHEA, BGP4HEH, BGP4HEHH BGP4PDLH, BGP5C1, BGP5C11, BGP5C13, BGP5C2, BGP5C3, BGP5CTL BGP5CTS, BGP5D5, BGP5D7, BGP5D72, BGP5DH, BGP5ELH, BGP5ELHH BGP5ETL, BGP5ETS, BGP5HEH BGP5HEHH, BGP5PDLH, BGP5TXV, BGP5TXV1, BGP6C1, BGP6C11 BGP6C12, BGP6C13, BGP6C2, BGP6C3, BGP6C31, BGP6CTL, BGP6CTS BGP6D5, BGP6D7, BGP6D72, BGP6DCHW, BGP6DH, BGP6ELH, BGP6ELHH BGP6ETL, BGP6HEH, BGP6HEHH, BGP6PDLH, BGP6XTV, BGP6XTV1 BGP7C1, BGP7C11, BGP7C12 BGP7C13, BGP7C2, BGP7C3, BGP7C31, BGP7CTL, BGP7CTS, BGP7D5 BGP7D7, BGP7D72, BGP7D721, BGP7D722, BGP7D731, BGP7D732
BGP7D75, BGP7DCHW, BGP7DH, BGP7ELH, BGP7ELHH, BGP7ETL, BGP7HEH	
\$\$VALI^XBIDIQ1	BGP7HEHH, BGP7PDLH, BGP7XTV, BGP7XTV1, BGPD30, BGPDAREA BGP3D3, BGP4D7, BGP5D7, BGP5EL4, BGP6C3, BGP6D7, BGP6EL4 BGP7C3, BGP7D7, BGP7EL4, BGPD22
^XBFMK	BGP3DA, BGP3DAR, BGP3DARA, BGP3DARH, BGP3DG, BGP3DL, BGP3ULF BGP3UTL, BGP4DA, BGP4DAP, BGP4DAR, BGP4DARA, BGP4DG, BGP4DHE
BGP4DHEA, BGP4DL, BGP4DPP, BGP4HUTL, BGP4ULF, BGP4UTL, BGP5CL	BGP5CTL, BGP5DAP, BGP5DAR, BGP5DEL, BGP5DELA, BGP5DGPA BGP5DGPU, BGP5DHE BGP5DHEA, BGP5DL, BGP5DNG, BGP5DPP, BGP5ETL, BGP5EUTL BGP5GAEL, BGP5GAGP, BGP5GAHE, BGP5GALL, BGP5GAPU, BGP5GCMS BGP5GCOM, BGP5GELD, BGP5GGPU, BGP5GHED, BGP5GLST, BGP5GNPL BGP5GNL, BGP5GPP, BGP5GUPL, BGP5HUTL, BGP5NPL, BGP5TXL BGP5ULF, BGP5UTL, BGP6AHW
BGP6CL, BGP6CTL, BGP6DAP, BGP6DAR, BGP6DCHA, BGP6DCL, BGP6DEL	
BGP6DELA, BGP6DGPA, BGP6DGPU, BGP6DHE, BGP6DHEA, BGP6DL	BGP6DLT, BGP6DMT, BGP6DNG, BGP6DPP, BGP6DSTM, BGP6ETL BGP6EUTL, BGP6GACW, BGP6GAEL, BGP6GAGP, BGP6GAHE, BGP6GALL BGP6GAPU, BGP6GCMP, BGP6GCMS BGP6GCOM, BGP6GELD, BGP6GGPU, BGP6GHED, BGP6GLST, BGP6GLTX BGP6GMTX, BGP6GNPL, BGP6GNST, BGP6GNL, BGP6GPP, BGP6GUPL

BGP6HUTL , BGP6NPL , BGP6ULF , BGP6UTL , BGP6XTL , BGP7CL , BGP7CTL
BGP7DAP , BGP7DAR , BGP7DCHA , BGP7DCL , BGP7DEL , BGP7DELA
BGP7DGPA , BGP7DGPU

BGP7DHE , BGP7DHEA , BGP7DL , BGP7DLT , BGP7DMT , BGP7DNE , BGP7DNG
BGP7DPE , BGP7DPEA , BGP7DPP , BGP7DSTM , BGP7ETL , BGP7EUTL
BGP7GACW , BGP7GAEL , BGP7GAGP , BGP7GAHE , BGP7GALL , BGP7GAPE
BGP7GAPU , BGP7GCMP , BGP7GCMS , BGP7GCOM , BGP7GELD , BGP7GGPU
BGP7GHED , BGP7GLST

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 475

Routine	is Invoked by:
^XBFBK	BGP7GLTX, BGP7GMTX, BGP7GNPL, BGP7GNST, BGP7GNTL, BGP7GNXP BGP7GPED, BGP7GPP, BGP7GUPL, BGP7HUTL, BGP7NPL, BGP7POS1 BGP7PUTL, BGP7ULF, BGP7UTL, BGP7XTL, BGPD, BGPDAREA, BGPDARP BGPDFTA, BGPDUP, BGPPPOST, BGPPPOST1
\$\$EXTSET^XBFFUNC	BGP7D74, BGP7DPE2
^XBGSAVE	
BGP3PDL, BGP3UTL, BGP4HEL, BGP4HUTL, BGP4PDL, BGP4UTL, BGP5ELL	BGP5EUTL, BGP5HEL, BGP5HUTL, BGP5PDL, BGP5UTL, BGP6DCHA BGP6ELL, BGP6EUTL, BGP6HEL, BGP6HUTL, BGP6PDL, BGP6UTL BGP7DNE, BGP7DPED, BGP7ELL, BGP7EUTL, BGP7HEL, BGP7HUTL BGP7PDL, BGP7PUTL, BGP7UTL BGPDARP, BGPDFTA
^XBKVAR	BGP5GAEL, BGP5GAGP, BGP5GAHE, BGP5GALL, BGP5GAPU, BGP5GCMS BGP5GCOM, BGP5GELD, BGP5GGPU, BGP5GHED, BGP5GLST, BGP5GNPL BGP5GNTL, BGP5GPP, BGP5GUPL, BGP6GACW, BGP6GAEL, BGP6GAGP BGP6GAHE, BGP6GALL, BGP6GAPU, BGP6GCMP, BGP6GCMS, BGP6GCOM BGP6GELD, BGP6GGPU BGP6GHED, BGP6GLST, BGP6GLTX, BGP6GMTX, BGP6GNPL, BGP6GNST BGP6GNTL, BGP6GPP, BGP6GUPL, BGP7GACW, BGP7GAEL, BGP7GAGP BGP7GAHE, BGP7GALL, BGP7GAPE, BGP7GAPU, BGP7GCMP, BGP7GCMS BGP7GCOM, BGP7GELD, BGP7GGPU, BGP7GHED, BGP7GLST, BGP7GLTX BGP7GMTX, BGP7GNPL BGP7GNST, BGP7GNTL, BGP7GNXP, BGP7GPED, BGP7GPP, BGP7GUPL
GUIR^XBLM	BGP5GAEL, BGP5GAGP, BGP5GAHE, BGP5GALL, BGP5GAPU, BGP5GCMS BGP5GCOM, BGP5GELD, BGP5GGPU, BGP5GHED, BGP5GNPL, BGP5GNTL BGP5GPP, BGP5GU, BGP6GAEL, BGP6GAGP, BGP6GAHE, BGP6GALL BGP6GAPU, BGP6GCMP, BGP6GCMS, BGP6GCOM, BGP6GELD, BGP6GGPU BGP6GHED, BGP6GLTX BGP6GMTX, BGP6GNPL, BGP6GNST, BGP6GNTL, BGP6GPP, BGP6GU BGP7GAEL, BGP7GAGP, BGP7GAHE, BGP7GALL, BGP7GAPE, BGP7GAPU BGP7GCMP, BGP7GCMS, BGP7GCOM, BGP7GELD, BGP7GGPU, BGP7GHED BGP7GLTX, BGP7GMTX, BGP7GNPL, BGP7GNST, BGP7GNTL, BGP7GNXP BGP7GPED, BGP7GPP, BGP7GU BGP7GNTL, BGP7GNXP, BGP7GPED, BGP7GPP, BGP7GUPL
EN^XBVK	BGP3DA, BGP3DAR, BGP3DARA, BGP3DARH, BGP3DG, BGP3DL, BGP3ULF BGP4DA, BGP4DAP, BGP4DAR, BGP4DARA, BGP4DG, BGP4DHE, BGP4DHEA BGP4DL, BGP4DPP, BGP4LSTF, BGP4ULF, BGP5CL, BGP5CTL, BGP5CTS BGP5DAP, BGP5DAR, BGP5DEL, BGP5DELA, BGP5DGPA, BGP5DGPU BGP5DHE, BGP5DHEA
BGP5DL, BGP5DNG, BGP5DPP, BGP5ETL, BGP5ETS, BGP5GAEL, BGP5GAGP	BGP5GAHE, BGP5GALL, BGP5GAPU, BGP5GCMS, BGP5GCOM, BGP5GELD BGP5GGPU, BGP5GHED, BGP5GLST, BGP5GNPL, BGP5GNTL, BGP5GPP
BGP5GR, BGP5GRA, BGP5GU, BGP5GUPL, BGP5LSTF, BGP5NPL, BGP5TXL	BGP5TXS, BGP5TXV
BGP5TXV1, BGP6AHW, BGP6CL, BGP6CTL, BGP6CTS, BGP6DAP, BGP6DAR	BGP6DCHA, BGP6DCL, BGP6DEL, BGP6DELA, BGP6DGPA, BGP6DGPU
BGP6DHE, BGP6DHEA, BGP6DL, BGP6DLT, BGP6DMT, BGP6DNG, BGP6DPP	

BGP6DSTM, BGP6ETL, BGP6GACW, BGP6GAEL, BGP6GAGP, BGP6GAHE
BGP6GALL, BGP6GAPU
BGP6GCMP, BGP6GCMS, BGP6GCOM, BGP6GELD, BGP6GGPU, BGP6GHED
BGP6GLST, BGP6GLTX, BGP6GMTX, BGP6GNPL, BGP6GNST, BGP6GNTL

BGP6GPP, BGP6GR, BGP6GR1, BGP6GRA, BGP6GU, BGP6GUPL, BGP6LSTF

BGP6NPL, BGP6XTL, BGP6XTS, BGP6XTV, BGP6XTV1, BGP7CL, BGP7CTL
BGP7CTS, BGP7DAP

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 476

Routine	is Invoked by:
EN^XBVK	BGP7DAR, BGP7DCHA, BGP7DCL, BGP7DEL, BGP7DELA, BGP7DGPA
BGP7DGPU, BGP7DHE, BGP7DHEA, BGP7DL, BGP7DLT, BGP7DMT, BGP7DNE	BGP7DNG, BGP7DPE, BGP7DPEA, BGP7DPP, BGP7DSTM, BGP7ETL BGP7GACW, BGP7GAEL, BGP7GAGP, BGP7GAHE, BGP7GALL, BGP7GAPE BGP7GAPU, BGP7GCMP, BGP7GCMS BGP7GCOM, BGP7GELD, BGP7GGPU, BGP7GHED, BGP7GLST, BGP7GLTX BGP7GMTX, BGP7GNPL, BGP7GNST, BGP7GNTL, BGP7GNXP, BGP7GPED BGP7GPP, BGP7GR, BGP7GR1, BGP7GR2, BGP7GRA, BGP7GU, BGP7GUPL BGP7LSTF, BGP7NPL, BGP7POS, BGP7XTL, BGP7XTS, BGP7XTV BGP7XTV1, BGP7D, BGP7DAREA
\$\$DT^XLFTD	BGP7DARP, BGP7DFTA, BGP7DUP, BGP7G, BGP7GRA, BGP7GU, BGP7POST BGP5GAEL, BGP5GAGP, BGP5GAHE, BGP5GALL, BGP5GAPU, BGP5GCMS BGP5GCOM, BGP5GELD, BGP5GGPU, BGP5GHED, BGP5GLST, BGP5GNPL BGP5GNTL, BGP5GPP, BGP5GUPL, BGP6GACW, BGP6GAEL, BGP6GAGP BGP6GAHE, BGP6GALL, BGP6GAPU, BGP6GCMP, BGP6GCMS, BGP6GCOM BGP6GELD, BGP6GGPU BGP6GHED, BGP6GLST, BGP6GLTX, BGP6GMTX, BGP6GNPL, BGP6GNST BGP6GNTL, BGP6GPP, BGP6GUPL, BGP7GACW, BGP7GAEL, BGP7GAGP BGP7GAHE, BGP7GALL, BGP7GAPE, BGP7GAPU, BGP7GCMP, BGP7GCMS BGP7GCOM, BGP7GELD, BGP7GGPU, BGP7GHED, BGP7GLST, BGP7GLTX BGP7GMTX, BGP7GNPL
\$\$FMADD^XLFTD	BGP7GNST, BGP7GNTL, BGP7GNXP, BGP7GPED, BGP7GPP, BGP7GUPL BGP3BAN, BGP3D1, BGP3D10, BGP3D3, BGP3D4, BGP3D5, BGP3D6 BGP3D7, BGP3D8, BGP3D9, BGP4BAN, BGP4D1, BGP4D10, BGP4D21 BGP4D3, BGP4D31, BGP4D4, BGP4D41, BGP4D5, BGP4D6, BGP4D61 BGP4D7, BGP4D8, BGP4D9, BGP5ASL, BGP5BAN, BGP5C1, BGP5C11 BGP5C12, BGP5C13, BGP5C3, BGP5CL BGP5D1, BGP5D10, BGP5D22, BGP5D3, BGP5D31, BGP5D4, BGP5D41 BGP5D5, BGP5D6, BGP5D61, BGP5D7, BGP5D71, BGP5D72, BGP5D73 BGP5D8, BGP5D9, BGP5DEL, BGP5DELA, BGP5DGPA, BGP5DGPU BGP5DHE1, BGP5EL3, BGP5GAEL, BGP5GAHE, BGP5GAPU, BGP5GCMS BGP5GELD, BGP5GGPU, BGP5GNPL
BGP5NPL, BGP5UTL1, BGP6AHW, BGP6ASL, BGP6BAN, BGP6C1, BGP6C11	BGP6C12, BGP6C13, BGP6C2, BGP6C3, BGP6C31, BGP6D1, BGP6D10 BGP6D21, BGP6D22, BGP6D23, BGP6D3, BGP6D31, BGP6D4, BGP6D41 BGP6D42, BGP6D5, BGP6D51, BGP6D52, BGP6D6, BGP6D61, BGP6D7 BGP6D71, BGP6D72
BGP6D73, BGP6D8, BGP6D82, BGP6D83, BGP6D84, BGP6D86, BGP6D863	BGP6D9, BGP6DCL, BGP6DEL, BGP6DELA, BGP6DGPA, BGP6DGPU BGP6DHE1, BGP6DSTM, BGP6DU, BGP6EL3, BGP6EL4, BGP6GAEL BGP6GAHE, BGP6GAPU, BGP6GCMP, BGP6GELD, BGP6GGPU, BGP6GNPL BGP6GNST, BGP6NPL, BGP6UTL
BGP6UTL1, BGP7ASL, BGP7BAN, BGP7C1, BGP7C11, BGP7C12, BGP7C13	BGP7C2, BGP7C3, BGP7C31, BGP7D1, BGP7D10, BGP7D21, BGP7D22 BGP7D23, BGP7D3, BGP7D31, BGP7D4, BGP7D41, BGP7D42, BGP7D5 BGP7D51, BGP7D52, BGP7D6, BGP7D61, BGP7D7, BGP7D71, BGP7D72 BGP7D721, BGP7D722

BGP7D73 ,BGP7D731 ,BGP7D732 ,BGP7D74 ,BGP7D75 ,BGP7D8 ,BGP7D82

BGP7D83 ,BGP7D84 ,BGP7D86 ,BGP7D863 ,BGP7D9 ,BGP7DAP ,BGP7DCL

BGP7DCLP ,BGP7DEL ,BGP7DELA ,BGP7DGPA ,BGP7DGPU ,BGP7DHE

BGP7DHE1 ,BGP7DHEA ,BGP7DL ,BGP7DNE ,BGP7DPE ,BGP7DPE1

BGP7DPEA ,BGP7DPP

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 477

```

Routine          is Invoked by:
  $$FMADD^XLFDT  BGP7DSTM,BGP7DU,BGP7EL3,BGP7EL4,BGP7GAEL,BGP7GAHE
                  BGP7GAPE,BGP7GAPU,BGP7GCMP,BGP7GELD,BGP7GGPU,BGP7GNPL
                  BGP7GNST,BGP7GPED,BGP7NPL,BGP7UTL,BGP7UTL1,BGPD,BGPD1
                  BGPD12,BGPD14,BGPD1B,BGPD2,BGPD23,BGPD24,BGPD29,BGPD3
                  BGPD4,BGPD5,BGPD6,BGPD7
                  BGPD8,BGPDA,BGPDB,BGPDC,BGPDD,BGPDU
  $$FMDIFF^XLFDT
BGP3D1,BGP3D3,BGP3D9,BGP4D1,BGP4D10,BGP4D7,BGP4D9,BGP5D1
                  BGP5D10,BGP5D32,BGP5D33,BGP5D34,BGP5D35,BGP5D7,BGP5D72
                  BGP5D9,BGP6D1,BGP6D10,BGP6D32,BGP6D33,BGP6D34,BGP6D35
                  BGP6D52,BGP6D7,BGP6D72,BGP6D82,BGP6D83,BGP6D84,BGP6D85
                  BGP6D861,BGP6D863
                  BGP6D9,BGP6GU,BGP7D1,BGP7D10,BGP7D32,BGP7D33,BGP7D34
                  BGP7D35,BGP7D52,BGP7D7,BGP7D72,BGP7D721,BGP7D722
                  BGP7D731,BGP7D732,BGP7D74,BGP7D75,BGP7D82,BGP7D83
                  BGP7D84,BGP7D85,BGP7D861,BGP7D863,BGP7D9,BGP7GU,BGPD2
                  BGPD22
  $$FMTE^XLFDT
BGP3D1,BGP3D10,BGP3D2,BGP3D3,BGP3D4,BGP3D5,BGP3D6,BGP3D7
                  BGP3D8,BGP3D9,BGP3DA,BGP3DAR,BGP3DARA,BGP3DARH,BGP3DG
                  BGP3DH,BGP3DL,BGP3DPH,BGP3DS,BGP3PDL,BGP3PDLs,BGP4D1
                  BGP4D10,BGP4D2,BGP4D21,BGP4D3,BGP4D31,BGP4D4,BGP4D41
                  BGP4D5,BGP4D6
                  BGP4D61,BGP4D7,BGP4D71,BGP4D8,BGP4D9,BGP4DA,BGP4DAP
                  BGP4DAR,BGP4DARA,BGP4DG,BGP4DH,BGP4DHE,BGP4DHEA,BGP4DL
BGP4DPH,BGP4DPP,BGP4DS,BGP4HEH,BGP4HEHH,BGP4HEL,BGP4HES
                  BGP4HESL,BGP4PDL,BGP4PDLH,BGP4PDLs,BGP5C1,BGP5C11
                  BGP5C13,BGP5C2,BGP5C3
                  BGP5CL,BGP5CP,BGP5D1,BGP5D10,BGP5D2,BGP5D21,BGP5D22
                  BGP5D3,BGP5D31,BGP5D4,BGP5D41,BGP5D5,BGP5D6,BGP5D61
                  BGP5D7,BGP5D71,BGP5D72,BGP5D8,BGP5D9,BGP5DAP,BGP5DAR
                  BGP5DEL,BGP5DELA,BGP5DGPA,BGP5DGPU,BGP5DH,BGP5DHE
                  BGP5DHE1,BGP5DHEA,BGP5DL
                  BGP5DNG,BGP5DPH,BGP5DPP,BGP5DS,BGP5EL2,BGP5EL3,BGP5EL4
BGP5ELH,BGP5ELHH,BGP5ELL,BGP5ELS,BGP5ELSL,BGP5GU,BGP5HEH
                  BGP5HEHH,BGP5HEL,BGP5HES,BGP5HESL,BGP5NPL,BGP5NPLD
                  BGP5NPLP,BGP5PDL,BGP5PDLH,BGP5PDLs,BGP6AHW,BGP6C1
                  BGP6C11,BGP6C12,BGP6C13
                  BGP6C2,BGP6C3,BGP6C31,BGP6CL,BGP6CP,BGP6D1,BGP6D10
                  BGP6D2,BGP6D21,BGP6D22,BGP6D23,BGP6D3,BGP6D31,BGP6D4
                  BGP6D41,BGP6D5,BGP6D51,BGP6D52,BGP6D6,BGP6D61,BGP6D7
                  BGP6D71,BGP6D72,BGP6D8,BGP6D81,BGP6D82,BGP6D83,BGP6D84
                  BGP6D861,BGP6D863,BGP6D9
                  BGP6DAP,BGP6DAR,BGP6DCHA,BGP6DCHW,BGP6DCL,BGP6DCLD
                  BGP6DCLP,BGP6DEL,BGP6DELA,BGP6DGPA,BGP6DGPU,BGP6DH
BGP6DHE,BGP6DHE1,BGP6DHEA,BGP6DL,BGP6DLT,BGP6DMT,BGP6DNG
                  BGP6DPH,BGP6DPP,BGP6DS,BGP6DSTM,BGP6DU,BGP6EL2,BGP6EL3
                  BGP6EL4,BGP6ELH,BGP6ELHH

```

BGP6ELL, BGP6ELS, BGP6ELSL, BGP6GU, BGP6HEH, BGP6HEHH, BGP6HEL
BGP6HES, BGP6HESL, BGP6NPL, BGP6NPLD, BGP6NPLP, BGP6PDL

BGP6PDLH, BGP6PDLs, BGP7C1, BGP7C11, BGP7C12, BGP7C13, BGP7C2
BGP7C3, BGP7C31, BGP7CL, BGP7CP, BGP7D1, BGP7D10, BGP7D2
BGP7D21, BGP7D211

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 478

Routine	is Invoked by:
\$\$FMTE^XLFDT	BGP7D21A, BGP7D22, BGP7D23, BGP7D3, BGP7D31, BGP7D4, BGP7D41 BGP7D5, BGP7D51, BGP7D52, BGP7D6, BGP7D61, BGP7D7, BGP7D71 BGP7D72, BGP7D721, BGP7D722, BGP7D731, BGP7D732, BGP7D74 BGP7D75, BGP7D8, BGP7D81, BGP7D9, BGP7DAP, BGP7DAR, BGP7DCHA BGP7DCHW, BGP7DCL BGP7DCLD, BGP7DCLP, BGP7DEL, BGP7DELA, BGP7DGPA, BGP7DGPU
BGP7DH, BGP7DHE, BGP7DHE1, BGP7DHE2, BGP7DHEA, BGP7DL, BGP7DLT	BGP7DMT, BGP7DNE, BGP7DNG, BGP7DPE, BGP7DPE1, BGP7DPE2 BGP7DPEA, BGP7DPEL, BGP7DPEP, BGP7DPH, BGP7DPP, BGP7DS BGP7DSTM, BGP7DU, BGP7EL2 BGP7EL3, BGP7EL4, BGP7ELH, BGP7ELHH, BGP7ELL, BGP7ELS BGP7ELSL, BGP7GU, BGP7HEH, BGP7HEHH, BGP7HEL, BGP7HES BGP7HESL, BGP7NPL, BGP7NPLD, BGP7NPLP, BGP7PDL, BGP7PDLH BGP7PDLs, BGPAPH, BGPD, BGPD1, BGPD12, BGPD13, BGPD14, BGPD1B BGPD2, BGPD22, BGPD24, BGPD29 BGPD3, BGPD31, BGPD4, BGPD5, BGPD6, BGPD7, BGPDA, BGPDAREA BGPDARP, BGPDB, BGPDD, BGPDFTA, BGPDH, BGPDPH, BGPDS, BGPDU BPGU
\$\$NOW^XLFDT	BGP3UTL, BGP4DA, BGP4DAR, BGP4DARA, BGP4DG, BGP5DAR, BGP5DNG BGP5GAEL, BGP5GAGP, BGP5GAHE, BGP5GALL, BGP5GAPU, BGP5GCMS BGP5GCOM, BGP5GELD, BGP5GGPU, BGP5GHED, BGP5GLST, BGP5GNPL BGP5GNTL, BGP5GPP, BGP5GUPL, BGP6DAR, BGP6DCHA, BGP6DNG BGP6GACW, BGP6GAEL BGP6GAGP, BGP6GAHE, BGP6GALL, BGP6GAPU, BGP6GCMP, BGP6GCMS BGP6GCOM, BGP6GELD, BGP6GGPU, BGP6GHED, BGP6GLST, BGP6GLTX BGP6GMTX, BGP6GNPL, BGP6GNST, BGP6GNTL, BGP6GPP, BGP6GUPL BGP7DAR, BGP7DCHA, BGP7DNE, BGP7DNG, BGP7GACW, BGP7GAEL BGP7GAGP, BGP7GAHE, BGP7GALL BGP7GAPE, BGP7GAPU, BGP7GCMP, BGP7GCMS, BGP7GCOM, BGP7GELD BGP7GGPU, BGP7GHED, BGP7GLST, BGP7GLTX, BGP7GMTX, BGP7GNPL BGP7GNST, BGP7GNTL, BGP7GNXP, BGP7GPED, BGP7GPP, BGP7GUPL BGP7D74
\$\$ABS^XLFMTH	BGP7D74
\$\$CJ^XLFSTR	BGP7POS, BGPPOST
\$\$REPEAT^XLFSTR	BGP6DLT, BGP6DMT, BGP6DSTM, BGP7DLT, BGP7DMT, BGP7DSTM
\$\$STRIP^XLFSTR	BGP4DP1, BGP4HEL1, BGP4HEP1, BGP4PDL3, BGP5DP1, BGP5ELL1 BGP5ELLN, BGP5ELP1, BGP5HEL1, BGP5HEP1, BGP5PDL3, BGP6D23 BGP6DCHW, BGP6DP1C, BGP6EL2, BGP6ELL1, BGP6ELLN, BGP6ELP1 BGP6HEL1, BGP6HEP1, BGP6PDL3, BGP7D23, BGP7DCHW, BGP7DP1C BGP7DPED, BGP7DPEP, BGP7EL2 BGP7ELL1, BGP7ELLN, BGP7ELP1, BGP7HEL1, BGP7HEP1, BGP7PDL3
\$\$SUP^XLFSTR	BGP3D2, BGP3ULF, BGP4D2, BGP4ULF, BGP5C1, BGP5C11, BGP5C12 BGP5C13, BGP5C2, BGP5C3, BGP5D2, BGP5D72, BGP5EL2, BGP5ULF BGP6C1, BGP6C11, BGP6C12, BGP6C13, BGP6C2, BGP6C3, BGP6D2 BGP6D72, BGP6EL2, BGP6ULF, BGP7C1, BGP7C11, BGP7C12, BGP7C13 BGP7C2, BGP7C3, BGP7D2 BGP7D72, BGP7D722, BGP7D731, BGP7D732, BGP7D75, BGP7EL2 BGP7ULF, BGPD2, BGPD4, dd90374.01, dd90530.01
^XMD	BGP50P1
<<<< 2.2.2.12 NOT MARKED AS ENTRY POINT >>>>	
MES^XPDUTL	BGP7POS, BGPPOST

***** Cross Reference of all Routines ***** DEC 14, 2006@06:51:55 page 479

Routine is Invoked by:
 DISP^XQORM1
 BGP3DSI, BGP3DSL, BGP4DHSL, BGP4DSI, BGP4DSL, BGP5ASL, BGP5CSI
 BGP5CTL, BGP5CTS, BGP5DESI, BGP5DESL, BGP5DHSL, BGP5DSI
 BGP5DSL, BGP5ETL, BGP5ETS, BGP5TXL, BGP5TXS, BGP5TXV, BGP5TXV1
 BGP6ASL, BGP6CSI, BGP6CTL, BGP6CTS, BGP6DESI, BGP6DESL
 BGP6DHSL, BGP6DSI
 BGP6DSL, BGP6ETL, BGP6XTL, BGP6XTS, BGP6XTV, BGP6XTV1, BGP7ASL
 BGP7CSI, BGP7CTL, BGP7CTS, BGP7DESI, BGP7DESL, BGP7DHSL
 BGP7DPSI, BGP7DSI, BGP7DSL, BGP7ETL, BGP7PESI, BGP7PESL
 BGP7XTL, BGP7XTS, BGP7XTV, BGP7XTV1, BGP7D0, BGP7D1
 \$\$\$NOJOURN^ZIBGCHAR BGP3D1, BGP4HUTL, BGP4UTL, BGP5EUTL, BGP5HUTL, BGP5UTL
 BGP6EUTL, BGP6HUTL, BGP6UTL, BGP7DPE1, BGP7EUTL, BGP7HUTL
 BGP7PUTL, BGP7UTL, BGP7D1

6.2 Callable Routines

There are no published entry points in this package.

6.3 Exported Options

Option Name	Description
BGP 02 MENU	GPRA+ FY02
BGP 03 AREA AREA	Run AREA Annual Clinical Report (Area Only)
BGP 03 AREA GPRA	Run AREA GPRA Report (to be used at Area only)
BGP 03 AREA MENU	Area Options
BGP 03 EDIT SITE PARAMETERS	Site Parameters
BGP 03 MAIN MENU	GPRA+ FY03
BGP 03 MENU REPORTS	Reports
BGP 03 MENU SETUP	System Setup
BGP 03 RPT AREA PERF	Run Area Annual Performance Report for Area Export
BGP 03 RPT GPRA AREA EXPT	Run GPRA Report for Area Export
BGP 03 RPT LOCAL USE	Run Indicator Reports for Local Use
BGP 03 TAXONOMY CHECK	Check for Taxonomies Required by the GPRA Report
BGP 03 TAXONOMY SETUP	Taxonomy Setup
BGP 03 UPLOAD TO AREA	Upload Report Files from Site
BGP 04 AREA AREA	Run AREA Area Director Report (Area Only)
BGP 04 AREA DIR REPORT	Area Director's Performance Report

Option Name	Description
BGP 04 AREA GPRA	Run AREA GPRA Report
BGP 04 AREA HEDIS	Run Area HEDIS Report
BGP 04 AREA MENU	Area Office Options
BGP 04 HEDIS AREA REPORT	Run AREA HEDIS Report
BGP 04 HEDIS REPORT	HEDIS Performance Report
BGP 04 LIST FILES	List files in a directory
BGP 04 MAIN MENU	GPRA+ FY04
BGP 04 REPORT GPRA	GPRA Performance Report
BGP 04 REPORTS MENU	Reports
BGP 04 SEL ALL PATS	Selected Indicators with No Community Defined
BGP 04 SEL PATIENT PANEL	Selected Indicators w/Patient Panel Population
BGP 04 SELECTED IND REPORT	Selected Indicators w/Community Taxonomy Defined
BGP 04 SITE PARAMETERS	Site Parameters
BGP 04 SYSTEM SETUP	System Setup
BGP 04 TAXONOMY CHECK	Check for Taxonomies Required by FY 04 GPRA Report
BGP 04 TAXONOMY SETUP	Taxonomy Setup
BGP 04 UPLOAD FILES	Upload Report Files from Site
BGP 05 AREA ELDER REPORT	Run AREA Elder Care Report
BGP 05 AREA GPRA	Run AREA National GPRA Report
BGP 05 AREA GPU	Run AREA GPRA Performance Report
BGP 05 AREA HEDIS	Run Area HEDIS Report
BGP 05 AREA MENU	Area Options
BGP 05 CMS PERF RPT	CMS Performance Report
BGP 05 CMS TAXONOMY CHECK	Taxonomy Check-CMS Report
BGP 05 CMS TAXONOMY SETUP	Taxonomy Setup-CMS Report
BGP 05 CRS TAXONOMY SETUP	Taxonomy Setup-All CRS Reports
BGP 05 ELDER REPORT	Elder Care Report
BGP 05 ELDER TAXONOMY CHECK	Taxonomy Check-Elder Care Report
BGP 05 GPU GPRA PERFORMANCE	GPRA Performance Report
BGP 05 HEDIS AREA REPORT	Run AREA HEDIS Report
BGP 05 HEDIS REPORT	HEDIS Performance Report

Option Name	Description
BGP 05 HEDIS TAXONOMY CHECK	Taxonomy Check-HEDIS Report
BGP 05 HOSPITAL REPORTS	Hospital Reports
BGP 05 LIST FILES	List files in a directory
BGP 05 LOCAL REPORTS MENU	Reports for Local Use: IHS Clinical Indicators
BGP 05 LOCAL TAXONOMY CHECK	Taxonomy Check-Local CRS Reports
BGP 05 MAIN MENU	CRS 2005
BGP 05 MENU NATIONAL	National GPRA Reports
BGP 05 NATIONAL GPRA REPORT	National GPRA Report
BGP 05 NATIONAL PAT LISTS	National GPRA Report Patient List
BGP 05 NGR TAXONOMY CHECK	Taxonomy Check-National GPRA/GPRA Performance Rpts
BGP 05 NGR TAXONOMY SETUP	Taxonomy Setup-National GPRA/GPRA Performance Rpts
BGP 05 OTHER CMS REPORTS	Other CMS Reports
BGP 05 OTHER NTL RPTS MENU	Other National Reports
BGP 05 REPORTS MENU	Reports
BGP 05 SEL ALL PATS	Selected Indicators w/All Communities
BGP 05 SEL PATIENT PANEL	Selected Indicators w/Patient Panel Population
BGP 05 SELECTED IND REPORT	Selected Indicators w/Community Specified
BGP 05 SITE PARAMETERS	Site Parameters
BGP 05 SYSTEM SETUP	System Setup
BGP 05 TAXONOMY CHECK	Taxonomy Check
BGP 05 TAXONOMY SETUP	Taxonomy Setup
BGP 05 UPLOAD FILES	Upload Report Files from Site
BGP 05 VIEW ANY TAXONOMY	View All CRS Taxonomies (Site-Defined&Hard Coded)
BGP 06 AREA CHILDHOOD WT/HT	Run AREA Childhood Height and Weight Data File
BGP 06 AREA ELDER REPORT	Run AREA Elder Care Report
BGP 06 AREA GPRA	Run AREA National GPRA Report
BGP 06 AREA GPU	Run AREA GPRA Performance Report
BGP 06 AREA HEDIS	Run Area HEDIS Report

Option Name	Description
BGP 06 AREA MENU	Area Options
BGP 06 CMS PERF RPT	CMS Performance Report
BGP 06 CMS TAXONOMY CHECK	Taxonomy Check-CMS Report
BGP 06 CMS TAXONOMY SETUP	Taxonomy Setup-CMS Report
BGP 06 COMPREHENSIVE PAT LIST	Comprehensive National GPRA Patient List
BGP 06 CRS TAXONOMY SETUP	Taxonomy Setup-All CRS Reports
BGP 06 ELDER REPORT	Elder Care Report
BGP 06 ELDER TAXONOMY CHECK	Taxonomy Check-Elder Care Report
BGP 06 GPU GPRA PERFORMANCE	GPRA Performance Report
BGP 06 HEDIS AREA REPORT	Run AREA HEDIS Report
BGP 06 HEDIS REPORT	HEDIS Performance Report
BGP 06 HEDIS TAXONOMY CHECK	Taxonomy Check-HEDIS Report
BGP 06 HOSPITAL REPORTS	Hospital Reports
BGP 06 LAB TAX REPORT	Lab Taxonomy Report
BGP 06 LIST FILES	List files in a directory
BGP 06 LOCAL REPORTS MENU	Reports for Local Use: IHS Clinical Measures
BGP 06 LOCAL TAXONOMY CHECK	Taxonomy Check-Local CRS Reports
BGP 06 MAIN MENU	CRS 2006
BGP 06 MED TAX REPORT	Medication Taxonomy Report
BGP 06 MENU NATIONAL	National GPRA Reports
BGP 06 NATIONAL GPRA REPORT	National GPRA Report
BGP 06 NPL SEARCH TEMPLATE	Create Search Template for National Patient Lists
BGP 06 NATIONAL PAT LISTS	National GPRA Report Patient List
BGP 06 NGR TAXONOMY CHECK	Taxonomy Check-National GPRA/GPRA Performance Rpts
BGP 06 NGR TAXONOMY SETUP	Taxonomy Setup-National GPRA/GPRA Performance Rpts
BGP 06 OTHER CMS REPORTS	Other CMS Reports

Option Name	Description
BGP 06 OTHER NTL RPTS MENU	Other National Reports
BGP 06 REPORTS MENU	Reports
BGP 06 SEL ALL PATS	Selected Measures w/All Communities
BGP 06 SEL PATIENT PANEL	Selected Measures w/Patient Panel Population
BGP 06 SELECTED IND REPORT	Selected Measures w/Community Specified
BGP 06 SITE PARAMETERS	Site Parameters
BGP 06 SYSTEM SETUP	System Setup
BGP 06 TAXONOMY CHECK	Taxonomy Check
BGP 06 TAXONOMY REPORTS	Taxonomy Reports
BGP 06 TAXONOMY SETUP	Taxonomy Setup
BGP 06 UPLOAD FILES	Upload Report Files from Site
BGP 06 VIEW ANY TAXONOMY	View All CRS Taxonomies (Site-Defined&Hard Coded)
BGP FAC TO AREA RPT PERIOD	Run Area Annual Report for Export to Area
BGP GPRA AREA REPORT	Run AREA GPRA Report (to be used at Area only)
BGP GPRA FAC TO AREA RPT	Run FY Report for Export to Area
BGP GPRA LOCAL REPORT	Run GPRA Report for Local Use
BGP GPRA UPLOAD FILE	Upload GPRA Data File from Site
BGP INDICATOR 32 TESTS	Enter/Edit Indicator 32 Lab Tests
BGP TAXONOMY CHECK	Check for Taxonomies Required by the GRPA Report
BGPGRPC	BGPG Application Context
BGP 07 AREA CHILDHOOD WT/HT	Run AREA Height and Weight Data File
BGP 07 AREA ELDER REPORT	Run AREA Elder Care Report
BGP 07 AREA GPRA	Run AREA National GPRA Report
BGP 07 AREA GPU	Run AREA GPRA Performance Report
BGP 07 AREA HEDIS	Run Area HEDIS Report
BGP 07 AREA MENU	Area Options
BGP 07 AREA PAT ED RPT	Run AREA Patient Education Report
BGP 07 CMS PERF RPT	CMS Performance Report
BGP 07 CMS TAXONOMY CHECK	Taxonomy Check-CMS Report
BGP 07 CMS TAXONOMY SETUP	Taxonomy Setup-CMS Report

Option Name	Description
BGP 07 COMP NATL EXPORT	Comprehensive National GPRA Export
BGP 07 COMPREHENSIVE PAT LIST	Comprehensive National GPRA Patient List
BGP 07 CRS TAXONOMY SETUP	Taxonomy Setup-All CRS Reports
BGP 07 ELDER REPORT	Elder Care Report
BGP 07 ELDER TAXONOMY CHECK	Taxonomy Check-Elder Care Report
BGP 07 GPU GPRA PERFORMANCE	GPRA Performance Report
BGP 07 HEDIS AREA REPORT	Run AREA HEDIS Report
BGP 07 HEDIS REPORT	HEDIS Performance Report
BGP 07 HEDIS TAXONOMY CHECK	Taxonomy Check-HEDIS Report
BGP 07 HOSPITAL REPORTS	Hospital Reports
BGP 07 LAB TAX REPORT	Lab Taxonomy Report
BGP 07 LIST FILES	List files in a directory
BGP 07 LOCAL REPORTS MENU	Reports for Local Use: IHS Clinical Measures
BGP MENU	IHS Clinical Reporting System (CRS) Main Menu

7.0 Internal Relations

All users should be given the access to the appropriate options and keys to them, as needed. All of the options in this system stand alone.

8.0 How to Generate Online Documentation

The file number range for this package is 90240 and 90370. The namespace is BGP. All templates, routines, screen forms, etc. begin with BGP.

This section describes some of the methods by which users can generate IHS RPMS CLINICAL Reporting system technical documentation. Online technical documentation pertaining to the IHS RPMS CLINICAL REPORTING reporting software, in addition to that which is located in the help prompts and on the help screens throughout the IHS RPMS CLINICAL Reporting package, can be generated through the use of several Kernel options. These include, but are not limited to, the following:

- %INDEX
- Menu Management
- Inquire Option
- Print Option File
- VA FileMan
- Data Dictionary Utilities
- List File Attributes

Entering question marks at the “Select...Option” prompts can also provide users with valuable technical information. For example, a single question mark (?) lists all options that can be accessed from the current option. Entering two question marks (??) lists all options accessible from the current one, showing the formal name and lock for each. Three question marks (???) displays a brief description for each option in a menu, whereas an option name preceded by a question mark (?OPTION) shows extended help, if available, for that option.

For a more exhaustive option listing and further information about other utilities that supply online technical information, please consult the DHCP Kernel Reference manual.

8.1 %INDEX

This option analyzes the structure of a routine to determine in part if the routine adheres to RPMS Programming Standards. The %INDEX output can include the following components:

- Compiled list of errors and warnings
- Routine listing
- Local variables
- Global variables
- Naked globals

- Label references
- External references

By running %INDEX for a specified set of routines, you are afforded the opportunity to discover any deviations from RPMS Programming Standards that exist in the selected routines and to see how routines interact with one another (i.e., which routines call or are called by other routines).

To run %INDEX for the IHS RPMS CLINICAL Reporting package, specify the BGP namespace at the Routine(s)?> prompt.

8.2 Inquire Option

This menu management option provides the following information about a specified option:

- Option name
- Menu text
- Option description
- Type of option
- Lock (if any)

In addition, all items on the menu are listed for each menu option. To secure information about IHS RPMS CLINICAL REPORTING Reporting options, you must specify the BGP namespace.

8.3 Print Option File

This utility generates a listing of options from the Option file (#19). You can choose to print all of the entries in this file or you can specify a single option or range of options. For a list of IHS RPMS CLINICAL Reporting options, please refer to the Exported Options section of this manual.

8.4 List File Attributes

This VA FileMan option allows you to generate documentation pertaining to files and file structure. Using the Standard format of this option yields the following data dictionary information for a specified file:

- File name and description
- Identifiers
- Cross-references
- Files pointed to by the file specified
- Files that point to the file specified input, print, and sort templates

In addition, the following applicable data is supplied for each field in the file:

- Field name, number, title, and description
- Global location
- Help prompt
- Cross-references
- Input transform
- Date last edited
- Notes

Using the Global Map format of this option generates an output that lists the following information:

- All cross-references for the file selected
- Global location of each field in the file
- Input, print, and sort templates

For a comprehensive listing of CLINICAL package files, please refer to the files section of this manual.

9.0 SAC Exemptions

This program uses a standard global for export of data to the Area office. The global is called ^BGPDATA(and is killed at the root level in various routines.

The kill of these unsubscripted globals has been approved by the Standards and Conventions committee.

10.0 Glossary

Term	Definition
Archiving	The storing of historical or little-used data off-line (often on tape).
Banner	A line of text with a user's name and domain.
Browser	An interactive application that displays ASCII text on a terminal that supports a scroll region. The text can be in the form of a word-processing field or sequential local or global array. The user is allowed to navigate freely within the document.
Callable Entry Points	Places in a routine that can be called from an application program.
Cross-reference	An indexing method whereby files can include pre-sorted lists of entries as part of the stored database. Cross-references (x-refs) facilitate look-up and reporting.
Entry Point	Entry point within a routine that is referenced by a "DO" or "GOTO" command from a routine internal to a package.
File	A set of related records or entries treated as a single unit.
FileMan	The database management system for RPMS.
Global	In MUMPS, global refers to a variable stored on disk (global variable) or the array to which the global variable may belong (global array).
ICD	International Classification of Diseases.
INDEX (%INDEX)	A Kernel utility used to verify routines and other MUMPS code associated with a package. Checking is done according to current ANSI MUMPS standards and RPMS programming standards. This tool can be invoked through an option or from direct mode (>D ^%INDEX).
IRM	Information Resource Management. The IHS personnel responsible for information systems management and security.
Kernel	The set of MUMPS software utilities that function as an intermediary between the host operating system and application packages, such as Laboratory and Pharmacy. The Kernel provides a standard and consistent user and programmer interface between application packages and the underlying MUMPS implementation. These utilities provide the foundation for RPMS.

Menu	A list of choices for computing activity. A menu is a type of option designed to identify a series of items (other options) for presentation to the user for selection. When displayed, menu-type options are preceded by the word “Select” and followed by the word “option” as in Select Menu Management option: (the menu’s select prompt).
Namespace	A unique set of 2 to 4 alpha characters that are assigned by the database administrator to a software application.
Option	An entry in the Option file. As an item on a menu, an option provides an opportunity for users to select it, thereby invoking the associated computing activity. Options may also be scheduled to run in the background, non-interactively, by TaskMan.
Patient Care Component (PCC)	The central repository for data in the Resource and Patient Management System (RPMS).
Queuing	Requesting that a job be processed at a later time rather than within the current session.
Routine	A program or sequence of instructions called by a program that may have some general or frequent use. MUMPS routines are groups of program lines that are saved, loaded, and called as a single unit via a specific name.
UCI	User Class Identification: a computing area.
Up-Hat (^)	A circumflex, also know as a “hat” or “caret,” that is used as a piece delimiter in a global. The up-hat is denoted as “^” and is typed by pressing Shift+6 on the keyboard.
Utility	A callable routine line tag or function. A universal routine usable by anyone.
Variable	A character or group of characters that refers to a value. MUMPS recognizes 3 types of variables: local variables, global variables, and special variables. Local variables exist in a partition of the main memory and disappear at sign-off. A global variable is stored on disk, potentially available to any user. Global variables usually exist as parts of global arrays.

11.0 Contact Information

If you have any questions or comments regarding this distribution, please contact the OIT Help Desk by:

Phone: (505) 248-4371 or
(888) 830-7280

Fax: (505) 248-4363

Web: <http://www.ihs.gov/GeneralWeb/HelpCenter/Helpdesk/index.cfm>

Email: ITSCHelp@mail.ihs.gov