



RESOURCE AND PATIENT MANAGEMENT SYSTEM

HIV Management System (HMS) (BKM)

Installation Guide and Release Notes

Version 1.0
September 2006

Office of Information Technology (OIT)
Division of Information Resources
Albuquerque, New Mexico

PREFACE

This manual provides information regarding only the release and installation of the HIV Management System (HMS) version 1.0 package.

TABLE OF CONTENTS

1.0	RELEASE NOTES	3
1.1	Product Overview	3
2.0	INSTALLATION NOTES.....	4
2.1	General Information	4
2.2	Contents of Distribution.....	4
2.3	Required Resources	4
3.0	INSTALLATION ISSUES	5
4.0	INSTALLATION INSTRUCTIONS	6
4.1	Pre – Installation	6
4.2	Installation Steps	6
5.0	SAMPLE INSTALL	8
6.0	INSTALLATION CONFIGURATION.....	13
7.0	POST-INSTALLATION INSTRUCTIONS	14
8.0	CONTACT INFORMATION.....	18
	APPENDIX A – ROUTINE CHECKSUMS	19

1.0 Release Notes

Version 1.0 is the first release of this new application. There are no modifications or enhancements to any existing applications. The intended audiences for this document are IHS Site Managers and other technical staff. Technical documentation accompanies this installation document in this release. Please refer to the HMS Technical Manual v1.0 for additional information.

This application introduces a new term for the application owner - Case File Manager, which is essentially the same as the ADPAC. The Site Manager should coordinate this installation with the Case File Manager regarding user roles and security keys, and scheduled tasks, which will be discussed later in this document.

This application is similar to existing case management modules in RPMS, such as Diabetes Management. However, due to the highly sensitive nature of healthcare information regarding HIV and related illnesses, HMS has additional functionality to protect the security of that information.

1.1 Product Overview

The HIV Management System software will be a component of the IHS Resource and Patient Management System (RPMS) that provides case management functions. Case management provides a mechanism for healthcare providers to be informed about and to follow up on clinical interventions and information on single or multiple disease states and conditions. These applications can assist healthcare providers to identify high-risk patients, proactively track care reminders and the health status of individuals or populations, provide appropriate care by embedding evidence-based guidelines, and report outcomes. Appropriate case management improves the delivery of care, as well as outcomes. In addition, there are documented fiscal benefits when appropriate case management is introduced into the health care business process model.

HMS is a key tool for accomplishing the IHS goal to improve compliance with clinical practice guidelines for five chronic diseases that are significant to the American Indian/Alaska Native (AI/AN) population: diabetes, HIV/AIDS, asthma, obesity, and cardiovascular disease. It provides the capability to establish and maintain a patient register for managing the health care of any select groupings of HIV/AIDS and/or At Risk patients defined by a facility's providers. It also provides the ability to collect data not currently available within RPMS that is specific to the health care and management of HIV/AIDS patient and provide relevant reminders and alerts to providers. Finally, HMS provides the capability to analyze and report patient data for individuals and populations.

2.0 Installation Notes

Prefix: BKM

Current Version: 1.0

Warning: Read the entire notes file (if included) prior to attempting any installation.

2.1 General Information

- Make a copy of this distribution for offline storage.
- Print all notes/readme files.
- It is recommended that the terminal output during the installation be captured using an auxport printer attached to the terminal at which you are performing the software installation. This will insure a printed audit trail if any problems should arise.

2.2 Contents of Distribution

- bkm_0100.k
- bkm_010i.pdf
- bkm_010t.pdf
- bkm_010u.pdf

2.3 Required Resources

Package	Minimum Version
VA FileMan (DI)	v22.0
VA Kernel (XU)	v8.0
VA List Manager (VALM)	v1.0 through patch 2
IHS/VA Utilities (XB)	v3.0 through patch 11
Adverse Reaction Tracking (GMRA)	v4.0
IHS Clinical Reporting (BGP)	v6.0
Patient Information Management System (PIMS)	v5.3 through patch 1004
PCC Data Entry (APCD)	v2.0 through patch 8
PCC Health Summary (APCH)	v2.0 through patch 13
Q-Man (AMQQ)	v2.0 through patch 18
Taxonomy (ATX)	v5.1 through patch 8

3.0 Installation Issues

Before installing this package, your system should be backed up.

Since this is a first release of an entirely new component, users may be logged into RPMS when this installation is run. The actual software installation should take about fifteen minutes.

4.0 Installation Instructions

4.1 Pre – Installation

Initially load the software into a test account, then into the production account. There are no routines or globals to remove after the installation is complete.

Refer to the Installation Steps (Section 4.2) and the Sample Installation (Section 5.0) for the applicable prompts and user responses. Prompts requiring responses appear in **bold** type.

Download the **bkm_0100.k** file from the anonymous directory.

4.2 Installation Steps

1. From the Kernel Installation and Distribution System Menu, select the Installation Menu.
2. From the Installation Menu select option **1**, Load a Distribution.
3. Enter **bkm_0100.k** at the “Enter a Host File:” prompt. You may need to precede the file name with the appropriate host path. This is the file that you retrieve from the appropriate software directory.
4. Enter **YES** (or press Enter to accept the default) at the “Want to Continue with Load? YES//” prompt.
5. OPTIONAL - At this time you may elect to exercise one of the following installation options. Please note these are OPTIONAL. If you elect to use one of these, use HIV MANAGEMENT SYSTEM 1.0 as the INSTALL NAME.
 - a. Backup a Transport Global – this option creates a MailMan message that will back up all current routines on your system that would be replaced by this release. (Since this is a new release, no files will be replaced.)
 - b. Compare Transport Global – this option allows you to view all changes that will result from the installation of this patch and compares these with the values currently loaded on your system (routines, DD’s, templates, etc.).
 - c. Verify Checksums in Transport Global – This option provides you with a mechanism for ensuring the integrity of your routines. It verifies checksums for the components of the Transport Global and reports any errors uncovered.
6. From the Installation Menu, select the **Install Packages(s)** option.

7. At the “Select Install Name:” prompt, enter **HIV MANAGEMENT SYSTEM 1.0**.
8. Enter the username associated with the Coordinator for Mail Group BKMI Alert.
9. Respond to the following prompts as indicated:
 - a. “Want KIDS to Rebuild Menu Trees Upon Completion of Install? YES//,” answer **YES**.
 - b. “Select the User who will be initialized as the HMS Case File Manager. This is the user who will set up all the other HMS users through the User Account Setup [BKMV USER MAINTENANCE] option.” Answer with the name of the user being designated as the HMS Case File Manager. This is required. Failure to complete will cause the installation to be terminated.
 - c. “Want KIDS to INHIBIT LOGONs during the install?//YES,” answer **NO**.
 - d. “Want to DISABLE Scheduled Options, Menu Options and Protocols?// YES,” answer **NO**.
10. Enter the Device you want to print the Install messages. You can queue the install by entering a “**Q**” at the device prompt. Enter a “**^**” to abort the install.
11. Exit the Installation Menu.
12. Proceed to Section 6.0 System Configuration, for further information and setup instructions.

5.0 Sample Install

Figure 5-1, which follows, depicts a sample installation.

Note: Prompts requiring user response appear in bold type.

```
Select OPTION NAME: XPD MAIN           Kernel Installation & Distribution System

      Edits and Distribution ...
      Utilities ...
      Installation ...

Select Kernel Installation & Distribution System Option: Installation

1      Load a Distribution
2      Verify Checksums in Transport Global
3      Print Transport Global
4      Compare Transport Global to Current System
5      Backup a Transport Global
6      Install Package(s)
      Restart Install of Package(s)
      Unload a Distribution

Select Installation Option: Load a Distribution

Enter a Host File: d:\rpms_tst\bkm_0100.k (Note: Enter the appropriate path
for your system.)

KIDS Distribution saved on Dec 08, 2005@17:35:08
Comment: HIV MANAGEMENT SYSTEM TEST V1.0

This Distribution contains Transport Globals for the following Package(s):
  HIV MANAGEMENT SYSTEM 1.0
Distribution OK!

Want to Continue with Load? YES// YES
Loading Distribution...

  HIV MANAGEMENT SYSTEM 1.0

Use INSTALL NAME: HIV MANAGEMENT SYSTEM 1.0 to install this Distribution.

1      Load a Distribution
2      Verify Checksums in Transport Global
3      Print Transport Global
4      Compare Transport Global to Current System
5      Backup a Transport Global
6      Install Package(s)
      Restart Install of Package(s)
      Unload a Distribution
```

```
Select Installation Option: Install Package(s)
Select INSTALL NAME: HIV MANAGEMENT SYSTEM 1.0           Loaded from Distribution
Loaded from Distribution 4/25/06@14:54:37
=> HMS Version 1.0 Build#101 ;Created on Apr 25, 2006@14:52:07

This Distribution was loaded on Apr 25, 2006@14:54:37 with header of
HMS Version 1.0 Build#101 ;Created on Apr 25, 2006@14:52:07
It consisted of the following Install(s):
HIV MANAGEMENT SYSTEM 1.0
Checking Install for Package HIV MANAGEMENT SYSTEM 1.0

Install Questions for HIV MANAGEMENT SYSTEM 1.0

Incoming Files:

90450      HMS BASE
90451      HMS REGISTRY
90451.1    HMS PCC ENTRY UPDATE
90451.2    HMS CANDIDATE FILE
90451.5    HMS ETIOLOGY (including data)
90451.7    HMS CLINICAL CLASSIFICATION (including data)
90452      HMS DOCUMENTS
90452.5    HMS ORDER FREQUENCY (including data)
90453      HMS PROVIDER GUIDELINES (including data)
90455      HMS AUDIT
90456      HMS HIV/AIDS STATE REPORT (including data)
90459      HMS PCC BUFFER

Incoming Mail Groups:

Enter the Coordinator for Mail Group 'BKMI ALERT': Enter Name Here (choose the
appropriate name for your site)

Want KIDS to Rebuild Menu Trees Upon Completion of Install? YES// YES

Select the User who will be initialized as the HMS Case File Manager.
This is the user who will set up all the other HMS users through the
User Account Setup [BKMV USER MAINTENANCE] option.
Select User Name: Enter Appropriate User Name Here. Failure to identify an HMS
Case File Manager will terminate the installation.

Want KIDS to INHIBIT LOGONs during the install? YES// NO
Want to DISABLE Scheduled Options, Menu Options, and Protocols? YES// NO
```

Enter the Device you want to print the Install messages.
You can queue the install by enter a 'Q' at the device prompt.
Enter a '^' to abort the install.

DEVICE: HOME// **VIRTUAL**

Install Started for HIV MANAGEMENT SYSTEM 1.0 :
Apr 25, 2006@14:55:32

Build Distribution Date: Apr 25, 2006

Installing Routines:
Apr 25, 2006@14:55:34

Running Pre-Install Routine: ^BKM1PRE

Installing Data Dictionaries:
Apr 25, 2006@14:55:35

Installing Data:
Apr 25, 2006@14:55:35

Installing PACKAGE COMPONENTS:

Installing BULLETIN

Installing SECURITY KEY

Installing PRINT TEMPLATE

Installing SORT TEMPLATE

Installing INPUT TEMPLATE

Installing FORM

Installing MAIL GROUP

Installing PROTOCOL

Installing LIST TEMPLATE

Installing OPTION
Apr 25, 2006@14:55:41

Running Post-Install Routine: ^BKM1POST

Creating [BKM TB MEDS] taxonomy...
Creating [BKMV FI MEDS] taxonomy...
Creating [BKMV MAC PROPH MEDS] taxonomy...
Creating [BKMV NNRTI MEDS] taxonomy...
Creating [BKMV NRTI MEDS] taxonomy...
Creating [BKMV PCP PROPH MEDS] taxonomy...
Creating [BKMV PI MEDS] taxonomy...
Updating [BGP CHLAMYDIA CPTS] taxonomy...
Updating [BGP CPT HIV TESTS] taxonomy...
Creating [BGP ALCOHOL HLTH FACTOR] taxonomy...
Creating [BGP ALCOHOL SCREEN ICDS] taxonomy...
Creating [BGP CD4 CPTS] taxonomy...

```

Creating [BGP CHLAMYDIA TEST PROCEDURES] taxonomy...
Creating [BGP DENTAL EXAM DENTAL CODE] taxonomy...
Creating [BGP EYE EXAM CPTS] taxonomy...
Creating [BGP FLU IZ CVX CODES] taxonomy...
Creating [BGP FLU IZ DXS] taxonomy...
Creating [BGP FLU IZ PROCEDURES] taxonomy...
Creating [BGP HIV VIRAL LOAD CPTS] taxonomy...
Creating [BGP HYSTERECTOMY PROCEDURES] taxonomy...
Creating [BGP IPV/DV COUNSELING ICDS] taxonomy...
Creating [BGP MAMMOGRAM ICDS] taxonomy...
Creating [BGP MAMMOGRAM PROCEDURES] taxonomy...
Creating [BGP NON TOBACCO USER FACTORS] taxonomy...
Creating [BGP PAP ICDS] taxonomy...
Creating [BGP PAP PROCEDURES] taxonomy...
Creating [BGP PNEUMO IZ CPTS] taxonomy...
Creating [BGP PNEUMO IZ DXS] taxonomy...
Creating [BGP PNEUMO IZ PROCEDURES] taxonomy...
Creating [BGP TOBACCO CESS DENTAL CODE] taxonomy...
Creating [BGP TOBACCO SCREEN HLTH FACTOR] taxonomy...
Creating [BGP TOBACCO USER HLTH FACTORS] taxonomy...
Creating [BKM BCG IZ CPTS] taxonomy...
Creating [BKM BCG IZ CVX CODES] taxonomy...
Creating [BKM BCG IZ PROCEDURE] taxonomy...
Creating [BKM CHLAMYDIA DXS] taxonomy...
Creating [BKM CMV IZ CPTS] taxonomy...
Creating [BKM CMV IZ CVX CODES] taxonomy...
Creating [BKM CMV LOINC CODES] taxonomy...
Creating [BKM CMV TEST CPTS] taxonomy...
Creating [BKM COCCI ANTIBODY CPTS] taxonomy...
Creating [BKM COCCI ANTIBODY LOINC CODES] taxonomy...
Creating [BKM DENTAL EXAMINATION] taxonomy...
Creating [BKM FAMILY PLANNING POV] taxonomy...
Creating [BKM FTA-ABS CPTS] taxonomy...
Creating [BKM FTA-ABS LOINC CODES] taxonomy...
Creating [BKM GONORRHEA DXS] taxonomy...
Creating [BKM GONORRHEA LOINC CODES] taxonomy...
Creating [BKM GONORRHEA TESTS CPTS] taxonomy...
Creating [BKM HEP A DXS] taxonomy...
Creating [BKM HEP A IZ CPTS] taxonomy...
Creating [BKM HEP A IZ CVX CODES] taxonomy...
Creating [BKM HEP B DXS] taxonomy...
Creating [BKM HEP B IZ CPTS] taxonomy...
Creating [BKM HEP B IZ CVX CODES] taxonomy...
Creating [BKM HEP B LOINC CODES] taxonomy...
Creating [BKM HEP B TESTS CPTS] taxonomy...
Creating [BKM HEP C DXS] taxonomy...
Creating [BKM HEP C EIA LOINC CODES] taxonomy...
Creating [BKM HEP C EIA TESTS CPTS] taxonomy...
Creating [BKM HEP PANEL LOINC CODES] taxonomy...
Creating [BKM HEP C RIBA LOINC CODES] taxonomy...
Creating [BKM HEP C RIBA TESTS CPTS] taxonomy...
Creating [BKM HEPATITIS PANEL CPTS] taxonomy...
Creating [BKM OTHER STD DXS] taxonomy...
Creating [BKM OTHER SUBSTANCE ABUSE DXS] taxonomy...
Creating [BKM PNEUMO IZ CVX CODES] taxonomy...
Creating [BKM PPD CPTS] taxonomy...
Creating [BKM PPD CVX CODES] taxonomy...
Creating [BKM PPD ICDS] taxonomy...
Creating [BKM PPD LOINC CODES] taxonomy...

```

```
Creating [BKM RPR CPTS] taxonomy...
Creating [BKM RPR LOINC CODES] taxonomy...
Creating [BKM SAFE SEX ED CODES] taxonomy...
Creating [BKM STATE REPORTABLE DXS] taxonomy...
Creating [BKM STD ED CODES] taxonomy...
Creating [BKM SYPHILIS DXS] taxonomy...
Creating [BKM TETANUS IZ CPTS] taxonomy...
Creating [BKM TETANUS IZ CVX CODES] taxonomy...
Creating [BKM TETANUS IZ DXS] taxonomy...
Creating [BKM TETANUS IZ PROCEDURES] taxonomy...
Creating [BKM TOXOPLASMOSIS CPTS] taxonomy...
Creating [BKM TOXOPLASMOSIS LOINC CODES] taxonomy...
Creating [BKM TRICHOMONIASIS DXS] taxonomy...
Creating [BKM TB MED NDCS] taxonomy...
Creating [BKMV AIDS DEF ILL DXS] taxonomy...
Creating [BKMV CD4 ABS CPTS] taxonomy...
Creating [BKMV CD4 ABS LOINC CODES] taxonomy...
Creating [BKMV FI MED NDCS] taxonomy...
Creating [BKMV HIV ED DXS] taxonomy...
Creating [BKMV HIV GENOTYPE CPTS] taxonomy...
Creating [BKMV HIV OPP INF DXS] taxonomy...
Creating [BKMV HIV PHENOTYPE CPTS] taxonomy...
Creating [BKMV MAC PROPH MED NDCS] taxonomy...
Creating [BKMV NNRTI MED NDCS] taxonomy...
Creating [BKMV NRTI MED NDCS] taxonomy...
Creating [BKMV PCP PROPH MED NDCS] taxonomy...
Creating [BKMV PI MED NDCS] taxonomy...
Creating [BKMV STD ED DXS] taxonomy...

The user [selected CFM] was added with an access role of CASE FILE MANAGER.

Updating Routine file...

The following Routines were created during this install:

Updating KIDS files...

HIV MANAGEMENT SYSTEM 1.0 Installed.
      Apr 25, 2006@14:55:50

Not a VA primary domain

Install Completed
```

Figure 5-1: Sample of installation screens

Exit the installation menu.

6.0 Installation Configuration

Two new globals will be created for this application during the installation process. Global ^BKM contains the patient register information. Under normal circumstances the growth of this global is expected to be relatively slow.

Global ^BKMV contains the audit data, which tracks user activity and changes to patient register data. This global will grow more rapidly than ^BKM, dependent on the number of patients assigned to the HIV Management Register, their activity at IHS healthcare sites, and the volume of user activity surrounding patient care and case management.

Globals for this application may be placed at the Site Manager's discretion.

Enable journaling for these two globals.

7.0 Post-Installation Instructions

The new HMS menu was automatically added to the XUCORE menu tree during installation. If a menu tree other than XUCORE is utilized at the site, the BKMVMENU option will need to be added to that menu manually using the System Manager Menu.

An HMS Key Delegate/Case File Manager was defined for the HIV Management System during the installation process and delegated the HMS security keys. This designated user can now add and set up any number of users in HMS using the User Maintenance [BKMV USER MAINTENANCE] menu option, (as shown below), including additional Case File Managers.

For additional information on user account setup, please refer to the HMS User Manual.

```
Setting up programmer environment
Terminal Type set to: C-VT320

Select OPTION NAME: XUCORE          Core Applications

    GIS    GIS Interface Menu ...
    HMS    HIV Management System ...
    QMAN   Q-Man (PCC Query Utility)
    QMGR   Q-Man Site Manager's Utilities

Select Core Applications Option: HMS  HIV Management System

                                     ***WARNING***

                                     ***RESTRICTED INFORMATION***
* * * * *
* All information in this system is protected by the Privacy Act of      *
* 1974. If you elect to proceed, you will be required to prove you have  *
* a need to know. Access of data in this system is tracked, and your    *
* station Security Officer may contact you for your justification.       *
* * * * *

Enter Yes or No: N// YES

    PM     Patient Management ...
    RM     Register Management ...
    MU     Manager Utilities ...

Select HIV Management System Option: MU  Manager Utilities

    UM     User Account Setup
```

```
TXCH  Taxonomy Check
SPS   Site Parameter Setup
RP    Register Population ...
EUP   Edit Unreviewed Patients
DOC   Generate Documents ...
AM    Register Audit Menu ...
TS    Taxonomy Setup
```

Select Manager Utilities Option: **UM** User Account Setup

Maintain HMS Users

Jul 18, 2005 15:48:15
[DEMO INDIAN HOSPITAL]

Page: 1 of 1

User Name: USER,DELEGATE

Access: CASE FILE MANAGER

Enter ?? for more actions

UM User Maintenance

Select Item(s): Quit// UM User Maintenance

(A)dd, (E)dit or (D)eleate Another User?: (A/E/D): **ADD**

Select User Name: USER

- 1 USER,CFM
- 2 USER,DELEGATE
- 3 USER,EDITOR
- 4 USER,IRM
- 5 USER,READER

CHOOSE 1-5: 1 USER,CFM

Add USER,CFM as an HMS user?: (Y/N): YES

Select one of the following:

- R READER
- E EDITOR
- M CASE FILE MANAGER

Occupation: **M** CASE FILE MANAGER

The user USER,CFM was added with an occupation of CASE FILE MANAGER.

(A)dd, (E)dit or (D)eleate Another User?: (A/E/D): **ADD**

Select User Name: USER

- 1 USER,CFM
- 2 USER,DELEGATE
- 3 USER,EDITOR
- 4 USER,IRM
- 5 USER,READER

CHOOSE 1-5: 3 USER,EDITOR

Add USER,EDITOR as an HMS user?: (Y/N): **YES**

Select one of the following:

- R READER


```

      E      EDITOR
      M      CASE FILE MANAGER

Occupation: EDITOR

The user USER,EDITOR was added with an occupation of EDITOR.

(A)dd, (E)dit or (D)delete Another User?: (A/E/D): ADD
Select User Name: USER
  1  USER,CFM
  2  USER,DELEGATE
  3  USER,EDITOR
  4  USER,IRM
  5  USER,READER
CHOOSE 1-5: 5  USER,READER
Add USER,READER as an HMS user?: (Y/N): YES

    Select one of the following:

      R      READER
      E      EDITOR
      M      CASE FILE MANAGER

Occupation: READER

The user USER,READER was added with an occupation of READER.

Maintain HMS Users          Jul 18, 2005 15:48:50          Page:    1 of    1
                           [DEMO INDIAN HOSPITAL]

      User Name: USER,DELEGATE          Access: CASE FILE MANAGER
      User Name: USER,CFM              Access: CASE FILE MANAGER
      User Name: USER,READER           Access: READER
      User Name: USER,EDITOR           Access: EDITOR

      Enter ?? for more actions

UM  User Maintenance

Select Item(s): Quit//

```

Figure 7-1: Example of the user maintenance option

The Site Manager must set the State HIV Reporting Required field to *No* if the local state does *not* require reporting of HIV cases (all states require AIDS case reporting, but may require HIV reporting as well). Contact your local Case File Manager for this information.

Consult with the Case File Manager to determine if the Autopopulate Register function should be a scheduled task, and schedule appropriately (TaskMan).

Add users to mail group BKMI Alert (new) for Register Alerts as designated by the Case File Manager. Alerts are sent by HMS when:

- A patient in the HMS register has a status change
- A patient in the HMS register has a register diagnosis change.

An automated post install routine will add taxonomies to your existing taxonomies. Taxonomies are in the Standard Data Dictionary, file #9002226.

Note: The default register (Human Immunodeficiency Virus) should *never* be deleted from the ICM base file (this cannot be done by a user).

8.0 Contact Information

If you have any questions or comments regarding this distribution, please contact the OIT Help Desk by:

Phone: (505) 248-4371 or
(888) 830-7280

Fax: (505) 248-4363

Web: <http://www.rpms.ihs.gov/TechSupp.asp>

Email: ITSCHelp@mail.ihs.gov

Appendix A – Routine Checksums

ROUTINE	CHECKSUM	
	OLD	NEW
BKM1POST	n/a	10005902
BKM1PRE	n/a	1848107
BKMALG	n/a	5931649
BKMAPG	n/a	4858957
BKMATX	n/a	1174218
BKMATXA	n/a	1303148
BKMATXB	n/a	1011658
BKMATXC	n/a	1030316
BKMATXD	n/a	1064397
BKMATXE	n/a	1023685
BKMATXF	n/a	1184115
BKMATXG	n/a	1077917
BKMATXH	n/a	1048206
BKMATXI	n/a	938840
BKMATXJ	n/a	1012445
BKMATXK	n/a	971570
BKMATXL	n/a	970317
BKMATXM	n/a	999617
BKMATXN	n/a	954499
BKMATXO	n/a	1442260
BKMATXP	n/a	1118768
BKMATXQ	n/a	924358
BKMATXR	n/a	989231
BKMATXS	n/a	1065613
BKMATXT	n/a	955828
BKMATXU	n/a	969867
BKMATXV	n/a	3559245
BKMATXW	n/a	1930654
BKMBTX	n/a	1086727
BKMBTXA	n/a	1071759
BKMBTXB	n/a	921335
BKMBTXC	n/a	932321
BKMBTXD	n/a	1280306
BKMBTXE	n/a	910435
BKMBTXF	n/a	920863
BKMBTXG	n/a	3084448
BKMBTXGB	n/a	4082038
BKMBTXGC	n/a	919806
BKMBTXH	n/a	1184294
BKMBTXI	n/a	947615
BKMBTXJ	n/a	2400253
BKMBTXK	n/a	945809
BKMBTXL	n/a	1039635
BKMBTXM	n/a	914976
BKMBTXN	n/a	3043833
BKMBTXO	n/a	917426
BKMBTXP	n/a	3534165
BKMBTXPB	n/a	631545
BKMBTXQ	n/a	1015553
BKMBTXR	n/a	966434
BKMBTXS	n/a	1074696
BKMBTXT	n/a	1190782
BKMBTXU	n/a	981148
BKMBTXV	n/a	1203502

BKMBTXW	n/a	1400338
BKMBTXX	n/a	3344683
BKMBTXXB	n/a	2427770
BKMBTXY	n/a	1059628
BKMCC1	n/a	2586373
BKMCTX	n/a	1089208
BKMCTXA	n/a	1185795
BKMCTXB	n/a	1034671
BKMCTXC	n/a	1014448
BKMCTXD	n/a	1042824
BKMCTXE	n/a	3423261
BKMCTXF	n/a	1091884
BKMCTXG	n/a	953044
BKMCTXH	n/a	1189843
BKMCTXI	n/a	1256096
BKMCTXJ	n/a	1044854
BKMCTXK	n/a	1060314
BKMCTXL	n/a	1062861
BKMCTXM	n/a	901292
BKMCTXN	n/a	936359
BKMCTXO	n/a	967756
BKMCTXP	n/a	1087259
BKMCTXQ	n/a	2850604
BKMCTXR	n/a	3263548
BKMCTXS	n/a	5999615
BKMCTXSB	n/a	3846999
BKMCTXT	n/a	914263
BKMCTXU	n/a	1316141
BKMCTXV	n/a	1349822
BKMCTXW	n/a	1186660
BKMCXXX	n/a	1026243
BKMDCMT	n/a	12607655
BKMDCMT1	n/a	6571519
BKMDCMT2	n/a	6573258
BKMDCMT3	n/a	6227844
BKMDOCE	n/a	536276
BKMDOCP	n/a	7503158
BKMDTX	n/a	803820
BKMDTXA	n/a	1012447
BKMDTXB	n/a	3243020
BKMDTXBB	n/a	3635439
BKMDTXC	n/a	943202
BKMDTXD	n/a	2760586
BKMDTXDB	n/a	4315572
BKMDTXDC	n/a	4794722
BKMDTXDD	n/a	4808532
BKMDTXDE	n/a	4790404
BKMDTXDF	n/a	1231192
BKMDTXE	n/a	3400968
BKMDTXEB	n/a	567570
BKMDTXF	n/a	917431
BKMDTXG	n/a	1053587
BKMDTXH	n/a	957459
BKMDTXI	n/a	913441
BKMDTXJ	n/a	940866
BKMDTXK	n/a	3315844
BKMDTXKB	n/a	1335851
BKMDTXL	n/a	946934
BKMDTXM	n/a	3572853
BKMDTXN	n/a	1700220

BKMDTX0	n/a	4063523
BKMDTX0B	n/a	816960
BKMDTXP	n/a	4168804
BKMDTXPB	n/a	823058
BKMDTXQ	n/a	3399307
BKMDTXR	n/a	913842
BKMED	n/a	450104
BKMIDTF	n/a	50821
BKMIHSM	n/a	1843458
BKMILK	n/a	4888059
BKMIMRP	n/a	24181430
BKMIMRP1	n/a	15539786
BKMIMRP2	n/a	15358774
BKMIMRP3	n/a	15343856
BKMIPDR	n/a	23465391
BKMIPRIV	n/a	1609943
BKMIRML	n/a	19897321
BKMIRML1	n/a	4278960
BKMIXX	n/a	8240270
BKMIXX1	n/a	8203930
BKMIXX2	n/a	9540935
BKMIXX3	n/a	10077313
BKMIXX4	n/a	9118523
BKMIXX5	n/a	6802845
BKMPHPR	n/a	17424558
BKMPLK	n/a	8939663
BKMPLK1	n/a	13764421
BKMPLK1I	n/a	3125484
BKMPLKB	n/a	3767036
BKMPLKH	n/a	2243795
BKMPLKI	n/a	6770701
BKMPLKID	n/a	6655732
BKMPLKP	n/a	1170199
BKMPLKUT	n/a	5636823
BKMPLKZ	n/a	185807
BKMSUPP	n/a	21193885
BKMSUPP1	n/a	18311728
BKMSUPP2	n/a	23154794
BKMSUPP3	n/a	13691751
BKMSUPP4	n/a	13548649
BKMSUPP5	n/a	13028407
BKMTX	n/a	1717937
BKMTXA	n/a	836264
BKMTXB	n/a	1029764
BKMTXC	n/a	866123
BKMTXD	n/a	1047922
BKMTXE	n/a	1176952
BKMTXF	n/a	867299
BKMTXG	n/a	1037406
BKMTXH	n/a	1129359
BKMTXI	n/a	1128323
BKMTXJ	n/a	2383002
BKMTXK	n/a	3420515
BKMTXKB	n/a	2040878
BKMTXL	n/a	1003385
BKMTXM	n/a	1132592
BKMTXN	n/a	1050939
BKMTX0	n/a	1042326
BKMTXP	n/a	1115654
BKMTXQ	n/a	1149421

BKMTXR	n/a	2915889
BKMTXRB	n/a	4106276
BKMTXRC	n/a	2392942
BKMTXS	n/a	1055035
BKMTXT	n/a	1483993
BKMTXU	n/a	1175507
BKMTXV	n/a	1373734
BKMTXW	n/a	2105326
BKMTXX	n/a	1585327
BKMTXY	n/a	1332906
BKMTXZ	n/a	1195525
BKMVA1	n/a	21677704
BKMVA1A	n/a	18614934
BKMVA1B	n/a	9205100
BKMVA1C	n/a	11946974
BKMVA1U	n/a	6063127
BKMVA2	n/a	26419597
BKMVA4	n/a	7456110
BKMVA41	n/a	5497418
BKMVA5	n/a	7314080
BKMVA51	n/a	3780908
BKMVA6	n/a	7787184
BKMVA61	n/a	2522994
BKMVA9	n/a	22250228
BKMVAUC	n/a	10907093
BKMVAUD	n/a	13963522
BKMVAUD1	n/a	5506849
BKMVAUD2	n/a	5314382
BKMVAUDN	n/a	1048189
BKMVAUP	n/a	4321428
BKMVAUP1	n/a	4158
BKMVAUP0	n/a	787256
BKMVAUTR	n/a	2430471
BKMVB0	n/a	12905121
BKMVB1	n/a	9943517
BKMVB2	n/a	5583174
BKMVC0	n/a	5382428
BKMVC1	n/a	3110911
BKMVC10	n/a	5954277
BKMVC11	n/a	5293433
BKMVC12	n/a	4315037
BKMVC2	n/a	28629229
BKMVC3	n/a	13376091
BKMVC4	n/a	2162255
BKMVC5	n/a	2526905
BKMVC6	n/a	7460079
BKMVC7	n/a	6074384
BKMVC8	n/a	2285548
BKMVC9	n/a	5146451
BKMVCD	n/a	22346328
BKMVD1	n/a	14133690
BKMVD11	n/a	10433070
BKMVD2	n/a	7016028
BKMVD21	n/a	3947965
BKMVD3	n/a	6634001
BKMVD31	n/a	2334497
BKMVD4	n/a	6896466
BKMVD41	n/a	3991524
BKMVD5	n/a	7023256
BKMVD51	n/a	3993793

BKMVD6	n/a	7337138
BKMVD61	n/a	4545377
BKMVD7	n/a	9090549
BKMVD71	n/a	4125216
BKMVD8	n/a	7349841
BKMVD81	n/a	4343544
BKMVDIC	n/a	2427873
BKMVDOD	n/a	9373709
BKMVDOD1	n/a	851240
BKMVDOD4	n/a	6616223
BKMVDOD8	n/a	10266493
BKMVF1	n/a	9163660
BKMVF3	n/a	14540517
BKMVF31	n/a	16124228
BKMVF32	n/a	8189174
BKMVF4	n/a	5189281
BKMVF5	n/a	1029133
BKMVFAP1	n/a	19940168
BKMVFAP1	n/a	7137330
BKMVFLD	n/a	2962916
BKMVFMP1	n/a	5312525
BKMVIST3	n/a	10856340
BKMVIST4	n/a	12246054
BKMVIST5	n/a	13420243
BKMVLC	n/a	14897631
BKMVLN	n/a	11530759
BKMVLP	n/a	8488880
BKMVNAR	n/a	1423711
BKMVND	n/a	19617356
BKMVPCM1	n/a	12332740
BKMVPCML	n/a	3460045
BKMVPG	n/a	6019274
BKMVPGP	n/a	7673671
BKMVPGP1	n/a	2339483
BKMVQCR	n/a	13121342
BKMVQCR1	n/a	5139792
BKMVQCR2	n/a	9205199
BKMVQCR3	n/a	7917168
BKMVQCR4	n/a	7745635
BKMVQCR5	n/a	5523688
BKMVQCR6	n/a	9008921
BKMVQCR7	n/a	13670317
BKMVQCR8	n/a	3734371
BKMVQCR9	n/a	7876259
BKMVQCRA	n/a	12089886
BKMVSCD	n/a	1233791
BKMVSRP	n/a	4746748
BKMVSRP1	n/a	6373098
BKMVSSR	n/a	7153773
BKMVSSR1	n/a	1815930
BKMVSSR2	n/a	1302723
BKMVSSR3	n/a	2784580
BKMVSSR4	n/a	5251490
BKMVSSR5	n/a	4417258
BKMVSSR6	n/a	9985099
BKMVTAX	n/a	13676954
BKMVTAX0	n/a	18314631
BKMVTAX1	n/a	9207407
BKMVTAX2	n/a	1612437
BKMVTAX4	n/a	5991306

BKMOVUTL	n/a	560229
----------	-----	--------