

RPMS Omnicell® & Pyxis® Configuration

```
Select Omnicell/ILC User Menu Option: ??  
TRB   Omnicell Interface Trouble shooting ... [VEFSO O TROUBLESHOOT MENU]  
DIS   Display Queue for a Patient [VEFSO O PATIENT DISPLAY]  
MSD   set quantity multiple for missed doses [VEFSO O MISSED DOSE QTY MULT]  
PICK  ILC Omnicell Pick List Re-Print [VEFSO O PICK LIST REPRINT]  
PUSH  Push a Single Drug [VEFSO OMNICELL PUSH A DRUG]  
SEND  Send one patient adt/orders to interface [VEFSO O TRANSMIT ONE PATIENT]  
SH    Show Ready Queue [VEFSO OMNICELL SHOW]  
SITE  Site Parameter Edit [VEFSO OMNICELL SITE]  
SUR   Send Surgery Data [VEFSO O SEND SURGERY]
```

Notes:

- The bag labels (v label) will no longer work after a change from Observation to Full Admit. In RPMS, type a caret followed by **send** (^**send**) to send the orders that are dropped from Omnicell back to Omnicell.
- Active medications are restored to the patient on the Full admission. However, though the IVPB or IV infusion orders appear to be active, the Pharmacy generated IV labels will not work for BCMA. Pharmacy must generate new labels and re-label the bags.