



Abbott

Important Product Notice

The SPECIMEN COLLECTION AND HANDLING for the ID NOW™ COVID-19 test has changed. Swabs eluted in VTM are no longer an appropriate sample type. Please refer to the updated product insert included in this kit.

ID NOW COVID-19 is intended for testing a swab directly without elution in viral transport media as dilution will result in decreased detection of low positive samples that are near the limit of detection of the test.

Due to this change, disposable transfer pipettes are no longer included in the kit.

© 2020 Abbott. All rights reserved.

All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

TB000041 Rev. 1 2020/04



Abbott

Important Product Notice

The SPECIMEN COLLECTION AND HANDLING for the ID NOW™ COVID-19 test has changed. Swabs eluted in VTM are no longer an appropriate sample type. Please refer to the updated product insert included in this kit.

ID NOW COVID-19 is intended for testing a swab directly without elution in viral transport media as dilution will result in decreased detection of low positive samples that are near the limit of detection of the test.

Due to this change, disposable transfer pipettes are no longer included in the kit.

© 2020 Abbott. All rights reserved.

All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

TB000041 Rev. 1 2020/04

Abbott
ID NOW
Tech Bulletin

Size:
8.5 in x 11 in



Black

PN: TB000041
Rev: 1

Revision Date:
1.2 2020/04/15